


[HPLC COLUMNS CAPABILITIES]



NEW INNOVATIONS ARE ELEVATING HPLC CAPABILITIES.
IS YOUR LABORATORY READY?

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



CREATING EXCEPTIONAL CHROMATOGRAPHY

Waters reputation is based on chromatography, but we do not create chromatography — You do. Innovative thinking within your laboratory creates the chromatographic methods and the assays that sustain your business. The metric of your success is monitored by the methods and the results that you produce and the HPLC column that you choose today needs to support your success for the future. Waters full line of state-of-the-art reversed-phase and HILIC HPLC columns are chosen by scientists who understand that performance and innovation are linked and their success depends on them.



XBRIDGE HPLC COLUMNS

XBridge™ HPLC columns are designed for one purpose – to maximize your productivity. Whether your goal is to create a quality control method or to develop a leading edge LC/MS assay, XBridge columns are designed to help you by:

Improving pH Stability

- increased column lifetime

Improving Column Reliability

- assay ruggedness

Maximizing Particle Efficiency


- unmatched peak shape and peak capacity

With a selection of 10 general purpose and application specific sorbents in the widest range of particle sizes available, no other HPLC column family gives you the tools you need for the most demanding chromatographic challenges. Whether you require robust HPLC methods, seamless UPLC® transferability, or preparative scaling for product isolation, count on the versatility of an XBridge column.

pH Stability


XBridge columns have been specifically designed to contain the most chemically stable chromatographic sorbent available allowing you to explore the full benefits of a wide pH (1-12) mobile-phase range. Chemical stability, especially for the extremes of pH, is built into the particle during the synthesis process and it cannot be duplicated using a conventional silica-based bonding process. No other column can match the chemical stability of an XBridge column.

Accelerated High pH Stability Test of Competitive Columns



Chromatograms, run at regular intervals during the high-pH lifetime study, verify that 86% of the original XBridge column efficiency remains after 300 hours at pH 10 and elevated temperature, with little change in peak shape or retention time.

Comparative separations may not be representative of all applications.



Ligand Type	Trifunctional C ₁₈	Trifunctional C ₈
Particle Size	2.5, 3.5, 5, 10 μm	2.5, 3.5, 5, 10 μm
Ligand Density*	3.1 μmol/m ²	3.2 μmol/m ²
Carbon Load*	18%	13%
Endcap Style	Proprietary	Proprietary
pH Range	1-12	1-12
Low pH Temp. Limit	80 °C	60 °C
High pH Temp. Limit	60 °C	60 °C
Pore Diameter*	130Å	130Å
Surface Area*	185 m ² /g	185 m ² /g

* Expected or approximate value.


Column Reliability

Much of the cost when developing a chromatographic method is associated with the rigorous testing and validation of the final method. We understand that revalidation of your method is not an option, so we thoroughly test each batch of sorbent and final column product to ensure that you get the most reproducible columns available. With an XBridge column, you have the confidence that the method you develop today will be repeatable for the lifetime of your assay.

Particle Efficiency

The BEH (Ethylene-Bridged Hybrid) particle offers many advantages over conventional silica-based particles, including the ability to control the silanol activity with great precision. By controlling the silanol activity, you control and reduce unwanted silanol interactions that increase peak tailing.

XBridge Family USP Tailing Factors




The combination of excellent particle and ligand stability as well as high chromatographic efficiencies makes XBridge columns an ideal choice for low and intermediate pH methods.

Comparative separations may not be representative of all applications.

Methods Transfer Using XP 2.5 µm Columns

All XBridge and XSelect HPLC columns are offered in eXtended Performance (XP) 2.5 µm column formats to help you transfer methods from HPLC to UPLC instrumentation. The XP 2.5 µm columns improve the performance of your current HPLC instrumentation and provide you with a pathway to gain maximum separation efficiency using sub-2-µm ACQUITY UPLC[®] technology.



Columns of different lengths and particle sizes were used to successfully reduce run times and maintain resolution.

LC Conditions

Columns: XBridge C₁₈ 2.1 x 150 mm, 5.0 µm
 XBridge C₁₈ 2.1 x 100 mm, 3.5 µm
 XBridge C₁₈ 2.1 x 75 mm, 2.5 µm
 ACQUITY UPLC BEH C₁₈ 2.1 x 50 mm, 1.7 µm

Sample Diluent: 3% acetonitrile in water with 0.1% formic acid

Sample Conc.: 25 µg/mL

Column Temp.: 38 °C

Detection: UV @ 280 nm

Sampling Rate: 40 pts/sec

Time Constant: 0.05

Instrument: ACQUITY UPLC with TUV detector


Mobile Phase A: 0.1% formic acid in water

Mobile Phase B: 0.1% formic acid in acetonitrile

Flow Rate: see chromatogram above

Injection Volume: see chromatogram above

Isocratic: 95% A:5% B



Monofunctional Embedded Polar Group (Shield RP18)	Trifunctional C ₆ Phenyl	HILIC	Amide	BEH130 C ₁₈ (Peptide Separation Technology)	BEH300 C ₁₈ (Peptide Separation Technology)	BEH300 C ₄ (Protein Separation Technology)	BEH C ₁₈ (Oligonucleotide Separation Technology)
2.5, 3.5, 5, 10 µm	2.5, 3.5, 5 µm	2.5, 3.5, 5 µm	2.5, 3.5 µm	3.5, 5, 10 µm	3.5, 5, 10 µm	3.5 µm	2.5 µm
3.3 µmol/m ²	3.0 µmol/m ²	N/A	7.5 µmol/m ²	3.1 µmol/m ²	3.1 µmol/m ²	2.4 µmol/m ²	3.1 µmol/m ²
17%	15%	unbonded	12%	18%	12%	8%	18%
TMS	Proprietary	N/A	None	Proprietary	Proprietary	None	Proprietary
2-11	1-12	1-9	2-11	1-12	1-12	1-10	1-12
50 °C	80 °C	45 °C	90 °C	80 °C	80 °C	80 °C	80 °C
45 °C	60 °C	45 °C	90 °C	60 °C	60 °C	50 °C	60 °C
130 Å	130 Å	130 Å	130 Å	130 Å	300 Å	300 Å	130 Å
185 m ² /g	185 m ² /g	185 m ² /g	185 m ² /g	185 m ² /g	90 m ² /g	90 m ² /g	185 m ² /g

XSELECT HPLC COLUMNS

XSelect™ HPLC columns are designed for the method development scientist. Scientists who develop chromatographic methods demand the most diverse selection of sorbents to easily separate the most difficult analyte co-elutions. XSelect columns are tools that are:

Designed for Selectivity

- improve your ability to separate closely eluting peaks

Intended for Isolation & Purification

- highest analyte mass loading available


Ideal for Rapid Method Development

- reduce the time and cost spent developing methods

The base particle or substrate critically influences analyte selectivity; the bonded ligand has a secondary influence. Each, when used alone does not create dramatic selectivity changes; however, in combination, the substrate and ligand create the ultimate tools for enhancing analyte selectivity. It is for this reason that the XSelect family of columns contains both High Strength Silica (HSS) and Charged Surface Hybrid (CSH™) technologies with a unique optimization of bonded ligands. The result is reproducibility and high selectivity.

Enhanced Selectivity

Selectivity and retentivity are the most powerful tools a method developer has to influence chromatographic behavior. The XSelect family offers a diverse range of reversed-phase C₁₈ columns (e.g., CSH C₁₈, HSS C₁₈, HSS C₁₈ SB) for general purpose separations; as well as, columns that offer improved polar retention (T3) and greater selectivity options (Phenyl-Hexyl, Fluoro-Phenyl and Cyano) for method development.



Observed selectivity differences for a mixture of basic analytes. Compounds: [1] aminopyrazine, [2] pindolol, [3] quinine, [4] labetalol, [5] verapamil, [6] diltiazem, [7] amitriptyline.


LC Conditions

Columns: 2.1 x 50 mm
 Mobile Phase A: 10 mM ammonium formate, pH 3.0
 Mobile Phase B: methanol
 Flow Rate: 0.4 mL/min
 Injection Volume: 1 µL
 Sample Diluent: water
 Column Temp.: 30 °C

Gradient:

Time (min)	%A	%B
0.00	70	30
3.00	15	85
3.50	15	85
3.51	70	30
4.50	70	30

Detection: UV @ 260 nm
 Sampling Rate: 20 pts/sec
 Filter Response: normal
 System: ACQUITY UPLC with ACQUITY UPLC PDA detector




Ligand and Particle Type	CSH C ₁₈	CSH Fluoro-Phenyl	CSH Phenyl-Hexyl
Particle Size	2.5, 3.5, 5 µm	2.5, 3.5, 5 µm	2.5, 3.5, 5 µm
Ligand Density*	2.3 µmol/m ²	2.3 µmol/m ²	2.3 µmol/m ²
Carbon Load*	15%	10%	14%
Endcap Style	Proprietary	None	Proprietary
pH Range	1-11	1-8	1-11
Low pH Temp. Limit	80 °C	60 °C	80 °C
High pH Temp. Limit	45 °C	45 °C	45 °C
Pore Diameter*	130Å	130Å	130Å
Surface Area*	185 m ² /g	185 m ² /g	185 m ² /g

* Expected or approximate value.

Isolation and Purification

High mass loading applications like compound purification and dissolution testing demand superior column performance. For these types of applications, column loading is limited by its inability to maintain symmetrical peak shape. This manifests itself as severe tailing of the main compound peak, which often overwhelms the trace impurities that you are trying to remove in the purification. XSelect CSH columns consistently provide narrow peaks under high loading conditions that allow the chromatographer the ability to separate trace-level impurities or degradants giving you more loading capacity with less time and solvent.



The improved mass loading of XSelect columns permits the separation, identification and quantification of closely eluting impurities or degradants.


Comparative separations may not be representative of all applications.

LC Conditions

Columns:	2.1 x 50 mm	Injection Volume:	5 μ L
Mobile Phase A:	water	Sample Diluent:	water
Mobile Phase B:	acetonitrile	Sample Conc.:	imipramine: 0.5 mg/mL; amitriptyline: as indicated (% of imipramine)
Mobile Phase C:	2% formic acid in water		
Gradient:		Column Temp.:	40 $^{\circ}$ C
	Time (min)	Flow (mL/min)	%A
	Initial	0.6	70
	2.0	0.6	60
	3.0	0.6	0
	3.1	0.6	70
	5.0	0.6	70
		%B	25
		%C	5
		Curve	Initial
			6
			6
			6
			6
		Detection:	UV @ 254 nm
		Wash Solvent:	50/50 acetonitrile/water
		Purge Solvent:	50/50 acetonitrile/water
		System:	ACQUITY UPLC H-Class with ACQUITY UPLC PDA detector

Method Development and Transfer


When developing methods, skilled chromatographers realize that any method developed using uniquely selective columns must be easily transferable across laboratories, independent of the LC system platform used. XSelect columns are engineered for method development and are fully compatible with all modern detection modes.



Reproducibility and scalability for gradient separations on 2.1 x 50 mm columns containing nine different batches of CSH Fluoro-Phenyl representing three (1.7, 3.5, and 5 μ m) particle sizes. Compounds: [1] thiourea, [2] resorcinol, [3] metoprolol, [4] 3-nitrophenol, [5] 2-chlorobenzoic acid, [6] amitriptyline, [7] diethylphthalate, [8] fenopropfen, [9] dipropylphthalate, [10] pyrenesulfonic acid.

LC Conditions

Columns:	2.1 x 50 mm	Column Temp.:	30 $^{\circ}$ C
Flow Rate:	0.5 mL/min	Detection:	UV @ 254 nm
Mobile Phase A:	15.4 mM ammonium formate, pH 3.0	Sampling Rate:	20 pts/sec
Mobile Phase B:	acetonitrile	Filter Response:	normal
Gradient:	5 to 90% B linear in 5 minutes	System:	ACQUITY UPLC with ACQUITY UPLC PDA detector
Injection Volume:	5 μ L		
Sample Diluent:	water		



HSS C ₁₈	HSS C ₁₈ SB	HSS T3	HSS PFP	HSS CN
2.5, 3.5, 5 μ m				
3.2 μ mol/m ²	1.6 μ mol/m ²	1.6 μ mol/m ²	3.2 μ mol/m ²	2.0 μ mol/m ²
15%	8%	11%	7%	5%
Proprietary	None	Proprietary	None	None
1-8	2-8	2-8	2-8	2-8
45 $^{\circ}$ C				
45 $^{\circ}$ C				
100 \AA				
230 m ² /g				

ATLANTIS HPLC COLUMNS

Atlantis® HPLC columns provide exceptional performance, versatility and retention for polar compounds, while also affording balanced retention for broad analyte mixtures. Built for polar compound retention, Atlantis HPLC columns provide you with:

Polar Compound Retention without Ion-Pairing Reagents

Eliminating ion-pairing reagents improves detection limits, method reproducibility and robustness, while reducing instrument maintenance due to harsh mobile-phase environments.

Compatibility with 100% Aqueous Mobile Phases


To maximize polar compound retention in reversed-phase methods, it is possible to use Atlantis reversed-phase HPLC columns with highly aqueous mobile phases and buffers without the risk of pore dewetting and hydrophobic collapse of the stationary phase.

Long Column Lifetimes Using Low-pH Mobile Phases


Atlantis columns resist ligand hydrolysis when using strongly acidic mobile phases, thus maintaining method efficiency, compound retention, and critical analyte selectivity.



20 hr Exposure to 0.5% TFA at 60 °C



Polar Compound Retention



Separating highly polar analytes on the Atlantis T3 column compared to competitive brands. Scientists rely on the uncompromised peak shape and retention that only Atlantis columns can provide. Compounds: [1] thiourea, [2] 5-fluorocystine, [3] adenine, [4] guanosine-5'-monophosphate, [5] thymine.

Comparative separations may not be representative of all applications.

LC Conditions

Column: 4.6 x 150 mm
 Mobile Phase: 10 mM ammonium formate, pH 3.0
 Flow Rate: 1.3 mL/min for 3 µm
 Injection Volume: 2.0 µL
 Column Temp.: 30 °C

Detection: 254 nm
 Sampling Rate: 10 pts/sec
 System: Alliance® 2695 with 2487 Dual-Wavelength Absorbance detector

Stationary Phase Characteristics

Intended Use	HPLC
Particle Type	Atlantis Silica
Available Chemistries	T3, dC ₁₈ , HILIC Silica
pH Range	T3: 2-8; dC ₁₈ : 3-7; HILIC Silica: 1-5
Maximum Rated Pressure	6000 psi [-400 bar]
Particle Size	3, 5, 10 µm
Pore Diameter / Volume	100Å/1.0 mL/g
Surface Area	330 m ² /g

SUNFIRE HPLC COLUMNS

SunFire™ columns set the standard for state-of-the-art bonded C₁₈- and C₈- silica HPLC columns. Benefiting from years of research and product development, SunFire columns represent the best in particle and bonding expertise and deliver industry-leading levels of chromatographic performance. Since their development in 2004, SunFire HPLC columns have gained the reputation for:

Excellent Low-pH Stability

Under low-pH mobile-phase conditions, SunFire columns exhibit superior column lifetimes that exceed many silica-based HPLC column brands.

High Efficiency


A synergistic combination of particle synthesis, packing technology, and hardware engineering is required for high efficiency. SunFire Intelligent Speed™ (IS™) and Optimum Bed Density (OBD™) columns were developed specifically from this knowledge.

Superior Peak Shape

For years, SunFire columns provide symmetrical peaks for improved resolution of acidic, neutral and basic compounds at low and moderate pH ranges (2-8).

Batch-to-Batch Reproducibility

Waters is dedicated to maintain the tightest specifications in the HPLC industry. Controlled manufacturing processes and column packing procedures ensure that you receive the best, most reproducible HPLC column available.



This excellent reproducibility is a result of our commitment to maintaining the tightest specifications in the HPLC column industry. SunFire columns start with high purity raw materials, and are produced using controlled manufacturing processes and column packing procedures that provide today's scientists with the best, most reproducible HPLC columns available.

Stationary Phase Characteristics


Intended Use	HPLC
Particle Type	SunFire Silica
Available Chemistries	C ₈ , C ₁₈ , Silica
pH Range	2-8
Maximum Rated Pressure	6000 psi [~400 bar]
Particle Size	2.5, 3.5, 5, 10 μm
Pore Diameter / Volume	100Å/1.0 mL/g
Surface Area	340 m ² /g

HILIC HPLC COLUMNS


Hydrophilic Interaction Chromatography (HILIC) is a versatile and effective alternative to reversed-phase retention for highly polar compounds. HILIC outperforms reversed-phase techniques when you are primarily interested in retaining small organic acids and bases as well as neutral compounds, like water-soluble vitamins and sugars.

HILIC for Orthogonal Selectivity

HILIC techniques are gaining popularity as part of HPLC method development strategies. As a complementary technique to traditional reversed-phase chromatography, HILIC separations not only retain highly polar compounds but provide, in some cases, large selectivity differences where there can be a complete reversal in elution order compared to a reversed-phase chromatogram.



Chemical structures of water soluble vitamins.



Column selectivity at low pH. Compounds: [1] nicotinamide, [2] pyridoxine, [3] riboflavin, [4] nicotinic acid, [5] thiamine, [6] ascorbic acid, [7] B12, [8] folic acid.


LC Conditions

Columns: 2.1 x 50 mm
 Flow Rate: 0.5 mL/min
 Mobile Phase A: 50/50 acetonitrile/water with 10 mM ammonium formate
 Mobile Phase B: 95/5 acetonitrile/water with 10 mM ammonium formate
 Gradient: 99.9% B to 0.1% B linear in 5 minutes

Injection Volume: 5 µL
 Sample Diluent: 75/25 acetonitrile/methanol
 Column Temp.: 30 °C
 Sampling Rate: 20 pts/sec
 System: ACQUITY UPLC with ACQUITY UPLC PDA detector

HILIC for Sample Throughput

Many common sample preparation techniques such as protein precipitation and elution from solid-phase extraction devices include the use of organic solvents such as acetonitrile. For reversed-phase LC, injection of a strong solvent from this type of extract may produce a chromatogram with peak distortion and smearing. This problem is eliminated when using a HILIC column. This not only saves time, but also improves variability in sampling because tedious sample evaporation and reconstitution steps can be eliminated.



XBridge HILIC provides direct injection of SPE eluents, thus increasing sample throughput.

Ligand Type	Atlantis HILIC Silica	XBridge HILIC	XBridge Amide
Particle Size	3, 5 μ m	2.5, 3.5, 5 μ m	2.5, 3.5 μ m
pH Range	1-5	1-9	2-11
Low pH Temp. Limit	45 °C	45 °C	90 °C
High pH Temp. Limit	45 °C	45 °C	90 °C
Pore Diameter*	100Å	130Å	130Å
Surface Area*	330 m ² /g	185 m ² /g	185 m ² /g

* Expected or approximate value.

TECHNICAL RESOURCES

Literature References

HPLC Troubleshooting Guide

Part Number: WA20769

HPLC Columns Theory, Technology, and Practice

Part Number: WAT038216

Comprehensive Guide to HILIC

Part Number: 715002531

Beginners Guide to Liquid Chromatography

Part Number: 715001531



Electronic Tools

Waters Reversed-Phase Column Selectivity Chart

www.waters.com/selectivitychart

Waters Column Advisor

www.waters.com/columnadvisor

Waters Part Selector & Selectivity Chart for iPad

www.waters.com/apps



SALES OFFICES:

**Austria and European Export
(Central South Eastern Europe, CIS
and Middle East)** 43 1 877 18 07

Australia 61 2 9933 1777

Belgium and Luxembourg 32 2 726 1000

Brazil 55 11 4134 3788

Canada 1 800 252 4752

China 86 21 6156 2666

Czech Republic 420 2 617 11384

Denmark 45 46 59 8080

Finland 358 9 5659 6288

France 33 1 30 48 72 00

Germany 49 6196 400 600

Hong Kong 852 2964 1800

Hungary 36 1 350 5086

India 91 80 2837 1900

Ireland 353 1 448 1500

Italy 39 02 265 0983

Japan 81 3 3471 7191

Korea 82 2 6300 4800

Mexico 52 55 52 00 1860

The Netherlands 31 76 508 7200

Norway 47 6 384 6050

Poland 48 22 101 5900

Puerto Rico 1 787 747 8445

Russia/CIS 7 495 727 4490/ 290 9737

Singapore 65 6593 7100

Spain 34 93 600 9300

Sweden 46 8 555 115 00

Switzerland 41 56 676 7000

Taiwan 886 2 2501 9928

United Kingdom 44 208 238 6100

All other countries:

Waters Corporation U.S.A.

1 508 478 2000

1 800 252 4752

www.waters.com

www.waters.com/hplccolumns

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



The quality management system of Waters' manufacturing facilities in Taunton, Massachusetts and Wexford, Ireland complies with the International Standard ISO 9001 Quality Management and Quality Assurance Standards. Waters' quality management system is periodically audited by the registering body to ensure compliance.

Waters, XBridge, ACQUITY UPLC, UPLC, XSelect, CSH, Atlantis, SunFire, Intelligent Speed, IS, OBD, Alliance, and The Science of What's Possible are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2012 Waters Corporation. Printed in the U.S.A.
October 2012 720003995EN VW-PDF