

## Silices en vrac

*im*Chem

164 Avenue Joseph Kessel - Bâtiment 7 - 78960 Voisins le Bretonneux - France

Tel. : +33 9 72 32 10 17 06 - Fax : +33 1 83 64 60 45 - [info@imchem.fr](mailto:info@imchem.fr) - [www.imchem.fr](http://www.imchem.fr)

# PHASES EN VRAC


## Choisir imChem

Une gamme complète de phases en vrac, de haute qualité, pour tous types de séparations chromatographiques.

<b>SILICE SURF - HPLC ANALYTIQUE</b>	<b>4</b>
<b>GAMME DE GREFFAGES</b>	<b>6</b>
<b>GAMME DE END-CAPPING</b>	<b>10</b>
<b>GAMME DE PARTICULES</b>	<b>13</b>
<b>SILICE CORAL - HPLC PREPARATIVE</b>	<b>14</b>
<b>POUR COMMANDER</b>	<b>16</b>

QUALITE

PERFORMANCE



# Silice en vrac pour HPLC

## La silice de la gamme Surf

Nous vous proposons des particules de silice naturelle, de haute pureté (99.99%), sphériques et très résistantes présentant d'excellentes caractéristiques de rétention grâce à une surface spécifique de 310 m<sup>2</sup>/g et un volume de pores optimaux 0,8 ml/g.

Disponibles de 1,5 à 10µm de diamètre, avec des porosités de 100 A à 300 A.

Disponibles vierges ou greffées.

Nous proposons les greffages suivants : C4, C8, C18, SO<sub>3</sub>H, DNH, PEI, ARG, EP2TNE, EP4TNE, Amino, CN, Phenyl, Phenyl-hexyl, Diol, ainsi que différents traitements de surface : Polar, Polar TriF, Wax, Extrême, AQ, AQ plus, TriF.

Notre silice ultra pure présente un taux de métaux inférieur à 10ppm, ce qui vous garantit une très bonne symétrie de pics sans interactions secondaires.

Métal	Na	Ti	Ca	Mg	Al	Fe	Zn
Taux en ppm	<10	<10	<10	<10	<10	<10	<10

Différents diamètres de particules proposées :

Diamètres des particules (en µm)	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	7.0	9.0	10.0





Diamètres des pores des particules proposées :

Diamètres des pores (en A)	100	150	200	300





# Silice en vrac pour HPLC

## Surf SILICA




microscope électronique



## Notre gamme de greffages

### Surf C18

Il s'agit de la phase la plus traditionnelle et la plus polyvalente en chromatographie, cette silice peu polaire est hydrophobe. Surf C18 permet une large gamme d'applications de séparations des composés acides, basiques et neutres en une seule analyse. Surf C18 a un bon mécanisme de rétention et peut être utilisée avec la plupart des éluants de phase inverse. Notre phase C18 existe avec ou sans end-capping



### Surf ARG

Notre silice Surf ARG a pour greffon l'acide aminé arginine. Notre silice ARG a une forte affinité pour les composés hydrophiles et possède une haute performance de séparation ainsi qu'une sélectivité différente par rapport aux autres silices, elle permet de séparer les composés hydrophiles tels que les acides aminés, les peptides, les vitamines et les acides nucléiques.

Différents diamètres de particules de silice Surf ARG sont disponibles pour l'analyse et la purification à grande échelle.

### Surf C8

Comme la Surf C18, cette silice Surf C8 est peu polaire et est hydrophobe, elle permet la séparation des composés acides, basiques et neutres en une seule analyse. Elle existe avec ou sans end-capping.



### Surf Diol

La silice Surf Diol est utilisable en phase normale ou en phase inverse. Elle permet la filtration de protéines et de peptides.




Surf C8 vous offre un temps de rétention plus court que Surf C18. La réduction de l'encombrement stérique des C8 comparée aux C18 vous apportera une sélectivité intéressante en fonction du poids de vos molécules.

En phase normale, sa sélectivité est intermédiaire entre celle de la silice nue et celle de la C18. Son temps de rétention plus long permet d'augmenter la performance de séparation et ainsi de séparer des isomères structuraux, ainsi que les oméga 3.


### Surf C4

Nos silices Surf C4 vous apportent une sélectivité complémentaire à celle des silices Surf C18 et C8 avec une rétention moindre permettant un temps d'analyse plus court.



### Surf NH2

Nos silices Surf NH<sub>2</sub> sont traitées en surface avec des groupements amine pour éviter les interactions ioniques d'adsorption entre le gel de silice et les composés basiques.



Cette différence de rétention est notable pour les molécules les plus hydrophobes sans pour autant affecter les molécules polaires. Avec leurs chaînes plus courtes et leur caractère plus hydrophile, les silices C4 se prêtent parfaitement à la rétention et à la séparation de nombreuses classes de molécules issues du domaine biologique, telles que des protéines, les bio-polymères et polypeptides.

Ce greffage provoque l'élution des composés basiques de la colonne, et est très appréciable pour la séparation des cyclodextrines et des oligo saccharides.


Nos Surf C4, C8 et C18 sont disponibles avec différentes tailles de pores, 100, 150, 200 et 300A pour permettre des séparations optimales d'une large gamme de masses moléculaires de composés. Elles existent aussi dans tous les diamètres de particules Surf que nous proposons.

### Surf SO<sub>3</sub>H

Les colonnes surf SO<sub>3</sub>H, greffées acide sulfonique font partie des colonnes de type SCX. Elles peuvent être utilisées en phase inverse ou pour de la chromatographie à échange d'ions. Le greffon SO<sub>3</sub>H est un greffon échangeur de cations fort.


### Surf Phenyl et Surf Phenyl-hexyl

Nos silices Surf Phenyl et Surf Phenyl-hexyl vous apportent une alternative à la rétention et à la sélectivité des phases C18 habituelles. Les interactions hydrophobes bien équilibrées et les interactions π-π uniques en font les candidates idéales pour les analytes composés de cycles aromatiques ou de composés ayant de longs systèmes conjugués.




Les silices Surf Phenyl et Phenyl-hexyl existent dans toute notre gamme de diamètres de particules Surf.

Choisissez votre silice selon l'acidité des composés que vous souhaitez analyser :



### Surf DNH



Surf DNH est un nouveau type de silice  $\text{NH}_2$ , idéale pour la séparation des composés basiques forts. Surf DNH ainsi que la Surf  $\text{NH}_2$  peuvent être utilisés pour les purifications.

La Surf DNH a une quantité de liaisons amine 2 fois plus élevée que celle de la Surf  $\text{NH}_2$  et a pour greffons une amine secondaire ainsi qu'une amine primaire.

Ces greffages rendent la liaison hydrogène et les interactions électrostatiques plus faibles en raison d'une chaîne carbonée plus longue que celle de la Surf  $\text{NH}_2$ . Il est donc possible d'utiliser un solvant de faible polarité car la vitesse d'éluion est plus rapide qu'avec la Surf  $\text{NH}_2$ .

Surf DNH est de ce fait utilisée pour la séparation de composés basiques à très basiques qui adsorbent fortement sur la silice nue et/ou se décomposent au contact de celle-ci.

### Surf PEI (Poly-Ethylène-Imine)



Surf PEI possède un polymère aminé macromoléculaire à sa surface. La silice est greffée de groupements amines secondaires et tertiaires, le groupe amine y est trois fois plus grand que pour la Surf  $\text{NH}$ . Pour cette raison, la basicité est plus forte que la Surf  $\text{NH}$  et la liaison hydrogène y est plus faible en raison de la chaîne carbonée plus longue.

Surf PEI convient à la séparation des composés basiques forts, elle est également utilisable pour la purification.

### Surf EP2TNE et Surf EP4TNE

Cette silice greffée 2 Ethyl Pyridine ou 4 Ethyl Pyridine est polaire et le plus souvent utilisée en SFC plutôt qu'en HPLC.

Elle permet de travailler dans des conditions de phase normale et en SFC (Chromatographie super critique).

Ces deux phases vous offrent une excellente reproductibilité. Leur utilisation vous offre pour les molécules achirales, rétention, finesse de pics et sélectivité.

Le EP2TNE est particulièrement bien adapté à la séparation de la famille des composés acides et présente également une grande sélectivité vis-à-vis de nombreux autres composés.

Le EP4TNE vient compléter de façon alternative la sélectivité de la phase 2EP.

### Surf CN

La silice Surf CN (Cyanopropyl) est légèrement polaire et offre une sélectivité complémentaire en phase normale ou en phase inverse en augmentant la reconnaissance hydrophile.



La porosité et l'homogénéité de notre silice assure une analyse de haute reproductibilité et une grande stabilité du lit de la colonne. En phase inverse, la phase CN permet de séparer des composés qui ne peuvent pas l'être sur des phases C18 ou C8.

(Le passage du mode phase normale à phase inverse nécessite un ré équilibrage de la colonne avant utilisation.)

Le mécanisme de séparation est basé sur des interactions  $\pi$ - $\pi$  ainsi que des interactions polaires et hydrophobes.

## Notre gamme de End-Capping (traitement de surface) pour les Surf C4, C8 et C18

### Surf

- End-capping traditionnel pour une séparation des composés acides, basiques et neutres en une seule analyse

### Surf Wax

- Séparation des composés acides, basiques et neutres en une seule analyse.
- end-capping haute densité pour notre gamme Surf Wax, désactivée pour les produits basiques, afin de réduire les traînées de pics et d'offrir une meilleure capacité de charge.

### Surf Polar

- Séparation des composés acides, basiques et neutres en une seule analyse,
- Sans end-capping pour des mélanges à polarités complexes

### Surf TriF

Silice née de notre dernière génération de technologie de greffage.

- Hydrophobicité contrôlée permettant des temps de rétention plus courts que sur un greffage traditionnel pour offrir une analyse plus rapide tout en consommant moins de phase mobile.
- Idéal pour une utilisation en spectrométrie de masse, greffons C18 Trifonctionnel.
- Grande robustesse en milieu extrême.
- Sélectivité stérique.
- end-capping normal pour la gamme Surf TriF.

### Surf Polar TriF

Silice née de notre dernière génération de technologie de greffage.

- Hydrophobicité contrôlée permettant des temps de rétention plus courts que sur un greffage traditionnel pour offrir une analyse plus rapide tout en consommant moins de phase mobile.
- Idéal pour une utilisation en spectrométrie de masse, Greffons C18 Trifonctionnel.
- Grande robustesse en milieu extrême.
- Sélectivité stérique.
- Pas de end-capping pour notre gamme Surf Polar TriF, idéale pour les mélanges à polarités complexes et la séparation des peptides.

### Surf EXTREME

Silice née de notre dernière génération de technologie de greffage.

- Hydrophobicité contrôlée permettant des temps de rétention plus courts que sur un greffage traditionnel pour offrir une analyse plus rapide tout en consommant moins de phase mobile.
- Idéal pour une utilisation en spectrométrie de masse, greffons C18 Trifonctionnel.
- Grande robustesse en milieu extrême.
- Sélectivité stérique end-capping haute densité pour notre gamme Surf Extreme permettant de travailler sur une large gamme de pH 1 à 10.

### Surf AQ

- Stable de 100% de phase organique à 100% de phase aqueuse.
- La phase Surf AQ possède une phase greffée extrêmement stable qui vous permet d'utiliser une seule colonne dans des conditions de phases mobiles diverses. Vous pouvez travailler aisément dans des conditions fortement organiques ou fortement aqueuses.
- Excellente rétention des analytes polaires avec des phases mobiles jusqu'à 100% aqueuses.
- Forte rétention des analytes hydrophobes avec des phases mobiles jusqu'à 100% organiques.



Phase	Mono-Fonctionnel (M) / Tri-fonctionnel (T)	End-capped	Hydrophobicité	Rétention hydrophobe	Hydrophile	Très hydrophile	Acides	Bases
Surf	M	OUI	••	•••	••	•	••	••
Surf Wax	M	HIGH DENSITY	•••	•••	••	•	•••	•••
Surf Polar	M	NON	••	••	••	•	•	•
Surf TriF	T	OUI	••	••	•••	•••	•	•
Surf Polar TriF	T	NON	•••	•••	••	••	••	••
Surf Extreme	T	HIGH DENSITY	•••	•••	•	•	•••	•••
Surf AQ	M	OUI	••	••	••	••	••	••

Grâce à leurs différentes sélectivités, nos gammes Surf permettent de relever les différents défis des chromatographistes :

Vous souhaitez	Optez pour
Analyser composés acides, basiques, neutres, en une seule analyse	Surf
Analyser des mélanges polaires complexes	Surf Polar
Analyser des composés basiques	Surf Wax
Plus de robustesse, raccourcir les temps d'analyses, économiser la phase mobile	Surf TriF
Analyser des mélanges à polarités complexes	Surf Polar
Utiliser une phase mobile à large pH : {1 à 10}	Surf Extreme
Travailler de 100% de phase aqueuse à 100% de phase organique	Surf AQ
Analyser des composés polaires dans phases 100% aqueuses	
Analyser des composés hydrophobes dans une de phase 100% organique	

## Notre gamme de End-Capping (traitement de surface) pour les Surf C4, C8 et C18

Plus précisément, selon le greffage :

Solubilité, poids moléculaire et difficulté de séparation...Utilisation de chaque phase :

Phase	Hydrophobe	Rétention hydrophobe	Hydrophile	Très hydrophile	Acide	Basique	Echange d'ions
C18	++	+++	++	+	++	++	
C18 Wax	+++	+++	++	+	+++	+++	
C18 Polar	++	++	++	++	+	+	+
C18 TriF	+++	+++	++	++	++	++	
C18 Extrême	+++	+++	+	+	+++	+++	
C18 Polar TriF	++	++	+++	+++	+	+	+
C18 AQ Plus	+	+	+++	+++	+	+	
C18 AQ	++	++	++	++	++	++	

C8	++	++	++	+	+	+	
C8 Wax	+++	+++	+	+	+	++	
C8 AQ	+	+	+++	+++	+	+	

C4	+	+	++	+	+	+	
C4 Wax	++	++	+	+	+	++	
C4 AQ	+	+	+++	+	+	+	

CN	+	+	++	+	+	+	
Diol	+		++	+	+	+	
NH <sub>2</sub>	+		+++	++		+++	+
DNH	+		+++	++		+++	++
PEI	+	+	+++	+++		+++	+++
ARG	+	+	+++	+++	+++	+	+++
SO <sub>3</sub> H	+		+	++	++		++

Taille de particules	Tailles des pores	Greffons apolaires possibles	End-Capping en phase RP	Greffons polaires possibles
Surf 1µm	100A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 2µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 2,5µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 3µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 3,5µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 4µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 4,5µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	ARG - EP2TNE - EP4TNE, NH <sub>2</sub> , CN - Phenyl - Phenyl-hexyl - Diol
Surf 5µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 5,5µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE, NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 7µm	100A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 9µm	100A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Surf 10µm	100A - 150A - 200A - 300A	C4 - C8 - C18	Wax - Polar - triF - Extreme - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol

## Notre gamme Coral pour du scale-up

Cette gamme est plus particulièrement dédiée à la HPLC préparative avec des avantages significatifs tels qu'une haute efficacité grâce à une grande surface spécifique, une haute capacité et une grande stabilité chimique compatible avec une grande quantité de solvants. De plus, notre silice Coral est extrêmement résistante mécaniquement parlant et facile à manipuler.

Nous vous proposons une grande variété de tailles de particules pour permettre une grande variété d'analyses (de 5µm à 140µm) avec un grand choix de tailles de pores (de 70A à 1000A selon les particules) pour optimiser la capacité.

Afin de vous permettre des performances reproductibles, les lots sont homogènes entre eux.

Le vrac vous ouvre des possibilités d'applications à grande échelle, sans limitation de mise à l'échelle, permettant ainsi des avantages économiques pour ces procédés.

Les silices Coral SO<sub>3</sub>H ; DNH ; PEI ; ARG ; EP2TNE ; EP4TNE ; Amino ; CN ; Phenyl ; Phenyl-hexyl et Diol ont les mêmes particularités que les Surf du même nom. Il en est de même pour les traitements de surface.

### Coral NH



Les groupes fonctionnels aminés réagissent avec la cétone et/ou le groupe carboxyle. Notre gel de silice Coral NH peut séparer plusieurs composés de sucre dans des conditions de solvant organique-H<sub>2</sub>O sans présenter de réaction chimique.

La silice Coral NH est utilisée dans des conditions normales de phase, en particulier en mode HILIC (Hydrophilic Interaction Chromatography mode), mode dans lequel la vitesse d'éluion peut être augmentée lorsque le pourcentage d'eau dans le solvant organique est augmenté permettant un gain de temps d'analyse.

Sont disponibles dans la gamme Coral :

Taille de particules	Tailles des pores	Greffons RP possibles	End-Capping en phase RP	Autres greffons possibles
Coral 5µm	70A - 100A - 150A - 200A - 300A - 500A - 800A	C4 - C8 - C18	Wax - Polar - triF - Extrême - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Coral 10µm	70A - 100A - 150A - 200A - 300A - 500A - 800A - 1000A	C4 - C8 - C18	Wax - Polar - triF - Extrême - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - EP2TNE - EP4TNE - NH <sub>2</sub> - CN - Phenyl - Phenyl-hexyl - Diol
Coral 15µm	70A - 100A - 150A - 200A - 300A - 500A - 800A - 1000A	C4 - C8 - C18	Wax - Polar - triF - Extrême - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - Amino - CN - Phenyl - Phenyl-hexyl - Diol
Coral 30µm	70A - 100A - 150A - 200A - 300A - 500A - 800A - 1000A	C4 - C8 - C18	Wax - Polar - triF - Extrême - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - Amino - CN - Phenyl - Phenyl-hexyl - Diol
Coral 60µm	55A - 70A - 100A - 300A - 500A - 800A - 1000A	C4 - C8 - C18	TriF - Extrême - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - Amino - CN - Diol
Coral 140µm	55A - 70A - 100A - 300A - 500A - 800A - 1000A	C4 - C8 - C18	TriF - Extrême - Polar - AQ plus - AQ - Polar triF	SO <sub>3</sub> H - DNH - PEI - ARG - Amino - CN - Diol



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 1,5µm</b>	
<b>Particule Surf 1,5µm, pore 100A</b>	
B-100A-1.5-SU-Si-100G	Bulk Surf 100A 1.5µm 100g
B-100A-1.5-SU-Si-1KG	Bulk Surf 100A 1.5µm 1kg
<b>C8</b>	
B-100A-1.5-SU-C8-100G	Bulk Surf 100A 1.5µm C8 100g
B-100A-1.5-SU-C8-1KG	Bulk Surf 100A 1.5µm C8 1kg
B-100A-1.5-SU-C8W-100G	Bulk Surf 100A 1.5µm C8 Wax 100g
B-100A-1.5-SU-C8W-1KG	Bulk Surf 100A 1.5µm C8 Wax 1kg
B-100A-1.5-SU-C8P-100G	Bulk Surf 100A 1.5µm C8 Polar 100g
B-100A-1.5-SU-C8P-1KG	Bulk Surf 100A 1.5µm C8 Polar 1kg
B-100A-1.5-SU-C8T-100G	Bulk Surf 100A 1.5µm C8 TriF 100g
B-100A-1.5-SU-C8T-1KG	Bulk Surf 100A 1.5µm C8 TriF 1kg
B-100A-1.5-SU-C8EX-100G	Bulk Surf 100A 1.5µm C8 Extreme 100g
B-100A-1.5-SU-C8EX-1KG	Bulk Surf 100A 1.5µm C8 Extreme 1kg
B-100A-1.5-SU-C8PT-100G	Bulk Surf 100A 1.5µm C8 Polar TriF 100g
B-100A-1.5-SU-C8PT-1KG	Bulk Surf 100A 1.5µm C8 Polar TriF 1kg
B-100A-1.5-SU-C8AQP-100G	Bulk Surf 100A 1.5µm C8 AQ Plus 100g
B-100A-1.5-SU-C8AQP-1KG	Bulk Surf 100A 1.5µm C8 AQ Plus 1kg
B-100A-1.5-SU-C8AQ-100G	Bulk Surf 100A 1.5µm C8 AQ 100g
B-100A-1.5-SU-C8AQ-1KG	Bulk Surf 100A 1.5µm C8 AQ 1kg
<b>C18</b>	
B-100A-1.5-SU-C18-100G	Bulk Surf 100A 1.5µm C18 100g
B-100A-1.5-SU-C18-1KG	Bulk Surf 100A 1.5µm C18 1kg
B-100A-1.5-SU-C18W-100G	Bulk Surf 100A 1.5µm C18 Wax 100g
B-100A-1.5-SU-C18W-1KG	Bulk Surf 100A 1.5µm C18 Wax 1kg
B-100A-1.5-SU-C18P-100G	Bulk Surf 100A 1.5µm C18 Polar 100g
B-100A-1.5-SU-C18P-1KG	Bulk Surf 100A 1.5µm C18 Polar 1kg
B-100A-1.5-SU-C18T-100G	Bulk Surf 100A 1.5µm C18 TriF 100g
B-100A-1.5-SU-C18T-1KG	Bulk Surf 100A 1.5µm C18 TriF 1kg
B-100A-1.5-SU-C18EX-100G	Bulk Surf 100A 1.5µm C18 Extreme 100g
B-100A-1.5-SU-C18EX-1KG	Bulk Surf 100A 1.5µm C18 Extreme 1kg
B-100A-1.5-SU-C18PT-100G	Bulk Surf 100A 1.5µm C18 Polar TriF 100g
B-100A-1.5-SU-C18PT-1KG	Bulk Surf 100A 1.5µm C18 Polar TriF 1kg
B-100A-1.5-SU-C18AQP-100G	Bulk Surf 100A 1.5µm C18 AQ Plus 100g
B-100A-1.5-SU-C18AQP-1KG	Bulk Surf 100A 1.5µm C18 AQ Plus 1kg
B-100A-1.5-SU-C18AQ-100G	Bulk Surf 100A 1.5µm C18 AQ 100g
B-100A-1.5-SU-C18AQ-1KG	Bulk Surf 100A 1.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-1.5-SU-C4-100G	Bulk Surf 100A 1.5µm C4 100g
B-100A-1.5-SU-C4-1KG	Bulk Surf 100A 1.5µm C4 1kg
B-100A-1.5-SU-C4W-100G	Bulk Surf 100A 1.5µm C4 Wax 100g
B-100A-1.5-SU-C4W-1KG	Bulk Surf 100A 1.5µm C4 Wax 1kg
B-100A-1.5-SU-C4P-100G	Bulk Surf 100A 1.5µm C4 Polar 100g
B-100A-1.5-SU-C4P-1KG	Bulk Surf 100A 1.5µm C4 Polar 1kg
B-100A-1.5-SU-C4T-100G	Bulk Surf 100A 1.5µm C4 TriF 100g
B-100A-1.5-SU-C4T-1KG	Bulk Surf 100A 1.5µm C4 TriF 1kg
B-100A-1.5-SU-C4EX-100G	Bulk Surf 100A 1.5µm C4 Extreme 100g
B-100A-1.5-SU-C4EX-1KG	Bulk Surf 100A 1.5µm C4 Extreme 1kg
B-100A-1.5-SU-C4PT-100G	Bulk Surf 100A 1.5µm C4 Polar TriF 100g
B-100A-1.5-SU-C4PT-1KG	Bulk Surf 100A 1.5µm C4 Polar TriF 1kg
B-100A-1.5-SU-C4AQP-100G	Bulk Surf 100A 1.5µm C4 AQ Plus 100g
B-100A-1.5-SU-C4AQP-1KG	Bulk Surf 100A 1.5µm C4 AQ Plus 1kg
B-100A-1.5-SU-C4AQ-100G	Bulk Surf 100A 1.5µm C4 AQ 100g
B-100A-1.5-SU-C4AQ-1KG	Bulk Surf 100A 1.5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-1.5-SU-ARG-100G	Bulk Surf 100A 1.5µm ARG 100g
B-100A-1.5-SU-ARG-1KG	Bulk Surf 100A 1.5µm ARG 1kg
B-100A-1.5-SU-EP2TNE-100G	Bulk Surf 100A 1.5µm EP2TNE 100g
B-100A-1.5-SU-EP2TNE-1KG	Bulk Surf 100A 1.5µm EP2TNE 1kg
B-100A-1.5-SU-EP4TNE-100G	Bulk Surf 100A 1.5µm EP4TNE 100g
B-100A-1.5-SU-EP4TNE-1KG	Bulk Surf 100A 1.5µm EP4TNE 1kg
B-100A-1.5-SU-NH2-100G	Bulk Surf 100A 1.5µm NH2 100g
B-100A-1.5-SU-NH2-1KG	Bulk Surf 100A 1.5µm NH2 1kg
B-100A-1.5-SU-CN-100G	Bulk Surf 100A 1.5µm CN 100g
B-100A-1.5-SU-CN-1KG	Bulk Surf 100A 1.5µm CN 1kg
B-100A-1.5-SU-Phenyl-100G	Bulk Surf 100A 1.5µm Phenyl 100g
B-100A-1.5-SU-Phenyl-1KG	Bulk Surf 100A 1.5µm Phenyl 1kg
B-100A-1.5-SU-Phenyl HE-100G	Bulk Surf 100A 1.5µm Phenyl HE 100g
B-100A-1.5-SU-Phenyl HE-1KG	Bulk Surf 100A 1.5µm Phenyl HE 1kg
B-100A-1.5-SU-Diol-100G	Bulk Surf 100A 1.5µm Diol 100g
B-100A-1.5-SU-Diol-1KG	Bulk Surf 100A 1.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 2 µm</b>	
<b>Particule Surf 2µm, pore 100A</b>	
Code article	Désignation
B-100A-2-SU-Si-100G	Bulk Surf Silice 100A 2µm 100g
B-100A-2-SU-Si-1KG	Bulk Surf Silice 100A 2µm 1kg
<b>C8</b>	
B-100A-2-SU-C8-100G	Bulk Surf 100A 2µm C8 100g
B-100A-2-SU-C8-1KG	Bulk Surf 100A 2µm C8 1kg
B-100A-2-SU-C8W-100G	Bulk Surf 100A 2µm C8 Wax 100g
B-100A-2-SU-C8W-1KG	Bulk Surf 100A 2µm C8 Wax 1kg
B-100A-2-SU-C8P-100G	Bulk Surf 100A 2µm C8 Polar 100g
B-100A-2-SU-C8P-1KG	Bulk Surf 100A 2µm C8 Polar 1kg
B-100A-2-SU-C8T-100G	Bulk Surf 100A 2µm C8 TriF 100g
B-100A-2-SU-C8T-1KG	Bulk Surf 100A 2µm C8 TriF 1kg
B-100A-2-SU-C8EX-100G	Bulk Surf 100A 2µm C8 Extreme 100g
B-100A-2-SU-C8EX-1KG	Bulk Surf 100A 2µm C8 Extreme 1kg
B-100A-2-SU-C8PT-100G	Bulk Surf 100A 2µm C8 Polar TriF 100g
B-100A-2-SU-C8PT-1KG	Bulk Surf 100A 2µm C8 Polar TriF 1kg
B-100A-2-SU-C8AQP-100G	Bulk Surf 100A 2µm C8 AQ Plus 100g
B-100A-2-SU-C8AQP-1KG	Bulk Surf 100A 2µm C8 AQ Plus 1kg
B-100A-2-SU-C8AQ-100G	Bulk Surf 100A 2µm C8 AQ 100g
B-100A-2-SU-C8AQ-1KG	Bulk Surf 100A 2µm C8 AQ 1kg
<b>C18</b>	
B-100A-2-SU-C18-100G	Bulk Surf 100A 2µm C18 100g
B-100A-2-SU-C18-1KG	Bulk Surf 100A 2µm C18 1kg
B-100A-2-SU-C18W-100G	Bulk Surf 100A 2µm C18 Wax 100g
B-100A-2-SU-C18W-1KG	Bulk Surf 100A 2µm C18 Wax 1kg
B-100A-2-SU-C18P-100G	Bulk Surf 100A 2µm C18 Polar 100g
B-100A-2-SU-C18P-1KG	Bulk Surf 100A 2µm C18 Polar 1kg
B-100A-2-SU-C18T-100G	Bulk Surf 100A 2µm C18 TriF 100g
B-100A-2-SU-C18T-1KG	Bulk Surf 100A 2µm C18 TriF 1kg
B-100A-2-SU-C18EX-100G	Bulk Surf 100A 2µm C18 Extreme 100g
B-100A-2-SU-C18EX-1KG	Bulk Surf 100A 2µm C18 Extreme 1kg
B-100A-2-SU-C18PT-100G	Bulk Surf 100A 2µm C18 Polar TriF 100g
B-100A-2-SU-C18PT-1KG	Bulk Surf 100A 2µm C18 Polar TriF 1kg
B-100A-2-SU-C18AQP-100G	Bulk Surf 100A 2µm C18 AQ Plus 100g
B-100A-2-SU-C18AQP-1KG	Bulk Surf 100A 2µm C18 AQ Plus 1kg
B-100A-2-SU-C18AQ-100G	Bulk Surf 100A 2µm C18 AQ 100g
B-100A-2-SU-C18AQ-1KG	Bulk Surf 100A 2µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-2-SU-C4-100G	Bulk Surf 100A 2µm C4 100g
B-100A-2-SU-C4-1KG	Bulk Surf 100A 2µm C4 1kg
B-100A-2-SU-C4W-100G	Bulk Surf 100A 2µm C4 Wax 100g
B-100A-2-SU-C4W-1KG	Bulk Surf 100A 2µm C4 Wax 1kg
B-100A-2-SU-C4P-100G	Bulk Surf 100A 2µm C4 Polar 100g
B-100A-2-SU-C4P-1KG	Bulk Surf 100A 2µm C4 Polar 1kg
B-100A-2-SU-C4T-100G	Bulk Surf 100A 2µm C4 TriF 100g
B-100A-2-SU-C4T-1KG	Bulk Surf 100A 2µm C4 TriF 1kg
B-100A-2-SU-C4EX-100G	Bulk Surf 100A 2µm C4 Extreme 100g
B-100A-2-SU-C4EX-1KG	Bulk Surf 100A 2µm C4 Extreme 1kg
B-100A-2-SU-C4PT-100G	Bulk Surf 100A 2µm C4 Polar TriF 100g
B-100A-2-SU-C4PT-1KG	Bulk Surf 100A 2µm C4 Polar TriF 1kg
B-100A-2-SU-C4AQP-100G	Bulk Surf 100A 2µm C4 AQ Plus 100g
B-100A-2-SU-C4AQP-1KG	Bulk Surf 100A 2µm C4 AQ Plus 1kg
B-100A-2-SU-C4AQ-100G	Bulk Surf 100A 2µm C4 AQ 100g
B-100A-2-SU-C4AQ-1KG	Bulk Surf 100A 2µm C4 AQ 1kg
<b>Autres</b>	
B-100A-2-SU-ARG-100G	Bulk Surf 100A 2µm ARG 100g
B-100A-2-SU-ARG-1KG	Bulk Surf 100A 2µm ARG 1kg
B-100A-2-SU-EP2TNE-100G	Bulk Surf 100A 2µm EP2TNE 100g
B-100A-2-SU-EP2TNE-1KG	Bulk Surf 100A 2µm EP2TNE 1kg
B-100A-2-SU-EP4TNE-100G	Bulk Surf 100A 2µm EP4TNE 100g
B-100A-2-SU-EP4TNE-1KG	Bulk Surf 100A 2µm EP4TNE 1kg
B-100A-2-SU-NH2-100G	Bulk Surf 100A 2µm NH2 100g
B-100A-2-SU-NH2-1KG	Bulk Surf 100A 2µm NH2 1kg
B-100A-2-SU-CN-100G	Bulk Surf 100A 2µm CN 100g
B-100A-2-SU-CN-1KG	Bulk Surf 100A 2µm CN 1kg
B-100A-2-SU-Phenyl-100G	Bulk Surf 100A 2µm Phenyl 100g
B-100A-2-SU-Phenyl-1KG	Bulk Surf 100A 2µm Phenyl 1kg
B-100A-2-SU-Phenyl HE-100G	Bulk Surf 100A 2µm Phenyl HE 100g
B-100A-2-SU-Phenyl HE-1KG	Bulk Surf 100A 2µm Phenyl HE 1kg
B-100A-2-SU-Diol-100G	Bulk Surf 100A 2µm Diol 100g
B-100A-2-SU-Diol-1KG	Bulk Surf 100A 2µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 2µm, pore 150A</b>	
B-150A-2-SU-Si-100G	Bulk Surf 150A 2µm 100g
B-150A-2-SU-Si-1KG	Bulk Surf 150A 2µm 1kg
<b>C8</b>	
B-150A-2-SU-C8-100G	Bulk Surf 150A 2µm C8 100g
B-150A-2-SU-C8-1KG	Bulk Surf 150A 2µm C8 1kg
B-150A-2-SU-C8W-100G	Bulk Surf 150A 2µm C8 Wax 100g
B-150A-2-SU-C8W-1KG	Bulk Surf 150A 2µm C8 Wax 1kg
B-150A-2-SU-C8P-100G	Bulk Surf 150A 2µm C8 Polar 100g
B-150A-2-SU-C8P-1KG	Bulk Surf 150A 2µm C8 Polar 1kg
B-150A-2-SU-C8T-100G	Bulk Surf 150A 2µm C8 TriF 100g
B-150A-2-SU-C8T-1KG	Bulk Surf 150A 2µm C8 TriF 1kg
B-150A-2-SU-C8EX-100G	Bulk Surf 150A 2µm C8 Extreme 100g
B-150A-2-SU-C8EX-1KG	Bulk Surf 150A 2µm C8 Extreme 1kg
B-150A-2-SU-C8PT-100G	Bulk Surf 150A 2µm C8 Polar TriF 100g
B-150A-2-SU-C8PT-1KG	Bulk Surf 150A 2µm C8 Polar TriF 1kg
B-150A-2-SU-C8AQP-100G	Bulk Surf 150A 2µm C8 AQ Plus 100g
B-150A-2-SU-C8AQP-1KG	Bulk Surf 150A 2µm C8 AQ Plus 1kg
B-150A-2-SU-C8AQ-100G	Bulk Surf 150A 2µm C8 AQ 100g
B-150A-2-SU-C8AQ-1KG	Bulk Surf 150A 2µm C8 AQ 1kg
<b>C18</b>	
B-150A-2-SU-C18-100G	Bulk Surf 150A 2µm C18 100g
B-150A-2-SU-C18-1KG	Bulk Surf 150A 2µm C18 1kg
B-150A-2-SU-C18W-100G	Bulk Surf 150A 2µm C18 Wax 100g
B-150A-2-SU-C18W-1KG	Bulk Surf 150A 2µm C18 Wax 1kg
B-150A-2-SU-C18P-100G	Bulk Surf 150A 2µm C18 Polar 100g
B-150A-2-SU-C18P-1KG	Bulk Surf 150A 2µm C18 Polar 1kg
B-150A-2-SU-C18T-100G	Bulk Surf 150A 2µm C18 TriF 100g
B-150A-2-SU-C18T-1KG	Bulk Surf 150A 2µm C18 TriF 1kg
B-150A-2-SU-C18EX-100G	Bulk Surf 150A 2µm C18 Extreme 100g
B-150A-2-SU-C18EX-1KG	Bulk Surf 150A 2µm C18 Extreme 1kg
B-150A-2-SU-C18PT-100G	Bulk Surf 150A 2µm C18 Polar TriF 100g
B-150A-2-SU-C18PT-1KG	Bulk Surf 150A 2µm C18 Polar TriF 1kg
B-150A-2-SU-C18AQP-100G	Bulk Surf 150A 2µm C18 AQ Plus 100g
B-150A-2-SU-C18AQP-1KG	Bulk Surf 150A 2µm C18 AQ Plus 1kg
B-150A-2-SU-C18AQ-100G	Bulk Surf 150A 2µm C18 AQ 100g
B-150A-2-SU-C18AQ-1KG	Bulk Surf 150A 2µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-150A-2-SU-C4-100G	Bulk Surf 150A 2µm C4 100g
B-150A-2-SU-C4-1KG	Bulk Surf 150A 2µm C4 1kg
B-150A-2-SU-C4W-100G	Bulk Surf 150A 2µm C4 Wax 100g
B-150A-2-SU-C4W-1KG	Bulk Surf 150A 2µm C4 Wax 1kg
B-150A-2-SU-C4P-100G	Bulk Surf 150A 2µm C4 Polar 100g
B-150A-2-SU-C4P-1KG	Bulk Surf 150A 2µm C4 Polar 1kg
B-150A-2-SU-C4T-100G	Bulk Surf 150A 2µm C4 TriF 100g
B-150A-2-SU-C4T-1KG	Bulk Surf 150A 2µm C4 TriF 1kg
B-150A-2-SU-C4EX-100G	Bulk Surf 150A 2µm C4 Extreme 100g
B-150A-2-SU-C4EX-1KG	Bulk Surf 150A 2µm C4 Extreme 1kg
B-150A-2-SU-C4PT-100G	Bulk Surf 150A 2µm C4 Polar TriF 100g
B-150A-2-SU-C4PT-1KG	Bulk Surf 150A 2µm C4 Polar TriF 1kg
B-150A-2-SU-C4AQP-100G	Bulk Surf 150A 2µm C4 AQ Plus 100g
B-150A-2-SU-C4AQP-1KG	Bulk Surf 150A 2µm C4 AQ Plus 1kg
B-150A-2-SU-C4AQ-100G	Bulk Surf 150A 2µm C4 AQ 100g
B-150A-2-SU-C4AQ-1KG	Bulk Surf 150A 2µm C4 AQ 1kg
<b>Autres</b>	
B-150A-2-SU-ARG-100G	Bulk Surf 150A 2µm ARG 100g
B-150A-2-SU-ARG-1KG	Bulk Surf 150A 2µm ARG 1kg
B-150A-2-SU-EP2TNE-100G	Bulk Surf 150A 2µm EP2TNE 100g
B-150A-2-SU-EP2TNE-1KG	Bulk Surf 150A 2µm EP2TNE 1kg
B-150A-2-SU-EP4TNE-100G	Bulk Surf 150A 2µm EP4TNE 100g
B-150A-2-SU-EP4TNE-1KG	Bulk Surf 150A 2µm EP4TNE 1kg
B-150A-2-SU-NH2-100G	Bulk Surf 150A 2µm NH2 100g
B-150A-2-SU-NH2-1KG	Bulk Surf 150A 2µm NH2 1kg
B-150A-2-SU-CN-100G	Bulk Surf 150A 2µm CN 100g
B-150A-2-SU-CN-1KG	Bulk Surf 150A 2µm CN 1kg
B-150A-2-SU-Phenyl-100G	Bulk Surf 150A 2µm Phenyl 100g
B-150A-2-SU-Phenyl-1KG	Bulk Surf 150A 2µm Phenyl 1kg
B-150A-2-SU-Phenyl HE-100G	Bulk Surf 150A 2µm Phenyl HE 100g
B-150A-2-SU-Phenyl HE-1KG	Bulk Surf 150A 2µm Phenyl HE 1kg
B-150A-2-SU-Diol-100G	Bulk Surf 150A 2µm Diol 100g
B-150A-2-SU-Diol-1KG	Bulk Surf 150A 2µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 2µm, pore 200A</b>	
B-200A-2-SU-Si-100G	Bulk Surf 200A 2µm 100g
B-200A-2-SU-Si-1KG	Bulk Surf 200A 2µm 1kg
<b>C8</b>	
B-200A-2-SU-C8-100G	Bulk Surf 200A 2µm C8 100g
B-200A-2-SU-C8-1KG	Bulk Surf 200A 2µm C8 1kg
B-200A-2-SU-C8W-100G	Bulk Surf 200A 2µm C8 Wax 100g
B-200A-2-SU-C8W-1KG	Bulk Surf 200A 2µm C8 Wax 1kg
B-200A-2-SU-C8P-100G	Bulk Surf 200A 2µm C8 Polar 100g
B-200A-2-SU-C8P-1KG	Bulk Surf 200A 2µm C8 Polar 1kg
B-200A-2-SU-C8T-100G	Bulk Surf 200A 2µm C8 TriF 100g
B-200A-2-SU-C8T-1KG	Bulk Surf 200A 2µm C8 TriF 1kg
B-200A-2-SU-C8EX-100G	Bulk Surf 200A 2µm C8 Extreme 100g
B-200A-2-SU-C8EX-1KG	Bulk Surf 200A 2µm C8 Extreme 1kg
B-200A-2-SU-C8PT-100G	Bulk Surf 200A 2µm C8 Polar TriF 100g
B-200A-2-SU-C8PT-1KG	Bulk Surf 200A 2µm C8 Polar TriF 1kg
B-200A-2-SU-C8AQP-100G	Bulk Surf 200A 2µm C8 AQ Plus 100g
B-200A-2-SU-C8AQP-1KG	Bulk Surf 200A 2µm C8 AQ Plus 1kg
B-200A-2-SU-C8AQ-100G	Bulk Surf 200A 2µm C8 AQ 100g
B-200A-2-SU-C8AQ-1KG	Bulk Surf 200A 2µm C8 AQ 1kg
<b>C18</b>	
B-200A-2-SU-C18-100G	Bulk Surf 200A 2µm C18 100g
B-200A-2-SU-C18-1KG	Bulk Surf 200A 2µm C18 1kg
B-200A-2-SU-C18W-100G	Bulk Surf 200A 2µm C18 Wax 100g
B-200A-2-SU-C18W-1KG	Bulk Surf 200A 2µm C18 Wax 1kg
B-200A-2-SU-C18P-100G	Bulk Surf 200A 2µm C18 Polar 100g
B-200A-2-SU-C18P-1KG	Bulk Surf 200A 2µm C18 Polar 1kg
B-200A-2-SU-C18T-100G	Bulk Surf 200A 2µm C18 TriF 100g
B-200A-2-SU-C18T-1KG	Bulk Surf 200A 2µm C18 TriF 1kg
B-200A-2-SU-C18EX-100G	Bulk Surf 200A 2µm C18 Extreme 100g
B-200A-2-SU-C18EX-1KG	Bulk Surf 200A 2µm C18 Extreme 1kg
B-200A-2-SU-C18PT-100G	Bulk Surf 200A 2µm C18 Polar TriF 100g
B-200A-2-SU-C18PT-1KG	Bulk Surf 200A 2µm C18 Polar TriF 1kg
B-200A-2-SU-C18AQP-100G	Bulk Surf 200A 2µm C18 AQ Plus 100g
B-200A-2-SU-C18AQP-1KG	Bulk Surf 200A 2µm C18 AQ Plus 1kg
B-200A-2-SU-C18AQ-100G	Bulk Surf 200A 2µm C18 AQ 100g
B-200A-2-SU-C18AQ-1KG	Bulk Surf 200A 2µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-200A-2-SU-C4-100G	Bulk Surf 200A 2µm C4 100g
B-200A-2-SU-C4-1KG	Bulk Surf 200A 2µm C4 1kg
B-200A-2-SU-C4W-100G	Bulk Surf 200A 2µm C4 Wax 100g
B-200A-2-SU-C4W-1KG	Bulk Surf 200A 2µm C4 Wax 1kg
B-200A-2-SU-C4P-100G	Bulk Surf 200A 2µm C4 Polar 100g
B-200A-2-SU-C4P-1KG	Bulk Surf 200A 2µm C4 Polar 1kg
B-200A-2-SU-C4T-100G	Bulk Surf 200A 2µm C4 TriF 100g
B-200A-2-SU-C4T-1KG	Bulk Surf 200A 2µm C4 TriF 1kg
B-200A-2-SU-C4EX-100G	Bulk Surf 200A 2µm C4 Extreme 100g
B-200A-2-SU-C4EX-1KG	Bulk Surf 200A 2µm C4 Extreme 1kg
B-200A-2-SU-C4PT-100G	Bulk Surf 200A 2µm C4 Polar TriF 100g
B-200A-2-SU-C4PT-1KG	Bulk Surf 200A 2µm C4 Polar TriF 1kg
B-200A-2-SU-C4AQP-100G	Bulk Surf 200A 2µm C4 AQ Plus 100g
B-200A-2-SU-C4AQP-1KG	Bulk Surf 200A 2µm C4 AQ Plus 1kg
B-200A-2-SU-C4AQ-100G	Bulk Surf 200A 2µm C4 AQ 100g
B-200A-2-SU-C4AQ-1KG	Bulk Surf 200A 2µm C4 AQ 1kg
<b>Autres</b>	
B-200A-2-SU-ARG-100G	Bulk Surf 200A 2µm ARG 100g
B-200A-2-SU-ARG-1KG	Bulk Surf 200A 2µm ARG 1kg
B-200A-2-SU-EP2TNE-100G	Bulk Surf 200A 2µm EP2TNE 100g
B-200A-2-SU-EP2TNE-1KG	Bulk Surf 200A 2µm EP2TNE 1kg
B-200A-2-SU-EP4TNE-100G	Bulk Surf 200A 2µm EP4TNE 100g
B-200A-2-SU-EP4TNE-1KG	Bulk Surf 200A 2µm EP4TNE 1kg
B-200A-2-SU-NH2-100G	Bulk Surf 200A 2µm NH2 100g
B-200A-2-SU-NH2-1KG	Bulk Surf 200A 2µm NH2 1kg
B-200A-2-SU-CN-100G	Bulk Surf 200A 2µm CN 100g
B-200A-2-SU-CN-1KG	Bulk Surf 200A 2µm CN 1kg
B-200A-2-SU-Phenyl-100G	Bulk Surf 200A 2µm Phenyl 100g
B-200A-2-SU-Phenyl-1KG	Bulk Surf 200A 2µm Phenyl 1kg
B-200A-2-SU-Phenyl HE-100G	Bulk Surf 200A 2µm Phenyl HE 100g
B-200A-2-SU-Phenyl HE-1KG	Bulk Surf 200A 2µm Phenyl HE 1kg
B-200A-2-SU-Diol-100G	Bulk Surf 200A 2µm Diol 100g
B-200A-2-SU-Diol-1KG	Bulk Surf 200A 2µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 2µm, pore 300A</b>	
B-300A-2-SU-Si-100G	Bulk Surf 300A 2µm 100g
B-300A-2-SU-Si-1KG	Bulk Surf 300A 2µm 1kg
<b>C8</b>	
B-300A-2-SU-C8-100G	Bulk Surf 300A 2µm C8 100g
B-300A-2-SU-C8-1KG	Bulk Surf 300A 2µm C8 1kg
B-300A-2-SU-C8W-100G	Bulk Surf 300A 2µm C8 Wax 100g
B-300A-2-SU-C8W-1KG	Bulk Surf 300A 2µm C8 Wax 1kg
B-300A-2-SU-C8P-100G	Bulk Surf 300A 2µm C8 Polar 100g
B-300A-2-SU-C8P-1KG	Bulk Surf 300A 2µm C8 Polar 1kg
B-300A-2-SU-C8T-100G	Bulk Surf 300A 2µm C8 TriF 100g
B-300A-2-SU-C8T-1KG	Bulk Surf 300A 2µm C8 TriF 1kg
B-300A-2-SU-C8EX-100G	Bulk Surf 300A 2µm C8 Extreme 100g
B-300A-2-SU-C8EX-1KG	Bulk Surf 300A 2µm C8 Extreme 1kg
B-300A-2-SU-C8PT-100G	Bulk Surf 300A 2µm C8 Polar TriF 100g
B-300A-2-SU-C8PT-1KG	Bulk Surf 300A 2µm C8 Polar TriF 1kg
B-300A-2-SU-C8AQP-100G	Bulk Surf 300A 2µm C8 AQ Plus 100g
B-300A-2-SU-C8AQP-1KG	Bulk Surf 300A 2µm C8 AQ Plus 1kg
B-300A-2-SU-C8AQ-100G	Bulk Surf 300A 2µm C8 AQ 100g
B-300A-2-SU-C8AQ-1KG	Bulk Surf 300A 2µm C8 AQ 1kg
<b>C18</b>	
B-300A-2-SU-C18-100G	Bulk Surf 300A 2µm C18 100g
B-300A-2-SU-C18-1KG	Bulk Surf 300A 2µm C18 1kg
B-300A-2-SU-C18W-100G	Bulk Surf 300A 2µm C18 Wax 100g
B-300A-2-SU-C18W-1KG	Bulk Surf 300A 2µm C18 Wax 1kg
B-300A-2-SU-C18P-100G	Bulk Surf 300A 2µm C18 Polar 100g
B-300A-2-SU-C18P-1KG	Bulk Surf 300A 2µm C18 Polar 1kg
B-300A-2-SU-C18T-100G	Bulk Surf 300A 2µm C18 TriF 100g
B-300A-2-SU-C18T-1KG	Bulk Surf 300A 2µm C18 TriF 1kg
B-300A-2-SU-C18EX-100G	Bulk Surf 300A 2µm C18 Extreme 100g
B-300A-2-SU-C18EX-1KG	Bulk Surf 300A 2µm C18 Extreme 1kg
B-300A-2-SU-C18PT-100G	Bulk Surf 300A 2µm C18 Polar TriF 100g
B-300A-2-SU-C18PT-1KG	Bulk Surf 300A 2µm C18 Polar TriF 1kg
B-300A-2-SU-C18AQP-100G	Bulk Surf 300A 2µm C18 AQ Plus 100g
B-300A-2-SU-C18AQP-1KG	Bulk Surf 300A 2µm C18 AQ Plus 1kg
B-300A-2-SU-C18AQ-100G	Bulk Surf 300A 2µm C18 AQ 100g
B-300A-2-SU-C18AQ-1KG	Bulk Surf 300A 2µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-2-SU-C4-100G	Bulk Surf 300A 2µm C4 100g
B-300A-2-SU-C4-1KG	Bulk Surf 300A 2µm C4 1kg
B-300A-2-SU-C4W-100G	Bulk Surf 300A 2µm C4 Wax 100g
B-300A-2-SU-C4W-1KG	Bulk Surf 300A 2µm C4 Wax 1kg
B-300A-2-SU-C4P-100G	Bulk Surf 300A 2µm C4 Polar 100g
B-300A-2-SU-C4P-1KG	Bulk Surf 300A 2µm C4 Polar 1kg
B-300A-2-SU-C4T-100G	Bulk Surf 300A 2µm C4 TriF 100g
B-300A-2-SU-C4T-1KG	Bulk Surf 300A 2µm C4 TriF 1kg
B-300A-2-SU-C4EX-100G	Bulk Surf 300A 2µm C4 Extreme 100g
B-300A-2-SU-C4EX-1KG	Bulk Surf 300A 2µm C4 Extreme 1kg
B-300A-2-SU-C4PT-100G	Bulk Surf 300A 2µm C4 Polar TriF 100g
B-300A-2-SU-C4PT-1KG	Bulk Surf 300A 2µm C4 Polar TriF 1kg
B-300A-2-SU-C4AQP-100G	Bulk Surf 300A 2µm C4 AQ Plus 100g
B-300A-2-SU-C4AQP-1KG	Bulk Surf 300A 2µm C4 AQ Plus 1kg
B-300A-2-SU-C4AQ-100G	Bulk Surf 300A 2µm C4 AQ 100g
B-300A-2-SU-C4AQ-1KG	Bulk Surf 300A 2µm C4 AQ 1kg
<b>Autres</b>	
B-300A-2-SU-ARG-100G	Bulk Surf 300A 2µm ARG 100g
B-300A-2-SU-ARG-1KG	Bulk Surf 300A 2µm ARG 1kg
B-300A-2-SU-EP2TNE-100G	Bulk Surf 300A 2µm EP2TNE 100g
B-300A-2-SU-EP2TNE-1KG	Bulk Surf 300A 2µm EP2TNE 1kg
B-300A-2-SU-EP4TNE-100G	Bulk Surf 300A 2µm EP4TNE 100g
B-300A-2-SU-EP4TNE-1KG	Bulk Surf 300A 2µm EP4TNE 1kg
B-300A-2-SU-NH2-100G	Bulk Surf 300A 2µm NH2 100g
B-300A-2-SU-NH2-1KG	Bulk Surf 300A 2µm NH2 1kg
B-300A-2-SU-CN-100G	Bulk Surf 300A 2µm CN 100g
B-300A-2-SU-CN-1KG	Bulk Surf 300A 2µm CN 1kg
B-300A-2-SU-Phenyl-100G	Bulk Surf 300A 2µm Phenyl 100g
B-300A-2-SU-Phenyl-1KG	Bulk Surf 300A 2µm Phenyl 1kg
B-300A-2-SU-Phenyl HE-100G	Bulk Surf 300A 2µm Phenyl HE 100g
B-300A-2-SU-Phenyl HE-1KG	Bulk Surf 300A 2µm Phenyl HE 1kg
B-300A-2-SU-Diol-100G	Bulk Surf 300A 2µm Diol 100g
B-300A-2-SU-Diol-1KG	Bulk Surf 300A 2µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 2,5 µm</b>	
<b>Particule Surf 2,5µm, pore 100A</b>	
B-100A-2.5-SU-Si-100G	Bulk Surf Silice 100A 2.5µm 100g
B-100A-2.5-SU-Si-1KG	Bulk Surf Silice 100A 2.5µm 1kg
<b>C8</b>	
B-100A-2.5-SU-C8-100G	Bulk Surf 100A 2.5µm C8 100g
B-100A-2.5-SU-C8-1KG	Bulk Surf 100A 2.5µm C8 1kg
B-100A-2.5-SU-C8W-100G	Bulk Surf 100A 2.5µm C8 Wax 100g
B-100A-2.5-SU-C8W-1KG	Bulk Surf 100A 2.5µm C8 Wax 1kg
B-100A-2.5-SU-C8P-100G	Bulk Surf 100A 2.5µm C8 Polar 100g
B-100A-2.5-SU-C8P-1KG	Bulk Surf 100A 2.5µm C8 Polar 1kg
B-100A-2.5-SU-C8T-100G	Bulk Surf 100A 2.5µm C8 TriF 100g
B-100A-2.5-SU-C8T-1KG	Bulk Surf 100A 2.5µm C8 TriF 1kg
B-100A-2.5-SU-C8EX-100G	Bulk Surf 100A 2.5µm C8 Extreme 100g
B-100A-2.5-SU-C8EX-1KG	Bulk Surf 100A 2.5µm C8 Extreme 1kg
B-100A-2.5-SU-C8PT-100G	Bulk Surf 100A 2.5µm C8 Polar TriF 100g
B-100A-2.5-SU-C8PT-1KG	Bulk Surf 100A 2.5µm C8 Polar TriF 1kg
B-100A-2.5-SU-C8AQP-100G	Bulk Surf 100A 2.5µm C8 AQ Plus 100g
B-100A-2.5-SU-C8AQP-1KG	Bulk Surf 100A 2.5µm C8 AQ Plus 1kg
B-100A-2.5-SU-C8AQ-100G	Bulk Surf 100A 2.5µm C8 AQ 100g
B-100A-2.5-SU-C8AQ-1KG	Bulk Surf 100A 2.5µm C8 AQ 1kg
<b>C18</b>	
B-100A-2.5-SU-C18-100G	Bulk Surf 100A 2.5µm C18 100g
B-100A-2.5-SU-C18-1KG	Bulk Surf 100A 2.5µm C18 1kg
B-100A-2.5-SU-C18W-100G	Bulk Surf 100A 2.5µm C18 Wax 100g
B-100A-2.5-SU-C18W-1KG	Bulk Surf 100A 2.5µm C18 Wax 1kg
B-100A-2.5-SU-C18P-100G	Bulk Surf 100A 2.5µm C18 Polar 100g
B-100A-2.5-SU-C18P-1KG	Bulk Surf 100A 2.5µm C18 Polar 1kg
B-100A-2.5-SU-C18T-100G	Bulk Surf 100A 2.5µm C18 TriF 100g
B-100A-2.5-SU-C18T-1KG	Bulk Surf 100A 2.5µm C18 TriF 1kg
B-100A-2.5-SU-C18EX-100G	Bulk Surf 100A 2.5µm C18 Extreme 100g
B-100A-2.5-SU-C18EX-1KG	Bulk Surf 100A 2.5µm C18 Extreme 1kg
B-100A-2.5-SU-C18PT-100G	Bulk Surf 100A 2.5µm C18 Polar TriF 100g
B-100A-2.5-SU-C18PT-1KG	Bulk Surf 100A 2.5µm C18 Polar TriF 1kg
B-100A-2.5-SU-C18AQP-100G	Bulk Surf 100A 2.5µm C18 AQ Plus 100g
B-100A-2.5-SU-C18AQP-1KG	Bulk Surf 100A 2.5µm C18 AQ Plus 1kg
B-100A-2.5-SU-C18AQ-100G	Bulk Surf 100A 2.5µm C18 AQ 100g
B-100A-2.5-SU-C18AQ-1KG	Bulk Surf 100A 2.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-2.5-SU-C4-100G	Bulk Surf 100A 2.5µm C4 100g
B-100A-2.5-SU-C4-1KG	Bulk Surf 100A 2.5µm C4 1kg
B-100A-2.5-SU-C4W-100G	Bulk Surf 100A 2.5µm C4 Wax 100g
B-100A-2.5-SU-C4W-1KG	Bulk Surf 100A 2.5µm C4 Wax 1kg
B-100A-2.5-SU-C4P-100G	Bulk Surf 100A 2.5µm C4 Polar 100g
B-100A-2.5-SU-C4P-1KG	Bulk Surf 100A 2.5µm C4 Polar 1kg
B-100A-2.5-SU-C4T-100G	Bulk Surf 100A 2.5µm C4 TriF 100g
B-100A-2.5-SU-C4T-1KG	Bulk Surf 100A 2.5µm C4 TriF 1kg
B-100A-2.5-SU-C4EX-100G	Bulk Surf 100A 2.5µm C4 Extreme 100g
B-100A-2.5-SU-C4EX-1KG	Bulk Surf 100A 2.5µm C4 Extreme 1kg
B-100A-2.5-SU-C4PT-100G	Bulk Surf 100A 2.5µm C4 Polar TriF 100g
B-100A-2.5-SU-C4PT-1KG	Bulk Surf 100A 2.5µm C4 Polar TriF 1kg
B-100A-2.5-SU-C4AQP-100G	Bulk Surf 100A 2.5µm C4 AQ Plus 100g
B-100A-2.5-SU-C4AQP-1KG	Bulk Surf 100A 2.5µm C4 AQ Plus 1kg
B-100A-2.5-SU-C4AQ-100G	Bulk Surf 100A 2.5µm C4 AQ 100g
B-100A-2.5-SU-C4AQ-1KG	Bulk Surf 100A 2.5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-2.5-SU-ARG-100G	Bulk Surf 100A 2.5µm ARG 100g
B-100A-2.5-SU-ARG-1KG	Bulk Surf 100A 2.5µm ARG 1kg
B-100A-2.5-SU-EP2TNE-100G	Bulk Surf 100A 2.5µm EP2TNE 100g
B-100A-2.5-SU-EP2TNE-1KG	Bulk Surf 100A 2.5µm EP2TNE 1kg
B-100A-2.5-SU-EP4TNE-100G	Bulk Surf 100A 2.5µm EP4TNE 100g
B-100A-2.5-SU-EP4TNE-1KG	Bulk Surf 100A 2.5µm EP4TNE 1kg
B-100A-2.5-SU-NH2-100G	Bulk Surf 100A 2.5µm NH2 100g
B-100A-2.5-SU-NH2-1KG	Bulk Surf 100A 2.5µm NH2 1kg
B-100A-2.5-SU-CN-100G	Bulk Surf 100A 2.5µm CN 100g
B-100A-2.5-SU-CN-1KG	Bulk Surf 100A 2.5µm CN 1kg
B-100A-2.5-SU-Phenyl-100G	Bulk Surf 100A 2.5µm Phenyl 100g
B-100A-2.5-SU-Phenyl-1KG	Bulk Surf 100A 2.5µm Phenyl 1kg
B-100A-2.5-SU-Phenyl HE-100G	Bulk Surf 100A 2.5µm Phenyl HE 100g
B-100A-2.5-SU-Phenyl HE-1KG	Bulk Surf 100A 2.5µm Phenyl HE 1kg
B-100A-2.5-SU-Diol-100G	Bulk Surf 100A 2.5µm Diol 100g
B-100A-2.5-SU-Diol-1KG	Bulk Surf 100A 2.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 2,5µm, pore 150A</b>	
B-150A-2.5-SU-Si-100G	Bulk Surf 150A 2.5µm 100g
B-150A-2.5-SU-Si-1KG	Bulk Surf 150A 2.5µm 1kg
<b>C8</b>	
B-150A-2.5-SU-C8-100G	Bulk Surf 150A 2.5µm C8 100g
B-150A-2.5-SU-C8-1KG	Bulk Surf 150A 2.5µm C8 1kg
B-150A-2.5-SU-C8W-100G	Bulk Surf 150A 2.5µm C8 Wax 100g
B-150A-2.5-SU-C8W-1KG	Bulk Surf 150A 2.5µm C8 Wax 1kg
B-150A-2.5-SU-C8P-100G	Bulk Surf 150A 2.5µm C8 Polar 100g
B-150A-2.5-SU-C8P-1KG	Bulk Surf 150A 2.5µm C8 Polar 1kg
B-150A-2.5-SU-C8T-100G	Bulk Surf 150A 2.5µm C8 TriF 100g
B-150A-2.5-SU-C8T-1KG	Bulk Surf 150A 2.5µm C8 TriF 1kg
B-150A-2.5-SU-C8EX-100G	Bulk Surf 150A 2.5µm C8 Extreme 100g
B-150A-2.5-SU-C8EX-1KG	Bulk Surf 150A 2.5µm C8 Extreme 1kg
B-150A-2.5-SU-C8PT-100G	Bulk Surf 150A 2.5µm C8 Polar TriF 100g
B-150A-2.5-SU-C8PT-1KG	Bulk Surf 150A 2.5µm C8 Polar TriF 1kg
B-150A-2.5-SU-C8AQP-100G	Bulk Surf 150A 2.5µm C8 AQ Plus 100g
B-150A-2.5-SU-C8AQP-1KG	Bulk Surf 150A 2.5µm C8 AQ Plus 1kg
B-150A-2.5-SU-C8AQ-100G	Bulk Surf 150A 2.5µm C8 AQ 100g
B-150A-2.5-SU-C8AQ-1KG	Bulk Surf 150A 2.5µm C8 AQ 1kg
<b>C18</b>	
B-150A-2.5-SU-C18-100G	Bulk Surf 150A 2.5µm C18 100g
B-150A-2.5-SU-C18-1KG	Bulk Surf 150A 2.5µm C18 1kg
B-150A-2.5-SU-C18W-100G	Bulk Surf 150A 2.5µm C18 Wax 100g
B-150A-2.5-SU-C18W-1KG	Bulk Surf 150A 2.5µm C18 Wax 1kg
B-150A-2.5-SU-C18P-100G	Bulk Surf 150A 2.5µm C18 Polar 100g
B-150A-2.5-SU-C18P-1KG	Bulk Surf 150A 2.5µm C18 Polar 1kg
B-150A-2.5-SU-C18T-100G	Bulk Surf 150A 2.5µm C18 TriF 100g
B-150A-2.5-SU-C18T-1KG	Bulk Surf 150A 2.5µm C18 TriF 1kg
B-150A-2.5-SU-C18EX-100G	Bulk Surf 150A 2.5µm C18 Extreme 100g
B-150A-2.5-SU-C18EX-1KG	Bulk Surf 150A 2.5µm C18 Extreme 1kg
B-150A-2.5-SU-C18PT-100G	Bulk Surf 150A 2.5µm C18 Polar TriF 100g
B-150A-2.5-SU-C18PT-1KG	Bulk Surf 150A 2.5µm C18 Polar TriF 1kg
B-150A-2.5-SU-C18AQP-100G	Bulk Surf 150A 2.5µm C18 AQ Plus 100g
B-150A-2.5-SU-C18AQP-1KG	Bulk Surf 150A 2.5µm C18 AQ Plus 1kg
B-150A-2.5-SU-C18AQ-100G	Bulk Surf 150A 2.5µm C18 AQ 100g
B-150A-2.5-SU-C18AQ-1KG	Bulk Surf 150A 2.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-150A-2.5-SU-C4-100G	Bulk Surf 150A 2.5µm C4 100g
B-150A-2.5-SU-C4-1KG	Bulk Surf 150A 2.5µm C4 1kg
B-150A-2.5-SU-C4W-100G	Bulk Surf 150A 2.5µm C4 Wax 100g
B-150A-2.5-SU-C4W-1KG	Bulk Surf 150A 2.5µm C4 Wax 1kg
B-150A-2.5-SU-C4P-100G	Bulk Surf 150A 2.5µm C4 Polar 100g
B-150A-2.5-SU-C4P-1KG	Bulk Surf 150A 2.5µm C4 Polar 1kg
B-150A-2.5-SU-C4T-100G	Bulk Surf 150A 2.5µm C4 TriF 100g
B-150A-2.5-SU-C4T-1KG	Bulk Surf 150A 2.5µm C4 TriF 1kg
B-150A-2.5-SU-C4EX-100G	Bulk Surf 150A 2.5µm C4 Extreme 100g
B-150A-2.5-SU-C4EX-1KG	Bulk Surf 150A 2.5µm C4 Extreme 1kg
B-150A-2.5-SU-C4PT-100G	Bulk Surf 150A 2.5µm C4 Polar TriF 100g
B-150A-2.5-SU-C4PT-1KG	Bulk Surf 150A 2.5µm C4 Polar TriF 1kg
B-150A-2.5-SU-C4AQP-100G	Bulk Surf 150A 2.5µm C4 AQ Plus 100g
B-150A-2.5-SU-C4AQP-1KG	Bulk Surf 150A 2.5µm C4 AQ Plus 1kg
B-150A-2.5-SU-C4AQ-100G	Bulk Surf 150A 2.5µm C4 AQ 100g
B-150A-2.5-SU-C4AQ-1KG	Bulk Surf 150A 2.5µm C4 AQ 1kg
<b>Autres</b>	
B-150A-2.5-SU-ARG-100G	Bulk Surf 150A 2.5µm ARG 100g
B-150A-2.5-SU-ARG-1KG	Bulk Surf 150A 2.5µm ARG 1kg
B-150A-2.5-SU-EP2TNE-100G	Bulk Surf 150A 2.5µm EP2TNE 100g
B-150A-2.5-SU-EP2TNE-1KG	Bulk Surf 150A 2.5µm EP2TNE 1kg
B-150A-2.5-SU-EP4TNE-100G	Bulk Surf 150A 2.5µm EP4TNE 100g
B-150A-2.5-SU-EP4TNE-1KG	Bulk Surf 150A 2.5µm EP4TNE 1kg
B-150A-2.5-SU-NH2-100G	Bulk Surf 150A 2.5µm NH2 100g
B-150A-2.5-SU-NH2-1KG	Bulk Surf 150A 2.5µm NH2 1kg
B-150A-2.5-SU-CN-100G	Bulk Surf 150A 2.5µm CN 100g
B-150A-2.5-SU-CN-1KG	Bulk Surf 150A 2.5µm CN 1kg
B-150A-2.5-SU-Phenyl-100G	Bulk Surf 150A 2.5µm Phenyl 100g
B-150A-2.5-SU-Phenyl-1KG	Bulk Surf 150A 2.5µm Phenyl 1kg
B-150A-2.5-SU-Phenyl HE-100G	Bulk Surf 150A 2.5µm Phenyl HE 100g
B-150A-2.5-SU-Phenyl HE-1KG	Bulk Surf 150A 2.5µm Phenyl HE 1kg
B-150A-2.5-SU-Diol-100G	Bulk Surf 150A 2.5µm Diol 100g
B-150A-2.5-SU-Diol-1KG	Bulk Surf 150A 2.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 2,5µm, pore 200A</b>	
B-200A-2.5-SU-Si-100G	Bulk Surf 200A 2.5µm 100g
B-200A-2.5-SU-Si-1KG	Bulk Surf 200A 2.5µm 1kg
<b>C8</b>	
B-200A-2.5-SU-C8-100G	Bulk Surf 200A 2.5µm C8 100g
B-200A-2.5-SU-C8-1KG	Bulk Surf 200A 2.5µm C8 1kg
B-200A-2.5-SU-C8W-100G	Bulk Surf 200A 2.5µm C8 Wax 100g
B-200A-2.5-SU-C8W-1KG	Bulk Surf 200A 2.5µm C8 Wax 1kg
B-200A-2.5-SU-C8P-100G	Bulk Surf 200A 2.5µm C8 Polar 100g
B-200A-2.5-SU-C8P-1KG	Bulk Surf 200A 2.5µm C8 Polar 1kg
B-200A-2.5-SU-C8T-100G	Bulk Surf 200A 2.5µm C8 TriF 100g
B-200A-2.5-SU-C8T-1KG	Bulk Surf 200A 2.5µm C8 TriF 1kg
B-200A-2.5-SU-C8EX-100G	Bulk Surf 200A 2.5µm C8 Extreme 100g
B-200A-2.5-SU-C8EX-1KG	Bulk Surf 200A 2.5µm C8 Extreme 1kg
B-200A-2.5-SU-C8PT-100G	Bulk Surf 200A 2.5µm C8 Polar TriF 100g
B-200A-2.5-SU-C8PT-1KG	Bulk Surf 200A 2.5µm C8 Polar TriF 1kg
B-200A-2.5-SU-C8AQP-100G	Bulk Surf 200A 2.5µm C8 AQ Plus 100g
B-200A-2.5-SU-C8AQP-1KG	Bulk Surf 200A 2.5µm C8 AQ Plus 1kg
B-200A-2.5-SU-C8AQ-100G	Bulk Surf 200A 2.5µm C8 AQ 100g
B-200A-2.5-SU-C8AQ-1KG	Bulk Surf 200A 2.5µm C8 AQ 1kg
<b>C18</b>	
B-200A-2.5-SU-C18-100G	Bulk Surf 200A 2.5µm C18 100g
B-200A-2.5-SU-C18-1KG	Bulk Surf 200A 2.5µm C18 1kg
B-200A-2.5-SU-C18W-100G	Bulk Surf 200A 2.5µm C18 Wax 100g
B-200A-2.5-SU-C18W-1KG	Bulk Surf 200A 2.5µm C18 Wax 1kg
B-200A-2.5-SU-C18P-100G	Bulk Surf 200A 2.5µm C18 Polar 100g
B-200A-2.5-SU-C18P-1KG	Bulk Surf 200A 2.5µm C18 Polar 1kg
B-200A-2.5-SU-C18T-100G	Bulk Surf 200A 2.5µm C18 TriF 100g
B-200A-2.5-SU-C18T-1KG	Bulk Surf 200A 2.5µm C18 TriF 1kg
B-200A-2.5-SU-C18EX-100G	Bulk Surf 200A 2.5µm C18 Extreme 100g
B-200A-2.5-SU-C18EX-1KG	Bulk Surf 200A 2.5µm C18 Extreme 1kg
B-200A-2.5-SU-C18PT-100G	Bulk Surf 200A 2.5µm C18 Polar TriF 100g
B-200A-2.5-SU-C18PT-1KG	Bulk Surf 200A 2.5µm C18 Polar TriF 1kg
B-200A-2.5-SU-C18AQP-100G	Bulk Surf 200A 2.5µm C18 AQ Plus 100g
B-200A-2.5-SU-C18AQP-1KG	Bulk Surf 200A 2.5µm C18 AQ Plus 1kg
B-200A-2.5-SU-C18AQ-100G	Bulk Surf 200A 2.5µm C18 AQ 100g
B-200A-2.5-SU-C18AQ-1KG	Bulk Surf 200A 2.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-200A-2.5-SU-C4-100G	Bulk Surf 200A 2.5µm C4 100g
B-200A-2.5-SU-C4-1KG	Bulk Surf 200A 2.5µm C4 1kg
B-200A-2.5-SU-C4W-100G	Bulk Surf 200A 2.5µm C4 Wax 100g
B-200A-2.5-SU-C4W-1KG	Bulk Surf 200A 2.5µm C4 Wax 1kg
B-200A-2.5-SU-C4P-100G	Bulk Surf 200A 2.5µm C4 Polar 100g
B-200A-2.5-SU-C4P-1KG	Bulk Surf 200A 2.5µm C4 Polar 1kg
B-200A-2.5-SU-C4T-100G	Bulk Surf 200A 2.5µm C4 TriF 100g
B-200A-2.5-SU-C4T-1KG	Bulk Surf 200A 2.5µm C4 TriF 1kg
B-200A-2.5-SU-C4EX-100G	Bulk Surf 200A 2.5µm C4 Extreme 100g
B-200A-2.5-SU-C4EX-1KG	Bulk Surf 200A 2.5µm C4 Extreme 1kg
B-200A-2.5-SU-C4PT-100G	Bulk Surf 200A 2.5µm C4 Polar TriF 100g
B-200A-2.5-SU-C4PT-1KG	Bulk Surf 200A 2.5µm C4 Polar TriF 1kg
B-200A-2.5-SU-C4AQP-100G	Bulk Surf 200A 2.5µm C4 AQ Plus 100g
B-200A-2.5-SU-C4AQP-1KG	Bulk Surf 200A 2.5µm C4 AQ Plus 1kg
B-200A-2.5-SU-C4AQ-100G	Bulk Surf 200A 2.5µm C4 AQ 100g
B-200A-2.5-SU-C4AQ-1KG	Bulk Surf 200A 2.5µm C4 AQ 1kg
<b>Autres</b>	
B-200A-2.5-SU-ARG-100G	Bulk Surf 200A 2.5µm ARG 100g
B-200A-2.5-SU-ARG-1KG	Bulk Surf 200A 2.5µm ARG 1kg
B-200A-2.5-SU-EP2TNE-100G	Bulk Surf 200A 2.5µm EP2TNE 100g
B-200A-2.5-SU-EP2TNE-1KG	Bulk Surf 200A 2.5µm EP2TNE 1kg
B-200A-2.5-SU-EP4TNE-100G	Bulk Surf 200A 2.5µm EP4TNE 100g
B-200A-2.5-SU-EP4TNE-1KG	Bulk Surf 200A 2.5µm EP4TNE 1kg
B-200A-2.5-SU-NH2-100G	Bulk Surf 200A 2.5µm NH2 100g
B-200A-2.5-SU-NH2-1KG	Bulk Surf 200A 2.5µm NH2 1kg
B-200A-2.5-SU-CN-100G	Bulk Surf 200A 2.5µm CN 100g
B-200A-2.5-SU-CN-1KG	Bulk Surf 200A 2.5µm CN 1kg
B-200A-2.5-SU-Phenyl-100G	Bulk Surf 200A 2.5µm Phenyl 100g
B-200A-2.5-SU-Phenyl-1KG	Bulk Surf 200A 2.5µm Phenyl 1kg
B-200A-2.5-SU-Phenyl HE-100G	Bulk Surf 200A 2.5µm Phenyl HE 100g
B-200A-2.5-SU-Phenyl HE-1KG	Bulk Surf 200A 2.5µm Phenyl HE 1kg
B-200A-2.5-SU-Diol-100G	Bulk Surf 200A 2.5µm Diol 100g
B-200A-2.5-SU-Diol-1KG	Bulk Surf 200A 2.5µm Diol 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 2,5µm, pore 200A</b>	
B-300A-2.5-SU-Si-100G	Bulk Surf 300A 2.5µm 100g
B-300A-2.5-SU-Si-1KG	Bulk Surf 300A 2.5µm 1kg
<b>C8</b>	
B-300A-2.5-SU-C8-100G	Bulk Surf 300A 2.5µm C8 100g
B-300A-2.5-SU-C8-1KG	Bulk Surf 300A 2.5µm C8 1kg
B-300A-2.5-SU-C8W-100G	Bulk Surf 300A 2.5µm C8 Wax 100g
B-300A-2.5-SU-C8W-1KG	Bulk Surf 300A 2.5µm C8 Wax 1kg
B-300A-2.5-SU-C8P-100G	Bulk Surf 300A 2.5µm C8 Polar 100g
B-300A-2.5-SU-C8P-1KG	Bulk Surf 300A 2.5µm C8 Polar 1kg
B-300A-2.5-SU-C8T-100G	Bulk Surf 300A 2.5µm C8 TriF 100g
B-300A-2.5-SU-C8T-1KG	Bulk Surf 300A 2.5µm C8 TriF 1kg
B-300A-2.5-SU-C8EX-100G	Bulk Surf 300A 2.5µm C8 Extreme 100g
B-300A-2.5-SU-C8EX-1KG	Bulk Surf 300A 2.5µm C8 Extreme 1kg
B-300A-2.5-SU-C8PT-100G	Bulk Surf 300A 2.5µm C8 Polar TriF 100g
B-300A-2.5-SU-C8PT-1KG	Bulk Surf 300A 2.5µm C8 Polar TriF 1kg
B-300A-2.5-SU-C8AQP-100G	Bulk Surf 300A 2.5µm C8 AQ Plus 100g
B-300A-2.5-SU-C8AQP-1KG	Bulk Surf 300A 2.5µm C8 AQ Plus 1kg
B-300A-2.5-SU-C8AQ-100G	Bulk Surf 300A 2.5µm C8 AQ 100g
B-300A-2.5-SU-C8AQ-1KG	Bulk Surf 300A 2.5µm C8 AQ 1kg
<b>C18</b>	
B-300A-2.5-SU-C18-100G	Bulk Surf 300A 2.5µm C18 100g
B-300A-2.5-SU-C18-1KG	Bulk Surf 300A 2.5µm C18 1kg
B-300A-2.5-SU-C18W-100G	Bulk Surf 300A 2.5µm C18 Wax 100g
B-300A-2.5-SU-C18W-1KG	Bulk Surf 300A 2.5µm C18 Wax 1kg
B-300A-2.5-SU-C18P-100G	Bulk Surf 300A 2.5µm C18 Polar 100g
B-300A-2.5-SU-C18P-1KG	Bulk Surf 300A 2.5µm C18 Polar 1kg
B-300A-2.5-SU-C18T-100G	Bulk Surf 300A 2.5µm C18 TriF 100g
B-300A-2.5-SU-C18T-1KG	Bulk Surf 300A 2.5µm C18 TriF 1kg
B-300A-2.5-SU-C18EX-100G	Bulk Surf 300A 2.5µm C18 Extreme 100g
B-300A-2.5-SU-C18EX-1KG	Bulk Surf 300A 2.5µm C18 Extreme 1kg
B-300A-2.5-SU-C18PT-100G	Bulk Surf 300A 2.5µm C18 Polar TriF 100g
B-300A-2.5-SU-C18PT-1KG	Bulk Surf 300A 2.5µm C18 Polar TriF 1kg
B-300A-2.5-SU-C18AQP-100G	Bulk Surf 300A 2.5µm C18 AQ Plus 100g
B-300A-2.5-SU-C18AQP-1KG	Bulk Surf 300A 2.5µm C18 AQ Plus 1kg
B-300A-2.5-SU-C18AQ-100G	Bulk Surf 300A 2.5µm C18 AQ 100g
B-300A-2.5-SU-C18AQ-1KG	Bulk Surf 300A 2.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-2.5-SU-C4-100G	Bulk Surf 300A 2.5µm C4 100g
B-300A-2.5-SU-C4-1KG	Bulk Surf 300A 2.5µm C4 1kg
B-300A-2.5-SU-C4W-100G	Bulk Surf 300A 2.5µm C4 Wax 100g
B-300A-2.5-SU-C4W-1KG	Bulk Surf 300A 2.5µm C4 Wax 1kg
B-300A-2.5-SU-C4P-100G	Bulk Surf 300A 2.5µm C4 Polar 100g
B-300A-2.5-SU-C4P-1KG	Bulk Surf 300A 2.5µm C4 Polar 1kg
B-300A-2.5-SU-C4T-100G	Bulk Surf 300A 2.5µm C4 TriF 100g
B-300A-2.5-SU-C4T-1KG	Bulk Surf 300A 2.5µm C4 TriF 1kg
B-300A-2.5-SU-C4EX-100G	Bulk Surf 300A 2.5µm C4 Extreme 100g
B-300A-2.5-SU-C4EX-1KG	Bulk Surf 300A 2.5µm C4 Extreme 1kg
B-300A-2.5-SU-C4PT-100G	Bulk Surf 300A 2.5µm C4 Polar TriF 100g
B-300A-2.5-SU-C4PT-1KG	Bulk Surf 300A 2.5µm C4 Polar TriF 1kg
B-300A-2.5-SU-C4AQP-100G	Bulk Surf 300A 2.5µm C4 AQ Plus 100g
B-300A-2.5-SU-C4AQP-1KG	Bulk Surf 300A 2.5µm C4 AQ Plus 1kg
B-300A-2.5-SU-C4AQ-100G	Bulk Surf 300A 2.5µm C4 AQ 100g
B-300A-2.5-SU-C4AQ-1KG	Bulk Surf 300A 2.5µm C4 AQ 1kg
<b>Autres</b>	
B-300A-2.5-SU-ARG-100G	Bulk Surf 300A 2.5µm ARG 100g
B-300A-2.5-SU-ARG-1KG	Bulk Surf 300A 2.5µm ARG 1kg
B-300A-2.5-SU-EP2TNE-100G	Bulk Surf 300A 2.5µm EP2TNE 100g
B-300A-2.5-SU-EP2TNE-1KG	Bulk Surf 300A 2.5µm EP2TNE 1kg
B-300A-2.5-SU-EP4TNE-100G	Bulk Surf 300A 2.5µm EP4TNE 100g
B-300A-2.5-SU-EP4TNE-1KG	Bulk Surf 300A 2.5µm EP4TNE 1kg
B-300A-2.5-SU-NH2-100G	Bulk Surf 300A 2.5µm NH2 100g
B-300A-2.5-SU-NH2-1KG	Bulk Surf 300A 2.5µm NH2 1kg
B-300A-2.5-SU-CN-100G	Bulk Surf 300A 2.5µm CN 100g
B-300A-2.5-SU-CN-1KG	Bulk Surf 300A 2.5µm CN 1kg
B-300A-2.5-SU-Phenyl-100G	Bulk Surf 300A 2.5µm Phenyl 100g
B-300A-2.5-SU-Phenyl-1KG	Bulk Surf 300A 2.5µm Phenyl 1kg
B-300A-2.5-SU-Phenyl HE-100G	Bulk Surf 300A 2.5µm Phenyl HE 100g
B-300A-2.5-SU-Phenyl HE-1KG	Bulk Surf 300A 2.5µm Phenyl HE 1kg
B-300A-2.5-SU-Diol-100G	Bulk Surf 300A 2.5µm Diol 100g
B-300A-2.5-SU-Diol-1KG	Bulk Surf 300A 2.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 3 µm</b>	
<b>Particule Surf 3µm, pore 100A</b>	
B-100A-3-SU-Si-100G	Bulk Surf Silice 100A 3µm 100g
B-100A-3-SU-Si-1KG	Bulk Surf Silice 100A 3µm 1kg
<b>C8</b>	
B-100A-3-SU-C8-100G	Bulk Surf 100A 3µm C8 100g
B-100A-3-SU-C8-1KG	Bulk Surf 100A 3µm C8 1kg
B-100A-3-SU-C8W-100G	Bulk Surf 100A 3µm C8 Wax 100g
B-100A-3-SU-C8W-1KG	Bulk Surf 100A 3µm C8 Wax 1kg
B-100A-3-SU-C8P-100G	Bulk Surf 100A 3µm C8 Polar 100g
B-100A-3-SU-C8P-1KG	Bulk Surf 100A 3µm C8 Polar 1kg
B-100A-3-SU-C8T-100G	Bulk Surf 100A 3µm C8 TriF 100g
B-100A-3-SU-C8T-1KG	Bulk Surf 100A 3µm C8 TriF 1kg
B-100A-3-SU-C8EX-100G	Bulk Surf 100A 3µm C8 Extreme 100g
B-100A-3-SU-C8EX-1KG	Bulk Surf 100A 3µm C8 Extreme 1kg
B-100A-3-SU-C8PT-100G	Bulk Surf 100A 3µm C8 Polar TriF 100g
B-100A-3-SU-C8PT-1KG	Bulk Surf 100A 3µm C8 Polar TriF 1kg
B-100A-3-SU-C8AQP-100G	Bulk Surf 100A 3µm C8 AQ Plus 100g
B-100A-3-SU-C8AQP-1KG	Bulk Surf 100A 3µm C8 AQ Plus 1kg
B-100A-3-SU-C8AQ-100G	Bulk Surf 100A 3µm C8 AQ 100g
B-100A-3-SU-C8AQ-1KG	Bulk Surf 100A 3µm C8 AQ 1kg
<b>C18</b>	
B-100A-3-SU-C18-100G	Bulk Surf 100A 3µm C18 100g
B-100A-3-SU-C18-1KG	Bulk Surf 100A 3µm C18 1kg
B-100A-3-SU-C18W-100G	Bulk Surf 100A 3µm C18 Wax 100g
B-100A-3-SU-C18W-1KG	Bulk Surf 100A 3µm C18 Wax 1kg
B-100A-3-SU-C18P-100G	Bulk Surf 100A 3µm C18 Polar 100g
B-100A-3-SU-C18P-1KG	Bulk Surf 100A 3µm C18 Polar 1kg
B-100A-3-SU-C18T-100G	Bulk Surf 100A 3µm C18 TriF 100g
B-100A-3-SU-C18T-1KG	Bulk Surf 100A 3µm C18 TriF 1kg
B-100A-3-SU-C18EX-100G	Bulk Surf 100A 3µm C18 Extreme 100g
B-100A-3-SU-C18EX-1KG	Bulk Surf 100A 3µm C18 Extreme 1kg
B-100A-3-SU-C18PT-100G	Bulk Surf 100A 3µm C18 Polar TriF 100g
B-100A-3-SU-C18PT-1KG	Bulk Surf 100A 3µm C18 Polar TriF 1kg
B-100A-3-SU-C18AQP-100G	Bulk Surf 100A 3µm C18 AQ Plus 100g
B-100A-3-SU-C18AQP-1KG	Bulk Surf 100A 3µm C18 AQ Plus 1kg
B-100A-3-SU-C18AQ-100G	Bulk Surf 100A 3µm C18 AQ 100g
B-100A-3-SU-C18AQ-1KG	Bulk Surf 100A 3µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-3-SU-C4-100G	Bulk Surf 100A 3µm C4 100g
B-100A-3-SU-C4-1KG	Bulk Surf 100A 3µm C4 1kg
B-100A-3-SU-C4W-100G	Bulk Surf 100A 3µm C4 Wax 100g
B-100A-3-SU-C4W-1KG	Bulk Surf 100A 3µm C4 Wax 1kg
B-100A-3-SU-C4P-100G	Bulk Surf 100A 3µm C4 Polar 100g
B-100A-3-SU-C4P-1KG	Bulk Surf 100A 3µm C4 Polar 1kg
B-100A-3-SU-C4T-100G	Bulk Surf 100A 3µm C4 TriF 100g
B-100A-3-SU-C4T-1KG	Bulk Surf 100A 3µm C4 TriF 1kg
B-100A-3-SU-C4EX-100G	Bulk Surf 100A 3µm C4 Extreme 100g
B-100A-3-SU-C4EX-1KG	Bulk Surf 100A 3µm C4 Extreme 1kg
B-100A-3-SU-C4PT-100G	Bulk Surf 100A 3µm C4 Polar TriF 100g
B-100A-3-SU-C4PT-1KG	Bulk Surf 100A 3µm C4 Polar TriF 1kg
B-100A-3-SU-C4AQP-100G	Bulk Surf 100A 3µm C4 AQ Plus 100g
B-100A-3-SU-C4AQP-1KG	Bulk Surf 100A 3µm C4 AQ Plus 1kg
B-100A-3-SU-C4AQ-100G	Bulk Surf 100A 3µm C4 AQ 100g
B-100A-3-SU-C4AQ-1KG	Bulk Surf 100A 3µm C4 AQ 1kg
<b>Autres</b>	
B-100A-3-SU-ARG-100G	Bulk Surf 100A 3µm ARG 100g
B-100A-3-SU-ARG-1KG	Bulk Surf 100A 3µm ARG 1kg
B-100A-3-SU-EP2TNE-100G	Bulk Surf 100A 3µm EP2TNE 100g
B-100A-3-SU-EP2TNE-1KG	Bulk Surf 100A 3µm EP2TNE 1kg
B-100A-3-SU-EP4TNE-100G	Bulk Surf 100A 3µm EP4TNE 100g
B-100A-3-SU-EP4TNE-1KG	Bulk Surf 100A 3µm EP4TNE 1kg
B-100A-3-SU-NH2-100G	Bulk Surf 100A 3µm NH2 100g
B-100A-3-SU-NH2-1KG	Bulk Surf 100A 3µm NH2 1kg
B-100A-3-SU-CN-100G	Bulk Surf 100A 3µm CN 100g
B-100A-3-SU-CN-1KG	Bulk Surf 100A 3µm CN 1kg
B-100A-3-SU-Phenyl-100G	Bulk Surf 100A 3µm Phenyl 100g
B-100A-3-SU-Phenyl-1KG	Bulk Surf 100A 3µm Phenyl 1kg
B-100A-3-SU-Phenyl HE-100G	Bulk Surf 100A 3µm Phenyl HE 100g
B-100A-3-SU-Phenyl HE-1KG	Bulk Surf 100A 3µm Phenyl HE 1kg
B-100A-3-SU-Diol-100G	Bulk Surf 100A 3µm Diol 100g
B-100A-3-SU-Diol-1KG	Bulk Surf 100A 3µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 3µm, pore 150A</b>	
B-150A-3-SU-Si-100G	Bulk Surf 150A 3µm 100g
B-150A-3-SU-Si-1KG	Bulk Surf 150A 3µm 1kg
<b>C8</b>	
B-150A-3-SU-C8-100G	Bulk Surf 150A 3µm C8 100g
B-150A-3-SU-C8-1KG	Bulk Surf 150A 3µm C8 1kg
B-150A-3-SU-C8W-100G	Bulk Surf 150A 3µm C8 Wax 100g
B-150A-3-SU-C8W-1KG	Bulk Surf 150A 3µm C8 Wax 1kg
B-150A-3-SU-C8P-100G	Bulk Surf 150A 3µm C8 Polar 100g
B-150A-3-SU-C8P-1KG	Bulk Surf 150A 3µm C8 Polar 1kg
B-150A-3-SU-C8T-100G	Bulk Surf 150A 3µm C8 TriF 100g
B-150A-3-SU-C8T-1KG	Bulk Surf 150A 3µm C8 TriF 1kg
B-150A-3-SU-C8EX-100G	Bulk Surf 150A 3µm C8 Extreme 100g
B-150A-3-SU-C8EX-1KG	Bulk Surf 150A 3µm C8 Extreme 1kg
B-150A-3-SU-C8PT-100G	Bulk Surf 150A 3µm C8 Polar TriF 100g
B-150A-3-SU-C8PT-1KG	Bulk Surf 150A 3µm C8 Polar TriF 1kg
B-150A-3-SU-C8AQP-100G	Bulk Surf 150A 3µm C8 AQ Plus 100g
B-150A-3-SU-C8AQP-1KG	Bulk Surf 150A 3µm C8 AQ Plus 1kg
B-150A-3-SU-C8AQ-100G	Bulk Surf 150A 3µm C8 AQ 100g
B-150A-3-SU-C8AQ-1KG	Bulk Surf 150A 3µm C8 AQ 1kg
<b>C18</b>	
B-150A-3-SU-C18-100G	Bulk Surf 150A 3µm C18 100g
B-150A-3-SU-C18-1KG	Bulk Surf 150A 3µm C18 1kg
B-150A-3-SU-C18W-100G	Bulk Surf 150A 3µm C18 Wax 100g
B-150A-3-SU-C18W-1KG	Bulk Surf 150A 3µm C18 Wax 1kg
B-150A-3-SU-C18P-100G	Bulk Surf 150A 3µm C18 Polar 100g
B-150A-3-SU-C18P-1KG	Bulk Surf 150A 3µm C18 Polar 1kg
B-150A-3-SU-C18T-100G	Bulk Surf 150A 3µm C18 TriF 100g
B-150A-3-SU-C18T-1KG	Bulk Surf 150A 3µm C18 TriF 1kg
B-150A-3-SU-C18EX-100G	Bulk Surf 150A 3µm C18 Extreme 100g
B-150A-3-SU-C18EX-1KG	Bulk Surf 150A 3µm C18 Extreme 1kg
B-150A-3-SU-C18PT-100G	Bulk Surf 150A 3µm C18 Polar TriF 100g
B-150A-3-SU-C18PT-1KG	Bulk Surf 150A 3µm C18 Polar TriF 1kg
B-150A-3-SU-C18AQP-100G	Bulk Surf 150A 3µm C18 AQ Plus 100g
B-150A-3-SU-C18AQP-1KG	Bulk Surf 150A 3µm C18 AQ Plus 1kg
B-150A-3-SU-C18AQ-100G	Bulk Surf 150A 3µm C18 AQ 100g
B-150A-3-SU-C18AQ-1KG	Bulk Surf 150A 3µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-150A-3-SU-C4-100G	Bulk Surf 150A 3µm C4 100g
B-150A-3-SU-C4-1KG	Bulk Surf 150A 3µm C4 1kg
B-150A-3-SU-C4W-100G	Bulk Surf 150A 3µm C4 Wax 100g
B-150A-3-SU-C4W-1KG	Bulk Surf 150A 3µm C4 Wax 1kg
B-150A-3-SU-C4P-100G	Bulk Surf 150A 3µm C4 Polar 100g
B-150A-3-SU-C4P-1KG	Bulk Surf 150A 3µm C4 Polar 1kg
B-150A-3-SU-C4T-100G	Bulk Surf 150A 3µm C4 TriF 100g
B-150A-3-SU-C4T-1KG	Bulk Surf 150A 3µm C4 TriF 1kg
B-150A-3-SU-C4EX-100G	Bulk Surf 150A 3µm C4 Extreme 100g
B-150A-3-SU-C4EX-1KG	Bulk Surf 150A 3µm C4 Extreme 1kg
B-150A-3-SU-C4PT-100G	Bulk Surf 150A 3µm C4 Polar TriF 100g
B-150A-3-SU-C4PT-1KG	Bulk Surf 150A 3µm C4 Polar TriF 1kg
B-150A-3-SU-C4AQP-100G	Bulk Surf 150A 3µm C4 AQ Plus 100g
B-150A-3-SU-C4AQP-1KG	Bulk Surf 150A 3µm C4 AQ Plus 1kg
B-150A-3-SU-C4AQ-100G	Bulk Surf 150A 3µm C4 AQ 100g
B-150A-3-SU-C4AQ-1KG	Bulk Surf 150A 3µm C4 AQ 1kg
<b>Autres</b>	
B-150A-3-SU-ARG-100G	Bulk Surf 150A 3µm ARG 100g
B-150A-3-SU-ARG-1KG	Bulk Surf 150A 3µm ARG 1kg
B-150A-3-SU-EP2TNE-100G	Bulk Surf 150A 3µm EP2TNE 100g
B-150A-3-SU-EP2TNE-1KG	Bulk Surf 150A 3µm EP2TNE 1kg
B-150A-3-SU-EP4TNE-100G	Bulk Surf 150A 3µm EP4TNE 100g
B-150A-3-SU-EP4TNE-1KG	Bulk Surf 150A 3µm EP4TNE 1kg
B-150A-3-SU-NH2-100G	Bulk Surf 150A 3µm NH2 100g
B-150A-3-SU-NH2-1KG	Bulk Surf 150A 3µm NH2 1kg
B-150A-3-SU-CN-100G	Bulk Surf 150A 3µm CN 100g
B-150A-3-SU-CN-1KG	Bulk Surf 150A 3µm CN 1kg
B-150A-3-SU-Phenyl-100G	Bulk Surf 150A 3µm Phenyl 100g
B-150A-3-SU-Phenyl-1KG	Bulk Surf 150A 3µm Phenyl 1kg
B-150A-3-SU-Phenyl HE-100G	Bulk Surf 150A 3µm Phenyl HE 100g
B-150A-3-SU-Phenyl HE-1KG	Bulk Surf 150A 3µm Phenyl HE 1kg
B-150A-3-SU-Diol-100G	Bulk Surf 150A 3µm Diol 100g
B-150A-3-SU-Diol-1KG	Bulk Surf 150A 3µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 3µm, pore 200A</b>	
B-200A-3-SU-Si-100G	Bulk Surf 200A 3µm 100g
B-200A-3-SU-Si-1KG	Bulk Surf 200A 3µm 1kg
<b>C8</b>	
B-200A-3-SU-C8-100G	Bulk Surf 200A 3µm C8 100g
B-200A-3-SU-C8-1KG	Bulk Surf 200A 3µm C8 1kg
B-200A-3-SU-C8W-100G	Bulk Surf 200A 3µm C8 Wax 100g
B-200A-3-SU-C8W-1KG	Bulk Surf 200A 3µm C8 Wax 1kg
B-200A-3-SU-C8P-100G	Bulk Surf 200A 3µm C8 Polar 100g
B-200A-3-SU-C8P-1KG	Bulk Surf 200A 3µm C8 Polar 1kg
B-200A-3-SU-C8T-100G	Bulk Surf 200A 3µm C8 TriF 100g
B-200A-3-SU-C8T-1KG	Bulk Surf 200A 3µm C8 TriF 1kg
B-200A-3-SU-C8EX-100G	Bulk Surf 200A 3µm C8 Extreme 100g
B-200A-3-SU-C8EX-1KG	Bulk Surf 200A 3µm C8 Extreme 1kg
B-200A-3-SU-C8PT-100G	Bulk Surf 200A 3µm C8 Polar TriF 100g
B-200A-3-SU-C8PT-1KG	Bulk Surf 200A 3µm C8 Polar TriF 1kg
B-200A-3-SU-C8AQP-100G	Bulk Surf 200A 3µm C8 AQ Plus 100g
B-200A-3-SU-C8AQP-1KG	Bulk Surf 200A 3µm C8 AQ Plus 1kg
B-200A-3-SU-C8AQ-100G	Bulk Surf 200A 3µm C8 AQ 100g
B-200A-3-SU-C8AQ-1KG	Bulk Surf 200A 3µm C8 AQ 1kg
<b>C18</b>	
B-200A-3-SU-C18-100G	Bulk Surf 200A 3µm C18 100g
B-200A-3-SU-C18-1KG	Bulk Surf 200A 3µm C18 1kg
B-200A-3-SU-C18W-100G	Bulk Surf 200A 3µm C18 Wax 100g
B-200A-3-SU-C18W-1KG	Bulk Surf 200A 3µm C18 Wax 1kg
B-200A-3-SU-C18P-100G	Bulk Surf 200A 3µm C18 Polar 100g
B-200A-3-SU-C18P-1KG	Bulk Surf 200A 3µm C18 Polar 1kg
B-200A-3-SU-C18T-100G	Bulk Surf 200A 3µm C18 TriF 100g
B-200A-3-SU-C18T-1KG	Bulk Surf 200A 3µm C18 TriF 1kg
B-200A-3-SU-C18EX-100G	Bulk Surf 200A 3µm C18 Extreme 100g
B-200A-3-SU-C18EX-1KG	Bulk Surf 200A 3µm C18 Extreme 1kg
B-200A-3-SU-C18PT-100G	Bulk Surf 200A 3µm C18 Polar TriF 100g
B-200A-3-SU-C18PT-1KG	Bulk Surf 200A 3µm C18 Polar TriF 1kg
B-200A-3-SU-C18AQP-100G	Bulk Surf 200A 3µm C18 AQ Plus 100g
B-200A-3-SU-C18AQP-1KG	Bulk Surf 200A 3µm C18 AQ Plus 1kg
B-200A-3-SU-C18AQ-100G	Bulk Surf 200A 3µm C18 AQ 100g
B-200A-3-SU-C18AQ-1KG	Bulk Surf 200A 3µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-200A-3-SU-C4-100G	Bulk Surf 200A 3µm C4 100g
B-200A-3-SU-C4-1KG	Bulk Surf 200A 3µm C4 1kg
B-200A-3-SU-C4W-100G	Bulk Surf 200A 3µm C4 Wax 100g
B-200A-3-SU-C4W-1KG	Bulk Surf 200A 3µm C4 Wax 1kg
B-200A-3-SU-C4P-100G	Bulk Surf 200A 3µm C4 Polar 100g
B-200A-3-SU-C4P-1KG	Bulk Surf 200A 3µm C4 Polar 1kg
B-200A-3-SU-C4T-100G	Bulk Surf 200A 3µm C4 TriF 100g
B-200A-3-SU-C4T-1KG	Bulk Surf 200A 3µm C4 TriF 1kg
B-200A-3-SU-C4EX-100G	Bulk Surf 200A 3µm C4 Extreme 100g
B-200A-3-SU-C4EX-1KG	Bulk Surf 200A 3µm C4 Extreme 1kg
B-200A-3-SU-C4PT-100G	Bulk Surf 200A 3µm C4 Polar TriF 100g
B-200A-3-SU-C4PT-1KG	Bulk Surf 200A 3µm C4 Polar TriF 1kg
B-200A-3-SU-C4AQP-100G	Bulk Surf 200A 3µm C4 AQ Plus 100g
B-200A-3-SU-C4AQP-1KG	Bulk Surf 200A 3µm C4 AQ Plus 1kg
B-200A-3-SU-C4AQ-100G	Bulk Surf 200A 3µm C4 AQ 100g
B-200A-3-SU-C4AQ-1KG	Bulk Surf 200A 3µm C4 AQ 1kg
<b>Autres</b>	
B-200A-3-SU-ARG-100G	Bulk Surf 200A 3µm ARG 100g
B-200A-3-SU-ARG-1KG	Bulk Surf 200A 3µm ARG 1kg
B-200A-3-SU-EP2TNE-100G	Bulk Surf 200A 3µm EP2TNE 100g
B-200A-3-SU-EP2TNE-1KG	Bulk Surf 200A 3µm EP2TNE 1kg
B-200A-3-SU-EP4TNE-100G	Bulk Surf 200A 3µm EP4TNE 100g
B-200A-3-SU-EP4TNE-1KG	Bulk Surf 200A 3µm EP4TNE 1kg
B-200A-3-SU-NH2-100G	Bulk Surf 200A 3µm NH2 100g
B-200A-3-SU-NH2-1KG	Bulk Surf 200A 3µm NH2 1kg
B-200A-3-SU-CN-100G	Bulk Surf 200A 3µm CN 100g
B-200A-3-SU-CN-1KG	Bulk Surf 200A 3µm CN 1kg
B-200A-3-SU-Phenyl-100G	Bulk Surf 200A 3µm Phenyl 100g
B-200A-3-SU-Phenyl-1KG	Bulk Surf 200A 3µm Phenyl 1kg
B-200A-3-SU-Phenyl HE-100G	Bulk Surf 200A 3µm Phenyl HE 100g
B-200A-3-SU-Phenyl HE-1KG	Bulk Surf 200A 3µm Phenyl HE 1kg
B-200A-3-SU-Diol-100G	Bulk Surf 200A 3µm Diol 100g
B-200A-3-SU-Diol-1KG	Bulk Surf 200A 3µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 3µm, pore 300A</b>	
B-300A-3-SU-Si-100G	Bulk Surf 300A 3µm 100g
B-300A-3-SU-Si-1KG	Bulk Surf 300A 3µm 1kg
<b>C8</b>	
B-300A-3-SU-C8-100G	Bulk Surf 300A 3µm C8 100g
B-300A-3-SU-C8-1KG	Bulk Surf 300A 3µm C8 1kg
B-300A-3-SU-C8W-100G	Bulk Surf 300A 3µm C8 Wax 100g
B-300A-3-SU-C8W-1KG	Bulk Surf 300A 3µm C8 Wax 1kg
B-300A-3-SU-C8P-100G	Bulk Surf 300A 3µm C8 Polar 100g
B-300A-3-SU-C8P-1KG	Bulk Surf 300A 3µm C8 Polar 1kg
B-300A-3-SU-C8T-100G	Bulk Surf 300A 3µm C8 TriF 100g
B-300A-3-SU-C8T-1KG	Bulk Surf 300A 3µm C8 TriF 1kg
B-300A-3-SU-C8EX-100G	Bulk Surf 300A 3µm C8 Extreme 100g
B-300A-3-SU-C8EX-1KG	Bulk Surf 300A 3µm C8 Extreme 1kg
B-300A-3-SU-C8PT-100G	Bulk Surf 300A 3µm C8 Polar TriF 100g
B-300A-3-SU-C8PT-1KG	Bulk Surf 300A 3µm C8 Polar TriF 1kg
B-300A-3-SU-C8AQP-100G	Bulk Surf 300A 3µm C8 AQ Plus 100g
B-300A-3-SU-C8AQP-1KG	Bulk Surf 300A 3µm C8 AQ Plus 1kg
B-300A-3-SU-C8AQ-100G	Bulk Surf 300A 3µm C8 AQ 100g
B-300A-3-SU-C8AQ-1KG	Bulk Surf 300A 3µm C8 AQ 1kg
<b>C18</b>	
B-300A-3-SU-C18-100G	Bulk Surf 300A 3µm C18 100g
B-300A-3-SU-C18-1KG	Bulk Surf 300A 3µm C18 1kg
B-300A-3-SU-C18W-100G	Bulk Surf 300A 3µm C18 Wax 100g
B-300A-3-SU-C18W-1KG	Bulk Surf 300A 3µm C18 Wax 1kg
B-300A-3-SU-C18P-100G	Bulk Surf 300A 3µm C18 Polar 100g
B-300A-3-SU-C18P-1KG	Bulk Surf 300A 3µm C18 Polar 1kg
B-300A-3-SU-C18T-100G	Bulk Surf 300A 3µm C18 TriF 100g
B-300A-3-SU-C18T-1KG	Bulk Surf 300A 3µm C18 TriF 1kg
B-300A-3-SU-C18EX-100G	Bulk Surf 300A 3µm C18 Extreme 100g
B-300A-3-SU-C18EX-1KG	Bulk Surf 300A 3µm C18 Extreme 1kg
B-300A-3-SU-C18PT-100G	Bulk Surf 300A 3µm C18 Polar TriF 100g
B-300A-3-SU-C18PT-1KG	Bulk Surf 300A 3µm C18 Polar TriF 1kg
B-300A-3-SU-C18AQP-100G	Bulk Surf 300A 3µm C18 AQ Plus 100g
B-300A-3-SU-C18AQP-1KG	Bulk Surf 300A 3µm C18 AQ Plus 1kg
B-300A-3-SU-C18AQ-100G	Bulk Surf 300A 3µm C18 AQ 100g
B-300A-3-SU-C18AQ-1KG	Bulk Surf 300A 3µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-3-SU-C4-100G	Bulk Surf 300A 3µm C4 100g
B-300A-3-SU-C4-1KG	Bulk Surf 300A 3µm C4 1kg
B-300A-3-SU-C4W-100G	Bulk Surf 300A 3µm C4 Wax 100g
B-300A-3-SU-C4W-1KG	Bulk Surf 300A 3µm C4 Wax 1kg
B-300A-3-SU-C4P-100G	Bulk Surf 300A 3µm C4 Polar 100g
B-300A-3-SU-C4P-1KG	Bulk Surf 300A 3µm C4 Polar 1kg
B-300A-3-SU-C4T-100G	Bulk Surf 300A 3µm C4 TriF 100g
B-300A-3-SU-C4T-1KG	Bulk Surf 300A 3µm C4 TriF 1kg
B-300A-3-SU-C4EX-100G	Bulk Surf 300A 3µm C4 Extreme 100g
B-300A-3-SU-C4EX-1KG	Bulk Surf 300A 3µm C4 Extreme 1kg
B-300A-3-SU-C4PT-100G	Bulk Surf 300A 3µm C4 Polar TriF 100g
B-300A-3-SU-C4PT-1KG	Bulk Surf 300A 3µm C4 Polar TriF 1kg
B-300A-3-SU-C4AQP-100G	Bulk Surf 300A 3µm C4 AQ Plus 100g
B-300A-3-SU-C4AQP-1KG	Bulk Surf 300A 3µm C4 AQ Plus 1kg
B-300A-3-SU-C4AQ-100G	Bulk Surf 300A 3µm C4 AQ 100g
B-300A-3-SU-C4AQ-1KG	Bulk Surf 300A 3µm C4 AQ 1kg
<b>Autres</b>	
B-300A-3-SU-ARG-100G	Bulk Surf 300A 3µm ARG 100g
B-300A-3-SU-ARG-1KG	Bulk Surf 300A 3µm ARG 1kg
B-300A-3-SU-EP2TNE-100G	Bulk Surf 300A 3µm EP2TNE 100g
B-300A-3-SU-EP2TNE-1KG	Bulk Surf 300A 3µm EP2TNE 1kg
B-300A-3-SU-EP4TNE-100G	Bulk Surf 300A 3µm EP4TNE 100g
B-300A-3-SU-EP4TNE-1KG	Bulk Surf 300A 3µm EP4TNE 1kg
B-300A-3-SU-NH2-100G	Bulk Surf 300A 3µm NH2 100g
B-300A-3-SU-NH2-1KG	Bulk Surf 300A 3µm NH2 1kg
B-300A-3-SU-CN-100G	Bulk Surf 300A 3µm CN 100g
B-300A-3-SU-CN-1KG	Bulk Surf 300A 3µm CN 1kg
B-300A-3-SU-Phenyl-100G	Bulk Surf 300A 3µm Phenyl 100g
B-300A-3-SU-Phenyl-1KG	Bulk Surf 300A 3µm Phenyl 1kg
B-300A-3-SU-Phenyl HE-100G	Bulk Surf 300A 3µm Phenyl HE 100g
B-300A-3-SU-Phenyl HE-1KG	Bulk Surf 300A 3µm Phenyl HE 1kg
B-300A-3-SU-Diol-100G	Bulk Surf 300A 3µm Diol 100g
B-300A-3-SU-Diol-1KG	Bulk Surf 300A 3µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 3,5 µm</b>	
<b>Particule Surf 3,5µm, pore 100A</b>	
B-100A-3.5-SU-Si-100G	Bulk Surf Silice 100A 3.5µm 100g
B-100A-3.5-SU-Si-1KG	Bulk Surf Silice 100A 3.5µm 1kg
<b>C8</b>	
B-100A-3.5-SU-C8-100G	Bulk Surf 100A 3.5µm C8 100g
B-100A-3.5-SU-C8-1KG	Bulk Surf 100A 3.5µm C8 1kg
B-100A-3.5-SU-C8W-100G	Bulk Surf 100A 3.5µm C8 Wax 100g
B-100A-3.5-SU-C8W-1KG	Bulk Surf 100A 3.5µm C8 Wax 1kg
B-100A-3.5-SU-C8P-100G	Bulk Surf 100A 3.5µm C8 Polar 100g
B-100A-3.5-SU-C8P-1KG	Bulk Surf 100A 3.5µm C8 Polar 1kg
B-100A-3.5-SU-C8T-100G	Bulk Surf 100A 3.5µm C8 TriF 100g
B-100A-3.5-SU-C8T-1KG	Bulk Surf 100A 3.5µm C8 TriF 1kg
B-100A-3.5-SU-C8EX-100G	Bulk Surf 100A 3.5µm C8 Extreme 100g
B-100A-3.5-SU-C8EX-1KG	Bulk Surf 100A 3.5µm C8 Extreme 1kg
B-100A-3.5-SU-C8PT-100G	Bulk Surf 100A 3.5µm C8 Polar TriF 100g
B-100A-3.5-SU-C8PT-1KG	Bulk Surf 100A 3.5µm C8 Polar TriF 1kg
B-100A-3.5-SU-C8AQP-100G	Bulk Surf 100A 3.5µm C8 AQ Plus 100g
B-100A-3.5-SU-C8AQP-1KG	Bulk Surf 100A 3.5µm C8 AQ Plus 1kg
B-100A-3.5-SU-C8AQ-100G	Bulk Surf 100A 3.5µm C8 AQ 100g
B-100A-3.5-SU-C8AQ-1KG	Bulk Surf 100A 3.5µm C8 AQ 1kg
<b>C18</b>	
B-100A-3.5-SU-C18-100G	Bulk Surf 100A 3.5µm C18 100g
B-100A-3.5-SU-C18-1KG	Bulk Surf 100A 3.5µm C18 1kg
B-100A-3.5-SU-C18W-100G	Bulk Surf 100A 3.5µm C18 Wax 100g
B-100A-3.5-SU-C18W-1KG	Bulk Surf 100A 3.5µm C18 Wax 1kg
B-100A-3.5-SU-C18P-100G	Bulk Surf 100A 3.5µm C18 Polar 100g
B-100A-3.5-SU-C18P-1KG	Bulk Surf 100A 3.5µm C18 Polar 1kg
B-100A-3.5-SU-C18T-100G	Bulk Surf 100A 3.5µm C18 TriF 100g
B-100A-3.5-SU-C18T-1KG	Bulk Surf 100A 3.5µm C18 TriF 1kg
B-100A-3.5-SU-C18EX-100G	Bulk Surf 100A 3.5µm C18 Extreme 100g
B-100A-3.5-SU-C18EX-1KG	Bulk Surf 100A 3.5µm C18 Extreme 1kg
B-100A-3.5-SU-C18PT-100G	Bulk Surf 100A 3.5µm C18 Polar TriF 100g
B-100A-3.5-SU-C18PT-1KG	Bulk Surf 100A 3.5µm C18 Polar TriF 1kg
B-100A-3.5-SU-C18AQP-100G	Bulk Surf 100A 3.5µm C18 AQ Plus 100g
B-100A-3.5-SU-C18AQP-1KG	Bulk Surf 100A 3.5µm C18 AQ Plus 1kg
B-100A-3.5-SU-C18AQ-100G	Bulk Surf 100A 3.5µm C18 AQ 100g
B-100A-3.5-SU-C18AQ-1KG	Bulk Surf 100A 3.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-3.5-SU-C4-100G	Bulk Surf 100A 3.5µm C4 100g
B-100A-3.5-SU-C4-1KG	Bulk Surf 100A 3.5µm C4 1kg
B-100A-3.5-SU-C4W-100G	Bulk Surf 100A 3.5µm C4 Wax 100g
B-100A-3.5-SU-C4W-1KG	Bulk Surf 100A 3.5µm C4 Wax 1kg
B-100A-3.5-SU-C4P-100G	Bulk Surf 100A 3.5µm C4 Polar 100g
B-100A-3.5-SU-C4P-1KG	Bulk Surf 100A 3.5µm C4 Polar 1kg
B-100A-3.5-SU-C4T-100G	Bulk Surf 100A 3.5µm C4 TriF 100g
B-100A-3.5-SU-C4T-1KG	Bulk Surf 100A 3.5µm C4 TriF 1kg
B-100A-3.5-SU-C4EX-100G	Bulk Surf 100A 3.5µm C4 Extreme 100g
B-100A-3.5-SU-C4EX-1KG	Bulk Surf 100A 3.5µm C4 Extreme 1kg
B-100A-3.5-SU-C4PT-100G	Bulk Surf 100A 3.5µm C4 Polar TriF 100g
B-100A-3.5-SU-C4PT-1KG	Bulk Surf 100A 3.5µm C4 Polar TriF 1kg
B-100A-3.5-SU-C4AQP-100G	Bulk Surf 100A 3.5µm C4 AQ Plus 100g
B-100A-3.5-SU-C4AQP-1KG	Bulk Surf 100A 3.5µm C4 AQ Plus 1kg
B-100A-3.5-SU-C4AQ-100G	Bulk Surf 100A 3.5µm C4 AQ 100g
B-100A-3.5-SU-C4AQ-1KG	Bulk Surf 100A 3.5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-3.5-SU-ARG-100G	Bulk Surf 100A 3.5µm ARG 100g
B-100A-3.5-SU-ARG-1KG	Bulk Surf 100A 3.5µm ARG 1kg
B-100A-3.5-SU-EP2TNE-100G	Bulk Surf 100A 3.5µm EP2TNE 100g
B-100A-3.5-SU-EP2TNE-1KG	Bulk Surf 100A 3.5µm EP2TNE 1kg
B-100A-3.5-SU-EP4TNE-100G	Bulk Surf 100A 3.5µm EP4TNE 100g
B-100A-3.5-SU-EP4TNE-1KG	Bulk Surf 100A 3.5µm EP4TNE 1kg
B-100A-3.5-SU-NH2-100G	Bulk Surf 100A 3.5µm NH2 100g
B-100A-3.5-SU-NH2-1KG	Bulk Surf 100A 3.5µm NH2 1kg
B-100A-3.5-SU-CN-100G	Bulk Surf 100A 3.5µm CN 100g
B-100A-3.5-SU-CN-1KG	Bulk Surf 100A 3.5µm CN 1kg
B-100A-3.5-SU-Phenyl-100G	Bulk Surf 100A 3.5µm Phenyl 100g
B-100A-3.5-SU-Phenyl-1KG	Bulk Surf 100A 3.5µm Phenyl 1kg
B-100A-3.5-SU-Phenyl HE-100G	Bulk Surf 100A 3.5µm Phenyl HE 100g
B-100A-3.5-SU-Phenyl HE-1KG	Bulk Surf 100A 3.5µm Phenyl HE 1kg
B-100A-3.5-SU-Diol-100G	Bulk Surf 100A 3.5µm Diol 100g
B-100A-3.5-SU-Diol-1KG	Bulk Surf 100A 3.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 3,5µm, pore 150A</b>	
B-150A-3.5-SU-Si-100G	Bulk Surf 150A 3.5µm 100g
B-150A-3.5-SU-Si-1KG	Bulk Surf 150A 3.5µm 1kg
<b>C8</b>	
B-150A-3.5-SU-C8-100G	Bulk Surf 150A 3.5µm C8 100g
B-150A-3.5-SU-C8-1KG	Bulk Surf 150A 3.5µm C8 1kg
B-150A-3.5-SU-C8W-100G	Bulk Surf 150A 3.5µm C8 Wax 100g
B-150A-3.5-SU-C8W-1KG	Bulk Surf 150A 3.5µm C8 Wax 1kg
B-150A-3.5-SU-C8P-100G	Bulk Surf 150A 3.5µm C8 Polar 100g
B-150A-3.5-SU-C8P-1KG	Bulk Surf 150A 3.5µm C8 Polar 1kg
B-150A-3.5-SU-C8T-100G	Bulk Surf 150A 3.5µm C8 TriF 100g
B-150A-3.5-SU-C8T-1KG	Bulk Surf 150A 3.5µm C8 TriF 1kg
B-150A-3.5-SU-C8EX-100G	Bulk Surf 150A 3.5µm C8 Extreme 100g
B-150A-3.5-SU-C8EX-1KG	Bulk Surf 150A 3.5µm C8 Extreme 1kg
B-150A-3.5-SU-C8PT-100G	Bulk Surf 150A 3.5µm C8 Polar TriF 100g
B-150A-3.5-SU-C8PT-1KG	Bulk Surf 150A 3.5µm C8 Polar TriF 1kg
B-150A-3.5-SU-C8AQP-100G	Bulk Surf 150A 3.5µm C8 AQ Plus 100g
B-150A-3.5-SU-C8AQP-1KG	Bulk Surf 150A 3.5µm C8 AQ Plus 1kg
B-150A-3.5-SU-C8AQ-100G	Bulk Surf 150A 3.5µm C8 AQ 100g
B-150A-3.5-SU-C8AQ-1KG	Bulk Surf 150A 3.5µm C8 AQ 1kg
<b>C18</b>	
B-150A-3.5-SU-C18-100G	Bulk Surf 150A 3.5µm C18 100g
B-150A-3.5-SU-C18-1KG	Bulk Surf 150A 3.5µm C18 1kg
B-150A-3.5-SU-C18W-100G	Bulk Surf 150A 3.5µm C18 Wax 100g
B-150A-3.5-SU-C18W-1KG	Bulk Surf 150A 3.5µm C18 Wax 1kg
B-150A-3.5-SU-C18P-100G	Bulk Surf 150A 3.5µm C18 Polar 100g
B-150A-3.5-SU-C18P-1KG	Bulk Surf 150A 3.5µm C18 Polar 1kg
B-150A-3.5-SU-C18T-100G	Bulk Surf 150A 3.5µm C18 TriF 100g
B-150A-3.5-SU-C18T-1KG	Bulk Surf 150A 3.5µm C18 TriF 1kg
B-150A-3.5-SU-C18EX-100G	Bulk Surf 150A 3.5µm C18 Extreme 100g
B-150A-3.5-SU-C18EX-1KG	Bulk Surf 150A 3.5µm C18 Extreme 1kg
B-150A-3.5-SU-C18PT-100G	Bulk Surf 150A 3.5µm C18 Polar TriF 100g
B-150A-3.5-SU-C18PT-1KG	Bulk Surf 150A 3.5µm C18 Polar TriF 1kg
B-150A-3.5-SU-C18AQP-100G	Bulk Surf 150A 3.5µm C18 AQ Plus 100g
B-150A-3.5-SU-C18AQP-1KG	Bulk Surf 150A 3.5µm C18 AQ Plus 1kg
B-150A-3.5-SU-C18AQ-100G	Bulk Surf 150A 3.5µm C18 AQ 100g
B-150A-3.5-SU-C18AQ-1KG	Bulk Surf 150A 3.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-150A-3.5-SU-C4-100G	Bulk Surf 150A 3.5µm C4 100g
B-150A-3.5-SU-C4-1KG	Bulk Surf 150A 3.5µm C4 1kg
B-150A-3.5-SU-C4W-100G	Bulk Surf 150A 3.5µm C4 Wax 100g
B-150A-3.5-SU-C4W-1KG	Bulk Surf 150A 3.5µm C4 Wax 1kg
B-150A-3.5-SU-C4P-100G	Bulk Surf 150A 3.5µm C4 Polar 100g
B-150A-3.5-SU-C4P-1KG	Bulk Surf 150A 3.5µm C4 Polar 1kg
B-150A-3.5-SU-C4T-100G	Bulk Surf 150A 3.5µm C4 TriF 100g
B-150A-3.5-SU-C4T-1KG	Bulk Surf 150A 3.5µm C4 TriF 1kg
B-150A-3.5-SU-C4EX-100G	Bulk Surf 150A 3.5µm C4 Extreme 100g
B-150A-3.5-SU-C4EX-1KG	Bulk Surf 150A 3.5µm C4 Extreme 1kg
B-150A-3.5-SU-C4PT-100G	Bulk Surf 150A 3.5µm C4 Polar TriF 100g
B-150A-3.5-SU-C4PT-1KG	Bulk Surf 150A 3.5µm C4 Polar TriF 1kg
B-150A-3.5-SU-C4AQP-100G	Bulk Surf 150A 3.5µm C4 AQ Plus 100g
B-150A-3.5-SU-C4AQP-1KG	Bulk Surf 150A 3.5µm C4 AQ Plus 1kg
B-150A-3.5-SU-C4AQ-100G	Bulk Surf 150A 3.5µm C4 AQ 100g
B-150A-3.5-SU-C4AQ-1KG	Bulk Surf 150A 3.5µm C4 AQ 1kg
<b>Autres</b>	
B-150A-3.5-SU-ARG-100G	Bulk Surf 150A 3.5µm ARG 100g
B-150A-3.5-SU-ARG-1KG	Bulk Surf 150A 3.5µm ARG 1kg
B-150A-3.5-SU-EP2TNE-100G	Bulk Surf 150A 3.5µm EP2TNE 100g
B-150A-3.5-SU-EP2TNE-1KG	Bulk Surf 150A 3.5µm EP2TNE 1kg
B-150A-3.5-SU-EP4TNE-100G	Bulk Surf 150A 3.5µm EP4TNE 100g
B-150A-3.5-SU-EP4TNE-1KG	Bulk Surf 150A 3.5µm EP4TNE 1kg
B-150A-3.5-SU-NH2-100G	Bulk Surf 150A 3.5µm NH2 100g
B-150A-3.5-SU-NH2-1KG	Bulk Surf 150A 3.5µm NH2 1kg
B-150A-3.5-SU-CN-100G	Bulk Surf 150A 3.5µm CN 100g
B-150A-3.5-SU-CN-1KG	Bulk Surf 150A 3.5µm CN 1kg
B-150A-3.5-SU-Phenyl-100G	Bulk Surf 150A 3.5µm Phenyl 100g
B-150A-3.5-SU-Phenyl-1KG	Bulk Surf 150A 3.5µm Phenyl 1kg
B-150A-3.5-SU-Phenyl HE-100G	Bulk Surf 150A 3.5µm Phenyl HE 100g
B-150A-3.5-SU-Phenyl HE-1KG	Bulk Surf 150A 3.5µm Phenyl HE 1kg
B-150A-3.5-SU-Diol-100G	Bulk Surf 150A 3.5µm Diol 100g
B-150A-3.5-SU-Diol-1KG	Bulk Surf 150A 3.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 3,5µm, pore 200A</b>	
B-200A-3.5-SU-Si-100G	Bulk Surf 200A 3.5µm 100g
B-200A-3.5-SU-Si-1KG	Bulk Surf 200A 3.5µm 1kg
<b>C8</b>	
B-200A-3.5-SU-C8-100G	Bulk Surf 200A 3.5µm C8 100g
B-200A-3.5-SU-C8-1KG	Bulk Surf 200A 3.5µm C8 1kg
B-200A-3.5-SU-C8W-100G	Bulk Surf 200A 3.5µm C8 Wax 100g
B-200A-3.5-SU-C8W-1KG	Bulk Surf 200A 3.5µm C8 Wax 1kg
B-200A-3.5-SU-C8P-100G	Bulk Surf 200A 3.5µm C8 Polar 100g
B-200A-3.5-SU-C8P-1KG	Bulk Surf 200A 3.5µm C8 Polar 1kg
B-200A-3.5-SU-C8T-100G	Bulk Surf 200A 3.5µm C8 TriF 100g
B-200A-3.5-SU-C8T-1KG	Bulk Surf 200A 3.5µm C8 TriF 1kg
B-200A-3.5-SU-C8EX-100G	Bulk Surf 200A 3.5µm C8 Extreme 100g
B-200A-3.5-SU-C8EX-1KG	Bulk Surf 200A 3.5µm C8 Extreme 1kg
B-200A-3.5-SU-C8PT-100G	Bulk Surf 200A 3.5µm C8 Polar TriF 100g
B-200A-3.5-SU-C8PT-1KG	Bulk Surf 200A 3.5µm C8 Polar TriF 1kg
B-200A-3.5-SU-C8AQP-100G	Bulk Surf 200A 3.5µm C8 AQ Plus 100g
B-200A-3.5-SU-C8AQP-1KG	Bulk Surf 200A 3.5µm C8 AQ Plus 1kg
B-200A-3.5-SU-C8AQ-100G	Bulk Surf 200A 3.5µm C8 AQ 100g
B-200A-3.5-SU-C8AQ-1KG	Bulk Surf 200A 3.5µm C8 AQ 1kg
<b>C18</b>	
B-200A-3.5-SU-C18-100G	Bulk Surf 200A 3.5µm C18 100g
B-200A-3.5-SU-C18-1KG	Bulk Surf 200A 3.5µm C18 1kg
B-200A-3.5-SU-C18W-100G	Bulk Surf 200A 3.5µm C18 Wax 100g
B-200A-3.5-SU-C18W-1KG	Bulk Surf 200A 3.5µm C18 Wax 1kg
B-200A-3.5-SU-C18P-100G	Bulk Surf 200A 3.5µm C18 Polar 100g
B-200A-3.5-SU-C18P-1KG	Bulk Surf 200A 3.5µm C18 Polar 1kg
B-200A-3.5-SU-C18T-100G	Bulk Surf 200A 3.5µm C18 TriF 100g
B-200A-3.5-SU-C18T-1KG	Bulk Surf 200A 3.5µm C18 TriF 1kg
B-200A-3.5-SU-C18EX-100G	Bulk Surf 200A 3.5µm C18 Extreme 100g
B-200A-3.5-SU-C18EX-1KG	Bulk Surf 200A 3.5µm C18 Extreme 1kg
B-200A-3.5-SU-C18PT-100G	Bulk Surf 200A 3.5µm C18 Polar TriF 100g
B-200A-3.5-SU-C18PT-1KG	Bulk Surf 200A 3.5µm C18 Polar TriF 1kg
B-200A-3.5-SU-C18AQP-100G	Bulk Surf 200A 3.5µm C18 AQ Plus 100g
B-200A-3.5-SU-C18AQP-1KG	Bulk Surf 200A 3.5µm C18 AQ Plus 1kg
B-200A-3.5-SU-C18AQ-100G	Bulk Surf 200A 3.5µm C18 AQ 100g
B-200A-3.5-SU-C18AQ-1KG	Bulk Surf 200A 3.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-200A-3.5-SU-C4-100G	Bulk Surf 200A 3.5µm C4 100g
B-200A-3.5-SU-C4-1KG	Bulk Surf 200A 3.5µm C4 1kg
B-200A-3.5-SU-C4W-100G	Bulk Surf 200A 3.5µm C4 Wax 100g
B-200A-3.5-SU-C4W-1KG	Bulk Surf 200A 3.5µm C4 Wax 1kg
B-200A-3.5-SU-C4P-100G	Bulk Surf 200A 3.5µm C4 Polar 100g
B-200A-3.5-SU-C4P-1KG	Bulk Surf 200A 3.5µm C4 Polar 1kg
B-200A-3.5-SU-C4T-100G	Bulk Surf 200A 3.5µm C4 TriF 100g
B-200A-3.5-SU-C4T-1KG	Bulk Surf 200A 3.5µm C4 TriF 1kg
B-200A-3.5-SU-C4EX-100G	Bulk Surf 200A 3.5µm C4 Extreme 100g
B-200A-3.5-SU-C4EX-1KG	Bulk Surf 200A 3.5µm C4 Extreme 1kg
B-200A-3.5-SU-C4PT-100G	Bulk Surf 200A 3.5µm C4 Polar TriF 100g
B-200A-3.5-SU-C4PT-1KG	Bulk Surf 200A 3.5µm C4 Polar TriF 1kg
B-200A-3.5-SU-C4AQP-100G	Bulk Surf 200A 3.5µm C4 AQ Plus 100g
B-200A-3.5-SU-C4AQP-1KG	Bulk Surf 200A 3.5µm C4 AQ Plus 1kg
B-200A-3.5-SU-C4AQ-100G	Bulk Surf 200A 3.5µm C4 AQ 100g
B-200A-3.5-SU-C4AQ-1KG	Bulk Surf 200A 3.5µm C4 AQ 1kg
<b>Autres</b>	
B-200A-3.5-SU-ARG-100G	Bulk Surf 200A 3.5µm ARG 100g
B-200A-3.5-SU-ARG-1KG	Bulk Surf 200A 3.5µm ARG 1kg
B-200A-3.5-SU-EP2TNE-100G	Bulk Surf 200A 3.5µm EP2TNE 100g
B-200A-3.5-SU-EP2TNE-1KG	Bulk Surf 200A 3.5µm EP2TNE 1kg
B-200A-3.5-SU-EP4TNE-100G	Bulk Surf 200A 3.5µm EP4TNE 100g
B-200A-3.5-SU-EP4TNE-1KG	Bulk Surf 200A 3.5µm EP4TNE 1kg
B-200A-3.5-SU-NH2-100G	Bulk Surf 200A 3.5µm NH2 100g
B-200A-3.5-SU-NH2-1KG	Bulk Surf 200A 3.5µm NH2 1kg
B-200A-3.5-SU-CN-100G	Bulk Surf 200A 3.5µm CN 100g
B-200A-3.5-SU-CN-1KG	Bulk Surf 200A 3.5µm CN 1kg
B-200A-3.5-SU-Phenyl-100G	Bulk Surf 200A 3.5µm Phenyl 100g
B-200A-3.5-SU-Phenyl-1KG	Bulk Surf 200A 3.5µm Phenyl 1kg
B-200A-3.5-SU-Phenyl HE-100G	Bulk Surf 200A 3.5µm Phenyl HE 100g
B-200A-3.5-SU-Phenyl HE-1KG	Bulk Surf 200A 3.5µm Phenyl HE 1kg
B-200A-3.5-SU-Diol-100G	Bulk Surf 200A 3.5µm Diol 100g
B-200A-3.5-SU-Diol-1KG	Bulk Surf 200A 3.5µm Diol 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 3,5µm, pore 300A</b>	
B-300A-3.5-SU-Si-100G	Bulk Surf 300A 3.5µm 100g
B-300A-3.5-SU-Si-1KG	Bulk Surf 300A 3.5µm 1kg
<b>C8</b>	
B-300A-3.5-SU-C8-100G	Bulk Surf 300A 3.5µm C8 100g
B-300A-3.5-SU-C8-1KG	Bulk Surf 300A 3.5µm C8 1kg
B-300A-3.5-SU-C8W-100G	Bulk Surf 300A 3.5µm C8 Wax 100g
B-300A-3.5-SU-C8W-1KG	Bulk Surf 300A 3.5µm C8 Wax 1kg
B-300A-3.5-SU-C8P-100G	Bulk Surf 300A 3.5µm C8 Polar 100g
B-300A-3.5-SU-C8P-1KG	Bulk Surf 300A 3.5µm C8 Polar 1kg
B-300A-3.5-SU-C8T-100G	Bulk Surf 300A 3.5µm C8 TriF 100g
B-300A-3.5-SU-C8T-1KG	Bulk Surf 300A 3.5µm C8 TriF 1kg
B-300A-3.5-SU-C8EX-100G	Bulk Surf 300A 3.5µm C8 Extreme 100g
B-300A-3.5-SU-C8EX-1KG	Bulk Surf 300A 3.5µm C8 Extreme 1kg
B-300A-3.5-SU-C8PT-100G	Bulk Surf 300A 3.5µm C8 Polar TriF 100g
B-300A-3.5-SU-C8PT-1KG	Bulk Surf 300A 3.5µm C8 Polar TriF 1kg
B-300A-3.5-SU-C8AQP-100G	Bulk Surf 300A 3.5µm C8 AQ Plus 100g
B-300A-3.5-SU-C8AQP-1KG	Bulk Surf 300A 3.5µm C8 AQ Plus 1kg
B-300A-3.5-SU-C8AQ-100G	Bulk Surf 300A 3.5µm C8 AQ 100g
B-300A-3.5-SU-C8AQ-1KG	Bulk Surf 300A 3.5µm C8 AQ 1kg
<b>C18</b>	
B-300A-3.5-SU-C18-100G	Bulk Surf 300A 3.5µm C18 100g
B-300A-3.5-SU-C18-1KG	Bulk Surf 300A 3.5µm C18 1kg
B-300A-3.5-SU-C18W-100G	Bulk Surf 300A 3.5µm C18 Wax 100g
B-300A-3.5-SU-C18W-1KG	Bulk Surf 300A 3.5µm C18 Wax 1kg
B-300A-3.5-SU-C18P-100G	Bulk Surf 300A 3.5µm C18 Polar 100g
B-300A-3.5-SU-C18P-1KG	Bulk Surf 300A 3.5µm C18 Polar 1kg
B-300A-3.5-SU-C18T-100G	Bulk Surf 300A 3.5µm C18 TriF 100g
B-300A-3.5-SU-C18T-1KG	Bulk Surf 300A 3.5µm C18 TriF 1kg
B-300A-3.5-SU-C18EX-100G	Bulk Surf 300A 3.5µm C18 Extreme 100g
B-300A-3.5-SU-C18EX-1KG	Bulk Surf 300A 3.5µm C18 Extreme 1kg
B-300A-3.5-SU-C18PT-100G	Bulk Surf 300A 3.5µm C18 Polar TriF 100g
B-300A-3.5-SU-C18PT-1KG	Bulk Surf 300A 3.5µm C18 Polar TriF 1kg
B-300A-3.5-SU-C18AQP-100G	Bulk Surf 300A 3.5µm C18 AQ Plus 100g
B-300A-3.5-SU-C18AQP-1KG	Bulk Surf 300A 3.5µm C18 AQ Plus 1kg
B-300A-3.5-SU-C18AQ-100G	Bulk Surf 300A 3.5µm C18 AQ 100g
B-300A-3.5-SU-C18AQ-1KG	Bulk Surf 300A 3.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-3.5-SU-C4-100G	Bulk Surf 300A 3.5µm C4 100g
B-300A-3.5-SU-C4-1KG	Bulk Surf 300A 3.5µm C4 1kg
B-300A-3.5-SU-C4W-100G	Bulk Surf 300A 3.5µm C4 Wax 100g
B-300A-3.5-SU-C4W-1KG	Bulk Surf 300A 3.5µm C4 Wax 1kg
B-300A-3.5-SU-C4P-100G	Bulk Surf 300A 3.5µm C4 Polar 100g
B-300A-3.5-SU-C4P-1KG	Bulk Surf 300A 3.5µm C4 Polar 1kg
B-300A-3.5-SU-C4T-100G	Bulk Surf 300A 3.5µm C4 TriF 100g
B-300A-3.5-SU-C4T-1KG	Bulk Surf 300A 3.5µm C4 TriF 1kg
B-300A-3.5-SU-C4EX-100G	Bulk Surf 300A 3.5µm C4 Extreme 100g
B-300A-3.5-SU-C4EX-1KG	Bulk Surf 300A 3.5µm C4 Extreme 1kg
B-300A-3.5-SU-C4PT-100G	Bulk Surf 300A 3.5µm C4 Polar TriF 100g
B-300A-3.5-SU-C4PT-1KG	Bulk Surf 300A 3.5µm C4 Polar TriF 1kg
B-300A-3.5-SU-C4AQP-100G	Bulk Surf 300A 3.5µm C4 AQ Plus 100g
B-300A-3.5-SU-C4AQP-1KG	Bulk Surf 300A 3.5µm C4 AQ Plus 1kg
B-300A-3.5-SU-C4AQ-100G	Bulk Surf 300A 3.5µm C4 AQ 100g
B-300A-3.5-SU-C4AQ-1KG	Bulk Surf 300A 3.5µm C4 AQ 1kg
<b>Autres</b>	
B-300A-3.5-SU-ARG-100G	Bulk Surf 300A 3.5µm ARG 100g
B-300A-3.5-SU-ARG-1KG	Bulk Surf 300A 3.5µm ARG 1kg
B-300A-3.5-SU-EP2TNE-100G	Bulk Surf 300A 3.5µm EP2TNE 100g
B-300A-3.5-SU-EP2TNE-1KG	Bulk Surf 300A 3.5µm EP2TNE 1kg
B-300A-3.5-SU-EP4TNE-100G	Bulk Surf 300A 3.5µm EP4TNE 100g
B-300A-3.5-SU-EP4TNE-1KG	Bulk Surf 300A 3.5µm EP4TNE 1kg
B-300A-3.5-SU-NH2-100G	Bulk Surf 300A 3.5µm NH2 100g
B-300A-3.5-SU-NH2-1KG	Bulk Surf 300A 3.5µm NH2 1kg
B-300A-3.5-SU-CN-100G	Bulk Surf 300A 3.5µm CN 100g
B-300A-3.5-SU-CN-1KG	Bulk Surf 300A 3.5µm CN 1kg
B-300A-3.5-SU-Phenyl-100G	Bulk Surf 300A 3.5µm Phenyl 100g
B-300A-3.5-SU-Phenyl-1KG	Bulk Surf 300A 3.5µm Phenyl 1kg
B-300A-3.5-SU-Phenyl HE-100G	Bulk Surf 300A 3.5µm Phenyl HE 100g
B-300A-3.5-SU-Phenyl HE-1KG	Bulk Surf 300A 3.5µm Phenyl HE 1kg
B-300A-3.5-SU-Diol-100G	Bulk Surf 300A 3.5µm Diol 100g
B-300A-3.5-SU-Diol-1KG	Bulk Surf 300A 3.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 4 µm</b>	
<b>Particule Surf 4µm, pore 100A</b>	
Code article ImChem	Désignation ImChem
B-100A-4-SU-Si-100G	Bulk Surf Silice 100A 4µm 100g
B-100A-4-SU-Si-1KG	Bulk Surf Silice 100A 4µm 1kg
<b>C8</b>	
B-100A-4-SU-C8-100G	Bulk Surf 100A 4µm C8 100g
B-100A-4-SU-C8-1KG	Bulk Surf 100A 4µm C8 1kg
B-100A-4-SU-C8W-100G	Bulk Surf 100A 4µm C8 Wax 100g
B-100A-4-SU-C8W-1KG	Bulk Surf 100A 4µm C8 Wax 1kg
B-100A-4-SU-C8P-100G	Bulk Surf 100A 4µm C8 Polar 100g
B-100A-4-SU-C8P-1KG	Bulk Surf 100A 4µm C8 Polar 1kg
B-100A-4-SU-C8T-100G	Bulk Surf 100A 4µm C8 TriF 100g
B-100A-4-SU-C8T-1KG	Bulk Surf 100A 4µm C8 TriF 1kg
B-100A-4-SU-C8EX-100G	Bulk Surf 100A 4µm C8 Extreme 100g
B-100A-4-SU-C8EX-1KG	Bulk Surf 100A 4µm C8 Extreme 1kg
B-100A-4-SU-C8PT-100G	Bulk Surf 100A 4µm C8 Polar TriF 100g
B-100A-4-SU-C8PT-1KG	Bulk Surf 100A 4µm C8 Polar TriF 1kg
B-100A-4-SU-C8AQP-100G	Bulk Surf 100A 4µm C8 AQ Plus 100g
B-100A-4-SU-C8AQP-1KG	Bulk Surf 100A 4µm C8 AQ Plus 1kg
B-100A-4-SU-C8AQ-100G	Bulk Surf 100A 4µm C8 AQ 100g
B-100A-4-SU-C8AQ-1KG	Bulk Surf 100A 4µm C8 AQ 1kg
<b>C18</b>	
B-100A-4-SU-C18-100G	Bulk Surf 100A 4µm C18 100g
B-100A-4-SU-C18-1KG	Bulk Surf 100A 4µm C18 1kg
B-100A-4-SU-C18W-100G	Bulk Surf 100A 4µm C18 Wax 100g
B-100A-4-SU-C18W-1KG	Bulk Surf 100A 4µm C18 Wax 1kg
B-100A-4-SU-C18P-100G	Bulk Surf 100A 4µm C18 Polar 100g
B-100A-4-SU-C18P-1KG	Bulk Surf 100A 4µm C18 Polar 1kg
B-100A-4-SU-C18T-100G	Bulk Surf 100A 4µm C18 TriF 100g
B-100A-4-SU-C18T-1KG	Bulk Surf 100A 4µm C18 TriF 1kg
B-100A-4-SU-C18EX-100G	Bulk Surf 100A 4µm C18 Extreme 100g
B-100A-4-SU-C18EX-1KG	Bulk Surf 100A 4µm C18 Extreme 1kg
B-100A-4-SU-C18PT-100G	Bulk Surf 100A 4µm C18 Polar TriF 100g
B-100A-4-SU-C18PT-1KG	Bulk Surf 100A 4µm C18 Polar TriF 1kg
B-100A-4-SU-C18AQP-100G	Bulk Surf 100A 4µm C18 AQ Plus 100g
B-100A-4-SU-C18AQP-1KG	Bulk Surf 100A 4µm C18 AQ Plus 1kg
B-100A-4-SU-C18AQ-100G	Bulk Surf 100A 4µm C18 AQ 100g
B-100A-4-SU-C18AQ-1KG	Bulk Surf 100A 4µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-4-SU-C4-100G	Bulk Surf 100A 4µm C4 100g
B-100A-4-SU-C4-1KG	Bulk Surf 100A 4µm C4 1kg
B-100A-4-SU-C4W-100G	Bulk Surf 100A 4µm C4 Wax 100g
B-100A-4-SU-C4W-1KG	Bulk Surf 100A 4µm C4 Wax 1kg
B-100A-4-SU-C4P-100G	Bulk Surf 100A 4µm C4 Polar 100g
B-100A-4-SU-C4P-1KG	Bulk Surf 100A 4µm C4 Polar 1kg
B-100A-4-SU-C4T-100G	Bulk Surf 100A 4µm C4 TriF 100g
B-100A-4-SU-C4T-1KG	Bulk Surf 100A 4µm C4 TriF 1kg
B-100A-4-SU-C4EX-100G	Bulk Surf 100A 4µm C4 Extreme 100g
B-100A-4-SU-C4EX-1KG	Bulk Surf 100A 4µm C4 Extreme 1kg
B-100A-4-SU-C4PT-100G	Bulk Surf 100A 4µm C4 Polar TriF 100g
B-100A-4-SU-C4PT-1KG	Bulk Surf 100A 4µm C4 Polar TriF 1kg
B-100A-4-SU-C4AQP-100G	Bulk Surf 100A 4µm C4 AQ Plus 100g
B-100A-4-SU-C4AQP-1KG	Bulk Surf 100A 4µm C4 AQ Plus 1kg
B-100A-4-SU-C4AQ-100G	Bulk Surf 100A 4µm C4 AQ 100g
B-100A-4-SU-C4AQ-1KG	Bulk Surf 100A 4µm C4 AQ 1kg
<b>Autres</b>	
B-100A-4-SU-ARG-100G	Bulk Surf 100A 4µm ARG 100g
B-100A-4-SU-ARG-1KG	Bulk Surf 100A 4µm ARG 1kg
B-100A-4-SU-EP2TNE-100G	Bulk Surf 100A 4µm EP2TNE 100g
B-100A-4-SU-EP2TNE-1KG	Bulk Surf 100A 4µm EP2TNE 1kg
B-100A-4-SU-EP4TNE-100G	Bulk Surf 100A 4µm EP4TNE 100g
B-100A-4-SU-EP4TNE-1KG	Bulk Surf 100A 4µm EP4TNE 1kg
B-100A-4-SU-NH2-100G	Bulk Surf 100A 4µm NH2 100g
B-100A-4-SU-NH2-1KG	Bulk Surf 100A 4µm NH2 1kg
B-100A-4-SU-CN-100G	Bulk Surf 100A 4µm CN 100g
B-100A-4-SU-CN-1KG	Bulk Surf 100A 4µm CN 1kg
B-100A-4-SU-Phenyl-100G	Bulk Surf 100A 4µm Phenyl 100g
B-100A-4-SU-Phenyl-1KG	Bulk Surf 100A 4µm Phenyl 1kg
B-100A-4-SU-Phenyl HE-100G	Bulk Surf 100A 4µm Phenyl HE 100g
B-100A-4-SU-Phenyl HE-1KG	Bulk Surf 100A 4µm Phenyl HE 1kg
B-100A-4-SU-Diol-100G	Bulk Surf 100A 4µm Diol 100g
B-100A-4-SU-Diol-1KG	Bulk Surf 100A 4µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 4µm, pore 150A</b>	
B-150A-4-SU-Si-100G	Bulk Surf 150A 4µm 100g
B-150A-4-SU-Si-1KG	Bulk Surf 150A 4µm 1kg
<b>C8</b>	
B-150A-4-SU-C8-100G	Bulk Surf 150A 4µm C8 100g
B-150A-4-SU-C8-1KG	Bulk Surf 150A 4µm C8 1kg
B-150A-4-SU-C8W-100G	Bulk Surf 150A 4µm C8 Wax 100g
B-150A-4-SU-C8W-1KG	Bulk Surf 150A 4µm C8 Wax 1kg
B-150A-4-SU-C8P-100G	Bulk Surf 150A 4µm C8 Polar 100g
B-150A-4-SU-C8P-1KG	Bulk Surf 150A 4µm C8 Polar 1kg
B-150A-4-SU-C8T-100G	Bulk Surf 150A 4µm C8 TriF 100g
B-150A-4-SU-C8T-1KG	Bulk Surf 150A 4µm C8 TriF 1kg
B-150A-4-SU-C8EX-100G	Bulk Surf 150A 4µm C8 Extreme 100g
B-150A-4-SU-C8EX-1KG	Bulk Surf 150A 4µm C8 Extreme 1kg
B-150A-4-SU-C8PT-100G	Bulk Surf 150A 4µm C8 Polar TriF 100g
B-150A-4-SU-C8PT-1KG	Bulk Surf 150A 4µm C8 Polar TriF 1kg
B-150A-4-SU-C8AQP-100G	Bulk Surf 150A 4µm C8 AQ Plus 100g
B-150A-4-SU-C8AQP-1KG	Bulk Surf 150A 4µm C8 AQ Plus 1kg
B-150A-4-SU-C8AQ-100G	Bulk Surf 150A 4µm C8 AQ 100g
B-150A-4-SU-C8AQ-1KG	Bulk Surf 150A 4µm C8 AQ 1kg
<b>C18</b>	
B-150A-4-SU-C18-100G	Bulk Surf 150A 4µm C18 100g
B-150A-4-SU-C18-1KG	Bulk Surf 150A 4µm C18 1kg
B-150A-4-SU-C18W-100G	Bulk Surf 150A 4µm C18 Wax 100g
B-150A-4-SU-C18W-1KG	Bulk Surf 150A 4µm C18 Wax 1kg
B-150A-4-SU-C18P-100G	Bulk Surf 150A 4µm C18 Polar 100g
B-150A-4-SU-C18P-1KG	Bulk Surf 150A 4µm C18 Polar 1kg
B-150A-4-SU-C18T-100G	Bulk Surf 150A 4µm C18 TriF 100g
B-150A-4-SU-C18T-1KG	Bulk Surf 150A 4µm C18 TriF 1kg
B-150A-4-SU-C18EX-100G	Bulk Surf 150A 4µm C18 Extreme 100g
B-150A-4-SU-C18EX-1KG	Bulk Surf 150A 4µm C18 Extreme 1kg
B-150A-4-SU-C18PT-100G	Bulk Surf 150A 4µm C18 Polar TriF 100g
B-150A-4-SU-C18PT-1KG	Bulk Surf 150A 4µm C18 Polar TriF 1kg
B-150A-4-SU-C18AQP-100G	Bulk Surf 150A 4µm C18 AQ Plus 100g
B-150A-4-SU-C18AQP-1KG	Bulk Surf 150A 4µm C18 AQ Plus 1kg
B-150A-4-SU-C18AQ-100G	Bulk Surf 150A 4µm C18 AQ 100g
B-150A-4-SU-C18AQ-1KG	Bulk Surf 150A 4µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-150A-4-SU-C4-100G	Bulk Surf 150A 4µm C4 100g
B-150A-4-SU-C4-1KG	Bulk Surf 150A 4µm C4 1kg
B-150A-4-SU-C4W-100G	Bulk Surf 150A 4µm C4 Wax 100g
B-150A-4-SU-C4W-1KG	Bulk Surf 150A 4µm C4 Wax 1kg
B-150A-4-SU-C4P-100G	Bulk Surf 150A 4µm C4 Polar 100g
B-150A-4-SU-C4P-1KG	Bulk Surf 150A 4µm C4 Polar 1kg
B-150A-4-SU-C4T-100G	Bulk Surf 150A 4µm C4 TriF 100g
B-150A-4-SU-C4T-1KG	Bulk Surf 150A 4µm C4 TriF 1kg
B-150A-4-SU-C4EX-100G	Bulk Surf 150A 4µm C4 Extreme 100g
B-150A-4-SU-C4EX-1KG	Bulk Surf 150A 4µm C4 Extreme 1kg
B-150A-4-SU-C4PT-100G	Bulk Surf 150A 4µm C4 Polar TriF 100g
B-150A-4-SU-C4PT-1KG	Bulk Surf 150A 4µm C4 Polar TriF 1kg
B-150A-4-SU-C4AQP-100G	Bulk Surf 150A 4µm C4 AQ Plus 100g
B-150A-4-SU-C4AQP-1KG	Bulk Surf 150A 4µm C4 AQ Plus 1kg
B-150A-4-SU-C4AQ-100G	Bulk Surf 150A 4µm C4 AQ 100g
B-150A-4-SU-C4AQ-1KG	Bulk Surf 150A 4µm C4 AQ 1kg
<b>Autres</b>	
B-150A-4-SU-ARG-100G	Bulk Surf 150A 4µm ARG 100g
B-150A-4-SU-ARG-1KG	Bulk Surf 150A 4µm ARG 1kg
B-150A-4-SU-EP2TNE-100G	Bulk Surf 150A 4µm EP2TNE 100g
B-150A-4-SU-EP2TNE-1KG	Bulk Surf 150A 4µm EP2TNE 1kg
B-150A-4-SU-EP4TNE-100G	Bulk Surf 150A 4µm EP4TNE 100g
B-150A-4-SU-EP4TNE-1KG	Bulk Surf 150A 4µm EP4TNE 1kg
B-150A-4-SU-NH2-100G	Bulk Surf 150A 4µm NH2 100g
B-150A-4-SU-NH2-1KG	Bulk Surf 150A 4µm NH2 1kg
B-150A-4-SU-CN-100G	Bulk Surf 150A 4µm CN 100g
B-150A-4-SU-CN-1KG	Bulk Surf 150A 4µm CN 1kg
B-150A-4-SU-Phenyl-100G	Bulk Surf 150A 4µm Phenyl 100g
B-150A-4-SU-Phenyl-1KG	Bulk Surf 150A 4µm Phenyl 1kg
B-150A-4-SU-Phenyl HE-100G	Bulk Surf 150A 4µm Phenyl HE 100g
B-150A-4-SU-Phenyl HE-1KG	Bulk Surf 150A 4µm Phenyl HE 1kg
B-150A-4-SU-Diol-100G	Bulk Surf 150A 4µm Diol 100g
B-150A-4-SU-Diol-1KG	Bulk Surf 150A 4µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 4µm, pore 200A</b>	
B-200A-4-SU-Si-100G	Bulk Surf 200A 4µm 100g
B-200A-4-SU-Si-1KG	Bulk Surf 200A 4µm 1kg
<b>C8</b>	
B-200A-4-SU-C8-100G	Bulk Surf 200A 4µm C8 100g
B-200A-4-SU-C8-1KG	Bulk Surf 200A 4µm C8 1kg
B-200A-4-SU-C8W-100G	Bulk Surf 200A 4µm C8 Wax 100g
B-200A-4-SU-C8W-1KG	Bulk Surf 200A 4µm C8 Wax 1kg
B-200A-4-SU-C8P-100G	Bulk Surf 200A 4µm C8 Polar 100g
B-200A-4-SU-C8P-1KG	Bulk Surf 200A 4µm C8 Polar 1kg
B-200A-4-SU-C8T-100G	Bulk Surf 200A 4µm C8 TriF 100g
B-200A-4-SU-C8T-1KG	Bulk Surf 200A 4µm C8 TriF 1kg
B-200A-4-SU-C8EX-100G	Bulk Surf 200A 4µm C8 Extreme 100g
B-200A-4-SU-C8EX-1KG	Bulk Surf 200A 4µm C8 Extreme 1kg
B-200A-4-SU-C8PT-100G	Bulk Surf 200A 4µm C8 Polar TriF 100g
B-200A-4-SU-C8PT-1KG	Bulk Surf 200A 4µm C8 Polar TriF 1kg
B-200A-4-SU-C8AQP-100G	Bulk Surf 200A 4µm C8 AQ Plus 100g
B-200A-4-SU-C8AQP-1KG	Bulk Surf 200A 4µm C8 AQ Plus 1kg
B-200A-4-SU-C8AQ-100G	Bulk Surf 200A 4µm C8 AQ 100g
B-200A-4-SU-C8AQ-1KG	Bulk Surf 200A 4µm C8 AQ 1kg
<b>C18</b>	
B-200A-4-SU-C18-100G	Bulk Surf 200A 4µm C18 100g
B-200A-4-SU-C18-1KG	Bulk Surf 200A 4µm C18 1kg
B-200A-4-SU-C18W-100G	Bulk Surf 200A 4µm C18 Wax 100g
B-200A-4-SU-C18W-1KG	Bulk Surf 200A 4µm C18 Wax 1kg
B-200A-4-SU-C18P-100G	Bulk Surf 200A 4µm C18 Polar 100g
B-200A-4-SU-C18P-1KG	Bulk Surf 200A 4µm C18 Polar 1kg
B-200A-4-SU-C18T-100G	Bulk Surf 200A 4µm C18 TriF 100g
B-200A-4-SU-C18T-1KG	Bulk Surf 200A 4µm C18 TriF 1kg
B-200A-4-SU-C18EX-100G	Bulk Surf 200A 4µm C18 Extreme 100g
B-200A-4-SU-C18EX-1KG	Bulk Surf 200A 4µm C18 Extreme 1kg
B-200A-4-SU-C18PT-100G	Bulk Surf 200A 4µm C18 Polar TriF 100g
B-200A-4-SU-C18PT-1KG	Bulk Surf 200A 4µm C18 Polar TriF 1kg
B-200A-4-SU-C18AQP-100G	Bulk Surf 200A 4µm C18 AQ Plus 100g
B-200A-4-SU-C18AQP-1KG	Bulk Surf 200A 4µm C18 AQ Plus 1kg
B-200A-4-SU-C18AQ-100G	Bulk Surf 200A 4µm C18 AQ 100g
B-200A-4-SU-C18AQ-1KG	Bulk Surf 200A 4µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-200A-4-SU-C4-100G	Bulk Surf 200A 4µm C4 100g
B-200A-4-SU-C4-1KG	Bulk Surf 200A 4µm C4 1kg
B-200A-4-SU-C4W-100G	Bulk Surf 200A 4µm C4 Wax 100g
B-200A-4-SU-C4W-1KG	Bulk Surf 200A 4µm C4 Wax 1kg
B-200A-4-SU-C4P-100G	Bulk Surf 200A 4µm C4 Polar 100g
B-200A-4-SU-C4P-1KG	Bulk Surf 200A 4µm C4 Polar 1kg
B-200A-4-SU-C4T-100G	Bulk Surf 200A 4µm C4 TriF 100g
B-200A-4-SU-C4T-1KG	Bulk Surf 200A 4µm C4 TriF 1kg
B-200A-4-SU-C4EX-100G	Bulk Surf 200A 4µm C4 Extreme 100g
B-200A-4-SU-C4EX-1KG	Bulk Surf 200A 4µm C4 Extreme 1kg
B-200A-4-SU-C4PT-100G	Bulk Surf 200A 4µm C4 Polar TriF 100g
B-200A-4-SU-C4PT-1KG	Bulk Surf 200A 4µm C4 Polar TriF 1kg
B-200A-4-SU-C4AQP-100G	Bulk Surf 200A 4µm C4 AQ Plus 100g
B-200A-4-SU-C4AQP-1KG	Bulk Surf 200A 4µm C4 AQ Plus 1kg
B-200A-4-SU-C4AQ-100G	Bulk Surf 200A 4µm C4 AQ 100g
B-200A-4-SU-C4AQ-1KG	Bulk Surf 200A 4µm C4 AQ 1kg
<b>Autres</b>	
B-200A-4-SU-ARG-100G	Bulk Surf 200A 4µm ARG 100g
B-200A-4-SU-ARG-1KG	Bulk Surf 200A 4µm ARG 1kg
B-200A-4-SU-EP2TNE-100G	Bulk Surf 200A 4µm EP2TNE 100g
B-200A-4-SU-EP2TNE-1KG	Bulk Surf 200A 4µm EP2TNE 1kg
B-200A-4-SU-EP4TNE-100G	Bulk Surf 200A 4µm EP4TNE 100g
B-200A-4-SU-EP4TNE-1KG	Bulk Surf 200A 4µm EP4TNE 1kg
B-200A-4-SU-NH2-100G	Bulk Surf 200A 4µm NH2 100g
B-200A-4-SU-NH2-1KG	Bulk Surf 200A 4µm NH2 1kg
B-200A-4-SU-CN-100G	Bulk Surf 200A 4µm CN 100g
B-200A-4-SU-CN-1KG	Bulk Surf 200A 4µm CN 1kg
B-200A-4-SU-Phenyl-100G	Bulk Surf 200A 4µm Phenyl 100g
B-200A-4-SU-Phenyl-1KG	Bulk Surf 200A 4µm Phenyl 1kg
B-200A-4-SU-Phenyl HE-100G	Bulk Surf 200A 4µm Phenyl HE 100g
B-200A-4-SU-Phenyl HE-1KG	Bulk Surf 200A 4µm Phenyl HE 1kg
B-200A-4-SU-Diol-100G	Bulk Surf 200A 4µm Diol 100g
B-200A-4-SU-Diol-1KG	Bulk Surf 200A 4µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 4µm, pore 300A</b>	
B-300A-4-SU-Si-100G	Bulk Surf 300A 4µm 100g
B-300A-4-SU-Si-1KG	Bulk Surf 300A 4µm 1kg
<b>C8</b>	
B-300A-4-SU-C8-100G	Bulk Surf 300A 4µm C8 100g
B-300A-4-SU-C8-1KG	Bulk Surf 300A 4µm C8 1kg
B-300A-4-SU-C8W-100G	Bulk Surf 300A 4µm C8 Wax 100g
B-300A-4-SU-C8W-1KG	Bulk Surf 300A 4µm C8 Wax 1kg
B-300A-4-SU-C8P-100G	Bulk Surf 300A 4µm C8 Polar 100g
B-300A-4-SU-C8P-1KG	Bulk Surf 300A 4µm C8 Polar 1kg
B-300A-4-SU-C8T-100G	Bulk Surf 300A 4µm C8 TriF 100g
B-300A-4-SU-C8T-1KG	Bulk Surf 300A 4µm C8 TriF 1kg
B-300A-4-SU-C8EX-100G	Bulk Surf 300A 4µm C8 Extreme 100g
B-300A-4-SU-C8EX-1KG	Bulk Surf 300A 4µm C8 Extreme 1kg
B-300A-4-SU-C8PT-100G	Bulk Surf 300A 4µm C8 Polar TriF 100g
B-300A-4-SU-C8PT-1KG	Bulk Surf 300A 4µm C8 Polar TriF 1kg
B-300A-4-SU-C8AQP-100G	Bulk Surf 300A 4µm C8 AQ Plus 100g
B-300A-4-SU-C8AQP-1KG	Bulk Surf 300A 4µm C8 AQ Plus 1kg
B-300A-4-SU-C8AQ-100G	Bulk Surf 300A 4µm C8 AQ 100g
B-300A-4-SU-C8AQ-1KG	Bulk Surf 300A 4µm C8 AQ 1kg
<b>C18</b>	
B-300A-4-SU-C18-100G	Bulk Surf 300A 4µm C18 100g
B-300A-4-SU-C18-1KG	Bulk Surf 300A 4µm C18 1kg
B-300A-4-SU-C18W-100G	Bulk Surf 300A 4µm C18 Wax 100g
B-300A-4-SU-C18W-1KG	Bulk Surf 300A 4µm C18 Wax 1kg
B-300A-4-SU-C18P-100G	Bulk Surf 300A 4µm C18 Polar 100g
B-300A-4-SU-C18P-1KG	Bulk Surf 300A 4µm C18 Polar 1kg
B-300A-4-SU-C18T-100G	Bulk Surf 300A 4µm C18 TriF 100g
B-300A-4-SU-C18T-1KG	Bulk Surf 300A 4µm C18 TriF 1kg
B-300A-4-SU-C18EX-100G	Bulk Surf 300A 4µm C18 Extreme 100g
B-300A-4-SU-C18EX-1KG	Bulk Surf 300A 4µm C18 Extreme 1kg
B-300A-4-SU-C18PT-100G	Bulk Surf 300A 4µm C18 Polar TriF 100g
B-300A-4-SU-C18PT-1KG	Bulk Surf 300A 4µm C18 Polar TriF 1kg
B-300A-4-SU-C18AQP-100G	Bulk Surf 300A 4µm C18 AQ Plus 100g
B-300A-4-SU-C18AQP-1KG	Bulk Surf 300A 4µm C18 AQ Plus 1kg
B-300A-4-SU-C18AQ-100G	Bulk Surf 300A 4µm C18 AQ 100g
B-300A-4-SU-C18AQ-1KG	Bulk Surf 300A 4µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-4-SU-C4-100G	Bulk Surf 300A 4µm C4 100g
B-300A-4-SU-C4-1KG	Bulk Surf 300A 4µm C4 1kg
B-300A-4-SU-C4W-100G	Bulk Surf 300A 4µm C4 Wax 100g
B-300A-4-SU-C4W-1KG	Bulk Surf 300A 4µm C4 Wax 1kg
B-300A-4-SU-C4P-100G	Bulk Surf 300A 4µm C4 Polar 100g
B-300A-4-SU-C4P-1KG	Bulk Surf 300A 4µm C4 Polar 1kg
B-300A-4-SU-C4T-100G	Bulk Surf 300A 4µm C4 TriF 100g
B-300A-4-SU-C4T-1KG	Bulk Surf 300A 4µm C4 TriF 1kg
B-300A-4-SU-C4EX-100G	Bulk Surf 300A 4µm C4 Extreme 100g
B-300A-4-SU-C4EX-1KG	Bulk Surf 300A 4µm C4 Extreme 1kg
B-300A-4-SU-C4PT-100G	Bulk Surf 300A 4µm C4 Polar TriF 100g
B-300A-4-SU-C4PT-1KG	Bulk Surf 300A 4µm C4 Polar TriF 1kg
B-300A-4-SU-C4AQP-100G	Bulk Surf 300A 4µm C4 AQ Plus 100g
B-300A-4-SU-C4AQP-1KG	Bulk Surf 300A 4µm C4 AQ Plus 1kg
B-300A-4-SU-C4AQ-100G	Bulk Surf 300A 4µm C4 AQ 100g
B-300A-4-SU-C4AQ-1KG	Bulk Surf 300A 4µm C4 AQ 1kg
<b>Autres</b>	
B-300A-4-SU-ARG-100G	Bulk Surf 300A 4µm ARG 100g
B-300A-4-SU-ARG-1KG	Bulk Surf 300A 4µm ARG 1kg
B-300A-4-SU-EP2TNE-100G	Bulk Surf 300A 4µm EP2TNE 100g
B-300A-4-SU-EP2TNE-1KG	Bulk Surf 300A 4µm EP2TNE 1kg
B-300A-4-SU-EP4TNE-100G	Bulk Surf 300A 4µm EP4TNE 100g
B-300A-4-SU-EP4TNE-1KG	Bulk Surf 300A 4µm EP4TNE 1kg
B-300A-4-SU-NH2-100G	Bulk Surf 300A 4µm NH2 100g
B-300A-4-SU-NH2-1KG	Bulk Surf 300A 4µm NH2 1kg
B-300A-4-SU-CN-100G	Bulk Surf 300A 4µm CN 100g
B-300A-4-SU-CN-1KG	Bulk Surf 300A 4µm CN 1kg
B-300A-4-SU-Phenyl-100G	Bulk Surf 300A 4µm Phenyl 100g
B-300A-4-SU-Phenyl-1KG	Bulk Surf 300A 4µm Phenyl 1kg
B-300A-4-SU-Phenyl HE-100G	Bulk Surf 300A 4µm Phenyl HE 100g
B-300A-4-SU-Phenyl HE-1KG	Bulk Surf 300A 4µm Phenyl HE 1kg
B-300A-4-SU-Diol-100G	Bulk Surf 300A 4µm Diol 100g
B-300A-4-SU-Diol-1KG	Bulk Surf 300A 4µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 4,5 µm</b>	
<b>Particule Surf 4,5µm, pore 100A</b>	
B-100A-4.5-SU-Si-100G	Bulk Surf Silice 100A 4.5µm 100g
B-100A-4.5-SU-Si-1KG	Bulk Surf Silice 100A 4.5µm 1kg
<b>C8</b>	
B-100A-4.5-SU-C8-100G	Bulk Surf 100A 4.5µm C8 100g
B-100A-4.5-SU-C8-1KG	Bulk Surf 100A 4.5µm C8 1kg
B-100A-4.5-SU-C8W-100G	Bulk Surf 100A 4.5µm C8 Wax 100g
B-100A-4.5-SU-C8W-1KG	Bulk Surf 100A 4.5µm C8 Wax 1kg
B-100A-4.5-SU-C8P-100G	Bulk Surf 100A 4.5µm C8 Polar 100g
B-100A-4.5-SU-C8P-1KG	Bulk Surf 100A 4.5µm C8 Polar 1kg
B-100A-4.5-SU-C8T-100G	Bulk Surf 100A 4.5µm C8 TriF 100g
B-100A-4.5-SU-C8T-1KG	Bulk Surf 100A 4.5µm C8 TriF 1kg
B-100A-4.5-SU-C8EX-100G	Bulk Surf 100A 4.5µm C8 Extreme 100g
B-100A-4.5-SU-C8EX-1KG	Bulk Surf 100A 4.5µm C8 Extreme 1kg
B-100A-4.5-SU-C8PT-100G	Bulk Surf 100A 4.5µm C8 Polar TriF 100g
B-100A-4.5-SU-C8PT-1KG	Bulk Surf 100A 4.5µm C8 Polar TriF 1kg
B-100A-4.5-SU-C8AQP-100G	Bulk Surf 100A 4.5µm C8 AQ Plus 100g
B-100A-4.5-SU-C8AQP-1KG	Bulk Surf 100A 4.5µm C8 AQ Plus 1kg
B-100A-4.5-SU-C8AQ-100G	Bulk Surf 100A 4.5µm C8 AQ 100g
B-100A-4.5-SU-C8AQ-1KG	Bulk Surf 100A 4.5µm C8 AQ 1kg
<b>C18</b>	
B-100A-4.5-SU-C18-100G	Bulk Surf 100A 4.5µm C18 100g
B-100A-4.5-SU-C18-1KG	Bulk Surf 100A 4.5µm C18 1kg
B-100A-4.5-SU-C18W-100G	Bulk Surf 100A 4.5µm C18 Wax 100g
B-100A-4.5-SU-C18W-1KG	Bulk Surf 100A 4.5µm C18 Wax 1kg
B-100A-4.5-SU-C18P-100G	Bulk Surf 100A 4.5µm C18 Polar 100g
B-100A-4.5-SU-C18P-1KG	Bulk Surf 100A 4.5µm C18 Polar 1kg
B-100A-4.5-SU-C18T-100G	Bulk Surf 100A 4.5µm C18 TriF 100g
B-100A-4.5-SU-C18T-1KG	Bulk Surf 100A 4.5µm C18 TriF 1kg
B-100A-4.5-SU-C18EX-100G	Bulk Surf 100A 4.5µm C18 Extreme 100g
B-100A-4.5-SU-C18EX-1KG	Bulk Surf 100A 4.5µm C18 Extreme 1kg
B-100A-4.5-SU-C18PT-100G	Bulk Surf 100A 4.5µm C18 Polar TriF 100g
B-100A-4.5-SU-C18PT-1KG	Bulk Surf 100A 4.5µm C18 Polar TriF 1kg
B-100A-4.5-SU-C18AQP-100G	Bulk Surf 100A 4.5µm C18 AQ Plus 100g
B-100A-4.5-SU-C18AQP-1KG	Bulk Surf 100A 4.5µm C18 AQ Plus 1kg
B-100A-4.5-SU-C18AQ-100G	Bulk Surf 100A 4.5µm C18 AQ 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-4.5-SU-C18AQ-1KG	Bulk Surf 100A 4.5µm C18 AQ 1kg
<b>C4</b>	
B-100A-4.5-SU-C4-100G	Bulk Surf 100A 4.5µm C4 100g
B-100A-4.5-SU-C4-1KG	Bulk Surf 100A 4.5µm C4 1kg
B-100A-4.5-SU-C4W-100G	Bulk Surf 100A 4.5µm C4 Wax 100g
B-100A-4.5-SU-C4W-1KG	Bulk Surf 100A 4.5µm C4 Wax 1kg
B-100A-4.5-SU-C4P-100G	Bulk Surf 100A 4.5µm C4 Polar 100g
B-100A-4.5-SU-C4P-1KG	Bulk Surf 100A 4.5µm C4 Polar 1kg
B-100A-4.5-SU-C4T-100G	Bulk Surf 100A 4.5µm C4 TriF 100g
B-100A-4.5-SU-C4T-1KG	Bulk Surf 100A 4.5µm C4 TriF 1kg
B-100A-4.5-SU-C4EX-100G	Bulk Surf 100A 4.5µm C4 Extreme 100g
B-100A-4.5-SU-C4EX-1KG	Bulk Surf 100A 4.5µm C4 Extreme 1kg
B-100A-4.5-SU-C4PT-100G	Bulk Surf 100A 4.5µm C4 Polar TriF 100g
B-100A-4.5-SU-C4PT-1KG	Bulk Surf 100A 4.5µm C4 Polar TriF 1kg
B-100A-4.5-SU-C4AQP-100G	Bulk Surf 100A 4.5µm C4 AQ Plus 100g
B-100A-4.5-SU-C4AQP-1KG	Bulk Surf 100A 4.5µm C4 AQ Plus 1kg
B-100A-4.5-SU-C4AQ-100G	Bulk Surf 100A 4.5µm C4 AQ 100g
B-100A-4.5-SU-C4AQ-1KG	Bulk Surf 100A 4.5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-4.5-SU-ARG-100G	Bulk Surf 100A 4.5µm ARG 100g
B-100A-4.5-SU-ARG-1KG	Bulk Surf 100A 4.5µm ARG 1kg
B-100A-4.5-SU-EP2TNE-100G	Bulk Surf 100A 4.5µm EP2TNE 100g
B-100A-4.5-SU-EP2TNE-1KG	Bulk Surf 100A 4.5µm EP2TNE 1kg
B-100A-4.5-SU-EP4TNE-100G	Bulk Surf 100A 4.5µm EP4TNE 100g
B-100A-4.5-SU-EP4TNE-1KG	Bulk Surf 100A 4.5µm EP4TNE 1kg
B-100A-4.5-SU-NH2-100G	Bulk Surf 100A 4.5µm NH2 100g
B-100A-4.5-SU-NH2-1KG	Bulk Surf 100A 4.5µm NH2 1kg
B-100A-4.5-SU-CN-100G	Bulk Surf 100A 4.5µm CN 100g
B-100A-4.5-SU-CN-1KG	Bulk Surf 100A 4.5µm CN 1kg
B-100A-4.5-SU-Phenyl-100G	Bulk Surf 100A 4.5µm Phenyl 100g
B-100A-4.5-SU-Phenyl-1KG	Bulk Surf 100A 4.5µm Phenyl 1kg
B-100A-4.5-SU-Phenyl HE-100G	Bulk Surf 100A 4.5µm Phenyl HE 100g
B-100A-4.5-SU-Phenyl HE-1KG	Bulk Surf 100A 4.5µm Phenyl HE 1kg
B-100A-4.5-SU-Diol-100G	Bulk Surf 100A 4.5µm Diol 100g
B-100A-4.5-SU-Diol-1KG	Bulk Surf 100A 4.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 4,5µm, pore 150A</b>	
B-150A-4.5-SU-Si-100G	Bulk Surf 150A 4.5µm 100g
B-150A-4.5-SU-Si-1KG	Bulk Surf 150A 4.5µm 1kg
<b>C8</b>	
B-150A-4.5-SU-C8-100G	Bulk Surf 150A 4.5µm C8 100g
B-150A-4.5-SU-C8-1KG	Bulk Surf 150A 4.5µm C8 1kg
B-150A-4.5-SU-C8W-100G	Bulk Surf 150A 4.5µm C8 Wax 100g
B-150A-4.5-SU-C8W-1KG	Bulk Surf 150A 4.5µm C8 Wax 1kg
B-150A-4.5-SU-C8P-100G	Bulk Surf 150A 4.5µm C8 Polar 100g
B-150A-4.5-SU-C8P-1KG	Bulk Surf 150A 4.5µm C8 Polar 1kg
B-150A-4.5-SU-C8T-100G	Bulk Surf 150A 4.5µm C8 TriF 100g
B-150A-4.5-SU-C8T-1KG	Bulk Surf 150A 4.5µm C8 TriF 1kg
B-150A-4.5-SU-C8EX-100G	Bulk Surf 150A 4.5µm C8 Extreme 100g
B-150A-4.5-SU-C8EX-1KG	Bulk Surf 150A 4.5µm C8 Extreme 1kg
B-150A-4.5-SU-C8PT-100G	Bulk Surf 150A 4.5µm C8 Polar TriF 100g
B-150A-4.5-SU-C8PT-1KG	Bulk Surf 150A 4.5µm C8 Polar TriF 1kg
B-150A-4.5-SU-C8AQP-100G	Bulk Surf 150A 4.5µm C8 AQ Plus 100g
B-150A-4.5-SU-C8AQP-1KG	Bulk Surf 150A 4.5µm C8 AQ Plus 1kg
B-150A-4.5-SU-C8AQ-100G	Bulk Surf 150A 4.5µm C8 AQ 100g
B-150A-4.5-SU-C8AQ-1KG	Bulk Surf 150A 4.5µm C8 AQ 1kg
<b>C18</b>	
B-150A-4.5-SU-C18-100G	Bulk Surf 150A 4.5µm C18 100g
B-150A-4.5-SU-C18-1KG	Bulk Surf 150A 4.5µm C18 1kg
B-150A-4.5-SU-C18W-100G	Bulk Surf 150A 4.5µm C18 Wax 100g
B-150A-4.5-SU-C18W-1KG	Bulk Surf 150A 4.5µm C18 Wax 1kg
B-150A-4.5-SU-C18P-100G	Bulk Surf 150A 4.5µm C18 Polar 100g
B-150A-4.5-SU-C18P-1KG	Bulk Surf 150A 4.5µm C18 Polar 1kg
B-150A-4.5-SU-C18T-100G	Bulk Surf 150A 4.5µm C18 TriF 100g
B-150A-4.5-SU-C18T-1KG	Bulk Surf 150A 4.5µm C18 TriF 1kg
B-150A-4.5-SU-C18EX-100G	Bulk Surf 150A 4.5µm C18 Extreme 100g
B-150A-4.5-SU-C18EX-1KG	Bulk Surf 150A 4.5µm C18 Extreme 1kg
B-150A-4.5-SU-C18PT-100G	Bulk Surf 150A 4.5µm C18 Polar TriF 100g
B-150A-4.5-SU-C18PT-1KG	Bulk Surf 150A 4.5µm C18 Polar TriF 1kg
B-150A-4.5-SU-C18AQP-100G	Bulk Surf 150A 4.5µm C18 AQ Plus 100g
B-150A-4.5-SU-C18AQP-1KG	Bulk Surf 150A 4.5µm C18 AQ Plus 1kg
B-150A-4.5-SU-C18AQ-100G	Bulk Surf 150A 4.5µm C18 AQ 100g
B-150A-4.5-SU-C18AQ-1KG	Bulk Surf 150A 4.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-150A-4.5-SU-C4-100G	Bulk Surf 150A 4.5µm C4 100g
B-150A-4.5-SU-C4-1KG	Bulk Surf 150A 4.5µm C4 1kg
B-150A-4.5-SU-C4W-100G	Bulk Surf 150A 4.5µm C4 Wax 100g
B-150A-4.5-SU-C4W-1KG	Bulk Surf 150A 4.5µm C4 Wax 1kg
B-150A-4.5-SU-C4P-100G	Bulk Surf 150A 4.5µm C4 Polar 100g
B-150A-4.5-SU-C4P-1KG	Bulk Surf 150A 4.5µm C4 Polar 1kg
B-150A-4.5-SU-C4T-100G	Bulk Surf 150A 4.5µm C4 TriF 100g
B-150A-4.5-SU-C4T-1KG	Bulk Surf 150A 4.5µm C4 TriF 1kg
B-150A-4.5-SU-C4EX-100G	Bulk Surf 150A 4.5µm C4 Extreme 100g
B-150A-4.5-SU-C4EX-1KG	Bulk Surf 150A 4.5µm C4 Extreme 1kg
B-150A-4.5-SU-C4PT-100G	Bulk Surf 150A 4.5µm C4 Polar TriF 100g
B-150A-4.5-SU-C4PT-1KG	Bulk Surf 150A 4.5µm C4 Polar TriF 1kg
B-150A-4.5-SU-C4AQP-100G	Bulk Surf 150A 4.5µm C4 AQ Plus 100g
B-150A-4.5-SU-C4AQP-1KG	Bulk Surf 150A 4.5µm C4 AQ Plus 1kg
B-150A-4.5-SU-C4AQ-100G	Bulk Surf 150A 4.5µm C4 AQ 100g
B-150A-4.5-SU-C4AQ-1KG	Bulk Surf 150A 4.5µm C4 AQ 1kg
<b>Autres</b>	
B-150A-4.5-SU-ARG-100G	Bulk Surf 150A 4.5µm ARG 100g
B-150A-4.5-SU-ARG-1KG	Bulk Surf 150A 4.5µm ARG 1kg
B-150A-4.5-SU-EP2TNE-100G	Bulk Surf 150A 4.5µm EP2TNE 100g
B-150A-4.5-SU-EP2TNE-1KG	Bulk Surf 150A 4.5µm EP2TNE 1kg
B-150A-4.5-SU-EP4TNE-100G	Bulk Surf 150A 4.5µm EP4TNE 100g
B-150A-4.5-SU-EP4TNE-1KG	Bulk Surf 150A 4.5µm EP4TNE 1kg
B-150A-4.5-SU-NH2-100G	Bulk Surf 150A 4.5µm NH2 100g
B-150A-4.5-SU-NH2-1KG	Bulk Surf 150A 4.5µm NH2 1kg
B-150A-4.5-SU-CN-100G	Bulk Surf 150A 4.5µm CN 100g
B-150A-4.5-SU-CN-1KG	Bulk Surf 150A 4.5µm CN 1kg
B-150A-4.5-SU-Phenyl-100G	Bulk Surf 150A 4.5µm Phenyl 100g
B-150A-4.5-SU-Phenyl-1KG	Bulk Surf 150A 4.5µm Phenyl 1kg
B-150A-4.5-SU-Phenyl HE-100G	Bulk Surf 150A 4.5µm Phenyl HE 100g
B-150A-4.5-SU-Phenyl HE-1KG	Bulk Surf 150A 4.5µm Phenyl HE 1kg
B-150A-4.5-SU-Diol-100G	Bulk Surf 150A 4.5µm Diol 100g
B-150A-4.5-SU-Diol-1KG	Bulk Surf 150A 4.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 4,5µm, pore 200A</b>	
B-200A-4.5-SU-Si-100G	Bulk Surf 200A 4.5µm 100g
B-200A-4.5-SU-Si-1KG	Bulk Surf 200A 4.5µm 1kg
<b>C8</b>	
B-200A-4.5-SU-C8-100G	Bulk Surf 200A 4.5µm C8 100g
B-200A-4.5-SU-C8-1KG	Bulk Surf 200A 4.5µm C8 1kg
B-200A-4.5-SU-C8W-100G	Bulk Surf 200A 4.5µm C8 Wax 100g
B-200A-4.5-SU-C8W-1KG	Bulk Surf 200A 4.5µm C8 Wax 1kg
B-200A-4.5-SU-C8P-100G	Bulk Surf 200A 4.5µm C8 Polar 100g
B-200A-4.5-SU-C8P-1KG	Bulk Surf 200A 4.5µm C8 Polar 1kg
B-200A-4.5-SU-C8T-100G	Bulk Surf 200A 4.5µm C8 TriF 100g
B-200A-4.5-SU-C8T-1KG	Bulk Surf 200A 4.5µm C8 TriF 1kg
B-200A-4.5-SU-C8EX-100G	Bulk Surf 200A 4.5µm C8 Extreme 100g
B-200A-4.5-SU-C8EX-1KG	Bulk Surf 200A 4.5µm C8 Extreme 1kg
B-200A-4.5-SU-C8PT-100G	Bulk Surf 200A 4.5µm C8 Polar TriF 100g
B-200A-4.5-SU-C8PT-1KG	Bulk Surf 200A 4.5µm C8 Polar TriF 1kg
B-200A-4.5-SU-C8AQP-100G	Bulk Surf 200A 4.5µm C8 AQ Plus 100g
B-200A-4.5-SU-C8AQP-1KG	Bulk Surf 200A 4.5µm C8 AQ Plus 1kg
B-200A-4.5-SU-C8AQ-100G	Bulk Surf 200A 4.5µm C8 AQ 100g
B-200A-4.5-SU-C8AQ-1KG	Bulk Surf 200A 4.5µm C8 AQ 1kg
<b>C18</b>	
B-200A-4.5-SU-C18-100G	Bulk Surf 200A 4.5µm C18 100g
B-200A-4.5-SU-C18-1KG	Bulk Surf 200A 4.5µm C18 1kg
B-200A-4.5-SU-C18W-100G	Bulk Surf 200A 4.5µm C18 Wax 100g
B-200A-4.5-SU-C18W-1KG	Bulk Surf 200A 4.5µm C18 Wax 1kg
B-200A-4.5-SU-C18P-100G	Bulk Surf 200A 4.5µm C18 Polar 100g
B-200A-4.5-SU-C18P-1KG	Bulk Surf 200A 4.5µm C18 Polar 1kg
B-200A-4.5-SU-C18T-100G	Bulk Surf 200A 4.5µm C18 TriF 100g
B-200A-4.5-SU-C18T-1KG	Bulk Surf 200A 4.5µm C18 TriF 1kg
B-200A-4.5-SU-C18EX-100G	Bulk Surf 200A 4.5µm C18 Extreme 100g
B-200A-4.5-SU-C18EX-1KG	Bulk Surf 200A 4.5µm C18 Extreme 1kg
B-200A-4.5-SU-C18PT-100G	Bulk Surf 200A 4.5µm C18 Polar TriF 100g
B-200A-4.5-SU-C18PT-1KG	Bulk Surf 200A 4.5µm C18 Polar TriF 1kg
B-200A-4.5-SU-C18AQP-100G	Bulk Surf 200A 4.5µm C18 AQ Plus 100g
B-200A-4.5-SU-C18AQP-1KG	Bulk Surf 200A 4.5µm C18 AQ Plus 1kg
B-200A-4.5-SU-C18AQ-100G	Bulk Surf 200A 4.5µm C18 AQ 100g
B-200A-4.5-SU-C18AQ-1KG	Bulk Surf 200A 4.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-200A-4.5-SU-C4-100G	Bulk Surf 200A 4.5µm C4 100g
B-200A-4.5-SU-C4-1KG	Bulk Surf 200A 4.5µm C4 1kg
B-200A-4.5-SU-C4W-100G	Bulk Surf 200A 4.5µm C4 Wax 100g
B-200A-4.5-SU-C4W-1KG	Bulk Surf 200A 4.5µm C4 Wax 1kg
B-200A-4.5-SU-C4P-100G	Bulk Surf 200A 4.5µm C4 Polar 100g
B-200A-4.5-SU-C4P-1KG	Bulk Surf 200A 4.5µm C4 Polar 1kg
B-200A-4.5-SU-C4T-100G	Bulk Surf 200A 4.5µm C4 TriF 100g
B-200A-4.5-SU-C4T-1KG	Bulk Surf 200A 4.5µm C4 TriF 1kg
B-200A-4.5-SU-C4EX-100G	Bulk Surf 200A 4.5µm C4 Extreme 100g
B-200A-4.5-SU-C4EX-1KG	Bulk Surf 200A 4.5µm C4 Extreme 1kg
B-200A-4.5-SU-C4PT-100G	Bulk Surf 200A 4.5µm C4 Polar TriF 100g
B-200A-4.5-SU-C4PT-1KG	Bulk Surf 200A 4.5µm C4 Polar TriF 1kg
B-200A-4.5-SU-C4AQP-100G	Bulk Surf 200A 4.5µm C4 AQ Plus 100g
B-200A-4.5-SU-C4AQP-1KG	Bulk Surf 200A 4.5µm C4 AQ Plus 1kg
B-200A-4.5-SU-C4AQ-100G	Bulk Surf 200A 4.5µm C4 AQ 100g
B-200A-4.5-SU-C4AQ-1KG	Bulk Surf 200A 4.5µm C4 AQ 1kg
<b>Autres</b>	
B-200A-4.5-SU-ARG-100G	Bulk Surf 200A 4.5µm ARG 100g
B-200A-4.5-SU-ARG-1KG	Bulk Surf 200A 4.5µm ARG 1kg
B-200A-4.5-SU-EP2TNE-100G	Bulk Surf 200A 4.5µm EP2TNE 100g
B-200A-4.5-SU-EP2TNE-1KG	Bulk Surf 200A 4.5µm EP2TNE 1kg
B-200A-4.5-SU-EP4TNE-100G	Bulk Surf 200A 4.5µm EP4TNE 100g
B-200A-4.5-SU-EP4TNE-1KG	Bulk Surf 200A 4.5µm EP4TNE 1kg
B-200A-4.5-SU-NH2-100G	Bulk Surf 200A 4.5µm NH2 100g
B-200A-4.5-SU-NH2-1KG	Bulk Surf 200A 4.5µm NH2 1kg
B-200A-4.5-SU-CN-100G	Bulk Surf 200A 4.5µm CN 100g
B-200A-4.5-SU-CN-1KG	Bulk Surf 200A 4.5µm CN 1kg
B-200A-4.5-SU-Phenyl-100G	Bulk Surf 200A 4.5µm Phenyl 100g
B-200A-4.5-SU-Phenyl-1KG	Bulk Surf 200A 4.5µm Phenyl 1kg
B-200A-4.5-SU-Phenyl HE-100G	Bulk Surf 200A 4.5µm Phenyl HE 100g
B-200A-4.5-SU-Phenyl HE-1KG	Bulk Surf 200A 4.5µm Phenyl HE 1kg
B-200A-4.5-SU-Diol-100G	Bulk Surf 200A 4.5µm Diol 100g
B-200A-4.5-SU-Diol-1KG	Bulk Surf 200A 4.5µm Diol 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 4,5µm, pore 300A</b>	
B-300A-4.5-SU-Si-100G	Bulk Surf 300A 4.5µm 100g
B-300A-4.5-SU-Si-1KG	Bulk Surf 300A 4.5µm 1kg
<b>C8</b>	
B-300A-4.5-SU-C8-100G	Bulk Surf 300A 4.5µm C8 100g
B-300A-4.5-SU-C8-1KG	Bulk Surf 300A 4.5µm C8 1kg
B-300A-4.5-SU-C8W-100G	Bulk Surf 300A 4.5µm C8 Wax 100g
B-300A-4.5-SU-C8W-1KG	Bulk Surf 300A 4.5µm C8 Wax 1kg
B-300A-4.5-SU-C8P-100G	Bulk Surf 300A 4.5µm C8 Polar 100g
B-300A-4.5-SU-C8P-1KG	Bulk Surf 300A 4.5µm C8 Polar 1kg
B-300A-4.5-SU-C8T-100G	Bulk Surf 300A 4.5µm C8 TriF 100g
B-300A-4.5-SU-C8T-1KG	Bulk Surf 300A 4.5µm C8 TriF 1kg
B-300A-4.5-SU-C8EX-100G	Bulk Surf 300A 4.5µm C8 Extreme 100g
B-300A-4.5-SU-C8EX-1KG	Bulk Surf 300A 4.5µm C8 Extreme 1kg
B-300A-4.5-SU-C8PT-100G	Bulk Surf 300A 4.5µm C8 Polar TriF 100g
B-300A-4.5-SU-C8PT-1KG	Bulk Surf 300A 4.5µm C8 Polar TriF 1kg
B-300A-4.5-SU-C8AQP-100G	Bulk Surf 300A 4.5µm C8 AQ Plus 100g
B-300A-4.5-SU-C8AQP-1KG	Bulk Surf 300A 4.5µm C8 AQ Plus 1kg
B-300A-4.5-SU-C8AQ-100G	Bulk Surf 300A 4.5µm C8 AQ 100g
B-300A-4.5-SU-C8AQ-1KG	Bulk Surf 300A 4.5µm C8 AQ 1kg
<b>C18</b>	
B-300A-4.5-SU-C18-100G	Bulk Surf 300A 4.5µm C18 100g
B-300A-4.5-SU-C18-1KG	Bulk Surf 300A 4.5µm C18 1kg
B-300A-4.5-SU-C18W-100G	Bulk Surf 300A 4.5µm C18 Wax 100g
B-300A-4.5-SU-C18W-1KG	Bulk Surf 300A 4.5µm C18 Wax 1kg
B-300A-4.5-SU-C18P-100G	Bulk Surf 300A 4.5µm C18 Polar 100g
B-300A-4.5-SU-C18P-1KG	Bulk Surf 300A 4.5µm C18 Polar 1kg
B-300A-4.5-SU-C18T-100G	Bulk Surf 300A 4.5µm C18 TriF 100g
B-300A-4.5-SU-C18T-1KG	Bulk Surf 300A 4.5µm C18 TriF 1kg
B-300A-4.5-SU-C18EX-100G	Bulk Surf 300A 4.5µm C18 Extreme 100g
B-300A-4.5-SU-C18EX-1KG	Bulk Surf 300A 4.5µm C18 Extreme 1kg
B-300A-4.5-SU-C18PT-100G	Bulk Surf 300A 4.5µm C18 Polar TriF 100g
B-300A-4.5-SU-C18PT-1KG	Bulk Surf 300A 4.5µm C18 Polar TriF 1kg
B-300A-4.5-SU-C18AQP-100G	Bulk Surf 300A 4.5µm C18 AQ Plus 100g
B-300A-4.5-SU-C18AQP-1KG	Bulk Surf 300A 4.5µm C18 AQ Plus 1kg
B-300A-4.5-SU-C18AQ-100G	Bulk Surf 300A 4.5µm C18 AQ 100g
B-300A-4.5-SU-C18AQ-1KG	Bulk Surf 300A 4.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-4.5-SU-C4-100G	Bulk Surf 300A 4.5µm C4 100g
B-300A-4.5-SU-C4-1KG	Bulk Surf 300A 4.5µm C4 1kg
B-300A-4.5-SU-C4W-100G	Bulk Surf 300A 4.5µm C4 Wax 100g
B-300A-4.5-SU-C4W-1KG	Bulk Surf 300A 4.5µm C4 Wax 1kg
B-300A-4.5-SU-C4P-100G	Bulk Surf 300A 4.5µm C4 Polar 100g
B-300A-4.5-SU-C4P-1KG	Bulk Surf 300A 4.5µm C4 Polar 1kg
B-300A-4.5-SU-C4T-100G	Bulk Surf 300A 4.5µm C4 TriF 100g
B-300A-4.5-SU-C4T-1KG	Bulk Surf 300A 4.5µm C4 TriF 1kg
B-300A-4.5-SU-C4EX-100G	Bulk Surf 300A 4.5µm C4 Extreme 100g
B-300A-4.5-SU-C4EX-1KG	Bulk Surf 300A 4.5µm C4 Extreme 1kg
B-300A-4.5-SU-C4PT-100G	Bulk Surf 300A 4.5µm C4 Polar TriF 100g
B-300A-4.5-SU-C4PT-1KG	Bulk Surf 300A 4.5µm C4 Polar TriF 1kg
B-300A-4.5-SU-C4AQP-100G	Bulk Surf 300A 4.5µm C4 AQ Plus 100g
B-300A-4.5-SU-C4AQP-1KG	Bulk Surf 300A 4.5µm C4 AQ Plus 1kg
B-300A-4.5-SU-C4AQ-100G	Bulk Surf 300A 4.5µm C4 AQ 100g
B-300A-4.5-SU-C4AQ-1KG	Bulk Surf 300A 4.5µm C4 AQ 1kg
<b>Autres</b>	
B-300A-4.5-SU-ARG-100G	Bulk Surf 300A 4.5µm ARG 100g
B-300A-4.5-SU-ARG-1KG	Bulk Surf 300A 4.5µm ARG 1kg
B-300A-4.5-SU-EP2TNE-100G	Bulk Surf 300A 4.5µm EP2TNE 100g
B-300A-4.5-SU-EP2TNE-1KG	Bulk Surf 300A 4.5µm EP2TNE 1kg
B-300A-4.5-SU-EP4TNE-100G	Bulk Surf 300A 4.5µm EP4TNE 100g
B-300A-4.5-SU-EP4TNE-1KG	Bulk Surf 300A 4.5µm EP4TNE 1kg
B-300A-4.5-SU-NH2-100G	Bulk Surf 300A 4.5µm NH2 100g
B-300A-4.5-SU-NH2-1KG	Bulk Surf 300A 4.5µm NH2 1kg
B-300A-4.5-SU-CN-100G	Bulk Surf 300A 4.5µm CN 100g
B-300A-4.5-SU-CN-1KG	Bulk Surf 300A 4.5µm CN 1kg
B-300A-4.5-SU-Phenyl-100G	Bulk Surf 300A 4.5µm Phenyl 100g
B-300A-4.5-SU-Phenyl-1KG	Bulk Surf 300A 4.5µm Phenyl 1kg
B-300A-4.5-SU-Phenyl HE-100G	Bulk Surf 300A 4.5µm Phenyl HE 100g
B-300A-4.5-SU-Phenyl HE-1KG	Bulk Surf 300A 4.5µm Phenyl HE 1kg
B-300A-4.5-SU-Diol-100G	Bulk Surf 300A 4.5µm Diol 100g
B-300A-4.5-SU-Diol-1KG	Bulk Surf 300A 4.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 5 µm</b>	
<b>Particule Surf 5µm, pore 100A</b>	
Code article ImChem	Désignation ImChem
B-100A-5-SU-Si-100G	Bulk Surf Silice 100A 5µm 100g
B-100A-5-SU-Si-1KG	Bulk Surf Silice 100A 5µm 1kg
<b>C8</b>	
B-100A-5-SU-C8-100G	Bulk Surf 100A 5µm C8 100g
B-100A-5-SU-C8-1KG	Bulk Surf 100A 5µm C8 1kg
B-100A-5-SU-C8W-100G	Bulk Surf 100A 5µm C8 Wax 100g
B-100A-5-SU-C8W-1KG	Bulk Surf 100A 5µm C8 Wax 1kg
B-100A-5-SU-C8P-100G	Bulk Surf 100A 5µm C8 Polar 100g
B-100A-5-SU-C8P-1KG	Bulk Surf 100A 5µm C8 Polar 1kg
B-100A-5-SU-C8T-100G	Bulk Surf 100A 5µm C8 TriF 100g
B-100A-5-SU-C8T-1KG	Bulk Surf 100A 5µm C8 TriF 1kg
B-100A-5-SU-C8EX-100G	Bulk Surf 100A 5µm C8 Extreme 100g
B-100A-5-SU-C8EX-1KG	Bulk Surf 100A 5µm C8 Extreme 1kg
B-100A-5-SU-C8PT-100G	Bulk Surf 100A 5µm C8 Polar TriF 100g
B-100A-5-SU-C8PT-1KG	Bulk Surf 100A 5µm C8 Polar TriF 1kg
B-100A-5-SU-C8AQP-100G	Bulk Surf 100A 5µm C8 AQ Plus 100g
B-100A-5-SU-C8AQP-1KG	Bulk Surf 100A 5µm C8 AQ Plus 1kg
B-100A-5-SU-C8AQ-100G	Bulk Surf 100A 5µm C8 AQ 100g
B-100A-5-SU-C8AQ-1KG	Bulk Surf 100A 5µm C8 AQ 1kg
<b>C18</b>	
B-100A-5-SU-C18-100G	Bulk Surf 100A 5µm C18 100g
B-100A-5-SU-C18-1KG	Bulk Surf 100A 5µm C18 1kg
B-100A-5-SU-C18W-100G	Bulk Surf 100A 5µm C18 Wax 100g
B-100A-5-SU-C18W-1KG	Bulk Surf 100A 5µm C18 Wax 1kg
B-100A-5-SU-C18P-100G	Bulk Surf 100A 5µm C18 Polar 100g
B-100A-5-SU-C18P-1KG	Bulk Surf 100A 5µm C18 Polar 1kg
B-100A-5-SU-C18T-100G	Bulk Surf 100A 5µm C18 TriF 100g
B-100A-5-SU-C18T-1KG	Bulk Surf 100A 5µm C18 TriF 1kg
B-100A-5-SU-C18EX-100G	Bulk Surf 100A 5µm C18 Extreme 100g
B-100A-5-SU-C18EX-1KG	Bulk Surf 100A 5µm C18 Extreme 1kg
B-100A-5-SU-C18PT-100G	Bulk Surf 100A 5µm C18 Polar TriF 100g
B-100A-5-SU-C18PT-1KG	Bulk Surf 100A 5µm C18 Polar TriF 1kg
B-100A-5-SU-C18AQP-100G	Bulk Surf 100A 5µm C18 AQ Plus 100g
B-100A-5-SU-C18AQP-1KG	Bulk Surf 100A 5µm C18 AQ Plus 1kg
B-100A-5-SU-C18AQ-100G	Bulk Surf 100A 5µm C18 AQ 100g
B-100A-5-SU-C18AQ-1KG	Bulk Surf 100A 5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-5-SU-C4-100G	Bulk Surf 100A 5µm C4 100g
B-100A-5-SU-C4-1KG	Bulk Surf 100A 5µm C4 1kg
B-100A-5-SU-C4W-100G	Bulk Surf 100A 5µm C4 Wax 100g
B-100A-5-SU-C4W-1KG	Bulk Surf 100A 5µm C4 Wax 1kg
B-100A-5-SU-C4P-100G	Bulk Surf 100A 5µm C4 Polar 100g
B-100A-5-SU-C4P-1KG	Bulk Surf 100A 5µm C4 Polar 1kg
B-100A-5-SU-C4T-100G	Bulk Surf 100A 5µm C4 TriF 100g
B-100A-5-SU-C4T-1KG	Bulk Surf 100A 5µm C4 TriF 1kg
B-100A-5-SU-C4EX-100G	Bulk Surf 100A 5µm C4 Extreme 100g
B-100A-5-SU-C4EX-1KG	Bulk Surf 100A 5µm C4 Extreme 1kg
B-100A-5-SU-C4PT-100G	Bulk Surf 100A 5µm C4 Polar TriF 100g
B-100A-5-SU-C4PT-1KG	Bulk Surf 100A 5µm C4 Polar TriF 1kg
B-100A-5-SU-C4AQP-100G	Bulk Surf 100A 5µm C4 AQ Plus 100g
B-100A-5-SU-C4AQP-1KG	Bulk Surf 100A 5µm C4 AQ Plus 1kg
B-100A-5-SU-C4AQ-100G	Bulk Surf 100A 5µm C4 AQ 100g
B-100A-5-SU-C4AQ-1KG	Bulk Surf 100A 5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-5-SU-COOH-100G	Bulk Surf 100A 5µm COOH 100g
B-100A-5-SU-COOH-1KG	Bulk Surf 100A 5µm COOH 1kg
B-100A-5-SU-SO3H-100G	Bulk Surf 100A 5µm SO3H 100g
B-100A-5-SU-SO3H-1KG	Bulk Surf 100A 5µm SO3H 1kg
B-100A-5-SU-DNH-100G	Bulk Surf 100A 5µm DNH 100g
B-100A-5-SU-DNH-1KG	Bulk Surf 100A 5µm DNH 1kg
B-100A-5-SU-PEI-100G	Bulk Surf 100A 5µm PEI 100g
B-100A-5-SU-PEI-1KG	Bulk Surf 100A 5µm PEI 1kg
B-100A-5-SU-ARG-100G	Bulk Surf 100A 5µm ARG 100g
B-100A-5-SU-ARG-1KG	Bulk Surf 100A 5µm ARG 1kg
B-100A-5-SU-EP2TNE-100G	Bulk Surf 100A 5µm EP2TNE 100g
B-100A-5-SU-EP2TNE-1KG	Bulk Surf 100A 5µm EP2TNE 1kg
B-100A-5-SU-EP4TNE-100G	Bulk Surf 100A 5µm EP4TNE 100g
B-100A-5-SU-EP4TNE-1KG	Bulk Surf 100A 5µm EP4TNE 1kg
B-100A-5-SU-NH2-100G	Bulk Surf 100A 5µm NH2 100g
B-100A-5-SU-NH2-1KG	Bulk Surf 100A 5µm NH2 1kg
B-100A-5-SU-CN-100G	Bulk Surf 100A 5µm CN 100g
B-100A-5-SU-CN-1KG	Bulk Surf 100A 5µm CN 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-5-SU-Phenyl-100G	Bulk Surf 100A 5µm Phenyl 100g
B-100A-5-SU-Phenyl-1KG	Bulk Surf 100A 5µm Phenyl 1kg
B-100A-5-SU-Phenyl HE-100G	Bulk Surf 100A 5µm Phenyl HE 100g
B-100A-5-SU-Phenyl HE-1KG	Bulk Surf 100A 5µm Phenyl HE 1kg
B-100A-5-SU-Diol-100G	Bulk Surf 100A 5µm Diol 100g
B-100A-5-SU-Diol-1KG	Bulk Surf 100A 5µm Diol 1kg
<b>Particule Surf 5µm, pore 150A</b>	
B-150A-5-SU-Si-100G	Bulk Surf Silice150A 5µm 100g
B-150A-5-SU-Si-1KG	Bulk Surf Silice 150A 5µm 1kg
<b>C8</b>	
B-150A-5-SU-C8-100G	Bulk Surf 150A 5µm C8 100g
B-150A-5-SU-C8-1KG	Bulk Surf 150A 5µm C8 1kg
B-150A-5-SU-C8W-100G	Bulk Surf 150A 5µm C8 Wax 100g
B-150A-5-SU-C8W-1KG	Bulk Surf 150A 5µm C8 Wax 1kg
B-150A-5-SU-C8P-100G	Bulk Surf 150A 5µm C8 Polar 100g
B-150A-5-SU-C8P-1KG	Bulk Surf 150A 5µm C8 Polar 1kg
B-150A-5-SU-C8T-100G	Bulk Surf 150A 5µm C8 TriF 100g
B-150A-5-SU-C8T-1KG	Bulk Surf 150A 5µm C8 TriF 1kg
B-150A-5-SU-C8EX-100G	Bulk Surf 150A 5µm C8 Extreme 100g
B-150A-5-SU-C8EX-1KG	Bulk Surf 150A 5µm C8 Extreme 1kg
B-150A-5-SU-C8PT-100G	Bulk Surf 150A 5µm C8 Polar TriF 100g
B-150A-5-SU-C8PT-1KG	Bulk Surf 150A 5µm C8 Polar TriF 1kg
B-150A-5-SU-C8AQP-100G	Bulk Surf 150A 5µm C8 AQ Plus 100g
B-150A-5-SU-C8AQP-1KG	Bulk Surf 150A 5µm C8 AQ Plus 1kg
B-150A-5-SU-C8AQ-100G	Bulk Surf 150A 5µm C8 AQ 100g
B-150A-5-SU-C8AQ-1KG	Bulk Surf 150A 5µm C8 AQ 1kg
<b>C18</b>	
B-150A-5-SU-C18-100G	Bulk Surf 150A 5µm C18 100g
B-150A-5-SU-C18-1KG	Bulk Surf 150A 5µm C18 1kg
B-150A-5-SU-C18W-100G	Bulk Surf 150A 5µm C18 Wax 100g
B-150A-5-SU-C18W-1KG	Bulk Surf 150A 5µm C18 Wax 1kg
B-150A-5-SU-C18P-100G	Bulk Surf 150A 5µm C18 Polar 100g
B-150A-5-SU-C18P-1KG	Bulk Surf 150A 5µm C18 Polar 1kg
B-150A-5-SU-C18T-100G	Bulk Surf 150A 5µm C18 TriF 100g
B-150A-5-SU-C18T-1KG	Bulk Surf 150A 5µm C18 TriF 1kg
B-150A-5-SU-C18EX-100G	Bulk Surf 150A 5µm C18 Extreme 100g
B-150A-5-SU-C18EX-1KG	Bulk Surf 150A 5µm C18 Extreme 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-5-SU-C18PT-100G	Bulk Surf 150A 5µm C18 Polar TriF 100g
B-150A-5-SU-C18PT-1KG	Bulk Surf 150A 5µm C18 Polar TriF 1kg
B-150A-5-SU-C18AQP-100G	Bulk Surf 150A 5µm C18 AQ Plus 100g
B-150A-5-SU-C18AQP-1KG	Bulk Surf 150A 5µm C18 AQ Plus 1kg
B-150A-5-SU-C18AQ-100G	Bulk Surf 150A 5µm C18 AQ 100g
B-150A-5-SU-C18AQ-1KG	Bulk Surf 150A 5µm C18 AQ 1kg
<b>C4</b>	
B-150A-5-SU-C4-100G	Bulk Surf 150A 5µm C4 100g
B-150A-5-SU-C4-1KG	Bulk Surf 150A 5µm C4 1kg
B-150A-5-SU-C4W-100G	Bulk Surf 150A 5µm C4 Wax 100g
B-150A-5-SU-C4W-1KG	Bulk Surf 150A 5µm C4 Wax 1kg
B-150A-5-SU-C4P-100G	Bulk Surf 150A 5µm C4 Polar 100g
B-150A-5-SU-C4P-1KG	Bulk Surf 150A 5µm C4 Polar 1kg
B-150A-5-SU-C4T-100G	Bulk Surf 150A 5µm C4 TriF 100g
B-150A-5-SU-C4T-1KG	Bulk Surf 150A 5µm C4 TriF 1kg
B-150A-5-SU-C4EX-100G	Bulk Surf 150A 5µm C4 Extreme 100g
B-150A-5-SU-C4EX-1KG	Bulk Surf 150A 5µm C4 Extreme 1kg
B-150A-5-SU-C4PT-100G	Bulk Surf 150A 5µm C4 Polar TriF 100g
B-150A-5-SU-C4PT-1KG	Bulk Surf 150A 5µm C4 Polar TriF 1kg
B-150A-5-SU-C4AQP-100G	Bulk Surf 150A 5µm C4 AQ Plus 100g
B-150A-5-SU-C4AQP-1KG	Bulk Surf 150A 5µm C4 AQ Plus 1kg
B-150A-5-SU-C4AQ-100G	Bulk Surf 150A 5µm C4 AQ 100g
B-150A-5-SU-C4AQ-1KG	Bulk Surf 150A 5µm C4 AQ 1kg
<b>Autres</b>	
B-150A-5-SU-COOH-100G	Bulk Surf 150A 5µm COOH 100g
B-150A-5-SU-COOH-1KG	Bulk Surf 150A 5µm COOH 1kg
B-150A-5-SU-SO3H-100G	Bulk Surf 150A 5µm SO3H 100g
B-150A-5-SU-SO3H-1KG	Bulk Surf 150A 5µm SO3H 1kg
B-150A-5-SU-DNH-100G	Bulk Surf 150A 5µm DNH 100g
B-150A-5-SU-DNH-1KG	Bulk Surf 150A 5µm DNH 1kg
B-150A-5-SU-PEI-100G	Bulk Surf 150A 5µm PEI 100g
B-150A-5-SU-PEI-1KG	Bulk Surf 150A 5µm PEI 1kg
B-150A-5-SU-ARG-100G	Bulk Surf 150A 5µm ARG 100g
B-150A-5-SU-ARG-1KG	Bulk Surf 150A 5µm ARG 1kg
B-150A-5-SU-EP2TNE-100G	Bulk Surf 150A 5µm EP2TNE 100g
B-150A-5-SU-EP2TNE-1KG	Bulk Surf 150A 5µm EP2TNE 1kg
B-150A-5-SU-EP4TNE-100G	Bulk Surf 150A 5µm EP4TNE 100g
B-150A-5-SU-EP4TNE-1KG	Bulk Surf 150A 5µm EP4TNE 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-5-SU-NH2-100G	Bulk Surf 150A 5µm NH2 100g
B-150A-5-SU-NH2-1KG	Bulk Surf 150A 5µm NH2 1kg
B-150A-5-SU-CN-100G	Bulk Surf 150A 5µm CN 100g
B-150A-5-SU-CN-1KG	Bulk Surf 150A 5µm CN 1kg
B-150A-5-SU-Phenyl-100G	Bulk Surf 150A 5µm Phenyl 100g
B-150A-5-SU-Phenyl-1KG	Bulk Surf 150A 5µm Phenyl 1kg
B-150A-5-SU-Phenyl HE-100G	Bulk Surf 150A 5µm Phenyl HE 100g
B-150A-5-SU-Phenyl HE-1KG	Bulk Surf 150A 5µm Phenyl HE 1kg
B-150A-5-SU-Diol-100G	Bulk Surf 150A 5µm Diol 100g
B-150A-5-SU-Diol-1KG	Bulk Surf 150A 5µm Diol 1kg
<b>Particule Surf 5µm, pore 200A</b>	
B-200A-5-SU-100G	Bulk Surf 200A 5µm 100g
B-200A-5-SU-1KG	Bulk Surf 200A 5µm 1kg
<b>C8</b>	
B-200A-5-SU-C8-100G	Bulk Surf 200A 5µm C8 100g
B-200A-5-SU-C8-1KG	Bulk Surf 200A 5µm C8 1kg
B-200A-5-SU-C8W-100G	Bulk Surf 200A 5µm C8 Wax 100g
B-200A-5-SU-C8W-1KG	Bulk Surf 200A 5µm C8 Wax 1kg
B-200A-5-SU-C8P-100G	Bulk Surf 200A 5µm C8 Polar 100g
B-200A-5-SU-C8P-1KG	Bulk Surf 200A 5µm C8 Polar 1kg
B-200A-5-SU-C8T-100G	Bulk Surf 200A 5µm C8 TriF 100g
B-200A-5-SU-C8T-1KG	Bulk Surf 200A 5µm C8 TriF 1kg
B-200A-5-SU-C8EX-100G	Bulk Surf 200A 5µm C8 Extreme 100g
B-200A-5-SU-C8EX-1KG	Bulk Surf 200A 5µm C8 Extreme 1kg
B-200A-5-SU-C8PT-100G	Bulk Surf 200A 5µm C8 Polar TriF 100g
B-200A-5-SU-C8PT-1KG	Bulk Surf 200A 5µm C8 Polar TriF 1kg
B-200A-5-SU-C8AQP-100G	Bulk Surf 200A 5µm C8 AQ Plus 100g
B-200A-5-SU-C8AQP-1KG	Bulk Surf 200A 5µm C8 AQ Plus 1kg
B-200A-5-SU-C8AQ-100G	Bulk Surf 200A 5µm C8 AQ 100g
B-200A-5-SU-C8AQ-1KG	Bulk Surf 200A 5µm C8 AQ 1kg
<b>C18</b>	
B-200A-5-SU-C18-100G	Bulk Surf 200A 5µm C18 100g
B-200A-5-SU-C18-1KG	Bulk Surf 200A 5µm C18 1kg
B-200A-5-SU-C18W-100G	Bulk Surf 200A 5µm C18 Wax 100g
B-200A-5-SU-C18W-1KG	Bulk Surf 200A 5µm C18 Wax 1kg
B-200A-5-SU-C18P-100G	Bulk Surf 200A 5µm C18 Polar 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-5-SU-C18P-1KG	Bulk Surf 200A 5µm C18 Polar 1kg
B-200A-5-SU-C18T-100G	Bulk Surf 200A 5µm C18 TriF 100g
B-200A-5-SU-C18T-1KG	Bulk Surf 200A 5µm C18 TriF 1kg
B-200A-5-SU-C18EX-100G	Bulk Surf 200A 5µm C18 Extreme 100g
B-200A-5-SU-C18EX-1KG	Bulk Surf 200A 5µm C18 Extreme 1kg
B-200A-5-SU-C18PT-100G	Bulk Surf 200A 5µm C18 Polar TriF 100g
B-200A-5-SU-C18PT-1KG	Bulk Surf 200A 5µm C18 Polar TriF 1kg
B-200A-5-SU-C18AQP-100G	Bulk Surf 200A 5µm C18 AQ Plus 100g
B-200A-5-SU-C18AQP-1KG	Bulk Surf 200A 5µm C18 AQ Plus 1kg
B-200A-5-SU-C18AQ-100G	Bulk Surf 200A 5µm C18 AQ 100g
B-200A-5-SU-C18AQ-1KG	Bulk Surf 200A 5µm C18 AQ 1kg
<b>C4</b>	
B-200A-5-SU-C4-100G	Bulk Surf 200A 5µm C4 100g
B-200A-5-SU-C4-1KG	Bulk Surf 200A 5µm C4 1kg
B-200A-5-SU-C4W-100G	Bulk Surf 200A 5µm C4 Wax 100g
B-200A-5-SU-C4W-1KG	Bulk Surf 200A 5µm C4 Wax 1kg
B-200A-5-SU-C4P-100G	Bulk Surf 200A 5µm C4 Polar 100g
B-200A-5-SU-C4P-1KG	Bulk Surf 200A 5µm C4 Polar 1kg
B-200A-5-SU-C4T-100G	Bulk Surf 200A 5µm C4 TriF 100g
B-200A-5-SU-C4T-1KG	Bulk Surf 200A 5µm C4 TriF 1kg
B-200A-5-SU-C4EX-100G	Bulk Surf 200A 5µm C4 Extreme 100g
B-200A-5-SU-C4EX-1KG	Bulk Surf 200A 5µm C4 Extreme 1kg
B-200A-5-SU-C4PT-100G	Bulk Surf 200A 5µm C4 Polar TriF 100g
B-200A-5-SU-C4PT-1KG	Bulk Surf 200A 5µm C4 Polar TriF 1kg
B-200A-5-SU-C4AQP-100G	Bulk Surf 200A 5µm C4 AQ Plus 100g
B-200A-5-SU-C4AQP-1KG	Bulk Surf 200A 5µm C4 AQ Plus 1kg
B-200A-5-SU-C4AQ-100G	Bulk Surf 200A 5µm C4 AQ 100g
B-200A-5-SU-C4AQ-1KG	Bulk Surf 200A 5µm C4 AQ 1kg
<b>Autres</b>	
B-200A-5-SU-COOH-100G	Bulk Surf 200A 5µm COOH 100g
B-200A-5-SU-COOH-1KG	Bulk Surf 200A 5µm COOH 1kg
B-200A-5-SU-SO3H-100G	Bulk Surf 200A 5µm SO3H 100g
B-200A-5-SU-SO3H-1KG	Bulk Surf 200A 5µm SO3H 1kg
B-200A-5-SU-DNH-100G	Bulk Surf 200A 5µm DNH 100g
B-200A-5-SU-DNH-1KG	Bulk Surf 200A 5µm DNH 1kg
B-200A-5-SU-PEI-100G	Bulk Surf 200A 5µm PEI 100g
B-200A-5-SU-PEI-1KG	Bulk Surf 200A 5µm PEI 1kg
B-200A-5-SU-ARG-100G	Bulk Surf 200A 5µm ARG 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-5-SU-ARG-1KG	Bulk Surf 200A 5µm ARG 1kg
B-200A-5-SU-EP2TNE-100G	Bulk Surf 200A 5µm EP2TNE 100g
B-200A-5-SU-EP2TNE-1KG	Bulk Surf 200A 5µm EP2TNE 1kg
B-200A-5-SU-EP4TNE-100G	Bulk Surf 200A 5µm EP4TNE 100g
B-200A-5-SU-EP4TNE-1KG	Bulk Surf 200A 5µm EP4TNE 1kg
B-200A-5-SU-NH2-100G	Bulk Surf 200A 5µm NH2 100g
B-200A-5-SU-NH2-1KG	Bulk Surf 200A 5µm NH2 1kg
B-200A-5-SU-CN-100G	Bulk Surf 200A 5µm CN 100g
B-200A-5-SU-CN-1KG	Bulk Surf 200A 5µm CN 1kg
B-200A-5-SU-Phenyl-100G	Bulk Surf 200A 5µm Phenyl 100g
B-200A-5-SU-Phenyl-1KG	Bulk Surf 200A 5µm Phenyl 1kg
B-200A-5-SU-Phenyl HE-100G	Bulk Surf 200A 5µm Phenyl HE 100g
B-200A-5-SU-Phenyl HE-1KG	Bulk Surf 200A 5µm Phenyl HE 1kg
B-200A-5-SU-Diol-100G	Bulk Surf 200A 5µm Diol 100g
B-200A-5-SU-Diol-1KG	Bulk Surf 200A 5µm Diol 1kg
<b>Particule Surf 5µm, pore 300A</b>	
B-300A-5-SU-100G	Bulk Surf 300A 5µm 100g
B-300A-5-SU-1KG	Bulk Surf 300A 5µm 1kg
<b>C8</b>	
B-300A-5-SU-C8-100G	Bulk Surf 300A 5µm C8 100g
B-300A-5-SU-C8-1KG	Bulk Surf 300A 5µm C8 1kg
B-300A-5-SU-C8W-100G	Bulk Surf 300A 5µm C8 Wax 100g
B-300A-5-SU-C8W-1KG	Bulk Surf 300A 5µm C8 Wax 1kg
B-300A-5-SU-C8P-100G	Bulk Surf 300A 5µm C8 Polar 100g
B-300A-5-SU-C8P-1KG	Bulk Surf 300A 5µm C8 Polar 1kg
B-300A-5-SU-C8T-100G	Bulk Surf 300A 5µm C8 TriF 100g
B-300A-5-SU-C8T-1KG	Bulk Surf 300A 5µm C8 TriF 1kg
B-300A-5-SU-C8EX-100G	Bulk Surf 300A 5µm C8 Extreme 100g
B-300A-5-SU-C8EX-1KG	Bulk Surf 300A 5µm C8 Extreme 1kg
B-300A-5-SU-C8PT-100G	Bulk Surf 300A 5µm C8 Polar TriF 100g
B-300A-5-SU-C8PT-1KG	Bulk Surf 300A 5µm C8 Polar TriF 1kg
B-300A-5-SU-C8AQP-100G	Bulk Surf 300A 5µm C8 AQ Plus 100g
B-300A-5-SU-C8AQP-1KG	Bulk Surf 300A 5µm C8 AQ Plus 1kg
B-300A-5-SU-C8AQ-100G	Bulk Surf 300A 5µm C8 AQ 100g
B-300A-5-SU-C8AQ-1KG	Bulk Surf 300A 5µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-300A-5-SU-C18-100G	Bulk Surf 300A 5µm C18 100g
B-300A-5-SU-C18-1KG	Bulk Surf 300A 5µm C18 1kg
B-300A-5-SU-C18W-100G	Bulk Surf 300A 5µm C18 Wax 100g
B-300A-5-SU-C18W-1KG	Bulk Surf 300A 5µm C18 Wax 1kg
B-300A-5-SU-C18P-100G	Bulk Surf 300A 5µm C18 Polar 100g
B-300A-5-SU-C18P-1KG	Bulk Surf 300A 5µm C18 Polar 1kg
B-300A-5-SU-C18T-100G	Bulk Surf 300A 5µm C18 TriF 100g
B-300A-5-SU-C18T-1KG	Bulk Surf 300A 5µm C18 TriF 1kg
B-300A-5-SU-C18EX-100G	Bulk Surf 300A 5µm C18 Extreme 100g
B-300A-5-SU-C18EX-1KG	Bulk Surf 300A 5µm C18 Extreme 1kg
B-300A-5-SU-C18PT-100G	Bulk Surf 300A 5µm C18 Polar TriF 100g
B-300A-5-SU-C18PT-1KG	Bulk Surf 300A 5µm C18 Polar TriF 1kg
B-300A-5-SU-C18AQP-100G	Bulk Surf 300A 5µm C18 AQ Plus 100g
B-300A-5-SU-C18AQP-1KG	Bulk Surf 300A 5µm C18 AQ Plus 1kg
B-300A-5-SU-C18AQ-100G	Bulk Surf 300A 5µm C18 AQ 100g
B-300A-5-SU-C18AQ-1KG	Bulk Surf 300A 5µm C18 AQ 1kg
<b>C4</b>	
B-300A-5-SU-C4-100G	Bulk Surf 300A 5µm C4 100g
B-300A-5-SU-C4-1KG	Bulk Surf 300A 5µm C4 1kg
B-300A-5-SU-C4W-100G	Bulk Surf 300A 5µm C4 Wax 100g
B-300A-5-SU-C4W-1KG	Bulk Surf 300A 5µm C4 Wax 1kg
B-300A-5-SU-C4P-100G	Bulk Surf 300A 5µm C4 Polar 100g
B-300A-5-SU-C4P-1KG	Bulk Surf 300A 5µm C4 Polar 1kg
B-300A-5-SU-C4T-100G	Bulk Surf 300A 5µm C4 TriF 100g
B-300A-5-SU-C4T-1KG	Bulk Surf 300A 5µm C4 TriF 1kg
B-300A-5-SU-C4EX-100G	Bulk Surf 300A 5µm C4 Extreme 100g
B-300A-5-SU-C4EX-1KG	Bulk Surf 300A 5µm C4 Extreme 1kg
B-300A-5-SU-C4PT-100G	Bulk Surf 300A 5µm C4 Polar TriF 100g
B-300A-5-SU-C4PT-1KG	Bulk Surf 300A 5µm C4 Polar TriF 1kg
B-300A-5-SU-C4AQP-100G	Bulk Surf 300A 5µm C4 AQ Plus 100g
B-300A-5-SU-C4AQP-1KG	Bulk Surf 300A 5µm C4 AQ Plus 1kg
B-300A-5-SU-C4AQ-100G	Bulk Surf 300A 5µm C4 AQ 100g
B-300A-5-SU-C4AQ-1KG	Bulk Surf 300A 5µm C4 AQ 1kg
<b>Autres</b>	
B-300A-5-SU-COOH-100G	Bulk Surf 300A 5µm COOH 100g
B-300A-5-SU-COOH-1KG	Bulk Surf 300A 5µm COOH 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-5-SU-S03H-100G	Bulk Surf 300A 5µm S03H 100g
B-300A-5-SU-S03H-1KG	Bulk Surf 300A 5µm S03H 1kg
B-300A-5-SU-DNH-100G	Bulk Surf 300A 5µm DNH 100g
B-300A-5-SU-DNH-1KG	Bulk Surf 300A 5µm DNH 1kg
B-300A-5-SU-PEI-100G	Bulk Surf 300A 5µm PEI 100g
B-300A-5-SU-PEI-1KG	Bulk Surf 300A 5µm PEI 1kg
B-300A-5-SU-ARG-100G	Bulk Surf 300A 5µm ARG 100g
B-300A-5-SU-ARG-1KG	Bulk Surf 300A 5µm ARG 1kg
B-300A-5-SU-EP2TNE-100G	Bulk Surf 300A 5µm EP2TNE 100g
B-300A-5-SU-EP2TNE-1KG	Bulk Surf 300A 5µm EP2TNE 1kg
B-300A-5-SU-EP4TNE-100G	Bulk Surf 300A 5µm EP4TNE 100g
B-300A-5-SU-EP4TNE-1KG	Bulk Surf 300A 5µm EP4TNE 1kg
B-300A-5-SU-NH2-100G	Bulk Surf 300A 5µm NH2 100g
B-300A-5-SU-NH2-1KG	Bulk Surf 300A 5µm NH2 1kg
B-300A-5-SU-CN-100G	Bulk Surf 300A 5µm CN 100g
B-300A-5-SU-CN-1KG	Bulk Surf 300A 5µm CN 1kg
B-300A-5-SU-Phenyl-100G	Bulk Surf 300A 5µm Phenyl 100g
B-300A-5-SU-Phenyl-1KG	Bulk Surf 300A 5µm Phenyl 1kg
B-300A-5-SU-Phenyl HE-100G	Bulk Surf 300A 5µm Phenyl HE 100g
B-300A-5-SU-Phenyl HE-1KG	Bulk Surf 300A 5µm Phenyl HE 1kg
B-300A-5-SU-Diol-100G	Bulk Surf 300A 5µm Diol 100g
B-300A-5-SU-Diol-1KG	Bulk Surf 300A 5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 5,5 µm</b>	
<b>Particule Surf 5,5µm, pore 100A</b>	
B-100A-5.5-SU-Si-100G	Bulk Surf Silice 100A 5.5µm 100g
B-100A-5.5-SU-Si-1KG	Bulk Surf Silice 100A 5.5µm 1kg
<b>C8</b>	
B-100A-5.5-SU-C8-100G	Bulk Surf 100A 5.5µm C8 100g
B-100A-5.5-SU-C8-1KG	Bulk Surf 100A 5.5µm C8 1kg
B-100A-5.5-SU-C8W-100G	Bulk Surf 100A 5.5µm C8 Wax 100g
B-100A-5.5-SU-C8W-1KG	Bulk Surf 100A 5.5µm C8 Wax 1kg
B-100A-5.5-SU-C8P-100G	Bulk Surf 100A 5.5µm C8 Polar 100g
B-100A-5.5-SU-C8P-1KG	Bulk Surf 100A 5.5µm C8 Polar 1kg
B-100A-5.5-SU-C8T-100G	Bulk Surf 100A 5.5µm C8 TriF 100g
B-100A-5.5-SU-C8T-1KG	Bulk Surf 100A 5.5µm C8 TriF 1kg
B-100A-5.5-SU-C8EX-100G	Bulk Surf 100A 5.5µm C8 Extreme 100g
B-100A-5.5-SU-C8EX-1KG	Bulk Surf 100A 5.5µm C8 Extreme 1kg
B-100A-5.5-SU-C8PT-100G	Bulk Surf 100A 5.5µm C8 Polar TriF 100g
B-100A-5.5-SU-C8PT-1KG	Bulk Surf 100A 5.5µm C8 Polar TriF 1kg
B-100A-5.5-SU-C8AQP-100G	Bulk Surf 100A 5.5µm C8 AQ Plus 100g
B-100A-5.5-SU-C8AQP-1KG	Bulk Surf 100A 5.5µm C8 AQ Plus 1kg
B-100A-5.5-SU-C8AQ-100G	Bulk Surf 100A 5.5µm C8 AQ 100g
B-100A-5.5-SU-C8AQ-1KG	Bulk Surf 100A 5.5µm C8 AQ 1kg
<b>C18</b>	
B-100A-5.5-SU-C18-100G	Bulk Surf 100A 5.5µm C18 100g
B-100A-5.5-SU-C18-1KG	Bulk Surf 100A 5.5µm C18 1kg
B-100A-5.5-SU-C18W-100G	Bulk Surf 100A 5.5µm C18 Wax 100g
B-100A-5.5-SU-C18W-1KG	Bulk Surf 100A 5.5µm C18 Wax 1kg
B-100A-5.5-SU-C18P-100G	Bulk Surf 100A 5.5µm C18 Polar 100g
B-100A-5.5-SU-C18P-1KG	Bulk Surf 100A 5.5µm C18 Polar 1kg
B-100A-5.5-SU-C18T-100G	Bulk Surf 100A 5.5µm C18 TriF 100g
B-100A-5.5-SU-C18T-1KG	Bulk Surf 100A 5.5µm C18 TriF 1kg
B-100A-5.5-SU-C18EX-100G	Bulk Surf 100A 5.5µm C18 Extreme 100g
B-100A-5.5-SU-C18EX-1KG	Bulk Surf 100A 5.5µm C18 Extreme 1kg
B-100A-5.5-SU-C18PT-100G	Bulk Surf 100A 5.5µm C18 Polar TriF 100g
B-100A-5.5-SU-C18PT-1KG	Bulk Surf 100A 5.5µm C18 Polar TriF 1kg
B-100A-5.5-SU-C18AQP-100G	Bulk Surf 100A 5.5µm C18 AQ Plus 100g
B-100A-5.5-SU-C18AQP-1KG	Bulk Surf 100A 5.5µm C18 AQ Plus 1kg
B-100A-5.5-SU-C18AQ-100G	Bulk Surf 100A 5.5µm C18 AQ 100g
B-100A-5.5-SU-C18AQ-1KG	Bulk Surf 100A 5.5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-5.5-SU-C4-100G	Bulk Surf 100A 5.5µm C4 100g
B-100A-5.5-SU-C4-1KG	Bulk Surf 100A 5.5µm C4 1kg
B-100A-5.5-SU-C4W-100G	Bulk Surf 100A 5.5µm C4 Wax 100g
B-100A-5.5-SU-C4W-1KG	Bulk Surf 100A 5.5µm C4 Wax 1kg
B-100A-5.5-SU-C4P-100G	Bulk Surf 100A 5.5µm C4 Polar 100g
B-100A-5.5-SU-C4P-1KG	Bulk Surf 100A 5.5µm C4 Polar 1kg
B-100A-5.5-SU-C4T-100G	Bulk Surf 100A 5.5µm C4 TriF 100g
B-100A-5.5-SU-C4T-1KG	Bulk Surf 100A 5.5µm C4 TriF 1kg
B-100A-5.5-SU-C4EX-100G	Bulk Surf 100A 5.5µm C4 Extreme 100g
B-100A-5.5-SU-C4EX-1KG	Bulk Surf 100A 5.5µm C4 Extreme 1kg
B-100A-5.5-SU-C4PT-100G	Bulk Surf 100A 5.5µm C4 Polar TriF 100g
B-100A-5.5-SU-C4PT-1KG	Bulk Surf 100A 5.5µm C4 Polar TriF 1kg
B-100A-5.5-SU-C4AQP-100G	Bulk Surf 100A 5.5µm C4 AQ Plus 100g
B-100A-5.5-SU-C4AQP-1KG	Bulk Surf 100A 5.5µm C4 AQ Plus 1kg
B-100A-5.5-SU-C4AQ-100G	Bulk Surf 100A 5.5µm C4 AQ 100g
B-100A-5.5-SU-C4AQ-1KG	Bulk Surf 100A 5.5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-5.5-SU-COOH-100G	Bulk Surf 100A 5.5µm COOH 100g
B-100A-5.5-SU-COOH-1KG	Bulk Surf 100A 5.5µm COOH 1kg
B-100A-5.5-SU-SO3H-100G	Bulk Surf 100A 5.5µm SO3H 100g
B-100A-5.5-SU-SO3H-1KG	Bulk Surf 100A 5.5µm SO3H 1kg
B-100A-5.5-SU-DNH-100G	Bulk Surf 100A 5.5µm DNH 100g
B-100A-5.5-SU-DNH-1KG	Bulk Surf 100A 5.5µm DNH 1kg
B-100A-5.5-SU-PEI-100G	Bulk Surf 100A 5.5µm PEI 100g
B-100A-5.5-SU-PEI-1KG	Bulk Surf 100A 5.5µm PEI 1kg
B-100A-5.5-SU-ARG-100G	Bulk Surf 100A 5.5µm ARG 100g
B-100A-5.5-SU-ARG-1KG	Bulk Surf 100A 5.5µm ARG 1kg
B-100A-5.5-SU-EP2TNE-100G	Bulk Surf 100A 5.5µm EP2TNE 100g
B-100A-5.5-SU-EP2TNE-1KG	Bulk Surf 100A 5.5µm EP2TNE 1kg
B-100A-5.5-SU-EP4TNE-100G	Bulk Surf 100A 5.5µm EP4TNE 100g
B-100A-5.5-SU-EP4TNE-1KG	Bulk Surf 100A 5.5µm EP4TNE 1kg
B-100A-5.5-SU-NH2-100G	Bulk Surf 100A 5.5µm NH2 100g
B-100A-5.5-SU-NH2-1KG	Bulk Surf 100A 5.5µm NH2 1kg
B-100A-5.5-SU-CN-100G	Bulk Surf 100A 5.5µm CN 100g
B-100A-5.5-SU-CN-1KG	Bulk Surf 100A 5.5µm CN 1kg
B-100A-5.5-SU-Phenyl-100G	Bulk Surf 100A 5.5µm Phenyl 100g
B-100A-5.5-SU-Phenyl-1KG	Bulk Surf 100A 5.5µm Phenyl 1kg
B-100A-5.5-SU-Phenyl HE-100G	Bulk Surf 100A 5.5µm Phenyl HE 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-5.5-SU-Phenyl HE-1KG	Bulk Surf 100A 5.5µm Phenyl HE 1kg
B-100A-5.5-SU-Diol-100G	Bulk Surf 100A 5.5µm Diol 100g
B-100A-5.5-SU-Diol-1KG	Bulk Surf 100A 5.5µm Diol 1kg
<b>Particule Surf 5,5µm, pore 150A</b>	
B-150A-5.5-SU-100G	Bulk Surf 150A 5.5µm 100g
B-150A-5.5-SU-1KG	Bulk Surf 150A 5.5µm 1kg
<b>C8</b>	
B-150A-5.5-SU-C8-100G	Bulk Surf 150A 5.5µm C8 100g
B-150A-5.5-SU-C8-1KG	Bulk Surf 150A 5.5µm C8 1kg
B-150A-5.5-SU-C8W-100G	Bulk Surf 150A 5.5µm C8 Wax 100g
B-150A-5.5-SU-C8W-1KG	Bulk Surf 150A 5.5µm C8 Wax 1kg
B-150A-5.5-SU-C8P-100G	Bulk Surf 150A 5.5µm C8 Polar 100g
B-150A-5.5-SU-C8P-1KG	Bulk Surf 150A 5.5µm C8 Polar 1kg
B-150A-5.5-SU-C8T-100G	Bulk Surf 150A 5.5µm C8 TriF 100g
B-150A-5.5-SU-C8T-1KG	Bulk Surf 150A 5.5µm C8 TriF 1kg
B-150A-5.5-SU-C8EX-100G	Bulk Surf 150A 5.5µm C8 Extreme 100g
B-150A-5.5-SU-C8EX-1KG	Bulk Surf 150A 5.5µm C8 Extreme 1kg
B-150A-5.5-SU-C8PT-100G	Bulk Surf 150A 5.5µm C8 Polar TriF 100g
B-150A-5.5-SU-C8PT-1KG	Bulk Surf 150A 5.5µm C8 Polar TriF 1kg
B-150A-5.5-SU-C8AQP-100G	Bulk Surf 150A 5.5µm C8 AQ Plus 100g
B-150A-5.5-SU-C8AQP-1KG	Bulk Surf 150A 5.5µm C8 AQ Plus 1kg
B-150A-5.5-SU-C8AQ-100G	Bulk Surf 150A 5.5µm C8 AQ 100g
B-150A-5.5-SU-C8AQ-1KG	Bulk Surf 150A 5.5µm C8 AQ 1kg
<b>C18</b>	
B-150A-5.5-SU-C18-100G	Bulk Surf 150A 5.5µm C18 100g
B-150A-5.5-SU-C18-1KG	Bulk Surf 150A 5.5µm C18 1kg
B-150A-5.5-SU-C18W-100G	Bulk Surf 150A 5.5µm C18 Wax 100g
B-150A-5.5-SU-C18W-1KG	Bulk Surf 150A 5.5µm C18 Wax 1kg
B-150A-5.5-SU-C18P-100G	Bulk Surf 150A 5.5µm C18 Polar 100g
B-150A-5.5-SU-C18P-1KG	Bulk Surf 150A 5.5µm C18 Polar 1kg
B-150A-5.5-SU-C18T-100G	Bulk Surf 150A 5.5µm C18 TriF 100g
B-150A-5.5-SU-C18T-1KG	Bulk Surf 150A 5.5µm C18 TriF 1kg
B-150A-5.5-SU-C18EX-100G	Bulk Surf 150A 5.5µm C18 Extreme 100g
B-150A-5.5-SU-C18EX-1KG	Bulk Surf 150A 5.5µm C18 Extreme 1kg
B-150A-5.5-SU-C18PT-100G	Bulk Surf 150A 5.5µm C18 Polar TriF 100g
B-150A-5.5-SU-C18PT-1KG	Bulk Surf 150A 5.5µm C18 Polar TriF 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-5.5-SU-C18AQP-100G	Bulk Surf 150A 5.5µm C18 AQ Plus 100g
B-150A-5.5-SU-C18AQP-1KG	Bulk Surf 150A 5.5µm C18 AQ Plus 1kg
B-150A-5.5-SU-C18AQ-100G	Bulk Surf 150A 5.5µm C18 AQ 100g
B-150A-5.5-SU-C18AQ-1KG	Bulk Surf 150A 5.5µm C18 AQ 1kg
<b>C4</b>	
B-150A-5.5-SU-C4-100G	Bulk Surf 150A 5.5µm C4 100g
B-150A-5.5-SU-C4-1KG	Bulk Surf 150A 5.5µm C4 1kg
B-150A-5.5-SU-C4W-100G	Bulk Surf 150A 5.5µm C4 Wax 100g
B-150A-5.5-SU-C4W-1KG	Bulk Surf 150A 5.5µm C4 Wax 1kg
B-150A-5.5-SU-C4P-100G	Bulk Surf 150A 5.5µm C4 Polar 100g
B-150A-5.5-SU-C4P-1KG	Bulk Surf 150A 5.5µm C4 Polar 1kg
B-150A-5.5-SU-C4T-100G	Bulk Surf 150A 5.5µm C4 TriF 100g
B-150A-5.5-SU-C4T-1KG	Bulk Surf 150A 5.5µm C4 TriF 1kg
B-150A-5.5-SU-C4EX-100G	Bulk Surf 150A 5.5µm C4 Extreme 100g
B-150A-5.5-SU-C4EX-1KG	Bulk Surf 150A 5.5µm C4 Extreme 1kg
B-150A-5.5-SU-C4PT-100G	Bulk Surf 150A 5.5µm C4 Polar TriF 100g
B-150A-5.5-SU-C4PT-1KG	Bulk Surf 150A 5.5µm C4 Polar TriF 1kg
B-150A-5.5-SU-C4AQP-100G	Bulk Surf 150A 5.5µm C4 AQ Plus 100g
B-150A-5.5-SU-C4AQP-1KG	Bulk Surf 150A 5.5µm C4 AQ Plus 1kg
B-150A-5.5-SU-C4AQ-100G	Bulk Surf 150A 5.5µm C4 AQ 100g
B-150A-5.5-SU-C4AQ-1KG	Bulk Surf 150A 5.5µm C4 AQ 1kg
<b>Autres</b>	
B-150A-5.5-SU-COOH-100G	Bulk Surf 150A 5.5µm COOH 100g
B-150A-5.5-SU-COOH-1KG	Bulk Surf 150A 5.5µm COOH 1kg
B-150A-5.5-SU-SO3H-100G	Bulk Surf 150A 5.5µm SO3H 100g
B-150A-5.5-SU-SO3H-1KG	Bulk Surf 150A 5.5µm SO3H 1kg
B-150A-5.5-SU-DNH-100G	Bulk Surf 150A 5.5µm DNH 100g
B-150A-5.5-SU-DNH-1KG	Bulk Surf 150A 5.5µm DNH 1kg
B-150A-5.5-SU-PEI-100G	Bulk Surf 150A 5.5µm PEI 100g
B-150A-5.5-SU-PEI-1KG	Bulk Surf 150A 5.5µm PEI 1kg
B-150A-5.5-SU-ARG-100G	Bulk Surf 150A 5.5µm ARG 100g
B-150A-5.5-SU-ARG-1KG	Bulk Surf 150A 5.5µm ARG 1kg
B-150A-5.5-SU-EP2TNE-100G	Bulk Surf 150A 5.5µm EP2TNE 100g
B-150A-5.5-SU-EP2TNE-1KG	Bulk Surf 150A 5.5µm EP2TNE 1kg
B-150A-5.5-SU-EP4TNE-100G	Bulk Surf 150A 5.5µm EP4TNE 100g
B-150A-5.5-SU-EP4TNE-1KG	Bulk Surf 150A 5.5µm EP4TNE 1kg
B-150A-5.5-SU-NH2-100G	Bulk Surf 150A 5.5µm NH2 100g
B-150A-5.5-SU-NH2-1KG	Bulk Surf 150A 5.5µm NH2 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-5.5-SU-CN-100G	Bulk Surf 150A 5.5µm CN 100g
B-150A-5.5-SU-CN-1KG	Bulk Surf 150A 5.5µm CN 1kg
B-150A-5.5-SU-Phenyl-100G	Bulk Surf 150A 5.5µm Phenyl 100g
B-150A-5.5-SU-Phenyl-1KG	Bulk Surf 150A 5.5µm Phenyl 1kg
B-150A-5.5-SU-Phenyl HE-100G	Bulk Surf 150A 5.5µm Phenyl HE 100g
B-150A-5.5-SU-Phenyl HE-1KG	Bulk Surf 150A 5.5µm Phenyl HE 1kg
B-150A-5.5-SU-Diol-100G	Bulk Surf 150A 5.5µm Diol 100g
B-150A-5.5-SU-Diol-1KG	Bulk Surf 150A 5.5µm Diol 1kg
<b>Particule Surf 5,5µm, pore 200A</b>	
B-200A-5.5-SU-100G	Bulk Surf 200A 5.5µm 100g
B-200A-5.5-SU-1KG	Bulk Surf 200A 5.5µm 1kg
<b>C8</b>	
B-200A-5.5-SU-C8-100G	Bulk Surf 200A 5.5µm C8 100g
B-200A-5.5-SU-C8-1KG	Bulk Surf 200A 5.5µm C8 1kg
B-200A-5.5-SU-C8W-100G	Bulk Surf 200A 5.5µm C8 Wax 100g
B-200A-5.5-SU-C8W-1KG	Bulk Surf 200A 5.5µm C8 Wax 1kg
B-200A-5.5-SU-C8P-100G	Bulk Surf 200A 5.5µm C8 Polar 100g
B-200A-5.5-SU-C8P-1KG	Bulk Surf 200A 5.5µm C8 Polar 1kg
B-200A-5.5-SU-C8T-100G	Bulk Surf 200A 5.5µm C8 TriF 100g
B-200A-5.5-SU-C8T-1KG	Bulk Surf 200A 5.5µm C8 TriF 1kg
B-200A-5.5-SU-C8EX-100G	Bulk Surf 200A 5.5µm C8 Extreme 100g
B-200A-5.5-SU-C8EX-1KG	Bulk Surf 200A 5.5µm C8 Extreme 1kg
B-200A-5.5-SU-C8PT-100G	Bulk Surf 200A 5.5µm C8 Polar TriF 100g
B-200A-5.5-SU-C8PT-1KG	Bulk Surf 200A 5.5µm C8 Polar TriF 1kg
B-200A-5.5-SU-C8AQP-100G	Bulk Surf 200A 5.5µm C8 AQ Plus 100g
B-200A-5.5-SU-C8AQP-1KG	Bulk Surf 200A 5.5µm C8 AQ Plus 1kg
B-200A-5.5-SU-C8AQ-100G	Bulk Surf 200A 5.5µm C8 AQ 100g
B-200A-5.5-SU-C8AQ-1KG	Bulk Surf 200A 5.5µm C8 AQ 1kg
<b>C18</b>	
B-200A-5.5-SU-C18-100G	Bulk Surf 200A 5.5µm C18 100g
B-200A-5.5-SU-C18-1KG	Bulk Surf 200A 5.5µm C18 1kg
B-200A-5.5-SU-C18W-100G	Bulk Surf 200A 5.5µm C18 Wax 100g
B-200A-5.5-SU-C18W-1KG	Bulk Surf 200A 5.5µm C18 Wax 1kg
B-200A-5.5-SU-C18P-100G	Bulk Surf 200A 5.5µm C18 Polar 100g
B-200A-5.5-SU-C18P-1KG	Bulk Surf 200A 5.5µm C18 Polar 1kg
B-200A-5.5-SU-C18T-100G	Bulk Surf 200A 5.5µm C18 TriF 100g



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-5.5-SU-C18T-1KG	Bulk Surf 200A 5.5µm C18 TriF 1kg
B-200A-5.5-SU-C18EX-100G	Bulk Surf 200A 5.5µm C18 Extreme 100g
B-200A-5.5-SU-C18EX-1KG	Bulk Surf 200A 5.5µm C18 Extreme 1kg
B-200A-5.5-SU-C18PT-100G	Bulk Surf 200A 5.5µm C18 Polar TriF 100g
B-200A-5.5-SU-C18PT-1KG	Bulk Surf 200A 5.5µm C18 Polar TriF 1kg
B-200A-5.5-SU-C18AQP-100G	Bulk Surf 200A 5.5µm C18 AQ Plus 100g
B-200A-5.5-SU-C18AQP-1KG	Bulk Surf 200A 5.5µm C18 AQ Plus 1kg
B-200A-5.5-SU-C18AQ-100G	Bulk Surf 200A 5.5µm C18 AQ 100g
B-200A-5.5-SU-C18AQ-1KG	Bulk Surf 200A 5.5µm C18 AQ 1kg
<b>C4</b>	
B-200A-5.5-SU-C4-100G	Bulk Surf 200A 5.5µm C4 100g
B-200A-5.5-SU-C4-1KG	Bulk Surf 200A 5.5µm C4 1kg
B-200A-5.5-SU-C4W-100G	Bulk Surf 200A 5.5µm C4 Wax 100g
B-200A-5.5-SU-C4W-1KG	Bulk Surf 200A 5.5µm C4 Wax 1kg
B-200A-5.5-SU-C4P-100G	Bulk Surf 200A 5.5µm C4 Polar 100g
B-200A-5.5-SU-C4P-1KG	Bulk Surf 200A 5.5µm C4 Polar 1kg
B-200A-5.5-SU-C4T-100G	Bulk Surf 200A 5.5µm C4 TriF 100g
B-200A-5.5-SU-C4T-1KG	Bulk Surf 200A 5.5µm C4 TriF 1kg
B-200A-5.5-SU-C4EX-100G	Bulk Surf 200A 5.5µm C4 Extreme 100g
B-200A-5.5-SU-C4EX-1KG	Bulk Surf 200A 5.5µm C4 Extreme 1kg
B-200A-5.5-SU-C4PT-100G	Bulk Surf 200A 5.5µm C4 Polar TriF 100g
B-200A-5.5-SU-C4PT-1KG	Bulk Surf 200A 5.5µm C4 Polar TriF 1kg
B-200A-5.5-SU-C4AQP-100G	Bulk Surf 200A 5.5µm C4 AQ Plus 100g
B-200A-5.5-SU-C4AQP-1KG	Bulk Surf 200A 5.5µm C4 AQ Plus 1kg
B-200A-5.5-SU-C4AQ-100G	Bulk Surf 200A 5.5µm C4 AQ 100g
B-200A-5.5-SU-C4AQ-1KG	Bulk Surf 200A 5.5µm C4 AQ 1kg
<b>Autres</b>	
B-200A-5.5-SU-COOH-100G	Bulk Surf 200A 5.5µm COOH 100g
B-200A-5.5-SU-COOH-1KG	Bulk Surf 200A 5.5µm COOH 1kg
B-200A-5.5-SU-SO3H-100G	Bulk Surf 200A 5.5µm SO3H 100g
B-200A-5.5-SU-SO3H-1KG	Bulk Surf 200A 5.5µm SO3H 1kg
B-200A-5.5-SU-DNH-100G	Bulk Surf 200A 5.5µm DNH 100g
B-200A-5.5-SU-DNH-1KG	Bulk Surf 200A 5.5µm DNH 1kg
B-200A-5.5-SU-PEI-100G	Bulk Surf 200A 5.5µm PEI 100g
B-200A-5.5-SU-PEI-1KG	Bulk Surf 200A 5.5µm PEI 1kg
B-200A-5.5-SU-ARG-100G	Bulk Surf 200A 5.5µm ARG 100g
B-200A-5.5-SU-ARG-1KG	Bulk Surf 200A 5.5µm ARG 1kg
B-200A-5.5-SU-EP2TNE-100G	Bulk Surf 200A 5.5µm EP2TNE 100g
B-200A-5.5-SU-EP2TNE-1KG	Bulk Surf 200A 5.5µm EP2TNE 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-5.5-SU-EP4TNE-100G	Bulk Surf 200A 5.5µm EP4TNE 100g
B-200A-5.5-SU-EP4TNE-1KG	Bulk Surf 200A 5.5µm EP4TNE 1kg
B-200A-5.5-SU-NH2-100G	Bulk Surf 200A 5.5µm NH2 100g
B-200A-5.5-SU-NH2-1KG	Bulk Surf 200A 5.5µm NH2 1kg
B-200A-5.5-SU-CN-100G	Bulk Surf 200A 5.5µm CN 100g
B-200A-5.5-SU-CN-1KG	Bulk Surf 200A 5.5µm CN 1kg
B-200A-5.5-SU-Phenyl-100G	Bulk Surf 200A 5.5µm Phenyl 100g
B-200A-5.5-SU-Phenyl-1KG	Bulk Surf 200A 5.5µm Phenyl 1kg
B-200A-5.5-SU-Phenyl HE-100G	Bulk Surf 200A 5.5µm Phenyl HE 100g
B-200A-5.5-SU-Phenyl HE-1KG	Bulk Surf 200A 5.5µm Phenyl HE 1kg
B-200A-5.5-SU-Diol-100G	Bulk Surf 200A 5.5µm Diol 100g
B-200A-5.5-SU-Diol-1KG	Bulk Surf 200A 5.5µm Diol 1kg
<b>Particule Surf 5,5µm, pore 300A</b>	
B-300A-5.5-SU-100G	Bulk Surf 300A 5.5µm 100g
B-300A-5.5-SU-1KG	Bulk Surf 300A 5.5µm 1kg
<b>C8</b>	
B-300A-5.5-SU-C8-100G	Bulk Surf 300A 5.5µm C8 100g
B-300A-5.5-SU-C8-1KG	Bulk Surf 300A 5.5µm C8 1kg
B-300A-5.5-SU-C8W-100G	Bulk Surf 300A 5.5µm C8 Wax 100g
B-300A-5.5-SU-C8W-1KG	Bulk Surf 300A 5.5µm C8 Wax 1kg
B-300A-5.5-SU-C8P-100G	Bulk Surf 300A 5.5µm C8 Polar 100g
B-300A-5.5-SU-C8P-1KG	Bulk Surf 300A 5.5µm C8 Polar 1kg
B-300A-5.5-SU-C8T-100G	Bulk Surf 300A 5.5µm C8 TriF 100g
B-300A-5.5-SU-C8T-1KG	Bulk Surf 300A 5.5µm C8 TriF 1kg
B-300A-5.5-SU-C8EX-100G	Bulk Surf 300A 5.5µm C8 Extreme 100g
B-300A-5.5-SU-C8EX-1KG	Bulk Surf 300A 5.5µm C8 Extreme 1kg
B-300A-5.5-SU-C8PT-100G	Bulk Surf 300A 5.5µm C8 Polar TriF 100g
B-300A-5.5-SU-C8PT-1KG	Bulk Surf 300A 5.5µm C8 Polar TriF 1kg
B-300A-5.5-SU-C8AQP-100G	Bulk Surf 300A 5.5µm C8 AQ Plus 100g
B-300A-5.5-SU-C8AQP-1KG	Bulk Surf 300A 5.5µm C8 AQ Plus 1kg
B-300A-5.5-SU-C8AQ-100G	Bulk Surf 300A 5.5µm C8 AQ 100g
B-300A-5.5-SU-C8AQ-1KG	Bulk Surf 300A 5.5µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-300A-5.5-SU-C18-100G	Bulk Surf 300A 5.5µm C18 100g
B-300A-5.5-SU-C18-1KG	Bulk Surf 300A 5.5µm C18 1kg
B-300A-5.5-SU-C18W-100G	Bulk Surf 300A 5.5µm C18 Wax 100g
B-300A-5.5-SU-C18W-1KG	Bulk Surf 300A 5.5µm C18 Wax 1kg
B-300A-5.5-SU-C18P-100G	Bulk Surf 300A 5.5µm C18 Polar 100g
B-300A-5.5-SU-C18P-1KG	Bulk Surf 300A 5.5µm C18 Polar 1kg
B-300A-5.5-SU-C18T-100G	Bulk Surf 300A 5.5µm C18 TriF 100g
B-300A-5.5-SU-C18T-1KG	Bulk Surf 300A 5.5µm C18 TriF 1kg
B-300A-5.5-SU-C18EX-100G	Bulk Surf 300A 5.5µm C18 Extreme 100g
B-300A-5.5-SU-C18EX-1KG	Bulk Surf 300A 5.5µm C18 Extreme 1kg
B-300A-5.5-SU-C18PT-100G	Bulk Surf 300A 5.5µm C18 Polar TriF 100g
B-300A-5.5-SU-C18PT-1KG	Bulk Surf 300A 5.5µm C18 Polar TriF 1kg
B-300A-5.5-SU-C18AQP-100G	Bulk Surf 300A 5.5µm C18 AQ Plus 100g
B-300A-5.5-SU-C18AQP-1KG	Bulk Surf 300A 5.5µm C18 AQ Plus 1kg
B-300A-5.5-SU-C18AQ-100G	Bulk Surf 300A 5.5µm C18 AQ 100g
B-300A-5.5-SU-C18AQ-1KG	Bulk Surf 300A 5.5µm C18 AQ 1kg
<b>C4</b>	
B-300A-5.5-SU-C4-100G	Bulk Surf 300A 5.5µm C4 100g
B-300A-5.5-SU-C4-1KG	Bulk Surf 300A 5.5µm C4 1kg
B-300A-5.5-SU-C4W-100G	Bulk Surf 300A 5.5µm C4 Wax 100g
B-300A-5.5-SU-C4W-1KG	Bulk Surf 300A 5.5µm C4 Wax 1kg
B-300A-5.5-SU-C4P-100G	Bulk Surf 300A 5.5µm C4 Polar 100g
B-300A-5.5-SU-C4P-1KG	Bulk Surf 300A 5.5µm C4 Polar 1kg
B-300A-5.5-SU-C4T-100G	Bulk Surf 300A 5.5µm C4 TriF 100g
B-300A-5.5-SU-C4T-1KG	Bulk Surf 300A 5.5µm C4 TriF 1kg
B-300A-5.5-SU-C4EX-100G	Bulk Surf 300A 5.5µm C4 Extreme 100g
B-300A-5.5-SU-C4EX-1KG	Bulk Surf 300A 5.5µm C4 Extreme 1kg
B-300A-5.5-SU-C4PT-100G	Bulk Surf 300A 5.5µm C4 Polar TriF 100g
B-300A-5.5-SU-C4PT-1KG	Bulk Surf 300A 5.5µm C4 Polar TriF 1kg
B-300A-5.5-SU-C4AQP-100G	Bulk Surf 300A 5.5µm C4 AQ Plus 100g
B-300A-5.5-SU-C4AQP-1KG	Bulk Surf 300A 5.5µm C4 AQ Plus 1kg
B-300A-5.5-SU-C4AQ-100G	Bulk Surf 300A 5.5µm C4 AQ 100g
B-300A-5.5-SU-C4AQ-1KG	Bulk Surf 300A 5.5µm C4 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Autres</b>	
B-300A-5.5-SU-COOH-100G	Bulk Surf 300A 5.5µm COOH 100g
B-300A-5.5-SU-COOH-1KG	Bulk Surf 300A 5.5µm COOH 1kg
B-300A-5.5-SU-SO3H-100G	Bulk Surf 300A 5.5µm SO3H 100g
B-300A-5.5-SU-SO3H-1KG	Bulk Surf 300A 5.5µm SO3H 1kg
B-300A-5.5-SU-DNH-100G	Bulk Surf 300A 5.5µm DNH 100g
B-300A-5.5-SU-DNH-1KG	Bulk Surf 300A 5.5µm DNH 1kg
B-300A-5.5-SU-PEI-100G	Bulk Surf 300A 5.5µm PEI 100g
B-300A-5.5-SU-PEI-1KG	Bulk Surf 300A 5.5µm PEI 1kg
B-300A-5.5-SU-ARG-100G	Bulk Surf 300A 5.5µm ARG 100g
B-300A-5.5-SU-ARG-1KG	Bulk Surf 300A 5.5µm ARG 1kg
B-300A-5.5-SU-EP2TNE-100G	Bulk Surf 300A 5.5µm EP2TNE 100g
B-300A-5.5-SU-EP2TNE-1KG	Bulk Surf 300A 5.5µm EP2TNE 1kg
B-300A-5.5-SU-EP4TNE-100G	Bulk Surf 300A 5.5µm EP4TNE 100g
B-300A-5.5-SU-EP4TNE-1KG	Bulk Surf 300A 5.5µm EP4TNE 1kg
B-300A-5.5-SU-NH2-100G	Bulk Surf 300A 5.5µm NH2 100g
B-300A-5.5-SU-NH2-1KG	Bulk Surf 300A 5.5µm NH2 1kg
B-300A-5.5-SU-CN-100G	Bulk Surf 300A 5.5µm CN 100g
B-300A-5.5-SU-CN-1KG	Bulk Surf 300A 5.5µm CN 1kg
B-300A-5.5-SU-Phenyl-100G	Bulk Surf 300A 5.5µm Phenyl 100g
B-300A-5.5-SU-Phenyl-1KG	Bulk Surf 300A 5.5µm Phenyl 1kg
B-300A-5.5-SU-Phenyl HE-100G	Bulk Surf 300A 5.5µm Phenyl HE 100g
B-300A-5.5-SU-Phenyl HE-1KG	Bulk Surf 300A 5.5µm Phenyl HE 1kg
B-300A-5.5-SU-Diol-100G	Bulk Surf 300A 5.5µm Diol 100g
B-300A-5.5-SU-Diol-1KG	Bulk Surf 300A 5.5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Surf 7 µm</b>	
<b>Particule Surf 7µm, pore 100A</b>	
B-100A-7-SU-Si-100G	Bulk Surf Silice 100A 7µm 100g
B-100A-7-SU-Si-1KG	Bulk Surf Silice 100A 7µm 1kg
<b>C8</b>	
B-100A-7-SU-C8-100G	Bulk Surf 100A 7µm C8 100g
B-100A-7-SU-C8-1KG	Bulk Surf 100A 7µm C8 1kg
B-100A-7-SU-C8W-100G	Bulk Surf 100A 7µm C8 Wax 100g
B-100A-7-SU-C8W-1KG	Bulk Surf 100A 7µm C8 Wax 1kg
B-100A-7-SU-C8P-100G	Bulk Surf 100A 7µm C8 Polar 100g
B-100A-7-SU-C8P-1KG	Bulk Surf 100A 7µm C8 Polar 1kg
B-100A-7-SU-C8T-100G	Bulk Surf 100A 7µm C8 TriF 100g
B-100A-7-SU-C8T-1KG	Bulk Surf 100A 7µm C8 TriF 1kg
B-100A-7-SU-C8EX-100G	Bulk Surf 100A 7µm C8 Extreme 100g
B-100A-7-SU-C8EX-1KG	Bulk Surf 100A 7µm C8 Extreme 1kg
B-100A-7-SU-C8PT-100G	Bulk Surf 100A 7µm C8 Polar TriF 100g
B-100A-7-SU-C8PT-1KG	Bulk Surf 100A 7µm C8 Polar TriF 1kg
B-100A-7-SU-C8AQP-100G	Bulk Surf 100A 7µm C8 AQ Plus 100g
B-100A-7-SU-C8AQP-1KG	Bulk Surf 100A 7µm C8 AQ Plus 1kg
B-100A-7-SU-C8AQ-100G	Bulk Surf 100A 7µm C8 AQ 100g
B-100A-7-SU-C8AQ-1KG	Bulk Surf 100A 7µm C8 AQ 1kg
<b>C18</b>	
B-100A-7-SU-C18-100G	Bulk Surf 100A 7µm C18 100g
B-100A-7-SU-C18-1KG	Bulk Surf 100A 7µm C18 1kg
B-100A-7-SU-C18W-100G	Bulk Surf 100A 7µm C18 Wax 100g
B-100A-7-SU-C18W-1KG	Bulk Surf 100A 7µm C18 Wax 1kg
B-100A-7-SU-C18P-100G	Bulk Surf 100A 7µm C18 Polar 100g
B-100A-7-SU-C18P-1KG	Bulk Surf 100A 7µm C18 Polar 1kg
B-100A-7-SU-C18T-100G	Bulk Surf 100A 7µm C18 TriF 100g
B-100A-7-SU-C18T-1KG	Bulk Surf 100A 7µm C18 TriF 1kg
B-100A-7-SU-C18EX-100G	Bulk Surf 100A 7µm C18 Extreme 100g
B-100A-7-SU-C18EX-1KG	Bulk Surf 100A 7µm C18 Extreme 1kg
B-100A-7-SU-C18PT-100G	Bulk Surf 100A 7µm C18 Polar TriF 100g
B-100A-7-SU-C18PT-1KG	Bulk Surf 100A 7µm C18 Polar TriF 1kg
B-100A-7-SU-C18AQP-100G	Bulk Surf 100A 7µm C18 AQ Plus 100g
B-100A-7-SU-C18AQP-1KG	Bulk Surf 100A 7µm C18 AQ Plus 1kg
B-100A-7-SU-C18AQ-100G	Bulk Surf 100A 7µm C18 AQ 100g
B-100A-7-SU-C18AQ-1KG	Bulk Surf 100A 7µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100A-7-SU-C4-100G	Bulk Surf 100A 7µm C4 100g
B-100A-7-SU-C4-1KG	Bulk Surf 100A 7µm C4 1kg
B-100A-7-SU-C4W-100G	Bulk Surf 100A 7µm C4 Wax 100g
B-100A-7-SU-C4W-1KG	Bulk Surf 100A 7µm C4 Wax 1kg
B-100A-7-SU-C4P-100G	Bulk Surf 100A 7µm C4 Polar 100g
B-100A-7-SU-C4P-1KG	Bulk Surf 100A 7µm C4 Polar 1kg
B-100A-7-SU-C4T-100G	Bulk Surf 100A 7µm C4 TriF 100g
B-100A-7-SU-C4T-1KG	Bulk Surf 100A 7µm C4 TriF 1kg
B-100A-7-SU-C4EX-100G	Bulk Surf 100A 7µm C4 Extreme 100g
B-100A-7-SU-C4EX-1KG	Bulk Surf 100A 7µm C4 Extreme 1kg
B-100A-7-SU-C4PT-100G	Bulk Surf 100A 7µm C4 Polar TriF 100g
B-100A-7-SU-C4PT-1KG	Bulk Surf 100A 7µm C4 Polar TriF 1kg
B-100A-7-SU-C4AQP-100G	Bulk Surf 100A 7µm C4 AQ Plus 100g
B-100A-7-SU-C4AQP-1KG	Bulk Surf 100A 7µm C4 AQ Plus 1kg
B-100A-7-SU-C4AQ-100G	Bulk Surf 100A 7µm C4 AQ 100g
B-100A-7-SU-C4AQ-1KG	Bulk Surf 100A 7µm C4 AQ 1kg
<b>Autres</b>	
B-100A-7-SU-COOH-100G	Bulk Surf 100A 7µm COOH 100g
B-100A-7-SU-COOH-1KG	Bulk Surf 100A 7µm COOH 1kg
B-100A-7-SU-SO3H-100G	Bulk Surf 100A 7µm SO3H 100g
B-100A-7-SU-SO3H-1KG	Bulk Surf 100A 7µm SO3H 1kg
B-100A-7-SU-DNH-100G	Bulk Surf 100A 7µm DNH 100g
B-100A-7-SU-DNH-1KG	Bulk Surf 100A 7µm DNH 1kg
B-100A-7-SU-PEI-100G	Bulk Surf 100A 7µm PEI 100g
B-100A-7-SU-PEI-1KG	Bulk Surf 100A 7µm PEI 1kg
B-100A-7-SU-ARG-100G	Bulk Surf 100A 7µm ARG 100g
B-100A-7-SU-ARG-1KG	Bulk Surf 100A 7µm ARG 1kg
B-100A-7-SU-EP2TNE-100G	Bulk Surf 100A 7µm EP2TNE 100g
B-100A-7-SU-EP2TNE-1KG	Bulk Surf 100A 7µm EP2TNE 1kg
B-100A-7-SU-EP4TNE-100G	Bulk Surf 100A 7µm EP4TNE 100g
B-100A-7-SU-EP4TNE-1KG	Bulk Surf 100A 7µm EP4TNE 1kg
B-100A-7-SU-NH2-100G	Bulk Surf 100A 7µm NH2 100g
B-100A-7-SU-NH2-1KG	Bulk Surf 100A 7µm NH2 1kg
B-100A-7-SU-CN-100G	Bulk Surf 100A 7µm CN 100g
B-100A-7-SU-CN-1KG	Bulk Surf 100A 7µm CN 1kg
B-100A-7-SU-Phenyl-100G	Bulk Surf 100A 7µm Phenyl 100g
B-100A-7-SU-Phenyl-1KG	Bulk Surf 100A 7µm Phenyl 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-7-SU-Phenyl HE-100G	Bulk Surf 100A 7µm Phenyl HE 100g
B-100A-7-SU-Phenyl HE-1KG	Bulk Surf 100A 7µm Phenyl HE 1kg
B-100A-7-SU-Diol-100G	Bulk Surf 100A 7µm Diol 100g
B-100A-7-SU-Diol-1KG	Bulk Surf 100A 7µm Diol 1kg
<b>Surf 9 µm</b>	
<b>Particule Surf 9µm, pore 100A</b>	
B-100A-9-SU-Si-100G	Bulk Surf Silice 100A 9µm 100g
B-100A-9-SU-Si-1KG	Bulk Surf Silice 100A 9µm 1kg
<b>C8</b>	
B-100A-9-SU-C8-100G	Bulk Surf 100A 9µm C8 100g
B-100A-9-SU-C8-1KG	Bulk Surf 100A 9µm C8 1kg
B-100A-9-SU-C8W-100G	Bulk Surf 100A 9µm C8 Wax 100g
B-100A-9-SU-C8W-1KG	Bulk Surf 100A 9µm C8 Wax 1kg
B-100A-9-SU-C8P-100G	Bulk Surf 100A 9µm C8 Polar 100g
B-100A-9-SU-C8P-1KG	Bulk Surf 100A 9µm C8 Polar 1kg
B-100A-9-SU-C8T-100G	Bulk Surf 100A 9µm C8 TriF 100g
B-100A-9-SU-C8T-1KG	Bulk Surf 100A 9µm C8 TriF 1kg
B-100A-9-SU-C8EX-100G	Bulk Surf 100A 9µm C8 Extreme 100g
B-100A-9-SU-C8EX-1KG	Bulk Surf 100A 9µm C8 Extreme 1kg
B-100A-9-SU-C8PT-100G	Bulk Surf 100A 9µm C8 Polar TriF 100g
B-100A-9-SU-C8PT-1KG	Bulk Surf 100A 9µm C8 Polar TriF 1kg
B-100A-9-SU-C8AQP-100G	Bulk Surf 100A 9µm C8 AQ Plus 100g
B-100A-9-SU-C8AQP-1KG	Bulk Surf 100A 9µm C8 AQ Plus 1kg
B-100A-9-SU-C8AQ-100G	Bulk Surf 100A 9µm C8 AQ 100g
B-100A-9-SU-C8AQ-1KG	Bulk Surf 100A 9µm C8 AQ 1kg
<b>C18</b>	
B-100A-9-SU-C18-100G	Bulk Surf 100A 9µm C18 100g
B-100A-9-SU-C18-1KG	Bulk Surf 100A 9µm C18 1kg
B-100A-9-SU-C18W-100G	Bulk Surf 100A 9µm C18 Wax 100g
B-100A-9-SU-C18W-1KG	Bulk Surf 100A 9µm C18 Wax 1kg
B-100A-9-SU-C18P-100G	Bulk Surf 100A 9µm C18 Polar 100g
B-100A-9-SU-C18P-1KG	Bulk Surf 100A 9µm C18 Polar 1kg
B-100A-9-SU-C18T-100G	Bulk Surf 100A 9µm C18 TriF 100g
B-100A-9-SU-C18T-1KG	Bulk Surf 100A 9µm C18 TriF 1kg
B-100A-9-SU-C18EX-100G	Bulk Surf 100A 9µm C18 Extreme 100g
B-100A-9-SU-C18EX-1KG	Bulk Surf 100A 9µm C18 Extreme 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-9-SU-C18PT-100G	Bulk Surf 100A 9µm C18 Polar TriF 100g
B-100A-9-SU-C18PT-1KG	Bulk Surf 100A 9µm C18 Polar TriF 1kg
B-100A-9-SU-C18AQP-100G	Bulk Surf 100A 9µm C18 AQ Plus 100g
B-100A-9-SU-C18AQP-1KG	Bulk Surf 100A 9µm C18 AQ Plus 1kg
B-100A-9-SU-C18AQ-100G	Bulk Surf 100A 9µm C18 AQ 100g
B-100A-9-SU-C18AQ-1KG	Bulk Surf 100A 9µm C18 AQ 1kg
<b>C4</b>	
B-100A-9-SU-C4-100G	Bulk Surf 100A 9µm C4 100g
B-100A-9-SU-C4-1KG	Bulk Surf 100A 9µm C4 1kg
B-100A-9-SU-C4W-100G	Bulk Surf 100A 9µm C4 Wax 100g
B-100A-9-SU-C4W-1KG	Bulk Surf 100A 9µm C4 Wax 1kg
B-100A-9-SU-C4P-100G	Bulk Surf 100A 9µm C4 Polar 100g
B-100A-9-SU-C4P-1KG	Bulk Surf 100A 9µm C4 Polar 1kg
B-100A-9-SU-C4T-100G	Bulk Surf 100A 9µm C4 TriF 100g
B-100A-9-SU-C4T-1KG	Bulk Surf 100A 9µm C4 TriF 1kg
B-100A-9-SU-C4EX-100G	Bulk Surf 100A 9µm C4 Extreme 100g
B-100A-9-SU-C4EX-1KG	Bulk Surf 100A 9µm C4 Extreme 1kg
B-100A-9-SU-C4PT-100G	Bulk Surf 100A 9µm C4 Polar TriF 100g
B-100A-9-SU-C4PT-1KG	Bulk Surf 100A 9µm C4 Polar TriF 1kg
B-100A-9-SU-C4AQP-100G	Bulk Surf 100A 9µm C4 AQ Plus 100g
B-100A-9-SU-C4AQP-1KG	Bulk Surf 100A 9µm C4 AQ Plus 1kg
B-100A-9-SU-C4AQ-100G	Bulk Surf 100A 9µm C4 AQ 100g
B-100A-9-SU-C4AQ-1KG	Bulk Surf 100A 9µm C4 AQ 1kg
<b>Autres</b>	
B-100A-9-SU-COOH-100G	Bulk Surf 100A 9µm COOH 100g
B-100A-9-SU-COOH-1KG	Bulk Surf 100A 9µm COOH 1kg
B-100A-9-SU-SO3H-100G	Bulk Surf 100A 9µm SO3H 100g
B-100A-9-SU-SO3H-1KG	Bulk Surf 100A 9µm SO3H 1kg
B-100A-9-SU-DNH-100G	Bulk Surf 100A 9µm DNH 100g
B-100A-9-SU-DNH-1KG	Bulk Surf 100A 9µm DNH 1kg
B-100A-9-SU-PEI-100G	Bulk Surf 100A 9µm PEI 100g
B-100A-9-SU-PEI-1KG	Bulk Surf 100A 9µm PEI 1kg
B-100A-9-SU-ARG-100G	Bulk Surf 100A 9µm ARG 100g
B-100A-9-SU-ARG-1KG	Bulk Surf 100A 9µm ARG 1kg
B-100A-9-SU-EP2TNE-100G	Bulk Surf 100A 9µm EP2TNE 100g
B-100A-9-SU-EP2TNE-1KG	Bulk Surf 100A 9µm EP2TNE 1kg
B-100A-9-SU-EP4TNE-100G	Bulk Surf 100A 9µm EP4TNE 100g
B-100A-9-SU-EP4TNE-1KG	Bulk Surf 100A 9µm EP4TNE 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-9-SU-NH2-100G	Bulk Surf 100A 9µm NH2 100g
B-100A-9-SU-NH2-1KG	Bulk Surf 100A 9µm NH2 1kg
B-100A-9-SU-CN-100G	Bulk Surf 100A 9µm CN 100g
B-100A-9-SU-CN-1KG	Bulk Surf 100A 9µm CN 1kg
B-100A-9-SU-Phenyl-100G	Bulk Surf 100A 9µm Phenyl 100g
B-100A-9-SU-Phenyl-1KG	Bulk Surf 100A 9µm Phenyl 1kg
B-100A-9-SU-Phenyl HE-100G	Bulk Surf 100A 9µm Phenyl HE 100g
B-100A-9-SU-Phenyl HE-1KG	Bulk Surf 100A 9µm Phenyl HE 1kg
B-100A-9-SU-Diol-100G	Bulk Surf 100A 9µm Diol 100g
B-100A-9-SU-Diol-1KG	Bulk Surf 100A 9µm Diol 1kg
<b>Surf 10 µm</b>	
<b>Particule Surf 10µm, pore 100A</b>	
B-100A-10-SU-Si-100G	Bulk Surf Silice 100A 10µm 100g
B-100A-10-SU-Si-1KG	Bulk Surf Silice 100A 10µm 1kg
<b>C8</b>	
B-100A-10-SU-C8-100G	Bulk Surf 100A 10µm C8 100g
B-100A-10-SU-C8-1KG	Bulk Surf 100A 10µm C8 1kg
B-100A-10-SU-C8W-100G	Bulk Surf 100A 10µm C8 Wax 100g
B-100A-10-SU-C8W-1KG	Bulk Surf 100A 10µm C8 Wax 1kg
B-100A-10-SU-C8P-100G	Bulk Surf 100A 10µm C8 Polar 100g
B-100A-10-SU-C8P-1KG	Bulk Surf 100A 10µm C8 Polar 1kg
B-100A-10-SU-C8T-100G	Bulk Surf 100A 10µm C8 TriF 100g
B-100A-10-SU-C8T-1KG	Bulk Surf 100A 10µm C8 TriF 1kg
B-100A-10-SU-C8EX-100G	Bulk Surf 100A 10µm C8 Extreme 100g
B-100A-10-SU-C8EX-1KG	Bulk Surf 100A 10µm C8 Extreme 1kg
B-100A-10-SU-C8PT-100G	Bulk Surf 100A 10µm C8 Polar TriF 100g
B-100A-10-SU-C8PT-1KG	Bulk Surf 100A 10µm C8 Polar TriF 1kg
B-100A-10-SU-C8AQP-100G	Bulk Surf 100A 10µm C8 AQ Plus 100g
B-100A-10-SU-C8AQP-1KG	Bulk Surf 100A 10µm C8 AQ Plus 1kg
B-100A-10-SU-C8AQ-100G	Bulk Surf 100A 10µm C8 AQ 100g
B-100A-10-SU-C8AQ-1KG	Bulk Surf 100A 10µm C8 AQ 1kg
<b>C18</b>	
B-100A-10-SU-C18-100G	Bulk Surf 100A 10µm C18 100g
B-100A-10-SU-C18-1KG	Bulk Surf 100A 10µm C18 1kg
B-100A-10-SU-C18W-100G	Bulk Surf 100A 10µm C18 Wax 100g
B-100A-10-SU-C18W-1KG	Bulk Surf 100A 10µm C18 Wax 1kg
B-100A-10-SU-C18P-100G	Bulk Surf 100A 10µm C18 Polar 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-10-SU-C18P-1KG	Bulk Surf 100A 10µm C18 Polar 1kg
B-100A-10-SU-C18T-100G	Bulk Surf 100A 10µm C18 TriF 100g
B-100A-10-SU-C18T-1KG	Bulk Surf 100A 10µm C18 TriF 1kg
B-100A-10-SU-C18EX-100G	Bulk Surf 100A 10µm C18 Extreme 100g
B-100A-10-SU-C18EX-1KG	Bulk Surf 100A 10µm C18 Extreme 1kg
B-100A-10-SU-C18PT-100G	Bulk Surf 100A 10µm C18 Polar TriF 100g
B-100A-10-SU-C18PT-1KG	Bulk Surf 100A 10µm C18 Polar TriF 1kg
B-100A-10-SU-C18AQP-100G	Bulk Surf 100A 10µm C18 AQ Plus 100g
B-100A-10-SU-C18AQP-1KG	Bulk Surf 100A 10µm C18 AQ Plus 1kg
B-100A-10-SU-C18AQ-100G	Bulk Surf 100A 10µm C18 AQ 100g
B-100A-10-SU-C18AQ-1KG	Bulk Surf 100A 10µm C18 AQ 1kg
<b>C4</b>	
B-100A-10-SU-C4-100G	Bulk Surf 100A 10µm C4 100g
B-100A-10-SU-C4-1KG	Bulk Surf 100A 10µm C4 1kg
B-100A-10-SU-C4W-100G	Bulk Surf 100A 10µm C4 Wax 100g
B-100A-10-SU-C4W-1KG	Bulk Surf 100A 10µm C4 Wax 1kg
B-100A-10-SU-C4P-100G	Bulk Surf 100A 10µm C4 Polar 100g
B-100A-10-SU-C4P-1KG	Bulk Surf 100A 10µm C4 Polar 1kg
B-100A-10-SU-C4T-100G	Bulk Surf 100A 10µm C4 TriF 100g
B-100A-10-SU-C4T-1KG	Bulk Surf 100A 10µm C4 TriF 1kg
B-100A-10-SU-C4EX-100G	Bulk Surf 100A 10µm C4 Extreme 100g
B-100A-10-SU-C4EX-1KG	Bulk Surf 100A 10µm C4 Extreme 1kg
B-100A-10-SU-C4PT-100G	Bulk Surf 100A 10µm C4 Polar TriF 100g
B-100A-10-SU-C4PT-1KG	Bulk Surf 100A 10µm C4 Polar TriF 1kg
B-100A-10-SU-C4AQP-100G	Bulk Surf 100A 10µm C4 AQ Plus 100g
B-100A-10-SU-C4AQP-1KG	Bulk Surf 100A 10µm C4 AQ Plus 1kg
B-100A-10-SU-C4AQ-100G	Bulk Surf 100A 10µm C4 AQ 100g
B-100A-10-SU-C4AQ-1KG	Bulk Surf 100A 10µm C4 AQ 1kg
<b>Autres</b>	
B-100A-10-SU-COOH-100G	Bulk Surf 100A 10µm COOH 100g
B-100A-10-SU-COOH-1KG	Bulk Surf 100A 10µm COOH 1kg
B-100A-10-SU-SO3H-100G	Bulk Surf 100A 10µm SO3H 100g
B-100A-10-SU-SO3H-1KG	Bulk Surf 100A 10µm SO3H 1kg
B-100A-10-SU-DNH-100G	Bulk Surf 100A 10µm DNH 100g
B-100A-10-SU-DNH-1KG	Bulk Surf 100A 10µm DNH 1kg
B-100A-10-SU-PEI-100G	Bulk Surf 100A 10µm PEI 100g
B-100A-10-SU-PEI-1KG	Bulk Surf 100A 10µm PEI 1kg
B-100A-10-SU-ARG-100G	Bulk Surf 100A 10µm ARG 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-10-SU-ARG-1KG	Bulk Surf 100A 10µm ARG 1kg
B-100A-10-SU-EP2TNE-100G	Bulk Surf 100A 10µm EP2TNE 100g
B-100A-10-SU-EP2TNE-1KG	Bulk Surf 100A 10µm EP2TNE 1kg
B-100A-10-SU-EP4TNE-100G	Bulk Surf 100A 10µm EP4TNE 100g
B-100A-10-SU-EP4TNE-1KG	Bulk Surf 100A 10µm EP4TNE 1kg
B-100A-10-SU-NH2-100G	Bulk Surf 100A 10µm NH2 100g
B-100A-10-SU-NH2-1KG	Bulk Surf 100A 10µm NH2 1kg
B-100A-10-SU-CN-100G	Bulk Surf 100A 10µm CN 100g
B-100A-10-SU-CN-1KG	Bulk Surf 100A 10µm CN 1kg
B-100A-10-SU-Phenyl-100G	Bulk Surf 100A 10µm Phenyl 100g
B-100A-10-SU-Phenyl-1KG	Bulk Surf 100A 10µm Phenyl 1kg
B-100A-10-SU-Phenyl HE-100G	Bulk Surf 100A 10µm Phenyl HE 100g
B-100A-10-SU-Phenyl HE-1KG	Bulk Surf 100A 10µm Phenyl HE 1kg
B-100A-10-SU-Diol-100G	Bulk Surf 100A 10µm Diol 100g
B-100A-10-SU-Diol-1KG	Bulk Surf 100A 10µm Diol 1kg
<b>Particule Surf 10µm, pore 150A</b>	
B-150A-10-SU-Si-100G	Bulk Surf 150A 10µm 100g
B-150A-10-SU-Si-1KG	Bulk Surf 150A 10µm 1kg
<b>C8</b>	
B-150A-10-SU-C8-100G	Bulk Surf 150A 10µm C8 100g
B-150A-10-SU-C8-1KG	Bulk Surf 150A 10µm C8 1kg
B-150A-10-SU-C8W-100G	Bulk Surf 150A 10µm C8 Wax 100g
B-150A-10-SU-C8W-1KG	Bulk Surf 150A 10µm C8 Wax 1kg
B-150A-10-SU-C8P-100G	Bulk Surf 150A 10µm C8 Polar 100g
B-150A-10-SU-C8P-1KG	Bulk Surf 150A 10µm C8 Polar 1kg
B-150A-10-SU-C8T-100G	Bulk Surf 150A 10µm C8 TriF 100g
B-150A-10-SU-C8T-1KG	Bulk Surf 150A 10µm C8 TriF 1kg
B-150A-10-SU-C8EX-100G	Bulk Surf 150A 10µm C8 Extreme 100g
B-150A-10-SU-C8EX-1KG	Bulk Surf 150A 10µm C8 Extreme 1kg
B-150A-10-SU-C8PT-100G	Bulk Surf 150A 10µm C8 Polar TriF 100g
B-150A-10-SU-C8PT-1KG	Bulk Surf 150A 10µm C8 Polar TriF 1kg
B-150A-10-SU-C8AQP-100G	Bulk Surf 150A 10µm C8 AQ Plus 100g
B-150A-10-SU-C8AQP-1KG	Bulk Surf 150A 10µm C8 AQ Plus 1kg
B-150A-10-SU-C8AQ-100G	Bulk Surf 150A 10µm C8 AQ 100g
B-150A-10-SU-C8AQ-1KG	Bulk Surf 150A 10µm C8 AQ 1kg
<b>C18</b>	

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-10-SU-C18-100G	Bulk Surf 150A 10µm C18 100g
B-150A-10-SU-C18-1KG	Bulk Surf 150A 10µm C18 1kg
B-150A-10-SU-C18W-100G	Bulk Surf 150A 10µm C18 Wax 100g
B-150A-10-SU-C18W-1KG	Bulk Surf 150A 10µm C18 Wax 1kg
B-150A-10-SU-C18P-100G	Bulk Surf 150A 10µm C18 Polar 100g
B-150A-10-SU-C18P-1KG	Bulk Surf 150A 10µm C18 Polar 1kg
B-150A-10-SU-C18T-100G	Bulk Surf 150A 10µm C18 TriF 100g
B-150A-10-SU-C18T-1KG	Bulk Surf 150A 10µm C18 TriF 1kg
B-150A-10-SU-C18EX-100G	Bulk Surf 150A 10µm C18 Extreme 100g
B-150A-10-SU-C18EX-1KG	Bulk Surf 150A 10µm C18 Extreme 1kg
B-150A-10-SU-C18PT-100G	Bulk Surf 150A 10µm C18 Polar TriF 100g
B-150A-10-SU-C18PT-1KG	Bulk Surf 150A 10µm C18 Polar TriF 1kg
B-150A-10-SU-C18AQP-100G	Bulk Surf 150A 10µm C18 AQ Plus 100g
B-150A-10-SU-C18AQP-1KG	Bulk Surf 150A 10µm C18 AQ Plus 1kg
B-150A-10-SU-C18AQ-100G	Bulk Surf 150A 10µm C18 AQ 100g
B-150A-10-SU-C18AQ-1KG	Bulk Surf 150A 10µm C18 AQ 1kg
<b>C4</b>	
B-150A-10-SU-C4-100G	Bulk Surf 150A 10µm C4 100g
B-150A-10-SU-C4-1KG	Bulk Surf 150A 10µm C4 1kg
B-150A-10-SU-C4W-100G	Bulk Surf 150A 10µm C4 Wax 100g
B-150A-10-SU-C4W-1KG	Bulk Surf 150A 10µm C4 Wax 1kg
B-150A-10-SU-C4P-100G	Bulk Surf 150A 10µm C4 Polar 100g
B-150A-10-SU-C4P-1KG	Bulk Surf 150A 10µm C4 Polar 1kg
B-150A-10-SU-C4T-100G	Bulk Surf 150A 10µm C4 TriF 100g
B-150A-10-SU-C4T-1KG	Bulk Surf 150A 10µm C4 TriF 1kg
B-150A-10-SU-C4EX-100G	Bulk Surf 150A 10µm C4 Extreme 100g
B-150A-10-SU-C4EX-1KG	Bulk Surf 150A 10µm C4 Extreme 1kg
B-150A-10-SU-C4PT-100G	Bulk Surf 150A 10µm C4 Polar TriF 100g
B-150A-10-SU-C4PT-1KG	Bulk Surf 150A 10µm C4 Polar TriF 1kg
B-150A-10-SU-C4AQP-100G	Bulk Surf 150A 10µm C4 AQ Plus 100g
B-150A-10-SU-C4AQP-1KG	Bulk Surf 150A 10µm C4 AQ Plus 1kg
B-150A-10-SU-C4AQ-100G	Bulk Surf 150A 10µm C4 AQ 100g
B-150A-10-SU-C4AQ-1KG	Bulk Surf 150A 10µm C4 AQ 1kg
Autres	
B-150A-10-SU-COOH-100G	Bulk Surf 150A 10µm COOH 100g
B-150A-10-SU-COOH-1KG	Bulk Surf 150A 10µm COOH 1kg
B-150A-10-SU-SO3H-100G	Bulk Surf 150A 10µm SO3H 100g
B-150A-10-SU-SO3H-1KG	Bulk Surf 150A 10µm SO3H 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-10-SU-DNH-100G	Bulk Surf 150A 10µm DNH 100g
B-150A-10-SU-DNH-1KG	Bulk Surf 150A 10µm DNH 1kg
B-150A-10-SU-PEI-100G	Bulk Surf 150A 10µm PEI 100g
B-150A-10-SU-PEI-1KG	Bulk Surf 150A 10µm PEI 1kg
B-150A-10-SU-ARG-100G	Bulk Surf 150A 10µm ARG 100g
B-150A-10-SU-ARG-1KG	Bulk Surf 150A 10µm ARG 1kg
B-150A-10-SU-EP2TNE-100G	Bulk Surf 150A 10µm EP2TNE 100g
B-150A-10-SU-EP2TNE-1KG	Bulk Surf 150A 10µm EP2TNE 1kg
B-150A-10-SU-EP4TNE-100G	Bulk Surf 150A 10µm EP4TNE 100g
B-150A-10-SU-EP4TNE-1KG	Bulk Surf 150A 10µm EP4TNE 1kg
B-150A-10-SU-NH2-100G	Bulk Surf 150A 10µm NH2 100g
B-150A-10-SU-NH2-1KG	Bulk Surf 150A 10µm NH2 1kg
B-150A-10-SU-CN-100G	Bulk Surf 150A 10µm CN 100g
B-150A-10-SU-CN-1KG	Bulk Surf 150A 10µm CN 1kg
B-150A-10-SU-Phenyl-100G	Bulk Surf 150A 10µm Phenyl 100g
B-150A-10-SU-Phenyl-1KG	Bulk Surf 150A 10µm Phenyl 1kg
B-150A-10-SU-Phenyl HE-100G	Bulk Surf 150A 10µm Phenyl HE 100g
B-150A-10-SU-Phenyl HE-1KG	Bulk Surf 150A 10µm Phenyl HE 1kg
B-150A-10-SU-Diol-100G	Bulk Surf 150A 10µm Diol 100g
B-150A-10-SU-Diol-1KG	Bulk Surf 150A 10µm Diol 1kg
<b>Particule Surf 10µm, pore 200A</b>	
B-200A-10-SU-Si-100G	Bulk Surf 200A 10µm 100g
B-200A-10-SU-Si-1KG	Bulk Surf 200A 10µm 1kg
<b>C8</b>	
B-200A-10-SU-C8-100G	Bulk Surf 200A 10µm C8 100g
B-200A-10-SU-C8-1KG	Bulk Surf 200A 10µm C8 1kg
B-200A-10-SU-C8W-100G	Bulk Surf 200A 10µm C8 Wax 100g
B-200A-10-SU-C8W-1KG	Bulk Surf 200A 10µm C8 Wax 1kg
B-200A-10-SU-C8P-100G	Bulk Surf 200A 10µm C8 Polar 100g
B-200A-10-SU-C8P-1KG	Bulk Surf 200A 10µm C8 Polar 1kg
B-200A-10-SU-C8T-100G	Bulk Surf 200A 10µm C8 TriF 100g
B-200A-10-SU-C8T-1KG	Bulk Surf 200A 10µm C8 TriF 1kg
B-200A-10-SU-C8EX-100G	Bulk Surf 200A 10µm C8 Extreme 100g
B-200A-10-SU-C8EX-1KG	Bulk Surf 200A 10µm C8 Extreme 1kg
B-200A-10-SU-C8PT-100G	Bulk Surf 200A 10µm C8 Polar TriF 100g
B-200A-10-SU-C8PT-1KG	Bulk Surf 200A 10µm C8 Polar TriF 1kg
B-200A-10-SU-C8AQP-100G	Bulk Surf 200A 10µm C8 AQ Plus 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-10-SU-C8AQP-1KG	Bulk Surf 200A 10µm C8 AQ Plus 1kg
B-200A-10-SU-C8AQ-100G	Bulk Surf 200A 10µm C8 AQ 100g
B-200A-10-SU-C8AQ-1KG	Bulk Surf 200A 10µm C8 AQ 1kg
<b>C18</b>	
B-200A-10-SU-C18-100G	Bulk Surf 200A 10µm C18 100g
B-200A-10-SU-C18-1KG	Bulk Surf 200A 10µm C18 1kg
B-200A-10-SU-C18W-100G	Bulk Surf 200A 10µm C18 Wax 100g
B-200A-10-SU-C18W-1KG	Bulk Surf 200A 10µm C18 Wax 1kg
B-200A-10-SU-C18P-100G	Bulk Surf 200A 10µm C18 Polar 100g
B-200A-10-SU-C18P-1KG	Bulk Surf 200A 10µm C18 Polar 1kg
B-200A-10-SU-C18T-100G	Bulk Surf 200A 10µm C18 TriF 100g
B-200A-10-SU-C18T-1KG	Bulk Surf 200A 10µm C18 TriF 1kg
B-200A-10-SU-C18EX-100G	Bulk Surf 200A 10µm C18 Extreme 100g
B-200A-10-SU-C18EX-1KG	Bulk Surf 200A 10µm C18 Extreme 1kg
B-200A-10-SU-C18PT-100G	Bulk Surf 200A 10µm C18 Polar TriF 100g
B-200A-10-SU-C18PT-1KG	Bulk Surf 200A 10µm C18 Polar TriF 1kg
B-200A-10-SU-C18AQP-100G	Bulk Surf 200A 10µm C18 AQ Plus 100g
B-200A-10-SU-C18AQP-1KG	Bulk Surf 200A 10µm C18 AQ Plus 1kg
B-200A-10-SU-C18AQ-100G	Bulk Surf 200A 10µm C18 AQ 100g
B-200A-10-SU-C18AQ-1KG	Bulk Surf 200A 10µm C18 AQ 1kg
<b>C4</b>	
B-200A-10-SU-C4-100G	Bulk Surf 200A 10µm C4 100g
B-200A-10-SU-C4-1KG	Bulk Surf 200A 10µm C4 1kg
B-200A-10-SU-C4W-100G	Bulk Surf 200A 10µm C4 Wax 100g
B-200A-10-SU-C4W-1KG	Bulk Surf 200A 10µm C4 Wax 1kg
B-200A-10-SU-C4P-100G	Bulk Surf 200A 10µm C4 Polar 100g
B-200A-10-SU-C4P-1KG	Bulk Surf 200A 10µm C4 Polar 1kg
B-200A-10-SU-C4T-100G	Bulk Surf 200A 10µm C4 TriF 100g
B-200A-10-SU-C4T-1KG	Bulk Surf 200A 10µm C4 TriF 1kg
B-200A-10-SU-C4EX-100G	Bulk Surf 200A 10µm C4 Extreme 100g
B-200A-10-SU-C4EX-1KG	Bulk Surf 200A 10µm C4 Extreme 1kg
B-200A-10-SU-C4PT-100G	Bulk Surf 200A 10µm C4 Polar TriF 100g
B-200A-10-SU-C4PT-1KG	Bulk Surf 200A 10µm C4 Polar TriF 1kg
B-200A-10-SU-C4AQP-100G	Bulk Surf 200A 10µm C4 AQ Plus 100g
B-200A-10-SU-C4AQP-1KG	Bulk Surf 200A 10µm C4 AQ Plus 1kg
B-200A-10-SU-C4AQ-100G	Bulk Surf 200A 10µm C4 AQ 100g
B-200A-10-SU-C4AQ-1KG	Bulk Surf 200A 10µm C4 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Autres</b>	
B-200A-10-SU-COOH-100G	Bulk Surf 200A 10µm COOH 100g
B-200A-10-SU-COOH-1KG	Bulk Surf 200A 10µm COOH 1kg
B-200A-10-SU-SO3H-100G	Bulk Surf 200A 10µm SO3H 100g
B-200A-10-SU-SO3H-1KG	Bulk Surf 200A 10µm SO3H 1kg
B-200A-10-SU-DNH-100G	Bulk Surf 200A 10µm DNH 100g
B-200A-10-SU-DNH-1KG	Bulk Surf 200A 10µm DNH 1kg
B-200A-10-SU-PEI-100G	Bulk Surf 200A 10µm PEI 100g
B-200A-10-SU-PEI-1KG	Bulk Surf 200A 10µm PEI 1kg
B-200A-10-SU-ARG-100G	Bulk Surf 200A 10µm ARG 100g
B-200A-10-SU-ARG-1KG	Bulk Surf 200A 10µm ARG 1kg
B-200A-10-SU-EP2TNE-100G	Bulk Surf 200A 10µm EP2TNE 100g
B-200A-10-SU-EP2TNE-1KG	Bulk Surf 200A 10µm EP2TNE 1kg
B-200A-10-SU-EP4TNE-100G	Bulk Surf 200A 10µm EP4TNE 100g
B-200A-10-SU-EP4TNE-1KG	Bulk Surf 200A 10µm EP4TNE 1kg
B-200A-10-SU-NH2-100G	Bulk Surf 200A 10µm NH2 100g
B-200A-10-SU-NH2-1KG	Bulk Surf 200A 10µm NH2 1kg
B-200A-10-SU-CN-100G	Bulk Surf 200A 10µm CN 100g
B-200A-10-SU-CN-1KG	Bulk Surf 200A 10µm CN 1kg
B-200A-10-SU-Phenyl-100G	Bulk Surf 200A 10µm Phenyl 100g
B-200A-10-SU-Phenyl-1KG	Bulk Surf 200A 10µm Phenyl 1kg
B-200A-10-SU-Phenyl HE-100G	Bulk Surf 200A 10µm Phenyl HE 100g
B-200A-10-SU-Phenyl HE-1KG	Bulk Surf 200A 10µm Phenyl HE 1kg
B-200A-10-SU-Diol-100G	Bulk Surf 200A 10µm Diol 100g
B-200A-10-SU-Diol-1KG	Bulk Surf 200A 10µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Surf 10µm, pore 300A</b>	
B-300A-10-SU-Si-100G	Bulk Surf 300A 10µm 100g
B-300A-10-SU-Si-1KG	Bulk Surf 300A 10µm 1kg
<b>C8</b>	
B-300A-10-SU-C8-100G	Bulk Surf 300A 10µm C8 100g
B-300A-10-SU-C8-1KG	Bulk Surf 300A 10µm C8 1kg
B-300A-10-SU-C8W-100G	Bulk Surf 300A 10µm C8 Wax 100g
B-300A-10-SU-C8W-1KG	Bulk Surf 300A 10µm C8 Wax 1kg
B-300A-10-SU-C8P-100G	Bulk Surf 300A 10µm C8 Polar 100g
B-300A-10-SU-C8P-1KG	Bulk Surf 300A 10µm C8 Polar 1kg
B-300A-10-SU-C8T-100G	Bulk Surf 300A 10µm C8 TriF 100g
B-300A-10-SU-C8T-1KG	Bulk Surf 300A 10µm C8 TriF 1kg
B-300A-10-SU-C8EX-100G	Bulk Surf 300A 10µm C8 Extreme 100g
B-300A-10-SU-C8EX-1KG	Bulk Surf 300A 10µm C8 Extreme 1kg
B-300A-10-SU-C8PT-100G	Bulk Surf 300A 10µm C8 Polar TriF 100g
B-300A-10-SU-C8PT-1KG	Bulk Surf 300A 10µm C8 Polar TriF 1kg
B-300A-10-SU-C8AQP-100G	Bulk Surf 300A 10µm C8 AQ Plus 100g
B-300A-10-SU-C8AQP-1KG	Bulk Surf 300A 10µm C8 AQ Plus 1kg
B-300A-10-SU-C8AQ-100G	Bulk Surf 300A 10µm C8 AQ 100g
B-300A-10-SU-C8AQ-1KG	Bulk Surf 300A 10µm C8 AQ 1kg
<b>C18</b>	
B-300A-10-SU-C18-100G	Bulk Surf 300A 10µm C18 100g
B-300A-10-SU-C18-1KG	Bulk Surf 300A 10µm C18 1kg
B-300A-10-SU-C18W-100G	Bulk Surf 300A 10µm C18 Wax 100g
B-300A-10-SU-C18W-1KG	Bulk Surf 300A 10µm C18 Wax 1kg
B-300A-10-SU-C18P-100G	Bulk Surf 300A 10µm C18 Polar 100g
B-300A-10-SU-C18P-1KG	Bulk Surf 300A 10µm C18 Polar 1kg
B-300A-10-SU-C18T-100G	Bulk Surf 300A 10µm C18 TriF 100g
B-300A-10-SU-C18T-1KG	Bulk Surf 300A 10µm C18 TriF 1kg
B-300A-10-SU-C18EX-100G	Bulk Surf 300A 10µm C18 Extreme 100g
B-300A-10-SU-C18EX-1KG	Bulk Surf 300A 10µm C18 Extreme 1kg
B-300A-10-SU-C18PT-100G	Bulk Surf 300A 10µm C18 Polar TriF 100g
B-300A-10-SU-C18PT-1KG	Bulk Surf 300A 10µm C18 Polar TriF 1kg
B-300A-10-SU-C18AQP-100G	Bulk Surf 300A 10µm C18 AQ Plus 100g
B-300A-10-SU-C18AQP-1KG	Bulk Surf 300A 10µm C18 AQ Plus 1kg
B-300A-10-SU-C18AQ-100G	Bulk Surf 300A 10µm C18 AQ 100g
B-300A-10-SU-C18AQ-1KG	Bulk Surf 300A 10µm C18 AQ 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-300A-10-SU-C4-100G	Bulk Surf 300A 10µm C4 100g
B-300A-10-SU-C4-1KG	Bulk Surf 300A 10µm C4 1kg
B-300A-10-SU-C4W-100G	Bulk Surf 300A 10µm C4 Wax 100g
B-300A-10-SU-C4W-1KG	Bulk Surf 300A 10µm C4 Wax 1kg
B-300A-10-SU-C4P-100G	Bulk Surf 300A 10µm C4 Polar 100g
B-300A-10-SU-C4P-1KG	Bulk Surf 300A 10µm C4 Polar 1kg
B-300A-10-SU-C4T-100G	Bulk Surf 300A 10µm C4 TriF 100g
B-300A-10-SU-C4T-1KG	Bulk Surf 300A 10µm C4 TriF 1kg
B-300A-10-SU-C4EX-100G	Bulk Surf 300A 10µm C4 Extreme 100g
B-300A-10-SU-C4EX-1KG	Bulk Surf 300A 10µm C4 Extreme 1kg
B-300A-10-SU-C4PT-100G	Bulk Surf 300A 10µm C4 Polar TriF 100g
B-300A-10-SU-C4PT-1KG	Bulk Surf 300A 10µm C4 Polar TriF 1kg
B-300A-10-SU-C4AQP-100G	Bulk Surf 300A 10µm C4 AQ Plus 100g
B-300A-10-SU-C4AQP-1KG	Bulk Surf 300A 10µm C4 AQ Plus 1kg
B-300A-10-SU-C4AQ-100G	Bulk Surf 300A 10µm C4 AQ 100g
B-300A-10-SU-C4AQ-1KG	Bulk Surf 300A 10µm C4 AQ 1kg
<b>Autres</b>	
B-300A-10-SU-COOH-100G	Bulk Surf 300A 10µm COOH 100g
B-300A-10-SU-COOH-1KG	Bulk Surf 300A 10µm COOH 1kg
B-300A-10-SU-SO3H-100G	Bulk Surf 300A 10µm SO3H 100g
B-300A-10-SU-SO3H-1KG	Bulk Surf 300A 10µm SO3H 1kg
B-300A-10-SU-DNH-100G	Bulk Surf 300A 10µm DNH 100g
B-300A-10-SU-DNH-1KG	Bulk Surf 300A 10µm DNH 1kg
B-300A-10-SU-PEI-100G	Bulk Surf 300A 10µm PEI 100g
B-300A-10-SU-PEI-1KG	Bulk Surf 300A 10µm PEI 1kg
B-300A-10-SU-ARG-100G	Bulk Surf 300A 10µm ARG 100g
B-300A-10-SU-ARG-1KG	Bulk Surf 300A 10µm ARG 1kg
B-300A-10-SU-EP2TNE-100G	Bulk Surf 300A 10µm EP2TNE 100g
B-300A-10-SU-EP2TNE-1KG	Bulk Surf 300A 10µm EP2TNE 1kg
B-300A-10-SU-EP4TNE-100G	Bulk Surf 300A 10µm EP4TNE 100g
B-300A-10-SU-EP4TNE-1KG	Bulk Surf 300A 10µm EP4TNE 1kg
B-300A-10-SU-NH2-100G	Bulk Surf 300A 10µm NH2 100g
B-300A-10-SU-NH2-1KG	Bulk Surf 300A 10µm NH2 1kg
B-300A-10-SU-CN-100G	Bulk Surf 300A 10µm CN 100g
B-300A-10-SU-CN-1KG	Bulk Surf 300A 10µm CN 1kg
B-300A-10-SU-Phenyl-100G	Bulk Surf 300A 10µm Phenyl 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-10-SU-Phenyl-1KG	Bulk Surf 300A 10µm Phenyl 1kg
B-300A-10-SU-Phenyl HE-100G	Bulk Surf 300A 10µm Phenyl HE 100g
B-300A-10-SU-Phenyl HE-1KG	Bulk Surf 300A 10µm Phenyl HE 1kg
B-300A-10-SU-Diol-100G	Bulk Surf 300A 10µm Diol 100g
B-300A-10-SU-Diol-1KG	Bulk Surf 300A 10µm Diol 1kg
<b>Coral 5µm</b>	
<b>Particule Coral 5, Pore 70A</b>	
B-70A-5-CO-Si-100G	Bulk Coral Silice 70A 5µm 100g
B-70A-5-CO-Si-1KG	Bulk Coral Silice 70A 5µm 1kg
<b>C8</b>	
B-70A-5-CO-C8-100G	Bulk Coral 70A 5µm C8 100g
B-70A-5-CO-C8-1KG	Bulk Coral 70A 5µm C8 1kg
B-70A-5-CO-C8W-100G	Bulk Coral 70A 5µm C8 Wax 100g
B-70A-5-CO-C8W-1KG	Bulk Coral 70A 5µm C8 Wax 1kg
B-70A-5-CO-C8P-100G	Bulk Coral 70A 5µm C8 Polar 100g
B-70A-5-CO-C8P-1KG	Bulk Coral 70A 5µm C8 Polar 1kg
B-70A-5-CO-C8T-100G	Bulk Coral 70A 5µm C8 TriF 100g
B-70A-5-CO-C8T-1KG	Bulk Coral 70A 5µm C8 TriF 1kg
B-70A-5-CO-C8EX-100G	Bulk Coral 70A 5µm C8 Extreme 100g
B-70A-5-CO-C8EX-1KG	Bulk Coral 70A 5µm C8 Extreme 1kg
B-70A-5-CO-C8PT-100G	Bulk Coral 70A 5µm C8 Polar TriF 100g
B-70A-5-CO-C8PT-1KG	Bulk Coral 70A 5µm C8 Polar TriF 1kg
B-70A-5-CO-C8AQP-100G	Bulk Coral 70A 5µm C8 AQ Plus 100g
B-70A-5-CO-C8AQP-1KG	Bulk Coral 70A 5µm C8 AQ Plus 1kg
B-70A-5-CO-C8AQ-100G	Bulk Coral 70A 5µm C8 AQ 100g
B-70A-5-CO-C8AQ-1KG	Bulk Coral 70A 5µm C8 AQ 1kg
<b>C18</b>	
B-70A-5-CO-C18-100G	Bulk Coral 70A 5µm C18 100g
B-70A-5-CO-C18-1KG	Bulk Coral 70A 5µm C18 1kg
B-70A-5-CO-C18W-100G	Bulk Coral 70A 5µm C18 Wax 100g
B-70A-5-CO-C18W-1KG	Bulk Coral 70A 5µm C18 Wax 1kg
B-70A-5-CO-C18P-100G	Bulk Coral 70A 5µm C18 Polar 100g
B-70A-5-CO-C18P-1KG	Bulk Coral 70A 5µm C18 Polar 1kg
B-70A-5-CO-C18T-100G	Bulk Coral 70A 5µm C18 TriF 100g
B-70A-5-CO-C18T-1KG	Bulk Coral 70A 5µm C18 TriF 1kg
B-70A-5-CO-C18EX-100G	Bulk Coral 70A 5µm C18 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-5-CO-C18EX-1KG	Bulk Coral 70A 5µm C18 Extreme 1kg
B-70A-5-CO-C18PT-100G	Bulk Coral 70A 5µm C18 Polar TriF 100g
B-70A-5-CO-C18PT-1KG	Bulk Coral 70A 5µm C18 Polar TriF 1kg
B-70A-5-CO-C18AQP-100G	Bulk Coral 70A 5µm C18 AQ Plus 100g
B-70A-5-CO-C18AQP-1KG	Bulk Coral 70A 5µm C18 AQ Plus 1kg
B-70A-5-CO-C18AQ-100G	Bulk Coral 70A 5µm C18 AQ 100g
B-70A-5-CO-C18AQ-1KG	Bulk Coral 70A 5µm C18 AQ 1kg
<b>C4</b>	
B-70A-5-CO-C4-100G	Bulk Coral 70A 5µm C4 100g
B-70A-5-CO-C4-1KG	Bulk Coral 70A 5µm C4 1kg
B-70A-5-CO-C4W-100G	Bulk Coral 70A 5µm C4 Wax 100g
B-70A-5-CO-C4W-1KG	Bulk Coral 70A 5µm C4 Wax 1kg
B-70A-5-CO-C4P-100G	Bulk Coral 70A 5µm C4 Polar 100g
B-70A-5-CO-C4P-1KG	Bulk Coral 70A 5µm C4 Polar 1kg
B-70A-5-CO-C4T-100G	Bulk Coral 70A 5µm C4 TriF 100g
B-70A-5-CO-C4T-1KG	Bulk Coral 70A 5µm C4 TriF 1kg
B-70A-5-CO-C4EX-100G	Bulk Coral 70A 5µm C4 Extreme 100g
B-70A-5-CO-C4EX-1KG	Bulk Coral 70A 5µm C4 Extreme 1kg
B-70A-5-CO-C4PT-100G	Bulk Coral 70A 5µm C4 Polar TriF 100g
B-70A-5-CO-C4PT-1KG	Bulk Coral 70A 5µm C4 Polar TriF 1kg
B-70A-5-CO-C4AQP-100G	Bulk Coral 70A 5µm C4 AQ Plus 100g
B-70A-5-CO-C4AQP-1KG	Bulk Coral 70A 5µm C4 AQ Plus 1kg
B-70A-5-CO-C4AQ-100G	Bulk Coral 70A 5µm C4 AQ 100g
B-70A-5-CO-C4AQ-1KG	Bulk Coral 70A 5µm C4 AQ 1kg
<b>Autres</b>	
B-70A-5-CO-COOH-100G	Bulk Coral 70A 5µm COOH 100g
B-70A-5-CO-COOH-1KG	Bulk Coral 70A 5µm COOH 1kg
B-70A-5-CO-SO3H-100G	Bulk Coral 70A 5µm SO3H 100g
B-70A-5-CO-SO3H-1KG	Bulk Coral 70A 5µm SO3H 1kg
B-70A-5-CO-DNH-100G	Bulk Coral 70A 5µm DNH 100g
B-70A-5-CO-DNH-1KG	Bulk Coral 70A 5µm DNH 1kg
B-70A-5-CO-PEI-100G	Bulk Coral 70A 5µm PEI 100g
B-70A-5-CO-PEI-1KG	Bulk Coral 70A 5µm PEI 1kg
B-70A-5-CO-ARG-100G	Bulk Coral 70A 5µm ARG 100g
B-70A-5-CO-ARG-1KG	Bulk Coral 70A 5µm ARG 1kg
B-70A-5-CO-EP2TNE-100G	Bulk Coral 70A 5µm EP2TNE 100g
B-70A-5-CO-EP2TNE-1KG	Bulk Coral 70A 5µm EP2TNE 1kg
B-70A-5-CO-EP4TNE-100G	Bulk Coral 70A 5µm EP4TNE 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-5-CO-EP4TNE-1KG	Bulk Coral 70A 5µm EP4TNE 1kg
B-70A-5-CO-NH2-100G	Bulk Coral 70A 5µm NH2 100g
B-70A-5-CO-NH2-1KG	Bulk Coral 70A 5µm NH2 1kg
B-70A-5-CO-CN-100G	Bulk Coral 70A 5µm CN 100g
B-70A-5-CO-CN-1KG	Bulk Coral 70A 5µm CN 1kg
B-70A-5-CO-Phenyl-100G	Bulk Coral 70A 5µm Phenyl 100g
B-70A-5-CO-Phenyl-1KG	Bulk Coral 70A 5µm Phenyl 1kg
B-70A-5-CO-Phenyl HE-100G	Bulk Coral 70A 5µm Phenyl HE 100g
B-70A-5-CO-Phenyl HE-1KG	Bulk Coral 70A 5µm Phenyl HE 1kg
B-70A-5-CO-Diol-100G	Bulk Coral 70A 5µm Diol 100g
B-70A-5-CO-Diol-1KG	Bulk Coral 70A 5µm Diol 1kg
<b>Particule Coral 5µ 100A</b>	
B-100A-5-CO-Si-100g	Bulk Coral 100A 5µm 100g
B-100A-5-CO-Si-1kg	Bulk Coral 100A 5µm 1kg
<b>C8</b>	
B-100A-5-CO-C8-100G	Bulk Coral 100A 5µm C8 100g
B-100A-5-CO-C8-1KG	Bulk Coral 100A 5µm C8 1kg
B-100A-5-CO-C8W-100G	Bulk Coral 100A 5µm C8 Wax 100g
B-100A-5-CO-C8W-1KG	Bulk Coral 100A 5µm C8 Wax 1kg
B-100A-5-CO-C8P-100G	Bulk Coral 100A 5µm C8 Polar 100g
B-100A-5-CO-C8P-1KG	Bulk Coral 100A 5µm C8 Polar 1kg
B-100A-5-CO-C8T-100G	Bulk Coral 100A 5µm C8 TriF 100g
B-100A-5-CO-C8T-1KG	Bulk Coral 100A 5µm C8 TriF 1kg
B-100A-5-CO-C8EX-100G	Bulk Coral 100A 5µm C8 Extreme 100g
B-100A-5-CO-C8EX-1KG	Bulk Coral 100A 5µm C8 Extreme 1kg
B-100A-5-CO-C8PT-100G	Bulk Coral 100A 5µm C8 Polar TriF 100g
B-100A-5-CO-C8PT-1KG	Bulk Coral 100A 5µm C8 Polar TriF 1kg
B-100A-5-CO-C8AQP-100G	Bulk Coral 100A 5µm C8 AQ Plus 100g
B-100A-5-CO-C8AQP-1KG	Bulk Coral 100A 5µm C8 AQ Plus 1kg
B-100A-5-CO-C8AQ-100G	Bulk Coral 100A 5µm C8 AQ 100g
B-100A-5-CO-C8AQ-1KG	Bulk Coral 100A 5µm C8 AQ 1kg
<b>C18</b>	
B-100A-5-CO-C18-100G	Bulk Coral 100A 5µm C18 100g
B-100A-5-CO-C18-1KG	Bulk Coral 100A 5µm C18 1kg
B-100A-5-CO-C18W-100G	Bulk Coral 100A 5µm C18 Wax 100g
B-100A-5-CO-C18W-1KG	Bulk Coral 100A 5µm C18 Wax 1kg
B-100A-5-CO-C18P-100G	Bulk Coral 100A 5µm C18 Polar 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-5-CO-C18P-1KG	Bulk Coral 100A 5µm C18 Polar 1kg
B-100A-5-CO-C18T-100G	Bulk Coral 100A 5µm C18 TriF 100g
B-100A-5-CO-C18T-1KG	Bulk Coral 100A 5µm C18 TriF 1kg
B-100A-5-CO-C18EX-100G	Bulk Coral 100A 5µm C18 Extreme 100g
B-100A-5-CO-C18EX-1KG	Bulk Coral 100A 5µm C18 Extreme 1kg
B-100A-5-CO-C18PT-100G	Bulk Coral 100A 5µm C18 Polar TriF 100g
B-100A-5-CO-C18PT-1KG	Bulk Coral 100A 5µm C18 Polar TriF 1kg
B-100A-5-CO-C18AQP-100G	Bulk Coral 100A 5µm C18 AQ Plus 100g
B-100A-5-CO-C18AQP-1KG	Bulk Coral 100A 5µm C18 AQ Plus 1kg
B-100A-5-CO-C18AQ-100G	Bulk Coral 100A 5µm C18 AQ 100g
B-100A-5-CO-C18AQ-1KG	Bulk Coral 100A 5µm C18 AQ 1kg
<b>C4</b>	
B-100A-5-CO-C4-100G	Bulk Coral 100A 5µm C4 100g
B-100A-5-CO-C4-1KG	Bulk Coral 100A 5µm C4 1kg
B-100A-5-CO-C4W-100G	Bulk Coral 100A 5µm C4 Wax 100g
B-100A-5-CO-C4W-1KG	Bulk Coral 100A 5µm C4 Wax 1kg
B-100A-5-CO-C4P-100G	Bulk Coral 100A 5µm C4 Polar 100g
B-100A-5-CO-C4P-1KG	Bulk Coral 100A 5µm C4 Polar 1kg
B-100A-5-CO-C4T-100G	Bulk Coral 100A 5µm C4 TriF 100g
B-100A-5-CO-C4T-1KG	Bulk Coral 100A 5µm C4 TriF 1kg
B-100A-5-CO-C4EX-100G	Bulk Coral 100A 5µm C4 Extreme 100g
B-100A-5-CO-C4EX-1KG	Bulk Coral 100A 5µm C4 Extreme 1kg
B-100A-5-CO-C4PT-100G	Bulk Coral 100A 5µm C4 Polar TriF 100g
B-100A-5-CO-C4PT-1KG	Bulk Coral 100A 5µm C4 Polar TriF 1kg
B-100A-5-CO-C4AQP-100G	Bulk Coral 100A 5µm C4 AQ Plus 100g
B-100A-5-CO-C4AQP-1KG	Bulk Coral 100A 5µm C4 AQ Plus 1kg
B-100A-5-CO-C4AQ-100G	Bulk Coral 100A 5µm C4 AQ 100g
B-100A-5-CO-C4AQ-1KG	Bulk Coral 100A 5µm C4 AQ 1kg
<b>Autres</b>	
B-100A-5-CO-COOH-100G	Bulk Coral 100A 5µm COOH 100g
B-100A-5-CO-COOH-1KG	Bulk Coral 100A 5µm COOH 1kg
B-100A-5-CO-SO3H-100G	Bulk Coral 100A 5µm SO3H 100g
B-100A-5-CO-SO3H-1KG	Bulk Coral 100A 5µm SO3H 1kg
B-100A-5-CO-DNH-100G	Bulk Coral 100A 5µm DNH 100g
B-100A-5-CO-DNH-1KG	Bulk Coral 100A 5µm DNH 1kg
B-100A-5-CO-PEI-100G	Bulk Coral 100A 5µm PEI 100g
B-100A-5-CO-PEI-1KG	Bulk Coral 100A 5µm PEI 1kg
B-100A-5-CO-ARG-100G	Bulk Coral 100A 5µm ARG 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-5-CO-ARG-1KG	Bulk Coral 100A 5µm ARG 1kg
B-100A-5-CO-EP2TNE-100G	Bulk Coral 100A 5µm EP2TNE 100g
B-100A-5-CO-EP2TNE-1KG	Bulk Coral 100A 5µm EP2TNE 1kg
B-100A-5-CO-EP4TNE-100G	Bulk Coral 100A 5µm EP4TNE 100g
B-100A-5-CO-EP4TNE-1KG	Bulk Coral 100A 5µm EP4TNE 1kg
B-100A-5-CO-NH2-100G	Bulk Coral 100A 5µm NH2 100g
B-100A-5-CO-NH2-1KG	Bulk Coral 100A 5µm NH2 1kg
B-100A-5-CO-CN-100G	Bulk Coral 100A 5µm CN 100g
B-100A-5-CO-CN-1KG	Bulk Coral 100A 5µm CN 1kg
B-100A-5-CO-Phenyl-100G	Bulk Coral 100A 5µm Phenyl 100g
B-100A-5-CO-Phenyl-1KG	Bulk Coral 100A 5µm Phenyl 1kg
B-100A-5-CO-Phenyl HE-100G	Bulk Coral 100A 5µm Phenyl HE 100g
B-100A-5-CO-Phenyl HE-1KG	Bulk Coral 100A 5µm Phenyl HE 1kg
B-100A-5-CO-Diol-100G	Bulk Coral 100A 5µm Diol 100g
B-100A-5-CO-Diol-1KG	Bulk Coral 100A 5µm Diol 1kg
<b>Particule Coral 5µ 150A</b>	
B-150A-5-CO-Si-100g	Bulk Coral 150A 5µm 100g
B-150A-5-CO-Si-1kg	Bulk Coral 150A 5µm 1kg
<b>C8</b>	
B-150A-5-CO-C8-100G	Bulk Coral 150A 5µm C8 100g
B-150A-5-CO-C8-1KG	Bulk Coral 150A 5µm C8 1kg
B-150A-5-CO-C8W-100G	Bulk Coral 150A 5µm C8 Wax 100g
B-150A-5-CO-C8W-1KG	Bulk Coral 150A 5µm C8 Wax 1kg
B-150A-5-CO-C8P-100G	Bulk Coral 150A 5µm C8 Polar 100g
B-150A-5-CO-C8P-1KG	Bulk Coral 150A 5µm C8 Polar 1kg
B-150A-5-CO-C8T-100G	Bulk Coral 150A 5µm C8 TriF 100g
B-150A-5-CO-C8T-1KG	Bulk Coral 150A 5µm C8 TriF 1kg
B-150A-5-CO-C8EX-100G	Bulk Coral 150A 5µm C8 Extreme 100g
B-150A-5-CO-C8EX-1KG	Bulk Coral 150A 5µm C8 Extreme 1kg
B-150A-5-CO-C8PT-100G	Bulk Coral 150A 5µm C8 Polar TriF 100g
B-150A-5-CO-C8PT-1KG	Bulk Coral 150A 5µm C8 Polar TriF 1kg
B-150A-5-CO-C8AQP-100G	Bulk Coral 150A 5µm C8 AQ Plus 100g
B-150A-5-CO-C8AQP-1KG	Bulk Coral 150A 5µm C8 AQ Plus 1kg
B-150A-5-CO-C8AQ-100G	Bulk Coral 150A 5µm C8 AQ 100g
B-150A-5-CO-C8AQ-1KG	Bulk Coral 150A 5µm C8 AQ 1kg
<b>C18</b>	
B-150A-5-CO-C18-100G	Bulk Coral 150A 5µm C18 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-5-CO-C18-1KG	Bulk Coral 150A 5µm C18 1kg
B-150A-5-CO-C18W-100G	Bulk Coral 150A 5µm C18 Wax 100g
B-150A-5-CO-C18W-1KG	Bulk Coral 150A 5µm C18 Wax 1kg
B-150A-5-CO-C18P-100G	Bulk Coral 150A 5µm C18 Polar 100g
B-150A-5-CO-C18P-1KG	Bulk Coral 150A 5µm C18 Polar 1kg
B-150A-5-CO-C18T-100G	Bulk Coral 150A 5µm C18 TriF 100g
B-150A-5-CO-C18T-1KG	Bulk Coral 150A 5µm C18 TriF 1kg
B-150A-5-CO-C18EX-100G	Bulk Coral 150A 5µm C18 Extreme 100g
B-150A-5-CO-C18EX-1KG	Bulk Coral 150A 5µm C18 Extreme 1kg
B-150A-5-CO-C18PT-100G	Bulk Coral 150A 5µm C18 Polar TriF 100g
B-150A-5-CO-C18PT-1KG	Bulk Coral 150A 5µm C18 Polar TriF 1kg
B-150A-5-CO-C18AQP-100G	Bulk Coral 150A 5µm C18 AQ Plus 100g
B-150A-5-CO-C18AQP-1KG	Bulk Coral 150A 5µm C18 AQ Plus 1kg
B-150A-5-CO-C18AQ-100G	Bulk Coral 150A 5µm C18 AQ 100g
B-150A-5-CO-C18AQ-1KG	Bulk Coral 150A 5µm C18 AQ 1kg
<b>C4</b>	
B-150A-5-CO-C4-100G	Bulk Coral 150A 5µm C4 100g
B-150A-5-CO-C4-1KG	Bulk Coral 150A 5µm C4 1kg
B-150A-5-CO-C4W-100G	Bulk Coral 150A 5µm C4 Wax 100g
B-150A-5-CO-C4W-1KG	Bulk Coral 150A 5µm C4 Wax 1kg
B-150A-5-CO-C4P-100G	Bulk Coral 150A 5µm C4 Polar 100g
B-150A-5-CO-C4P-1KG	Bulk Coral 150A 5µm C4 Polar 1kg
B-150A-5-CO-C4T-100G	Bulk Coral 150A 5µm C4 TriF 100g
B-150A-5-CO-C4T-1KG	Bulk Coral 150A 5µm C4 TriF 1kg
B-150A-5-CO-C4EX-100G	Bulk Coral 150A 5µm C4 Extreme 100g
B-150A-5-CO-C4EX-1KG	Bulk Coral 150A 5µm C4 Extreme 1kg
B-150A-5-CO-C4PT-100G	Bulk Coral 150A 5µm C4 Polar TriF 100g
B-150A-5-CO-C4PT-1KG	Bulk Coral 150A 5µm C4 Polar TriF 1kg
B-150A-5-CO-C4AQP-100G	Bulk Coral 150A 5µm C4 AQ Plus 100g
B-150A-5-CO-C4AQP-1KG	Bulk Coral 150A 5µm C4 AQ Plus 1kg
B-150A-5-CO-C4AQ-100G	Bulk Coral 150A 5µm C4 AQ 100g
B-150A-5-CO-C4AQ-1KG	Bulk Coral 150A 5µm C4 AQ 1kg
<b>Autres</b>	
B-150A-5-CO-COOH-100G	Bulk Coral 150A 5µm COOH 100g
B-150A-5-CO-COOH-1KG	Bulk Coral 150A 5µm COOH 1kg
B-150A-5-CO-SO3H-100G	Bulk Coral 150A 5µm SO3H 100g
B-150A-5-CO-SO3H-1KG	Bulk Coral 150A 5µm SO3H 1kg
B-150A-5-CO-DNH-100G	Bulk Coral 150A 5µm DNH 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-5-CO-DNH-1KG	Bulk Coral 150A 5µm DNH 1kg
B-150A-5-CO-PEI-100G	Bulk Coral 150A 5µm PEI 100g
B-150A-5-CO-PEI-1KG	Bulk Coral 150A 5µm PEI 1kg
B-150A-5-CO-ARG-100G	Bulk Coral 150A 5µm ARG 100g
B-150A-5-CO-ARG-1KG	Bulk Coral 150A 5µm ARG 1kg
B-150A-5-CO-EP2TNE-100G	Bulk Coral 150A 5µm EP2TNE 100g
B-150A-5-CO-EP2TNE-1KG	Bulk Coral 150A 5µm EP2TNE 1kg
B-150A-5-CO-EP4TNE-100G	Bulk Coral 150A 5µm EP4TNE 100g
B-150A-5-CO-EP4TNE-1KG	Bulk Coral 150A 5µm EP4TNE 1kg
B-150A-5-CO-NH2-100G	Bulk Coral 150A 5µm NH2 100g
B-150A-5-CO-NH2-1KG	Bulk Coral 150A 5µm NH2 1kg
B-150A-5-CO-CN-100G	Bulk Coral 150A 5µm CN 100g
B-150A-5-CO-CN-1KG	Bulk Coral 150A 5µm CN 1kg
B-150A-5-CO-Phenyl-100G	Bulk Coral 150A 5µm Phenyl 100g
B-150A-5-CO-Phenyl-1KG	Bulk Coral 150A 5µm Phenyl 1kg
B-150A-5-CO-Phenyl HE-100G	Bulk Coral 150A 5µm Phenyl HE 100g
B-150A-5-CO-Phenyl HE-1KG	Bulk Coral 150A 5µm Phenyl HE 1kg
B-150A-5-CO-Diol-100G	Bulk Coral 150A 5µm Diol 100g
B-150A-5-CO-Diol-1KG	Bulk Coral 150A 5µm Diol 1kg
<b>Particule Coral 5µ 200A</b>	
B-200A-5-CO-Si-100g	Bulk Coral 200A 5µm 100g
B-200A-5-CO-Si-1kg	Bulk Coral 200A 5µm 1kg
<b>C8</b>	
B-200A-5-CO-C8-100G	Bulk Coral 200A 5µm C8 100g
B-200A-5-CO-C8-1KG	Bulk Coral 200A 5µm C8 1kg
B-200A-5-CO-C8W-100G	Bulk Coral 200A 5µm C8 Wax 100g
B-200A-5-CO-C8W-1KG	Bulk Coral 200A 5µm C8 Wax 1kg
B-200A-5-CO-C8P-100G	Bulk Coral 200A 5µm C8 Polar 100g
B-200A-5-CO-C8P-1KG	Bulk Coral 200A 5µm C8 Polar 1kg
B-200A-5-CO-C8T-100G	Bulk Coral 200A 5µm C8 TriF 100g
B-200A-5-CO-C8T-1KG	Bulk Coral 200A 5µm C8 TriF 1kg
B-200A-5-CO-C8EX-100G	Bulk Coral 200A 5µm C8 Extreme 100g
B-200A-5-CO-C8EX-1KG	Bulk Coral 200A 5µm C8 Extreme 1kg
B-200A-5-CO-C8PT-100G	Bulk Coral 200A 5µm C8 Polar TriF 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-5-CO-C8PT-1KG	Bulk Coral 200A 5µm C8 Polar TriF 1kg
B-200A-5-CO-C8AQP-100G	Bulk Coral 200A 5µm C8 AQ Plus 100g
B-200A-5-CO-C8AQP-1KG	Bulk Coral 200A 5µm C8 AQ Plus 1kg
B-200A-5-CO-C8AQ-100G	Bulk Coral 200A 5µm C8 AQ 100g
B-200A-5-CO-C8AQ-1KG	Bulk Coral 200A 5µm C8 AQ 1kg
<b>C18</b>	
B-200A-5-CO-C18-100G	Bulk Coral 200A 5µm C18 100g
B-200A-5-CO-C18-1KG	Bulk Coral 200A 5µm C18 1kg
B-200A-5-CO-C18W-100G	Bulk Coral 200A 5µm C18 Wax 100g
B-200A-5-CO-C18W-1KG	Bulk Coral 200A 5µm C18 Wax 1kg
B-200A-5-CO-C18P-100G	Bulk Coral 200A 5µm C18 Polar 100g
B-200A-5-CO-C18P-1KG	Bulk Coral 200A 5µm C18 Polar 1kg
B-200A-5-CO-C18T-100G	Bulk Coral 200A 5µm C18 TriF 100g
B-200A-5-CO-C18T-1KG	Bulk Coral 200A 5µm C18 TriF 1kg
B-200A-5-CO-C18EX-100G	Bulk Coral 200A 5µm C18 Extreme 100g
B-200A-5-CO-C18EX-1KG	Bulk Coral 200A 5µm C18 Extreme 1kg
B-200A-5-CO-C18PT-100G	Bulk Coral 200A 5µm C18 Polar TriF 100g
B-200A-5-CO-C18PT-1KG	Bulk Coral 200A 5µm C18 Polar TriF 1kg
B-200A-5-CO-C18AQP-100G	Bulk Coral 200A 5µm C18 AQ Plus 100g
B-200A-5-CO-C18AQP-1KG	Bulk Coral 200A 5µm C18 AQ Plus 1kg
B-200A-5-CO-C18AQ-100G	Bulk Coral 200A 5µm C18 AQ 100g
B-200A-5-CO-C18AQ-1KG	Bulk Coral 200A 5µm C18 AQ 1kg
<b>C4</b>	
B-200A-5-CO-C4-100G	Bulk Coral 200A 5µm C4 100g
B-200A-5-CO-C4-1KG	Bulk Coral 200A 5µm C4 1kg
B-200A-5-CO-C4W-100G	Bulk Coral 200A 5µm C4 Wax 100g
B-200A-5-CO-C4W-1KG	Bulk Coral 200A 5µm C4 Wax 1kg
B-200A-5-CO-C4P-100G	Bulk Coral 200A 5µm C4 Polar 100g
B-200A-5-CO-C4P-1KG	Bulk Coral 200A 5µm C4 Polar 1kg
B-200A-5-CO-C4T-100G	Bulk Coral 200A 5µm C4 TriF 100g
B-200A-5-CO-C4T-1KG	Bulk Coral 200A 5µm C4 TriF 1kg
B-200A-5-CO-C4EX-100G	Bulk Coral 200A 5µm C4 Extreme 100g
B-200A-5-CO-C4EX-1KG	Bulk Coral 200A 5µm C4 Extreme 1kg
B-200A-5-CO-C4PT-100G	Bulk Coral 200A 5µm C4 Polar TriF 100g
B-200A-5-CO-C4PT-1KG	Bulk Coral 200A 5µm C4 Polar TriF 1kg
B-200A-5-CO-C4AQP-100G	Bulk Coral 200A 5µm C4 AQ Plus 100g
B-200A-5-CO-C4AQP-1KG	Bulk Coral 200A 5µm C4 AQ Plus 1kg
B-200A-5-CO-C4AQ-100G	Bulk Coral 200A 5µm C4 AQ 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-5-CO-C4AQ-1KG	Bulk Coral 200A 5µm C4 AQ 1kg
<b>Autres</b>	
B-200A-5-CO-COOH-100G	Bulk Coral 200A 5µm COOH 100g
B-200A-5-CO-COOH-1KG	Bulk Coral 200A 5µm COOH 1kg
B-200A-5-CO-SO3H-100G	Bulk Coral 200A 5µm SO3H 100g
B-200A-5-CO-SO3H-1KG	Bulk Coral 200A 5µm SO3H 1kg
B-200A-5-CO-DNH-100G	Bulk Coral 200A 5µm DNH 100g
B-200A-5-CO-DNH-1KG	Bulk Coral 200A 5µm DNH 1kg
B-200A-5-CO-PEI-100G	Bulk Coral 200A 5µm PEI 100g
B-200A-5-CO-PEI-1KG	Bulk Coral 200A 5µm PEI 1kg
B-200A-5-CO-ARG-100G	Bulk Coral 200A 5µm ARG 100g
B-200A-5-CO-ARG-1KG	Bulk Coral 200A 5µm ARG 1kg
B-200A-5-CO-EP2TNE-100G	Bulk Coral 200A 5µm EP2TNE 100g
B-200A-5-CO-EP2TNE-1KG	Bulk Coral 200A 5µm EP2TNE 1kg
B-200A-5-CO-EP4TNE-100G	Bulk Coral 200A 5µm EP4TNE 100g
B-200A-5-CO-EP4TNE-1KG	Bulk Coral 200A 5µm EP4TNE 1kg
B-200A-5-CO-NH2-100G	Bulk Coral 200A 5µm NH2 100g
B-200A-5-CO-NH2-1KG	Bulk Coral 200A 5µm NH2 1kg
B-200A-5-CO-CN-100G	Bulk Coral 200A 5µm CN 100g
B-200A-5-CO-CN-1KG	Bulk Coral 200A 5µm CN 1kg
B-200A-5-CO-Phenyl-100G	Bulk Coral 200A 5µm Phenyl 100g
B-200A-5-CO-Phenyl-1KG	Bulk Coral 200A 5µm Phenyl 1kg
B-200A-5-CO-Phenyl HE-100G	Bulk Coral 200A 5µm Phenyl HE 100g
B-200A-5-CO-Phenyl HE-1KG	Bulk Coral 200A 5µm Phenyl HE 1kg
B-200A-5-CO-Diol-100G	Bulk Coral 200A 5µm Diol 100g
B-200A-5-CO-Diol-1KG	Bulk Coral 200A 5µm Diol 1kg
<b>Particule Coral 5µ 300A</b>	
B-300A-5-CO-Si-100g	Bulk Coral 300A 5µm 100g
B-300A-5-CO-Si-1kg	Bulk Coral 300A 5µm 1kg
<b>C8</b>	
B-300A-5-CO-C8-100G	Bulk Coral 300A 5µm C8 100g
B-300A-5-CO-C8-1KG	Bulk Coral 300A 5µm C8 1kg
B-300A-5-CO-C8W-100G	Bulk Coral 300A 5µm C8 Wax 100g
B-300A-5-CO-C8W-1KG	Bulk Coral 300A 5µm C8 Wax 1kg
B-300A-5-CO-C8P-100G	Bulk Coral 300A 5µm C8 Polar 100g
B-300A-5-CO-C8P-1KG	Bulk Coral 300A 5µm C8 Polar 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-5-CO-C8T-100G	Bulk Coral 300A 5µm C8 TriF 100g
B-300A-5-CO-C8T-1KG	Bulk Coral 300A 5µm C8 TriF 1kg
B-300A-5-CO-C8EX-100G	Bulk Coral 300A 5µm C8 Extreme 100g
B-300A-5-CO-C8EX-1KG	Bulk Coral 300A 5µm C8 Extreme 1kg
B-300A-5-CO-C8PT-100G	Bulk Coral 300A 5µm C8 Polar TriF 100g
B-300A-5-CO-C8PT-1KG	Bulk Coral 300A 5µm C8 Polar TriF 1kg
B-300A-5-CO-C8AQP-100G	Bulk Coral 300A 5µm C8 AQ Plus 100g
B-300A-5-CO-C8AQP-1KG	Bulk Coral 300A 5µm C8 AQ Plus 1kg
B-300A-5-CO-C8AQ-100G	Bulk Coral 300A 5µm C8 AQ 100g
B-300A-5-CO-C8AQ-1KG	Bulk Coral 300A 5µm C8 AQ 1kg
<b>C18</b>	
B-300A-5-CO-C18-100G	Bulk Coral 300A 5µm C18 100g
B-300A-5-CO-C18-1KG	Bulk Coral 300A 5µm C18 1kg
B-300A-5-CO-C18W-100G	Bulk Coral 300A 5µm C18 Wax 100g
B-300A-5-CO-C18W-1KG	Bulk Coral 300A 5µm C18 Wax 1kg
B-300A-5-CO-C18P-100G	Bulk Coral 300A 5µm C18 Polar 100g
B-300A-5-CO-C18P-1KG	Bulk Coral 300A 5µm C18 Polar 1kg
B-300A-5-CO-C18T-100G	Bulk Coral 300A 5µm C18 TriF 100g
B-300A-5-CO-C18T-1KG	Bulk Coral 300A 5µm C18 TriF 1kg
B-300A-5-CO-C18EX-100G	Bulk Coral 300A 5µm C18 Extreme 100g
B-300A-5-CO-C18EX-1KG	Bulk Coral 300A 5µm C18 Extreme 1kg
B-300A-5-CO-C18PT-100G	Bulk Coral 300A 5µm C18 Polar TriF 100g
B-300A-5-CO-C18PT-1KG	Bulk Coral 300A 5µm C18 Polar TriF 1kg
B-300A-5-CO-C18AQP-100G	Bulk Coral 300A 5µm C18 AQ Plus 100g
B-300A-5-CO-C18AQP-1KG	Bulk Coral 300A 5µm C18 AQ Plus 1kg
B-300A-5-CO-C18AQ-100G	Bulk Coral 300A 5µm C18 AQ 100g
B-300A-5-CO-C18AQ-1KG	Bulk Coral 300A 5µm C18 AQ 1kg
<b>C4</b>	
B-300A-5-CO-C4-100G	Bulk Coral 300A 5µm C4 100g
B-300A-5-CO-C4-1KG	Bulk Coral 300A 5µm C4 1kg
B-300A-5-CO-C4W-100G	Bulk Coral 300A 5µm C4 Wax 100g
B-300A-5-CO-C4W-1KG	Bulk Coral 300A 5µm C4 Wax 1kg
B-300A-5-CO-C4P-100G	Bulk Coral 300A 5µm C4 Polar 100g
B-300A-5-CO-C4P-1KG	Bulk Coral 300A 5µm C4 Polar 1kg
B-300A-5-CO-C4T-100G	Bulk Coral 300A 5µm C4 TriF 100g
B-300A-5-CO-C4T-1KG	Bulk Coral 300A 5µm C4 TriF 1kg
B-300A-5-CO-C4EX-100G	Bulk Coral 300A 5µm C4 Extreme 100g
B-300A-5-CO-C4EX-1KG	Bulk Coral 300A 5µm C4 Extreme 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-5-CO-C4PT-100G	Bulk Coral 300A 5µm C4 Polar TriF 100g
B-300A-5-CO-C4PT-1KG	Bulk Coral 300A 5µm C4 Polar TriF 1kg
B-300A-5-CO-C4AQP-100G	Bulk Coral 300A 5µm C4 AQ Plus 100g
B-300A-5-CO-C4AQP-1KG	Bulk Coral 300A 5µm C4 AQ Plus 1kg
B-300A-5-CO-C4AQ-100G	Bulk Coral 300A 5µm C4 AQ 100g
B-300A-5-CO-C4AQ-1KG	Bulk Coral 300A 5µm C4 AQ 1kg
<b>Autres</b>	
B-300A-5-CO-COOH-100G	Bulk Coral 300A 5µm COOH 100g
B-300A-5-CO-COOH-1KG	Bulk Coral 300A 5µm COOH 1kg
B-300A-5-CO-SO3H-100G	Bulk Coral 300A 5µm SO3H 100g
B-300A-5-CO-SO3H-1KG	Bulk Coral 300A 5µm SO3H 1kg
B-300A-5-CO-DNH-100G	Bulk Coral 300A 5µm DNH 100g
B-300A-5-CO-DNH-1KG	Bulk Coral 300A 5µm DNH 1kg
B-300A-5-CO-PEI-100G	Bulk Coral 300A 5µm PEI 100g
B-300A-5-CO-PEI-1KG	Bulk Coral 300A 5µm PEI 1kg
B-300A-5-CO-ARG-100G	Bulk Coral 300A 5µm ARG 100g
B-300A-5-CO-ARG-1KG	Bulk Coral 300A 5µm ARG 1kg
B-300A-5-CO-EP2TNE-100G	Bulk Coral 300A 5µm EP2TNE 100g
B-300A-5-CO-EP2TNE-1KG	Bulk Coral 300A 5µm EP2TNE 1kg
B-300A-5-CO-EP4TNE-100G	Bulk Coral 300A 5µm EP4TNE 100g
B-300A-5-CO-EP4TNE-1KG	Bulk Coral 300A 5µm EP4TNE 1kg
B-300A-5-CO-NH2-100G	Bulk Coral 300A 5µm NH2 100g
B-300A-5-CO-NH2-1KG	Bulk Coral 300A 5µm NH2 1kg
B-300A-5-CO-CN-100G	Bulk Coral 300A 5µm CN 100g
B-300A-5-CO-CN-1KG	Bulk Coral 300A 5µm CN 1kg
B-300A-5-CO-Phenyl-100G	Bulk Coral 300A 5µm Phenyl 100g
B-300A-5-CO-Phenyl-1KG	Bulk Coral 300A 5µm Phenyl 1kg
B-300A-5-CO-Phenyl HE-100G	Bulk Coral 300A 5µm Phenyl HE 100g
B-300A-5-CO-Phenyl HE-1KG	Bulk Coral 300A 5µm Phenyl HE 1kg
B-300A-5-CO-Diol-100G	Bulk Coral 300A 5µm Diol 100g
B-300A-5-CO-Diol-1KG	Bulk Coral 300A 5µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Coral 5µ 500A</b>	
B-500A-5-CO-Si-100g	Bulk Coral 500A 5µm 100g
B-500A-5-CO-Si-1kg	Bulk Coral 500A 5µm 1kg
<b>C8</b>	
B-500A-5-CO-C8-100G	Bulk Coral 500A 5µm C8 100g
B-500A-5-CO-C8-1KG	Bulk Coral 500A 5µm C8 1kg
B-500A-5-CO-C8W-100G	Bulk Coral 500A 5µm C8 Wax 100g
B-500A-5-CO-C8W-1KG	Bulk Coral 500A 5µm C8 Wax 1kg
B-500A-5-CO-C8P-100G	Bulk Coral 500A 5µm C8 Polar 100g
B-500A-5-CO-C8P-1KG	Bulk Coral 500A 5µm C8 Polar 1kg
B-500A-5-CO-C8T-100G	Bulk Coral 500A 5µm C8 TriF 100g
B-500A-5-CO-C8T-1KG	Bulk Coral 500A 5µm C8 TriF 1kg
B-500A-5-CO-C8EX-100G	Bulk Coral 500A 5µm C8 Extreme 100g
B-500A-5-CO-C8EX-1KG	Bulk Coral 500A 5µm C8 Extreme 1kg
B-500A-5-CO-C8PT-100G	Bulk Coral 500A 5µm C8 Polar TriF 100g
B-500A-5-CO-C8PT-1KG	Bulk Coral 500A 5µm C8 Polar TriF 1kg
B-500A-5-CO-C8AQP-100G	Bulk Coral 500A 5µm C8 AQ Plus 100g
B-500A-5-CO-C8AQP-1KG	Bulk Coral 500A 5µm C8 AQ Plus 1kg
B-500A-5-CO-C8AQ-100G	Bulk Coral 500A 5µm C8 AQ 100g
B-500A-5-CO-C8AQ-1KG	Bulk Coral 500A 5µm C8 AQ 1kg
<b>C18</b>	
B-500A-5-CO-C18-100G	Bulk Coral 500A 5µm C18 100g
B-500A-5-CO-C18-1KG	Bulk Coral 500A 5µm C18 1kg
B-500A-5-CO-C18W-100G	Bulk Coral 500A 5µm C18 Wax 100g
B-500A-5-CO-C18W-1KG	Bulk Coral 500A 5µm C18 Wax 1kg
B-500A-5-CO-C18P-100G	Bulk Coral 500A 5µm C18 Polar 100g
B-500A-5-CO-C18P-1KG	Bulk Coral 500A 5µm C18 Polar 1kg
B-500A-5-CO-C18T-100G	Bulk Coral 500A 5µm C18 TriF 100g
B-500A-5-CO-C18T-1KG	Bulk Coral 500A 5µm C18 TriF 1kg
B-500A-5-CO-C18EX-100G	Bulk Coral 500A 5µm C18 Extreme 100g
B-500A-5-CO-C18EX-1KG	Bulk Coral 500A 5µm C18 Extreme 1kg
B-500A-5-CO-C18PT-100G	Bulk Coral 500A 5µm C18 Polar TriF 100g
B-500A-5-CO-C18PT-1KG	Bulk Coral 500A 5µm C18 Polar TriF 1kg
B-500A-5-CO-C18AQP-100G	Bulk Coral 500A 5µm C18 AQ Plus 100g
B-500A-5-CO-C18AQP-1KG	Bulk Coral 500A 5µm C18 AQ Plus 1kg
B-500A-5-CO-C18AQ-100G	Bulk Coral 500A 5µm C18 AQ 100g
B-500A-5-CO-C18AQ-1KG	Bulk Coral 500A 5µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-500A-5-CO-C4-100G	Bulk Coral 500A 5µm C4 100g
B-500A-5-CO-C4-1KG	Bulk Coral 500A 5µm C4 1kg
B-500A-5-CO-C4W-100G	Bulk Coral 500A 5µm C4 Wax 100g
B-500A-5-CO-C4W-1KG	Bulk Coral 500A 5µm C4 Wax 1kg
B-500A-5-CO-C4P-100G	Bulk Coral 500A 5µm C4 Polar 100g
B-500A-5-CO-C4P-1KG	Bulk Coral 500A 5µm C4 Polar 1kg
B-500A-5-CO-C4T-100G	Bulk Coral 500A 5µm C4 TriF 100g
B-500A-5-CO-C4T-1KG	Bulk Coral 500A 5µm C4 TriF 1kg
B-500A-5-CO-C4EX-100G	Bulk Coral 500A 5µm C4 Extreme 100g
B-500A-5-CO-C4EX-1KG	Bulk Coral 500A 5µm C4 Extreme 1kg
B-500A-5-CO-C4PT-100G	Bulk Coral 500A 5µm C4 Polar TriF 100g
B-500A-5-CO-C4PT-1KG	Bulk Coral 500A 5µm C4 Polar TriF 1kg
B-500A-5-CO-C4AQP-100G	Bulk Coral 500A 5µm C4 AQ Plus 100g
B-500A-5-CO-C4AQP-1KG	Bulk Coral 500A 5µm C4 AQ Plus 1kg
B-500A-5-CO-C4AQ-100G	Bulk Coral 500A 5µm C4 AQ 100g
B-500A-5-CO-C4AQ-1KG	Bulk Coral 500A 5µm C4 AQ 1kg
<b>Autres</b>	
B-500A-5-CO-COOH-100G	Bulk Coral 500A 5µm COOH 100g
B-500A-5-CO-COOH-1KG	Bulk Coral 500A 5µm COOH 1kg
B-500A-5-CO-SO3H-100G	Bulk Coral 500A 5µm SO3H 100g
B-500A-5-CO-SO3H-1KG	Bulk Coral 500A 5µm SO3H 1kg
B-500A-5-CO-DNH-100G	Bulk Coral 500A 5µm DNH 100g
B-500A-5-CO-DNH-1KG	Bulk Coral 500A 5µm DNH 1kg
B-500A-5-CO-PEI-100G	Bulk Coral 500A 5µm PEI 100g
B-500A-5-CO-PEI-1KG	Bulk Coral 500A 5µm PEI 1kg
B-500A-5-CO-ARG-100G	Bulk Coral 500A 5µm ARG 100g
B-500A-5-CO-ARG-1KG	Bulk Coral 500A 5µm ARG 1kg
B-500A-5-CO-EP2TNE-100G	Bulk Coral 500A 5µm EP2TNE 100g
B-500A-5-CO-EP2TNE-1KG	Bulk Coral 500A 5µm EP2TNE 1kg
B-500A-5-CO-EP4TNE-100G	Bulk Coral 500A 5µm EP4TNE 100g
B-500A-5-CO-EP4TNE-1KG	Bulk Coral 500A 5µm EP4TNE 1kg
B-500A-5-CO-NH2-100G	Bulk Coral 500A 5µm NH2 100g
B-500A-5-CO-NH2-1KG	Bulk Coral 500A 5µm NH2 1kg
B-500A-5-CO-CN-100G	Bulk Coral 500A 5µm CN 100g
B-500A-5-CO-CN-1KG	Bulk Coral 500A 5µm CN 1kg
B-500A-5-CO-Phenyl-100G	Bulk Coral 500A 5µm Phenyl 100g
B-500A-5-CO-Phenyl-1KG	Bulk Coral 500A 5µm Phenyl 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-500A-5-CO-Phenyl HE-100G	Bulk Coral 500A 5µm Phenyl HE 100g
B-500A-5-CO-Phenyl HE-1KG	Bulk Coral 500A 5µm Phenyl HE 1kg
B-500A-5-CO-Diol-100G	Bulk Coral 500A 5µm Diol 100g
B-500A-5-CO-Diol-1KG	Bulk Coral 500A 5µm Diol 1kg
<b>Particule Coral 5µ 800A</b>	
B-800A-5-CO-Si-100g	Bulk Coral 800A 5µm 100g
B-800A-5-CO-Si-1kg	Bulk Coral 800A 5µm 1kg
<b>C8</b>	
B-800A-5-CO-C8-100G	Bulk Coral 800A 5µm C8 100g
B-800A-5-CO-C8-1KG	Bulk Coral 800A 5µm C8 1kg
B-800A-5-CO-C8W-100G	Bulk Coral 800A 5µm C8 Wax 100g
B-800A-5-CO-C8W-1KG	Bulk Coral 800A 5µm C8 Wax 1kg
B-800A-5-CO-C8P-100G	Bulk Coral 800A 5µm C8 Polar 100g
B-800A-5-CO-C8P-1KG	Bulk Coral 800A 5µm C8 Polar 1kg
B-800A-5-CO-C8T-100G	Bulk Coral 800A 5µm C8 TriF 100g
B-800A-5-CO-C8T-1KG	Bulk Coral 800A 5µm C8 TriF 1kg
B-800A-5-CO-C8EX-100G	Bulk Coral 800A 5µm C8 Extreme 100g
B-800A-5-CO-C8EX-1KG	Bulk Coral 800A 5µm C8 Extreme 1kg
B-800A-5-CO-C8PT-100G	Bulk Coral 800A 5µm C8 Polar TriF 100g
B-800A-5-CO-C8PT-1KG	Bulk Coral 800A 5µm C8 Polar TriF 1kg
B-800A-5-CO-C8AQP-100G	Bulk Coral 800A 5µm C8 AQ Plus 100g
B-800A-5-CO-C8AQP-1KG	Bulk Coral 800A 5µm C8 AQ Plus 1kg
B-800A-5-CO-C8AQ-100G	Bulk Coral 800A 5µm C8 AQ 100g
B-800A-5-CO-C8AQ-1KG	Bulk Coral 800A 5µm C8 AQ 1kg
<b>C18</b>	
B-800A-5-CO-C18-100G	Bulk Coral 800A 5µm C18 100g
B-800A-5-CO-C18-1KG	Bulk Coral 800A 5µm C18 1kg
B-800A-5-CO-C18W-100G	Bulk Coral 800A 5µm C18 Wax 100g
B-800A-5-CO-C18W-1KG	Bulk Coral 800A 5µm C18 Wax 1kg
B-800A-5-CO-C18P-100G	Bulk Coral 800A 5µm C18 Polar 100g
B-800A-5-CO-C18P-1KG	Bulk Coral 800A 5µm C18 Polar 1kg
B-800A-5-CO-C18T-100G	Bulk Coral 800A 5µm C18 TriF 100g
B-800A-5-CO-C18T-1KG	Bulk Coral 800A 5µm C18 TriF 1kg
B-800A-5-CO-C18EX-100G	Bulk Coral 800A 5µm C18 Extreme 100g
B-800A-5-CO-C18EX-1KG	Bulk Coral 800A 5µm C18 Extreme 1kg
B-800A-5-CO-C18PT-100G	Bulk Coral 800A 5µm C18 Polar TriF 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-5-CO-C18PT-1KG	Bulk Coral 800A 5µm C18 Polar TriF 1kg
B-800A-5-CO-C18AQP-100G	Bulk Coral 800A 5µm C18 AQ Plus 100g
B-800A-5-CO-C18AQP-1KG	Bulk Coral 800A 5µm C18 AQ Plus 1kg
B-800A-5-CO-C18AQ-100G	Bulk Coral 800A 5µm C18 AQ 100g
B-800A-5-CO-C18AQ-1KG	Bulk Coral 800A 5µm C18 AQ 1kg
<b>C4</b>	
B-800A-5-CO-C4-100G	Bulk Coral 800A 5µm C4 100g
B-800A-5-CO-C4-1KG	Bulk Coral 800A 5µm C4 1kg
B-800A-5-CO-C4W-100G	Bulk Coral 800A 5µm C4 Wax 100g
B-800A-5-CO-C4W-1KG	Bulk Coral 800A 5µm C4 Wax 1kg
B-800A-5-CO-C4P-100G	Bulk Coral 800A 5µm C4 Polar 100g
B-800A-5-CO-C4P-1KG	Bulk Coral 800A 5µm C4 Polar 1kg
B-800A-5-CO-C4T-100G	Bulk Coral 800A 5µm C4 TriF 100g
B-800A-5-CO-C4T-1KG	Bulk Coral 800A 5µm C4 TriF 1kg
B-800A-5-CO-C4EX-100G	Bulk Coral 800A 5µm C4 Extreme 100g
B-800A-5-CO-C4EX-1KG	Bulk Coral 800A 5µm C4 Extreme 1kg
B-800A-5-CO-C4PT-100G	Bulk Coral 800A 5µm C4 Polar TriF 100g
B-800A-5-CO-C4PT-1KG	Bulk Coral 800A 5µm C4 Polar TriF 1kg
B-800A-5-CO-C4AQP-100G	Bulk Coral 800A 5µm C4 AQ Plus 100g
B-800A-5-CO-C4AQP-1KG	Bulk Coral 800A 5µm C4 AQ Plus 1kg
B-800A-5-CO-C4AQ-100G	Bulk Coral 800A 5µm C4 AQ 100g
B-800A-5-CO-C4AQ-1KG	Bulk Coral 800A 5µm C4 AQ 1kg
<b>Autres</b>	
B-800A-5-CO-COOH-100G	Bulk Coral 800A 5µm COOH 100g
B-800A-5-CO-COOH-1KG	Bulk Coral 800A 5µm COOH 1kg
B-800A-5-CO-SO3H-100G	Bulk Coral 800A 5µm SO3H 100g
B-800A-5-CO-SO3H-1KG	Bulk Coral 800A 5µm SO3H 1kg
B-800A-5-CO-DNH-100G	Bulk Coral 800A 5µm DNH 100g
B-800A-5-CO-DNH-1KG	Bulk Coral 800A 5µm DNH 1kg
B-800A-5-CO-PEI-100G	Bulk Coral 800A 5µm PEI 100g
B-800A-5-CO-PEI-1KG	Bulk Coral 800A 5µm PEI 1kg
B-800A-5-CO-ARG-100G	Bulk Coral 800A 5µm ARG 100g
B-800A-5-CO-ARG-1KG	Bulk Coral 800A 5µm ARG 1kg
B-800A-5-CO-EP2TNE-100G	Bulk Coral 800A 5µm EP2TNE 100g
B-800A-5-CO-EP2TNE-1KG	Bulk Coral 800A 5µm EP2TNE 1kg
B-800A-5-CO-EP4TNE-100G	Bulk Coral 800A 5µm EP4TNE 100g
B-800A-5-CO-EP4TNE-1KG	Bulk Coral 800A 5µm EP4TNE 1kg
B-800A-5-CO-NH2-100G	Bulk Coral 800A 5µm NH2 100g



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-5-CO-NH2-1KG	Bulk Coral 800A 5µm NH2 1kg
B-800A-5-CO-CN-100G	Bulk Coral 800A 5µm CN 100g
B-800A-5-CO-CN-1KG	Bulk Coral 800A 5µm CN 1kg
B-800A-5-CO-Phenyl-100G	Bulk Coral 800A 5µm Phenyl 100g
B-800A-5-CO-Phenyl-1KG	Bulk Coral 800A 5µm Phenyl 1kg
B-800A-5-CO-Phenyl HE-100G	Bulk Coral 800A 5µm Phenyl HE 100g
B-800A-5-CO-Phenyl HE-1KG	Bulk Coral 800A 5µm Phenyl HE 1kg
B-800A-5-CO-Diol-100G	Bulk Coral 800A 5µm Diol 100g
B-800A-5-CO-Diol-1KG	Bulk Coral 800A 5µm Diol 1kg
<b>Coral 10µm</b>	
<b>Particule Coral 10, Pore 70A</b>	
B-70A-10-CO-Si-100G	Bulk Coral Silice 70A 10µm 100g
B-70A-10-CO-Si-1KG	Bulk Coral Silice 70A 10µm 1kg
<b>C8</b>	
B-70A-10-CO-C8-100G	Bulk Coral 70A 10µm C8 100g
B-70A-10-CO-C8-1KG	Bulk Coral 70A 10µm C8 1kg
B-70A-10-CO-C8W-100G	Bulk Coral 70A 10µm C8 Wax 100g
B-70A-10-CO-C8W-1KG	Bulk Coral 70A 10µm C8 Wax 1kg
B-70A-10-CO-C8P-100G	Bulk Coral 70A 10µm C8 Polar 100g
B-70A-10-CO-C8P-1KG	Bulk Coral 70A 10µm C8 Polar 1kg
B-70A-10-CO-C8T-100G	Bulk Coral 70A 10µm C8 triF 100g
B-70A-10-CO-C8T-1KG	Bulk Coral 70A 10µm C8 triF 1kg
B-70A-10-CO-C8EX-100G	Bulk Coral 70A 10µm C8 extreme 100g
B-70A-10-CO-C8EX-1KG	Bulk Coral 70A 10µm C8 extreme 1kg
B-70A-10-CO-C8PT-100G	Bulk Coral 70A 10µm C8 Polar TriF 100g
B-70A-10-CO-C8PT-1KG	Bulk Coral 70A 10µm C8 Polar TriF 1kg
B-70A-10-CO-C8AQP-100G	Bulk Coral 70A 10µm C8 AQ Plus 100g
B-70A-10-CO-C8AQP-1KG	Bulk Coral 70A 10µm C8 AQ Plus 1kg
B-70A-10-CO-C8AQ-100G	Bulk Coral 70A 10µm C8 AQ 100g
B-70A-10-CO-C8AQ-1KG	Bulk Coral 70A 10µm C8 AQ 1kg
<b>C18</b>	
B-70A-10-CO-C18-100G	Bulk Coral 70A 10µm C18 100g
B-70A-10-CO-C18-1KG	Bulk Coral 70A 10µm C18 1kg
B-70A-10-CO-C18W-100G	Bulk Coral 70A 10µm C18 Wax 100g
B-70A-10-CO-C18W-1KG	Bulk Coral 70A 10µm C18 Wax 1kg
B-70A-10-CO-C18P-100G	Bulk Coral 70A 10µm C18 Polar 100g
B-70A-10-CO-C18P-1KG	Bulk Coral 70A 10µm C18 Polar 1kg
B-70A-10-CO-C18T-100G	Bulk Coral 70A 10µm C18 triF 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-10-CO-C18T-1KG	Bulk Coral 70A 10µm C18 triF 1kg
B-70A-10-CO-C18EX-100G	Bulk Coral 70A 10µm C18 extreme 100g
B-70A-10-CO-C18EX-1KG	Bulk Coral 70A 10µm C18 extreme 1kg
B-70A-10-CO-C18PT-100G	Bulk Coral 70A 10µm C18 Polar TriF 100g
B-70A-10-CO-C18PT-1KG	Bulk Coral 70A 10µm C18 Polar TriF 1kg
B-70A-10-CO-C18AQP-100G	Bulk Coral 70A 10µm C18 AQ Plus 100g
B-70A-10-CO-C18AQP-1KG	Bulk Coral 70A 10µm C18 AQ Plus 1kg
B-70A-10-CO-C18AQ-100G	Bulk Coral 70A 10µm C18 AQ 100g
B-70A-10-CO-C18AQ-1KG	Bulk Coral 70A 10µm C18 AQ 1kg
<b>C4</b>	
B-70A-10-CO-C4-100G	Bulk Coral 70A 10µm C4 100g
B-70A-10-CO-C4-1KG	Bulk Coral 70A 10µm C4 1kg
B-70A-10-CO-C4W-100G	Bulk Coral 70A 10µm C4 Wax 100g
B-70A-10-CO-C4W-1KG	Bulk Coral 70A 10µm C4 Wax 1kg
B-70A-10-CO-C4P-100G	Bulk Coral 70A 10µm C4 Polar 100g
B-70A-10-CO-C4P-1KG	Bulk Coral 70A 10µm C4 Polar 1kg
B-70A-10-CO-C4T-100G	Bulk Coral 70A 10µm C4 TriF100g
B-70A-10-CO-C4T-1KG	Bulk Coral 70A 10µm C4 TriF1kg
B-70A-10-CO-C4EX-100G	Bulk Coral 70A 10µm C4 Extreme 100g
B-70A-10-CO-C4EX-1KG	Bulk Coral 70A 10µm C4 Extreme 1kg
B-70A-10-CO-C4PT-100G	Bulk Coral 70A 10µm C4 Polar TriF 100g
B-70A-10-CO-C4PT-1KG	Bulk Coral 70A 10µm C4 Polar TriF 1kg
B-70A-10-CO-C4AQP-100G	Bulk Coral 70A 10µm C4 AQ Plus 100g
B-70A-10-CO-C4AQP-1KG	Bulk Coral 70A 10µm C4 AQ Plus 1kg
B-70A-10-CO-C4AQ-100G	Bulk Coral 70A 10µm C4 AQ 100g
B-70A-10-CO-C4AQ-1KG	Bulk Coral 70A 10µm C4 AQ 1kg
<b>Autres</b>	
B-70A-10-CO-COOH-100G	Bulk Coral 70A 10µm COOH 100g
B-70A-10-CO-COOH-1KG	Bulk Coral 70A 10µm COOH 1kg
B-70A-10-CO-SO3H-100G	Bulk Coral 70A 10µm SO3H 100g
B-70A-10-CO-SO3H-1KG	Bulk Coral 70A 10µm SO3H 1kg
B-70A-10-CO-DNH-100G	Bulk Coral 70A 10µm DNH 100g
B-70A-10-CO-DNH-1KG	Bulk Coral 70A 10µm DNH 1kg
B-70A-10-CO-PEI-100G	Bulk Coral 70A 10µm PEI 100g
B-70A-10-CO-PEI-1KG	Bulk Coral 70A 10µm PEI 1kg
B-70A-10-CO-ARG-100G	Bulk Coral 70A 10µm ARG 100g
B-70A-10-CO-ARG-1KG	Bulk Coral 70A 10µm ARG 1kg
B-70A-10-CO-EP2TNE-100G	Bulk Coral 70A 10µm EP2TNE 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-10-CO-EP2TNE-1KG	Bulk Coral 70A 10µm EP2TNE 1kg
B-70A-10-CO-EP4TNE-100G	Bulk Coral 70A 10µm EP4TNE 100g
B-70A-10-CO-EP4TNE-1KG	Bulk Coral 70A 10µm EP4TNE 1kg
B-70A-10-CO-NH2-100G	Bulk Coral 70A 10µm NH2 100g
B-70A-10-CO-NH2-1KG	Bulk Coral 70A 10µm NH2 1kg
B-70A-10-CO-CN-100G	Bulk Coral 70A 10µm CN 100g
B-70A-10-CO-CN-1KG	Bulk Coral 70A 10µm CN 1kg
B-70A-10-CO-Phenyl-100G	Bulk Coral 70A 10µm Phenyl 100g
B-70A-10-CO-Phenyl-1KG	Bulk Coral 70A 10µm Phenyl 1kg
B-70A-10-CO-Phenyl HE-100G	Bulk Coral 70A 10µm Phenyl HE 100g
B-70A-10-CO-Phenyl HE-1KG	Bulk Coral 70A 10µm Phenyl HE 1kg
B-70A-10-CO-Diol-100G	Bulk Coral 70A 10µm Diol 100g
B-70A-10-CO-Diol-1KG	Bulk Coral 70A 10µm Diol 1kg
<b>Particule Coral 10, Pore 100A</b>	
B-100A-10-CO-Si-100G	Bulk Coral 100A 10µm 100g
B-100A-10-CO-Si-1KG	Bulk Coral 100A 10µm 1kg
<b>C8</b>	
B-100A-10-CO-C8-100G	Bulk Coral 100A 10µm C8 100g
B-100A-10-CO-C8-1KG	Bulk Coral 100A 10µm C8 1kg
B-100A-10-CO-C8W-100G	Bulk Coral 100A 10µm C8 Wax 100g
B-100A-10-CO-C8W-1KG	Bulk Coral 100A 10µm C8 Wax 1kg
B-100A-10-CO-C8P-100G	Bulk Coral 100A 10µm C8 Polar 100g
B-100A-10-CO-C8P-1KG	Bulk Coral 100A 10µm C8 Polar 1kg
B-100A-10-CO-C8T-100G	Bulk Coral 100A 10µm C8 triF 100g
B-100A-10-CO-C8T-1KG	Bulk Coral 100A 10µm C8 triF 1kg
B-100A-10-CO-C8EX-100G	Bulk Coral 100A 10µm C8 extreme 100g
B-100A-10-CO-C8EX-1KG	Bulk Coral 100A 10µm C8 extreme 1kg
B-100A-10-CO-C8PT-100G	Bulk Coral 100A 10µm C8 Polar TriF 100g
B-100A-10-CO-C8PT-1KG	Bulk Coral 100A 10µm C8 Polar TriF 1kg
B-100A-10-CO-C8AQP-100G	Bulk Coral 100A 10µm C8 AQ Plus 100g
B-100A-10-CO-C8AQP-1KG	Bulk Coral 100A 10µm C8 AQ Plus 1kg
B-100A-10-CO-C8AQ-100G	Bulk Coral 100A 10µm C8 AQ 100g
B-100A-10-CO-C8AQ-1KG	Bulk Coral 100A 10µm C8 AQ 1kg
<b>C18</b>	
B-100A-10-CO-C18-100G	Bulk Coral 100A 10µm C18 100g
B-100A-10-CO-C18-1KG	Bulk Coral 100A 10µm C18 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-10-CO-C18W-100G	Bulk Coral 100A 10µm C18 Wax 100g
B-100A-10-CO-C18W-1KG	Bulk Coral 100A 10µm C18 Wax 1kg
B-100A-10-CO-C18P-100G	Bulk Coral 100A 10µm C18 Polar 100g
B-100A-10-CO-C18P-1KG	Bulk Coral 100A 10µm C18 Polar 1kg
B-100A-10-CO-C18T-100G	Bulk Coral 100A 10µm C18 triF 100g
B-100A-10-CO-C18T-1KG	Bulk Coral 100A 10µm C18 triF 1kg
B-100A-10-CO-C18EX-100G	Bulk Coral 100A 10µm C18 extreme 100g
B-100A-10-CO-C18EX-1KG	Bulk Coral 100A 10µm C18 extreme 1kg
B-100A-10-CO-C18PT-100G	Bulk Coral 100A 10µm C18 Polar TriF 100g
B-100A-10-CO-C18PT-1KG	Bulk Coral 100A 10µm C18 Polar TriF 1kg
B-100A-10-CO-C18AQP-100G	Bulk Coral 100A 10µm C18 AQ Plus 100g
B-100A-10-CO-C18AQP-1KG	Bulk Coral 100A 10µm C18 AQ Plus 1kg
B-100A-10-CO-C18AQ-100G	Bulk Coral 100A 10µm C18 AQ 100g
B-100A-10-CO-C18AQ-1KG	Bulk Coral 100A 10µm C18 AQ 1kg
<b>C4</b>	
B-100A-10-CO-C4-100G	Bulk Coral 100A 10µm C4 100g
B-100A-10-CO-C4-1KG	Bulk Coral 100A 10µm C4 1kg
B-100A-10-CO-C4W-100G	Bulk Coral 100A 10µm C4 Wax 100g
B-100A-10-CO-C4W-1KG	Bulk Coral 100A 10µm C4 Wax 1kg
B-100A-10-CO-C4P-100G	Bulk Coral 100A 10µm C4 Polar 100g
B-100A-10-CO-C4P-1KG	Bulk Coral 100A 10µm C4 Polar 1kg
B-100A-10-CO-C4T-100G	Bulk Coral 100A 10µm C4 TriF100g
B-100A-10-CO-C4T-1KG	Bulk Coral 100A 10µm C4 TriF1kg
B-100A-10-CO-C4EX-100G	Bulk Coral 100A 10µm C4 Extreme 100g
B-100A-10-CO-C4EX-1KG	Bulk Coral 100A 10µm C4 Extreme 1kg
B-100A-10-CO-C4PT-100G	Bulk Coral 100A 10µm C4 Polar TriF 100g
B-100A-10-CO-C4PT-1KG	Bulk Coral 100A 10µm C4 Polar TriF 1kg
B-100A-10-CO-C4AQP-100G	Bulk Coral 100A 10µm C4 AQ Plus 100g
B-100A-10-CO-C4AQP-1KG	Bulk Coral 100A 10µm C4 AQ Plus 1kg
B-100A-10-CO-C4AQ-100G	Bulk Coral 100A 10µm C4 AQ 100g
B-100A-10-CO-C4AQ-1KG	Bulk Coral 100A 10µm C4 AQ 1kg
<b>Autre</b>	
B-100A-10-CO-COOH-100G	Bulk Coral 100A 10µm COOH 100g
B-100A-10-CO-COOH-1KG	Bulk Coral 100A 10µm COOH 1kg
B-100A-10-CO-SO3H-100G	Bulk Coral 100A 10µm SO3H 100g
B-100A-10-CO-SO3H-1KG	Bulk Coral 100A 10µm SO3H 1kg
B-100A-10-CO-DNH-100G	Bulk Coral 100A 10µm DNH 100g
B-100A-10-CO-DNH-1KG	Bulk Coral 100A 10µm DNH 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-10-CO-PEI-100G	Bulk Coral 100A 10µm PEI 100g
B-100A-10-CO-PEI-1KG	Bulk Coral 100A 10µm PEI 1kg
B-100A-10-CO-ARG-100G	Bulk Coral 100A 10µm ARG 100g
B-100A-10-CO-ARG-1KG	Bulk Coral 100A 10µm ARG 1kg
B-100A-10-CO-EP2TNE-100G	Bulk Coral 100A 10µm EP2TNE 100g
B-100A-10-CO-EP2TNE-1KG	Bulk Coral 100A 10µm EP2TNE 1kg
B-100A-10-CO-EP4TNE-100G	Bulk Coral 100A 10µm EP4TNE 100g
B-100A-10-CO-EP4TNE-1KG	Bulk Coral 100A 10µm EP4TNE 1kg
B-100A-10-CO-NH2-100G	Bulk Coral 100A 10µm NH2 100g
B-100A-10-CO-NH2-1KG	Bulk Coral 100A 10µm NH2 1kg
B-100A-10-CO-CN-100G	Bulk Coral 100A 10µm CN 100g
B-100A-10-CO-CN-1KG	Bulk Coral 100A 10µm CN 1kg
B-100A-10-CO-Phenyl-100G	Bulk Coral 100A 10µm Phenyl 100g
B-100A-10-CO-Phenyl-1KG	Bulk Coral 100A 10µm Phenyl 1kg
B-100A-10-CO-Phenyl HE-100G	Bulk Coral 100A 10µm Phenyl HE 100g
B-100A-10-CO-Phenyl HE-1KG	Bulk Coral 100A 10µm Phenyl HE 1kg
B-100A-10-CO-Diol-100G	Bulk Coral 100A 10µm Diol 100g
B-100A-10-CO-Diol-1KG	Bulk Coral 100A 10µm Diol 1kg
<b>Particule Coral 10, Pore 150A</b>	
B-150A-10-CO-Si-100G	Bulk Coral 150A 10µm 100g
B-150A-10-CO-Si-1KG	Bulk Coral 150A 10µm 1kg
<b>C8</b>	
B-150A-10-CO-C8-100G	Bulk Coral 150A 10µm C8 100g
B-150A-10-CO-C8-1KG	Bulk Coral 150A 10µm C8 1kg
B-150A-10-CO-C8W-100G	Bulk Coral 150A 10µm C8 Wax 100g
B-150A-10-CO-C8W-1KG	Bulk Coral 150A 10µm C8 Wax 1kg
B-150A-10-CO-C8P-100G	Bulk Coral 150A 10µm C8 Polar 100g
B-150A-10-CO-C8P-1KG	Bulk Coral 150A 10µm C8 Polar 1kg
B-150A-10-CO-C8T-100G	Bulk Coral 150A 10µm C8 triF 100g
B-150A-10-CO-C8T-1KG	Bulk Coral 150A 10µm C8 triF 1kg
B-150A-10-CO-C8EX-100G	Bulk Coral 150A 10µm C8 extreme 100g
B-150A-10-CO-C8EX-1KG	Bulk Coral 150A 10µm C8 extreme 1kg
B-150A-10-CO-C8PT-100G	Bulk Coral 150A 10µm C8 Polar TriF 100g
B-150A-10-CO-C8PT-1KG	Bulk Coral 150A 10µm C8 Polar TriF 1kg
B-150A-10-CO-C8AQP-100G	Bulk Coral 150A 10µm C8 AQ Plus 100g
B-150A-10-CO-C8AQP-1KG	Bulk Coral 150A 10µm C8 AQ Plus 1kg
B-150A-10-CO-C8AQ-100G	Bulk Coral 150A 10µm C8 AQ 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-10-CO-C8AQ-1KG	Bulk Coral 150A 10µm C8 AQ 1kg
<b>C18</b>	
B-150A-10-CO-C18-100G	Bulk Coral 150A 10µm C18 100g
B-150A-10-CO-C18-1KG	Bulk Coral 150A 10µm C18 1kg
B-150A-10-CO-C18W-100G	Bulk Coral 150A 10µm C18 Wax 100g
B-150A-10-CO-C18W-1KG	Bulk Coral 150A 10µm C18 Wax 1kg
B-150A-10-CO-C18P-100G	Bulk Coral 150A 10µm C18 Polar 100g
B-150A-10-CO-C18P-1KG	Bulk Coral 150A 10µm C18 Polar 1kg
B-150A-10-CO-C18T-100G	Bulk Coral 150A 10µm C18 triF 100g
B-150A-10-CO-C18T-1KG	Bulk Coral 150A 10µm C18 triF 1kg
B-150A-10-CO-C18EX-100G	Bulk Coral 150A 10µm C18 extreme 100g
B-150A-10-CO-C18EX-1KG	Bulk Coral 150A 10µm C18 extreme 1kg
B-150A-10-CO-C18PT-100G	Bulk Coral 150A 10µm C18 Polar TriF 100g
B-150A-10-CO-C18PT-1KG	Bulk Coral 150A 10µm C18 Polar TriF 1kg
B-150A-10-CO-C18AQP-100G	Bulk Coral 150A 10µm C18 AQ Plus 100g
B-150A-10-CO-C18AQP-1KG	Bulk Coral 150A 10µm C18 AQ Plus 1kg
B-150A-10-CO-C18AQ-100G	Bulk Coral 150A 10µm C18 AQ 100g
B-150A-10-CO-C18AQ-1KG	Bulk Coral 150A 10µm C18 AQ 1kg
<b>C4</b>	
B-150A-10-CO-C4-100G	Bulk Coral 150A 10µm C4 100g
B-150A-10-CO-C4-1KG	Bulk Coral 150A 10µm C4 1kg
B-150A-10-CO-C4W-100G	Bulk Coral 150A 10µm C4 Wax 100g
B-150A-10-CO-C4W-1KG	Bulk Coral 150A 10µm C4 Wax 1kg
B-150A-10-CO-C4P-100G	Bulk Coral 150A 10µm C4 Polar 100g
B-150A-10-CO-C4P-1KG	Bulk Coral 150A 10µm C4 Polar 1kg
B-150A-10-CO-C4T-100G	Bulk Coral 150A 10µm C4 TriF100g
B-150A-10-CO-C4T-1KG	Bulk Coral 150A 10µm C4 TriF1kg
B-150A-10-CO-C4EX-100G	Bulk Coral 150A 10µm C4 Extreme 100g
B-150A-10-CO-C4EX-1KG	Bulk Coral 150A 10µm C4 Extreme 1kg
B-150A-10-CO-C4PT-100G	Bulk Coral 150A 10µm C4 Polar TriF 100g
B-150A-10-CO-C4PT-1KG	Bulk Coral 150A 10µm C4 Polar TriF 1kg
B-150A-10-CO-C4AQP-100G	Bulk Coral 150A 10µm C4 AQ Plus 100g
B-150A-10-CO-C4AQP-1KG	Bulk Coral 150A 10µm C4 AQ Plus 1kg
B-150A-10-CO-C4AQ-100G	Bulk Coral 150A 10µm C4 AQ 100g
B-150A-10-CO-C4AQ-1KG	Bulk Coral 150A 10µm C4 AQ 1kg
<b>Autres</b>	
B-150A-10-CO-COOH-100G	Bulk Coral 150A 10µm COOH 100g
B-150A-10-CO-COOH-1KG	Bulk Coral 150A 10µm COOH 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-10-CO-SO3H-100G	Bulk Coral 150A 10µm SO3H 100g
B-150A-10-CO-SO3H-1KG	Bulk Coral 150A 10µm SO3H 1kg
B-150A-10-CO-DNH-100G	Bulk Coral 150A 10µm DNH 100g
B-150A-10-CO-DNH-1KG	Bulk Coral 150A 10µm DNH 1kg
B-150A-10-CO-PEI-100G	Bulk Coral 150A 10µm PEI 100g
B-150A-10-CO-PEI-1KG	Bulk Coral 150A 10µm PEI 1kg
B-150A-10-CO-ARG-100G	Bulk Coral 150A 10µm ARG 100g
B-150A-10-CO-ARG-1KG	Bulk Coral 150A 10µm ARG 1kg
B-150A-10-CO-EP2TNE-100G	Bulk Coral 150A 10µm EP2TNE 100g
B-150A-10-CO-EP2TNE-1KG	Bulk Coral 150A 10µm EP2TNE 1kg
B-150A-10-CO-EP4TNE-100G	Bulk Coral 150A 10µm EP4TNE 100g
B-150A-10-CO-EP4TNE-1KG	Bulk Coral 150A 10µm EP4TNE 1kg
B-150A-10-CO-NH2-100G	Bulk Coral 150A 10µm NH2 100g
B-150A-10-CO-NH2-1KG	Bulk Coral 150A 10µm NH2 1kg
B-150A-10-CO-CN-100G	Bulk Coral 150A 10µm CN 100g
B-150A-10-CO-CN-1KG	Bulk Coral 150A 10µm CN 1kg
B-150A-10-CO-Phenyl-100G	Bulk Coral 150A 10µm Phenyl 100g
B-150A-10-CO-Phenyl-1KG	Bulk Coral 150A 10µm Phenyl 1kg
B-150A-10-CO-Phenyl HE-100G	Bulk Coral 150A 10µm Phenyl HE 100g
B-150A-10-CO-Phenyl HE-1KG	Bulk Coral 150A 10µm Phenyl HE 1kg
B-150A-10-CO-Diol-100G	Bulk Coral 150A 10µm Diol 100g
B-150A-10-CO-Diol-1KG	Bulk Coral 150A 10µm Diol 1kg
<b>Particule Coral 10, Pore 200A</b>	
B-200A-10-CO-Si-100G	Bulk Coral 200A 10µm Si 100g
B-200A-10-CO-Si-1KG	Bulk Coral 200A 10µm Si 1kg
<b>C8</b>	
B-200A-10-CO-C8-100G	Bulk Coral 200A 10µm C8 100g
B-200A-10-CO-C8-1KG	Bulk Coral 200A 10µm C8 1kg
B-200A-10-CO-C8W-100G	Bulk Coral 200A 10µm C8 Wax 100g
B-200A-10-CO-C8W-1KG	Bulk Coral 200A 10µm C8 Wax 1kg
B-200A-10-CO-C8P-100G	Bulk Coral 200A 10µm C8 Polar 100g
B-200A-10-CO-C8P-1KG	Bulk Coral 200A 10µm C8 Polar 1kg
B-200A-10-CO-C8T-100G	Bulk Coral 200A 10µm C8 triF 100g
B-200A-10-CO-C8T-1KG	Bulk Coral 200A 10µm C8 triF 1kg
B-200A-10-CO-C8EX-100G	Bulk Coral 200A 10µm C8 extreme 100g
B-200A-10-CO-C8EX-1KG	Bulk Coral 200A 10µm C8 extreme 1kg
B-200A-10-CO-C8PT-100G	Bulk Coral 200A 10µm C8 Polar TriF 100g
B-200A-10-CO-C8PT-1KG	Bulk Coral 200A 10µm C8 Polar TriF 1kg
B-200A-10-CO-C8AQP-100G	Bulk Coral 200A 10µm C8 AQ Plus 100g
B-200A-10-CO-C8AQP-1KG	Bulk Coral 200A 10µm C8 AQ Plus 1kg
B-200A-10-CO-C8AQ-100G	Bulk Coral 200A 10µm C8 AQ 100g
B-200A-10-CO-C8AQ-1KG	Bulk Coral 200A 10µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-10-CO-C8AQP-100G	Bulk Coral 200A 10µm C8 AQ Plus 100g
B-200A-10-CO-C8AQP-1KG	Bulk Coral 200A 10µm C8 AQ Plus 1kg
B-200A-10-CO-C8AQ-100G	Bulk Coral 200A 10µm C8 AQ 100g
B-200A-10-CO-C8AQ-1KG	Bulk Coral 200A 10µm C8 AQ 1kg
<b>C18</b>	
B-200A-10-CO-C18-100G	Bulk Coral 200A 10µm C18 100g
B-200A-10-CO-C18-1KG	Bulk Coral 200A 10µm C18 1kg
B-200A-10-CO-C18W-100G	Bulk Coral 200A 10µm C18 Wax 100g
B-200A-10-CO-C18W-1KG	Bulk Coral 200A 10µm C18 Wax 1kg
B-200A-10-CO-C18P-100G	Bulk Coral 200A 10µm C18 Polar 100g
B-200A-10-CO-C18P-1KG	Bulk Coral 200A 10µm C18 Polar 1kg
B-200A-10-CO-C18T-100G	Bulk Coral 200A 10µm C18 triF 100g
B-200A-10-CO-C18T-1KG	Bulk Coral 200A 10µm C18 triF 1kg
B-200A-10-CO-C18EX-100G	Bulk Coral 200A 10µm C18 extreme 100g
B-200A-10-CO-C18EX-1KG	Bulk Coral 200A 10µm C18 extreme 1kg
B-200A-10-CO-C18PT-100G	Bulk Coral 200A 10µm C18 Polar TriF 100g
B-200A-10-CO-C18PT-1KG	Bulk Coral 200A 10µm C18 Polar TriF 1kg
B-200A-10-CO-C18AQP-100G	Bulk Coral 200A 10µm C18 AQ Plus 100g
B-200A-10-CO-C18AQP-1KG	Bulk Coral 200A 10µm C18 AQ Plus 1kg
B-200A-10-CO-C18AQ-100G	Bulk Coral 200A 10µm C18 AQ 100g
B-200A-10-CO-C18AQ-1KG	Bulk Coral 200A 10µm C18 AQ 1kg
<b>C4</b>	
B-200A-10-CO-C4-100G	Bulk Coral 200A 10µm C4 100g
B-200A-10-CO-C4-1KG	Bulk Coral 200A 10µm C4 1kg
B-200A-10-CO-C4W-100G	Bulk Coral 200A 10µm C4 Wax 100g
B-200A-10-CO-C4W-1KG	Bulk Coral 200A 10µm C4 Wax 1kg
B-200A-10-CO-C4P-100G	Bulk Coral 200A 10µm C4 Polar 100g
B-200A-10-CO-C4P-1KG	Bulk Coral 200A 10µm C4 Polar 1kg
B-200A-10-CO-C4T-100G	Bulk Coral 200A 10µm C4 TriF100g
B-200A-10-CO-C4T-1KG	Bulk Coral 200A 10µm C4 TriF1kg
B-200A-10-CO-C4EX-100G	Bulk Coral 200A 10µm C4 Extreme 100g
B-200A-10-CO-C4EX-1KG	Bulk Coral 200A 10µm C4 Extreme 1kg
B-200A-10-CO-C4PT-100G	Bulk Coral 200A 10µm C4 Polar TriF 100g
B-200A-10-CO-C4PT-1KG	Bulk Coral 200A 10µm C4 Polar TriF 1kg
B-200A-10-CO-C4AQP-100G	Bulk Coral 200A 10µm C4 AQ Plus 100g
B-200A-10-CO-C4AQP-1KG	Bulk Coral 200A 10µm C4 AQ Plus 1kg
B-200A-10-CO-C4AQ-100G	Bulk Coral 200A 10µm C4 AQ 100g
B-200A-10-CO-C4AQ-1KG	Bulk Coral 200A 10µm C4 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Autres</b>	
B-200A-10-CO-COOH-100G	Bulk Coral 200A 10µm COOH 100g
B-200A-10-CO-COOH-1KG	Bulk Coral 200A 10µm COOH 1kg
B-200A-10-CO-S03H-100G	Bulk Coral 200A 10µm S03H 100g
B-200A-10-CO-S03H-1KG	Bulk Coral 200A 10µm S03H 1kg
B-200A-10-CO-DNH-100G	Bulk Coral 200A 10µm DNH 100g
B-200A-10-CO-DNH-1KG	Bulk Coral 200A 10µm DNH 1kg
B-200A-10-CO-PEI-100G	Bulk Coral 200A 10µm PEI 100g
B-200A-10-CO-PEI-1KG	Bulk Coral 200A 10µm PEI 1kg
B-200A-10-CO-ARG-100G	Bulk Coral 200A 10µm ARG 100g
B-200A-10-CO-ARG-1KG	Bulk Coral 200A 10µm ARG 1kg
B-200A-10-CO-EP2TNE-100G	Bulk Coral 200A 10µm EP2TNE 100g
B-200A-10-CO-EP2TNE-1KG	Bulk Coral 200A 10µm EP2TNE 1kg
B-200A-10-CO-EP4TNE-100G	Bulk Coral 200A 10µm EP4TNE 100g
B-200A-10-CO-EP4TNE-1KG	Bulk Coral 200A 10µm EP4TNE 1kg
B-200A-10-CO-NH2-100G	Bulk Coral 200A 10µm NH2 100g
B-200A-10-CO-NH2-1KG	Bulk Coral 200A 10µm NH2 1kg
B-200A-10-CO-CN-100G	Bulk Coral 200A 10µm CN 100g
B-200A-10-CO-CN-1KG	Bulk Coral 200A 10µm CN 1kg
B-200A-10-CO-Phenyl-100G	Bulk Coral 200A 10µm Phenyl 100g
B-200A-10-CO-Phenyl-1KG	Bulk Coral 200A 10µm Phenyl 1kg
B-200A-10-CO-Phenyl HE-100G	Bulk Coral 200A 10µm Phenyl HE 100g
B-200A-10-CO-Phenyl HE-1KG	Bulk Coral 200A 10µm Phenyl HE 1kg
B-200A-10-CO-Diol-100G	Bulk Coral 200A 10µm Diol 100g
B-200A-10-CO-Diol-1KG	Bulk Coral 200A 10µm Diol 1kg
<b>Particule Coral 10, Pore 300A</b>	
B-300A-10-CO-Si-100G	Bulk Coral 300A 10µm Si 100g
B-300A-10-CO-Si-1KG	Bulk Coral 300A 10µm Si 1kg
<b>C8</b>	
B-300A-10-CO-C8-100G	Bulk Coral 300A 10µm C8 100g
B-300A-10-CO-C8-1KG	Bulk Coral 300A 10µm C8 1kg
B-300A-10-CO-C8W-100G	Bulk Coral 300A 10µm C8 Wax 100g
B-300A-10-CO-C8W-1KG	Bulk Coral 300A 10µm C8 Wax 1kg
B-300A-10-CO-C8P-100G	Bulk Coral 300A 10µm C8 Polar 100g
B-300A-10-CO-C8P-1KG	Bulk Coral 300A 10µm C8 Polar 1kg
B-300A-10-CO-C8T-100G	Bulk Coral 300A 10µm C8 triF 100g
B-300A-10-CO-C8T-1KG	Bulk Coral 300A 10µm C8 triF 1kg
B-300A-10-CO-C8EX-100G	Bulk Coral 300A 10µm C8 extreme 100g
B-300A-10-CO-C8EX-1KG	Bulk Coral 300A 10µm C8 extreme 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-10-CO-C8EX-1KG	Bulk Coral 300A 10µm C8 extreme 1kg
B-300A-10-CO-C8PT-100G	Bulk Coral 300A 10µm C8 Polar TriF 100g
B-300A-10-CO-C8PT-1KG	Bulk Coral 300A 10µm C8 Polar TriF 1kg
B-300A-10-CO-C8AQP-100G	Bulk Coral 300A 10µm C8 AQ Plus 100g
B-300A-10-CO-C8AQP-1KG	Bulk Coral 300A 10µm C8 AQ Plus 1kg
B-300A-10-CO-C8AQ-100G	Bulk Coral 300A 10µm C8 AQ 100g
B-300A-10-CO-C8AQ-1KG	Bulk Coral 300A 10µm C8 AQ 1kg
<b>C18</b>	
B-300A-10-CO-C18-100G	Bulk Coral 300A 10µm C18 100g
B-300A-10-CO-C18-1KG	Bulk Coral 300A 10µm C18 1kg
B-300A-10-CO-C18W-100G	Bulk Coral 300A 10µm C18 Wax 100g
B-300A-10-CO-C18W-1KG	Bulk Coral 300A 10µm C18 Wax 1kg
B-300A-10-CO-C18P-100G	Bulk Coral 300A 10µm C18 Polar 100g
B-300A-10-CO-C18P-1KG	Bulk Coral 300A 10µm C18 Polar 1kg
B-300A-10-CO-C18T-100G	Bulk Coral 300A 10µm C18 triF 100g
B-300A-10-CO-C18T-1KG	Bulk Coral 300A 10µm C18 triF 1kg
B-300A-10-CO-C18EX-100G	Bulk Coral 300A 10µm C18 extreme 100g
B-300A-10-CO-C18EX-1KG	Bulk Coral 300A 10µm C18 extreme 1kg
B-300A-10-CO-C18PT-100G	Bulk Coral 300A 10µm C18 Polar TriF 100g
B-300A-10-CO-C18PT-1KG	Bulk Coral 300A 10µm C18 Polar TriF 1kg
B-300A-10-CO-C18AQP-100G	Bulk Coral 300A 10µm C18 AQ Plus 100g
B-300A-10-CO-C18AQP-1KG	Bulk Coral 300A 10µm C18 AQ Plus 1kg
B-300A-10-CO-C18AQ-100G	Bulk Coral 300A 10µm C18 AQ 100g
B-300A-10-CO-C18AQ-1KG	Bulk Coral 300A 10µm C18 AQ 1kg
<b>C4</b>	
B-300A-10-CO-C4-100G	Bulk Coral 300A 10µm C4 100g
B-300A-10-CO-C4-1KG	Bulk Coral 300A 10µm C4 1kg
B-300A-10-CO-C4W-100G	Bulk Coral 300A 10µm C4 Wax 100g
B-300A-10-CO-C4W-1KG	Bulk Coral 300A 10µm C4 Wax 1kg
B-300A-10-CO-C4P-100G	Bulk Coral 300A 10µm C4 Polar 100g
B-300A-10-CO-C4P-1KG	Bulk Coral 300A 10µm C4 Polar 1kg
B-300A-10-CO-C4T-100G	Bulk Coral 300A 10µm C4 TriF100g
B-300A-10-CO-C4T-1KG	Bulk Coral 300A 10µm C4 TriF1kg
B-300A-10-CO-C4EX-100G	Bulk Coral 300A 10µm C4 Extreme 100g
B-300A-10-CO-C4EX-1KG	Bulk Coral 300A 10µm C4 Extreme 1kg
B-300A-10-CO-C4PT-100G	Bulk Coral 300A 10µm C4 Polar TriF 100g
B-300A-10-CO-C4PT-1KG	Bulk Coral 300A 10µm C4 Polar TriF 1kg
B-300A-10-CO-C4AQP-100G	Bulk Coral 300A 10µm C4 AQ Plus 100g
B-300A-10-CO-C4AQP-1KG	Bulk Coral 300A 10µm C4 AQ Plus 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-10-CO-C4AQ-100G	Bulk Coral 300A 10µm C4 AQ 100g
B-300A-10-CO-C4AQ-1KG	Bulk Coral 300A 10µm C4 AQ 1kg
<b>Autres</b>	
B-300A-10-CO-COOH-100G	Bulk Coral 300A 10µm COOH 100g
B-300A-10-CO-COOH-1KG	Bulk Coral 300A 10µm COOH 1kg
B-300A-10-CO-SO3H-100G	Bulk Coral 300A 10µm SO3H 100g
B-300A-10-CO-SO3H-1KG	Bulk Coral 300A 10µm SO3H 1kg
B-300A-10-CO-DNH-100G	Bulk Coral 300A 10µm DNH 100g
B-300A-10-CO-DNH-1KG	Bulk Coral 300A 10µm DNH 1kg
B-300A-10-CO-PEI-100G	Bulk Coral 300A 10µm PEI 100g
B-300A-10-CO-PEI-1KG	Bulk Coral 300A 10µm PEI 1kg
B-300A-10-CO-ARG-100G	Bulk Coral 300A 10µm ARG 100g
B-300A-10-CO-ARG-1KG	Bulk Coral 300A 10µm ARG 1kg
B-300A-10-CO-EP2TNE-100G	Bulk Coral 300A 10µm EP2TNE 100g
B-300A-10-CO-EP2TNE-1KG	Bulk Coral 300A 10µm EP2TNE 1kg
B-300A-10-CO-EP4TNE-100G	Bulk Coral 300A 10µm EP4TNE 100g
B-300A-10-CO-EP4TNE-1KG	Bulk Coral 300A 10µm EP4TNE 1kg
B-300A-10-CO-NH2-100G	Bulk Coral 300A 10µm NH2 100g
B-300A-10-CO-NH2-1KG	Bulk Coral 300A 10µm NH2 1kg
B-300A-10-CO-CN-100G	Bulk Coral 300A 10µm CN 100g
B-300A-10-CO-CN-1KG	Bulk Coral 300A 10µm CN 1kg
B-300A-10-CO-Phenyl-100G	Bulk Coral 300A 10µm Phenyl 100g
B-300A-10-CO-Phenyl-1KG	Bulk Coral 300A 10µm Phenyl 1kg
B-300A-10-CO-Phenyl HE-100G	Bulk Coral 300A 10µm Phenyl HE 100g
B-300A-10-CO-Phenyl HE-1KG	Bulk Coral 300A 10µm Phenyl HE 1kg
B-300A-10-CO-Diol-100G	Bulk Coral 300A 10µm Diol 100g
B-300A-10-CO-Diol-1KG	Bulk Coral 300A 10µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Coral 10, Pore 500A</b>	
B-500A-10-CO-Si-100G	Bulk Coral 500A 10µm 100g
B-500A-10-CO-Si-1KG	Bulk Coral 500A 10µm 1kg
<b>C8</b>	
B-500A-10-CO-C8-100G	Bulk Coral 500A 10µm C8 100g
B-500A-10-CO-C8-1KG	Bulk Coral 500A 10µm C8 1kg
B-500A-10-CO-C8W-100G	Bulk Coral 500A 10µm C8 Wax 100g
B-500A-10-CO-C8W-1KG	Bulk Coral 500A 10µm C8 Wax 1kg
B-500A-10-CO-C8P-100G	Bulk Coral 500A 10µm C8 Polar 100g
B-500A-10-CO-C8P-1KG	Bulk Coral 500A 10µm C8 Polar 1kg
B-500A-10-CO-C8T-100G	Bulk Coral 500A 10µm C8 triF 100g
B-500A-10-CO-C8T-1KG	Bulk Coral 500A 10µm C8 triF 1kg
B-500A-10-CO-C8EX-100G	Bulk Coral 500A 10µm C8 extreme 100g
B-500A-10-CO-C8EX-1KG	Bulk Coral 500A 10µm C8 extreme 1kg
B-500A-10-CO-C8PT-100G	Bulk Coral 500A 10µm C8 Polar TriF 100g
B-500A-10-CO-C8PT-1KG	Bulk Coral 500A 10µm C8 Polar TriF 1kg
B-500A-10-CO-C8AQP-100G	Bulk Coral 500A 10µm C8 AQ Plus 100g
B-500A-10-CO-C8AQP-1KG	Bulk Coral 500A 10µm C8 AQ Plus 1kg
B-500A-10-CO-C8AQ-100G	Bulk Coral 500A 10µm C8 AQ 100g
B-500A-10-CO-C8AQ-1KG	Bulk Coral 500A 10µm C8 AQ 1kg
<b>C18</b>	
B-500A-10-CO-C18-100G	Bulk Coral 500A 10µm C18 100g
B-500A-10-CO-C18-1KG	Bulk Coral 500A 10µm C18 1kg
B-500A-10-CO-C18W-100G	Bulk Coral 500A 10µm C18 Wax 100g
B-500A-10-CO-C18W-1KG	Bulk Coral 500A 10µm C18 Wax 1kg
B-500A-10-CO-C18P-100G	Bulk Coral 500A 10µm C18 Polar 100g
B-500A-10-CO-C18P-1KG	Bulk Coral 500A 10µm C18 Polar 1kg
B-500A-10-CO-C18T-100G	Bulk Coral 500A 10µm C18 triF 100g
B-500A-10-CO-C18T-1KG	Bulk Coral 500A 10µm C18 triF 1kg
B-500A-10-CO-C18EX-100G	Bulk Coral 500A 10µm C18 extreme 100g
B-500A-10-CO-C18EX-1KG	Bulk Coral 500A 10µm C18 extreme 1kg
B-500A-10-CO-C18PT-100G	Bulk Coral 500A 10µm C18 Polar TriF 100g
B-500A-10-CO-C18PT-1KG	Bulk Coral 500A 10µm C18 Polar TriF 1kg
B-500A-10-CO-C18AQP-100G	Bulk Coral 500A 10µm C18 AQ Plus 100g
B-500A-10-CO-C18AQP-1KG	Bulk Coral 500A 10µm C18 AQ Plus 1kg
B-500A-10-CO-C18AQ-100G	Bulk Coral 500A 10µm C18 AQ 100g
B-500A-10-CO-C18AQ-1KG	Bulk Coral 500A 10µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-500A-10-CO-C4-100G	Bulk Coral 500A 10µm C4 100g
B-500A-10-CO-C4-1KG	Bulk Coral 500A 10µm C4 1kg
B-500A-10-CO-C4W-100G	Bulk Coral 500A 10µm C4 Wax 100g
B-500A-10-CO-C4W-1KG	Bulk Coral 500A 10µm C4 Wax 1kg
B-500A-10-CO-C4P-100G	Bulk Coral 500A 10µm C4 Polar 100g
B-500A-10-CO-C4P-1KG	Bulk Coral 500A 10µm C4 Polar 1kg
B-500A-10-CO-C4T-100G	Bulk Coral 500A 10µm C4 TriF100g
B-500A-10-CO-C4T-1KG	Bulk Coral 500A 10µm C4 TriF1kg
B-500A-10-CO-C4EX-100G	Bulk Coral 500A 10µm C4 Extreme 100g
B-500A-10-CO-C4EX-1KG	Bulk Coral 500A 10µm C4 Extreme 1kg
B-500A-10-CO-C4PT-100G	Bulk Coral 500A 10µm C4 Polar TriF 100g
B-500A-10-CO-C4PT-1KG	Bulk Coral 500A 10µm C4 Polar TriF 1kg
B-500A-10-CO-C4AQP-100G	Bulk Coral 500A 10µm C4 AQ Plus 100g
B-500A-10-CO-C4AQP-1KG	Bulk Coral 500A 10µm C4 AQ Plus 1kg
B-500A-10-CO-C4AQ-100G	Bulk Coral 500A 10µm C4 AQ 100g
B-500A-10-CO-C4AQ-1KG	Bulk Coral 500A 10µm C4 AQ 1kg
<b>Autres</b>	
B-500A-10-CO-COOH-100G	Bulk Coral 500A 10µm COOH 100g
B-500A-10-CO-COOH-1KG	Bulk Coral 500A 10µm COOH 1kg
B-500A-10-CO-SO3H-100G	Bulk Coral 500A 10µm SO3H 100g
B-500A-10-CO-SO3H-1KG	Bulk Coral 500A 10µm SO3H 1kg
B-500A-10-CO-DNH-100G	Bulk Coral 500A 10µm DNH 100g
B-500A-10-CO-DNH-1KG	Bulk Coral 500A 10µm DNH 1kg
B-500A-10-CO-PEI-100G	Bulk Coral 500A 10µm PEI 100g
B-500A-10-CO-PEI-1KG	Bulk Coral 500A 10µm PEI 1kg
B-500A-10-CO-ARG-100G	Bulk Coral 500A 10µm ARG 100g
B-500A-10-CO-ARG-1KG	Bulk Coral 500A 10µm ARG 1kg
B-500A-10-CO-EP2TNE-100G	Bulk Coral 500A 10µm EP2TNE 100g
B-500A-10-CO-EP2TNE-1KG	Bulk Coral 500A 10µm EP2TNE 1kg
B-500A-10-CO-EP4TNE-100G	Bulk Coral 500A 10µm EP4TNE 100g
B-500A-10-CO-EP4TNE-1KG	Bulk Coral 500A 10µm EP4TNE 1kg
B-500A-10-CO-NH2-100G	Bulk Coral 500A 10µm NH2 100g
B-500A-10-CO-NH2-1KG	Bulk Coral 500A 10µm NH2 1kg
B-500A-10-CO-CN-100G	Bulk Coral 500A 10µm CN 100g
B-500A-10-CO-CN-1KG	Bulk Coral 500A 10µm CN 1kg
B-500A-10-CO-Phenyl-100G	Bulk Coral 500A 10µm Phenyl 100g
B-500A-10-CO-Phenyl-1KG	Bulk Coral 500A 10µm Phenyl 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-500A-10-CO-Phenyl HE-100G	Bulk Coral 500A 10µm Phenyl HE 100g
B-500A-10-CO-Phenyl HE-1KG	Bulk Coral 500A 10µm Phenyl HE 1kg
B-500A-10-CO-Diol-100G	Bulk Coral 500A 10µm Diol 100g
B-500A-10-CO-Diol-1KG	Bulk Coral 500A 10µm Diol 1kg
<b>Particule Coral 10, Pore 800A</b>	
B-800A-10-CO-Si-100G	Bulk Coral 800A 10µm 100g
B-800A-10-CO-Si-1KG	Bulk Coral 800A 10µm 1kg
<b>C8</b>	
B-800A-10-CO-C8-100G	Bulk Coral 800A 10µm C8 100g
B-800A-10-CO-C8-1KG	Bulk Coral 800A 10µm C8 1kg
B-800A-10-CO-C8W-100G	Bulk Coral 800A 10µm C8 Wax 100g
B-800A-10-CO-C8W-1KG	Bulk Coral 800A 10µm C8 Wax 1kg
B-800A-10-CO-C8P-100G	Bulk Coral 800A 10µm C8 Polar 100g
B-800A-10-CO-C8P-1KG	Bulk Coral 800A 10µm C8 Polar 1kg
B-800A-10-CO-C8T-100G	Bulk Coral 800A 10µm C8 triF 100g
B-800A-10-CO-C8T-1KG	Bulk Coral 800A 10µm C8 triF 1kg
B-800A-10-CO-C8EX-100G	Bulk Coral 800A 10µm C8 extreme 100g
B-800A-10-CO-C8EX-1KG	Bulk Coral 800A 10µm C8 extreme 1kg
B-800A-10-CO-C8PT-100G	Bulk Coral 800A 10µm C8 Polar TriF 100g
B-800A-10-CO-C8PT-1KG	Bulk Coral 800A 10µm C8 Polar TriF 1kg
B-800A-10-CO-C8AQP-100G	Bulk Coral 800A 10µm C8 AQ Plus 100g
B-800A-10-CO-C8AQP-1KG	Bulk Coral 800A 10µm C8 AQ Plus 1kg
B-800A-10-CO-C8AQ-100G	Bulk Coral 800A 10µm C8 AQ 100g
B-800A-10-CO-C8AQ-1KG	Bulk Coral 800A 10µm C8 AQ 1kg
<b>C18</b>	
B-800A-10-CO-C18-100G	Bulk Coral 800A 10µm C18 100g
B-800A-10-CO-C18-1KG	Bulk Coral 800A 10µm C18 1kg
B-800A-10-CO-C18W-100G	Bulk Coral 800A 10µm C18 Wax 100g
B-800A-10-CO-C18W-1KG	Bulk Coral 800A 10µm C18 Wax 1kg
B-800A-10-CO-C18P-100G	Bulk Coral 800A 10µm C18 Polar 100g
B-800A-10-CO-C18P-1KG	Bulk Coral 800A 10µm C18 Polar 1kg
B-800A-10-CO-C18T-100G	Bulk Coral 800A 10µm C18 triF 100g
B-800A-10-CO-C18T-1KG	Bulk Coral 800A 10µm C18 triF 1kg
B-800A-10-CO-C18EX-100G	Bulk Coral 800A 10µm C18 extreme 100g
B-800A-10-CO-C18EX-1KG	Bulk Coral 800A 10µm C18 extreme 1kg
B-800A-10-CO-C18PT-100G	Bulk Coral 800A 10µm C18 Polar TriF 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-10-CO-C18PT-1KG	Bulk Coral 800A 10µm C18 Polar TriF 1kg
B-800A-10-CO-C18AQP-100G	Bulk Coral 800A 10µm C18 AQ Plus 100g
B-800A-10-CO-C18AQP-1KG	Bulk Coral 800A 10µm C18 AQ Plus 1kg
B-800A-10-CO-C18AQ-100G	Bulk Coral 800A 10µm C18 AQ 100g
B-800A-10-CO-C18AQ-1KG	Bulk Coral 800A 10µm C18 AQ 1kg
<b>C4</b>	
B-800A-10-CO-C4-100G	Bulk Coral 800A 10µm C4 100g
B-800A-10-CO-C4-1KG	Bulk Coral 800A 10µm C4 1kg
B-800A-10-CO-C4W-100G	Bulk Coral 800A 10µm C4 Wax 100g
B-800A-10-CO-C4W-1KG	Bulk Coral 800A 10µm C4 Wax 1kg
B-800A-10-CO-C4P-100G	Bulk Coral 800A 10µm C4 Polar 100g
B-800A-10-CO-C4P-1KG	Bulk Coral 800A 10µm C4 Polar 1kg
B-800A-10-CO-C4T-100G	Bulk Coral 800A 10µm C4 TriF100g
B-800A-10-CO-C4T-1KG	Bulk Coral 800A 10µm C4 TriF1kg
B-800A-10-CO-C4EX-100G	Bulk Coral 800A 10µm C4 Extreme 100g
B-800A-10-CO-C4EX-1KG	Bulk Coral 800A 10µm C4 Extreme 1kg
B-800A-10-CO-C4PT-100G	Bulk Coral 800A 10µm C4 Polar TriF 100g
B-800A-10-CO-C4PT-1KG	Bulk Coral 800A 10µm C4 Polar TriF 1kg
B-800A-10-CO-C4AQP-100G	Bulk Coral 800A 10µm C4 AQ Plus 100g
B-800A-10-CO-C4AQP-1KG	Bulk Coral 800A 10µm C4 AQ Plus 1kg
B-800A-10-CO-C4AQ-100G	Bulk Coral 800A 10µm C4 AQ 100g
B-800A-10-CO-C4AQ-1KG	Bulk Coral 800A 10µm C4 AQ 1kg
<b>Autres</b>	
B-800A-10-CO-COOH-100G	Bulk Coral 800A 10µm COOH 100g
B-800A-10-CO-COOH-1KG	Bulk Coral 800A 10µm COOH 1kg
B-800A-10-CO-SO3H-100G	Bulk Coral 800A 10µm SO3H 100g
B-800A-10-CO-SO3H-1KG	Bulk Coral 800A 10µm SO3H 1kg
B-800A-10-CO-DNH-100G	Bulk Coral 800A 10µm DNH 100g
B-800A-10-CO-DNH-1KG	Bulk Coral 800A 10µm DNH 1kg
B-800A-10-CO-PEI-100G	Bulk Coral 800A 10µm PEI 100g
B-800A-10-CO-PEI-1KG	Bulk Coral 800A 10µm PEI 1kg
B-800A-10-CO-ARG-100G	Bulk Coral 800A 10µm ARG 100g
B-800A-10-CO-ARG-1KG	Bulk Coral 800A 10µm ARG 1kg
B-800A-10-CO-EP2TNE-100G	Bulk Coral 800A 10µm EP2TNE 100g
B-800A-10-CO-EP2TNE-1KG	Bulk Coral 800A 10µm EP2TNE 1kg
B-800A-10-CO-EP4TNE-100G	Bulk Coral 800A 10µm EP4TNE 100g
B-800A-10-CO-EP4TNE-1KG	Bulk Coral 800A 10µm EP4TNE 1kg
B-800A-10-CO-NH2-100G	Bulk Coral 800A 10µm NH2 100g
B-800A-10-CO-NH2-1KG	Bulk Coral 800A 10µm NH2 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-10-CO-CN-100G	Bulk Coral 800A 10µm CN 100g
B-800A-10-CO-CN-1KG	Bulk Coral 800A 10µm CN 1kg
B-800A-10-CO-Phenyl-100G	Bulk Coral 800A 10µm Phenyl 100g
B-800A-10-CO-Phenyl-1KG	Bulk Coral 800A 10µm Phenyl 1kg
B-800A-10-CO-Phenyl HE-100G	Bulk Coral 800A 10µm Phenyl HE 100g
B-800A-10-CO-Phenyl HE-1KG	Bulk Coral 800A 10µm Phenyl HE 1kg
B-800A-10-CO-Diol-100G	Bulk Coral 800A 10µm Diol 100g
B-800A-10-CO-Diol-1KG	Bulk Coral 800A 10µm Diol 1kg
<b>Particule Coral 10, Pore 1000A</b>	
B-1000A-10-CO-Si-100G	Bulk Coral 1000A 10µm 100g
B-1000A-10-CO-Si-1KG	Bulk Coral 1000A 10µm 1kg
<b>C8</b>	
B-1000A-10-CO-C8-100G	Bulk Coral 1000A 10µm C8 100g
B-1000A-10-CO-C8-1KG	Bulk Coral 1000A 10µm C8 1kg
B-1000A-10-CO-C8W-100G	Bulk Coral 1000A 10µm C8 Wax 100g
B-1000A-10-CO-C8W-1KG	Bulk Coral 1000A 10µm C8 Wax 1kg
B-1000A-10-CO-C8P-100G	Bulk Coral 1000A 10µm C8 Polar 100g
B-1000A-10-CO-C8P-1KG	Bulk Coral 1000A 10µm C8 Polar 1kg
B-1000A-10-CO-C8T-100G	Bulk Coral 1000A 10µm C8 triF 100g
B-1000A-10-CO-C8T-1KG	Bulk Coral 1000A 10µm C8 triF 1kg
B-1000A-10-CO-C8EX-100G	Bulk Coral 1000A 10µm C8 extreme 100g
B-1000A-10-CO-C8EX-1KG	Bulk Coral 1000A 10µm C8 extreme 1kg
B-1000A-10-CO-C8PT-100G	Bulk Coral 1000A 10µm C8 Polar TriF 100g
B-1000A-10-CO-C8PT-1KG	Bulk Coral 1000A 10µm C8 Polar TriF 1kg
B-1000A-10-CO-C8AQP-100G	Bulk Coral 1000A 10µm C8 AQ Plus 100g
B-1000A-10-CO-C8AQP-1KG	Bulk Coral 1000A 10µm C8 AQ Plus 1kg
B-1000A-10-CO-C8AQ-100G	Bulk Coral 1000A 10µm C8 AQ 100g
B-1000A-10-CO-C8AQ-1KG	Bulk Coral 1000A 10µm C8 AQ 1kg
<b>C18</b>	
B-1000A-10-CO-C18-100G	Bulk Coral 1000A 10µm C18 100g
B-1000A-10-CO-C18-1KG	Bulk Coral 1000A 10µm C18 1kg
B-1000A-10-CO-C18W-100G	Bulk Coral 1000A 10µm C18 Wax 100g
B-1000A-10-CO-C18W-1KG	Bulk Coral 1000A 10µm C18 Wax 1kg
B-1000A-10-CO-C18P-100G	Bulk Coral 1000A 10µm C18 Polar 100g
B-1000A-10-CO-C18P-1KG	Bulk Coral 1000A 10µm C18 Polar 1kg
B-1000A-10-CO-C18T-100G	Bulk Coral 1000A 10µm C18 triF 100g
B-1000A-10-CO-C18T-1KG	Bulk Coral 1000A 10µm C18 triF 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-1000A-10-CO-C18EX-100G	Bulk Coral 1000A 10µm C18 extreme 100g
B-1000A-10-CO-C18EX-1KG	Bulk Coral 1000A 10µm C18 extreme 1kg
B-1000A-10-CO-C18PT-100G	Bulk Coral 1000A 10µm C18 Polar TriF 100g
B-1000A-10-CO-C18PT-1KG	Bulk Coral 1000A 10µm C18 Polar TriF 1kg
B-1000A-10-CO-C18AQP-100G	Bulk Coral 1000A 10µm C18 AQ Plus 100g
B-1000A-10-CO-C18AQP-1KG	Bulk Coral 1000A 10µm C18 AQ Plus 1kg
B-1000A-10-CO-C18AQ-100G	Bulk Coral 1000A 10µm C18 AQ 100g
B-1000A-10-CO-C18AQ-1KG	Bulk Coral 1000A 10µm C18 AQ 1kg
<b>C4</b>	
B-1000A-10-CO-C4-100G	Bulk Coral 1000A 10µm C4 100g
B-1000A-10-CO-C4-1KG	Bulk Coral 1000A 10µm C4 1kg
B-1000A-10-CO-C4W-100G	Bulk Coral 1000A 10µm C4 Wax 100g
B-1000A-10-CO-C4W-1KG	Bulk Coral 1000A 10µm C4 Wax 1kg
B-1000A-10-CO-C4P-100G	Bulk Coral 1000A 10µm C4 Polar 100g
B-1000A-10-CO-C4P-1KG	Bulk Coral 1000A 10µm C4 Polar 1kg
B-1000A-10-CO-C4T-100G	Bulk Coral 1000A 10µm C4 TriF100g
B-1000A-10-CO-C4T-1KG	Bulk Coral 1000A 10µm C4 TriF1kg
B-1000A-10-CO-C4EX-100G	Bulk Coral 1000A 10µm C4 Extreme 100g
B-1000A-10-CO-C4EX-1KG	Bulk Coral 1000A 10µm C4 Extreme 1kg
B-1000A-10-CO-C4PT-100G	Bulk Coral 1000A 10µm C4 Polar TriF 100g
B-1000A-10-CO-C4PT-1KG	Bulk Coral 1000A 10µm C4 Polar TriF 1kg
B-1000A-10-CO-C4AQP-100G	Bulk Coral 1000A 10µm C4 AQ Plus 100g
B-1000A-10-CO-C4AQP-1KG	Bulk Coral 1000A 10µm C4 AQ Plus 1kg
B-1000A-10-CO-C4AQ-100G	Bulk Coral 1000A 10µm C4 AQ 100g
B-1000A-10-CO-C4AQ-1KG	Bulk Coral 1000A 10µm C4 AQ 1kg
<b>Autres</b>	
B-1000A-10-CO-COOH-100G	Bulk Coral 1000A 10µm COOH 100g
B-1000A-10-CO-COOH-1KG	Bulk Coral 1000A 10µm COOH 1kg
B-1000A-10-CO-SO3H-100G	Bulk Coral 1000A 10µm SO3H 100g
B-1000A-10-CO-SO3H-1KG	Bulk Coral 1000A 10µm SO3H 1kg
B-1000A-10-CO-DNH-100G	Bulk Coral 1000A 10µm DNH 100g
B-1000A-10-CO-DNH-1KG	Bulk Coral 1000A 10µm DNH 1kg
B-1000A-10-CO-PEI-100G	Bulk Coral 1000A 10µm PEI 100g
B-1000A-10-CO-PEI-1KG	Bulk Coral 1000A 10µm PEI 1kg
B-1000A-10-CO-ARG-100G	Bulk Coral 1000A 10µm ARG 100g
B-1000A-10-CO-ARG-1KG	Bulk Coral 1000A 10µm ARG 1kg
B-1000A-10-CO-EP2TNE-100G	Bulk Coral 1000A 10µm EP2TNE 100g
B-1000A-10-CO-EP2TNE-1KG	Bulk Coral 1000A 10µm EP2TNE 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-1000A-10-CO-EP4TNE-100G	Bulk Coral 1000A 10µm EP4TNE 100g
B-1000A-10-CO-EP4TNE-1KG	Bulk Coral 1000A 10µm EP4TNE 1kg
B-1000A-10-CO-NH2-100G	Bulk Coral 1000A 10µm NH2 100g
B-1000A-10-CO-NH2-1KG	Bulk Coral 1000A 10µm NH2 1kg
B-1000A-10-CO-CN-100G	Bulk Coral 1000A 10µm CN 100g
B-1000A-10-CO-CN-1KG	Bulk Coral 1000A 10µm CN 1kg
B-1000A-10-CO-Phenyl-100G	Bulk Coral 1000A 10µm Phenyl 100g
B-1000A-10-CO-Phenyl-1KG	Bulk Coral 1000A 10µm Phenyl 1kg
B-1000A-10-CO-Phenyl HE-100G	Bulk Coral 1000A 10µm Phenyl HE 100g
B-1000A-10-CO-Phenyl HE-1KG	Bulk Coral 1000A 10µm Phenyl HE 1kg
B-1000A-10-CO-Diol-100G	Bulk Coral 1000A 10µm Diol 100g
B-1000A-10-CO-Diol-1KG	Bulk Coral 1000A 10µm Diol 1kg
<b>Coral 15µm</b>	
<b>Particule Coral 15, Pore 70A</b>	
B-70A-15-CO-Si-100G	Bulk Coral Silice 70A 15µm 100g
B-70A-15-CO-Si-1KG	Bulk Coral Silice 70A 15µm 1kg
<b>C8</b>	
B-70A-15-CO-C8-100G	Bulk Coral 70A 15µm C8 100g
B-70A-15-CO-C8-1KG	Bulk Coral 70A 15µm C8 1kg
B-70A-15-CO-C8W-100G	Bulk Coral 70A 15µm C8 Wax 100g
B-70A-15-CO-C8W-1KG	Bulk Coral 70A 15µm C8 Wax 1kg
B-70A-15-CO-C8P-100G	Bulk Coral 70A 15µm C8 Polar 100g
B-70A-15-CO-C8P-1KG	Bulk Coral 70A 15µm C8 Polar 1kg
B-70A-15-CO-C8T-100G	Bulk Coral 70A 15µm C8 TriF 100g
B-70A-15-CO-C8T-1KG	Bulk Coral 70A 15µm C8 TriF 1kg
B-70A-15-CO-C8EX-100G	Bulk Coral 70A 15µm C8 Extreme 100g
B-70A-15-CO-C8EX-1KG	Bulk Coral 70A 15µm C8 Extreme 1kg
B-70A-15-CO-C8PT-100G	Bulk Coral 70A 15µm C8 Polar TriF 100g
B-70A-15-CO-C8PT-1KG	Bulk Coral 70A 15µm C8 Polar TriF 1kg
B-70A-15-CO-C8AQP-100G	Bulk Coral 70A 15µm C8 AQ Plus 100g
B-70A-15-CO-C8AQP-1KG	Bulk Coral 70A 15µm C8 AQ Plus 1kg
B-70A-15-CO-C8AQ-100G	Bulk Coral 70A 15µm C8 AQ 100g
B-70A-15-CO-C8AQ-1KG	Bulk Coral 70A 15µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-70A-15-CO-C18-100G	Bulk Coral 70A 15µm C18 100g
B-70A-15-CO-C18-1KG	Bulk Coral 70A 15µm C18 1kg
B-70A-15-CO-C18W-100G	Bulk Coral 70A 15µm C18 Wax 100g
B-70A-15-CO-C18W-1KG	Bulk Coral 70A 15µm C18 Wax 1kg
B-70A-15-CO-C18P-100G	Bulk Coral 70A 15µm C18 Polar 100g
B-70A-15-CO-C18P-1KG	Bulk Coral 70A 15µm C18 Polar 1kg
B-70A-15-CO-C18T-100G	Bulk Coral 70A 15µm C18 TriF 100g
B-70A-15-CO-C18T-1KG	Bulk Coral 70A 15µm C18 TriF 1kg
B-70A-15-CO-C18EX-100G	Bulk Coral 70A 15µm C18 Extreme 100g
B-70A-15-CO-C18EX-1KG	Bulk Coral 70A 15µm C18 Extreme 1kg
B-70A-15-CO-C18PT-100G	Bulk Coral 70A 15µm C18 Polar TriF 100g
B-70A-15-CO-C18PT-1KG	Bulk Coral 70A 15µm C18 Polar TriF 1kg
B-70A-15-CO-C18AQP-100G	Bulk Coral 70A 15µm C18 AQ Plus 100g
B-70A-15-CO-C18AQP-1KG	Bulk Coral 70A 15µm C18 AQ Plus 1kg
B-70A-15-CO-C18AQ-100G	Bulk Coral 70A 15µm C18 AQ 100g
B-70A-15-CO-C18AQ-1KG	Bulk Coral 70A 15µm C18 AQ 1kg
<b>C4</b>	
B-70A-15-CO-C4-100G	Bulk Coral 70A 15µm C4 100g
B-70A-15-CO-C4-1KG	Bulk Coral 70A 15µm C4 1kg
B-70A-15-CO-C4W-100G	Bulk Coral 70A 15µm C4 Wax 100g
B-70A-15-CO-C4W-1KG	Bulk Coral 70A 15µm C4 Wax 1kg
B-70A-15-CO-C4P-100G	Bulk Coral 70A 15µm C4 Polar 100g
B-70A-15-CO-C4P-1KG	Bulk Coral 70A 15µm C4 Polar 1kg
B-70A-15-CO-C4T-100G	Bulk Coral 70A 15µm C4 TriF 100g
B-70A-15-CO-C4T-1KG	Bulk Coral 70A 15µm C4 TriF 1kg
B-70A-15-CO-C4EX-100G	Bulk Coral 70A 15µm C4 Extreme 100g
B-70A-15-CO-C4EX-1KG	Bulk Coral 70A 15µm C4 Extreme 1kg
B-70A-15-CO-C4PT-100G	Bulk Coral 70A 15µm C4 Polar TriF 100g
B-70A-15-CO-C4PT-1KG	Bulk Coral 70A 15µm C4 Polar TriF 1kg
B-70A-15-CO-C4AQP-100G	Bulk Coral 70A 15µm C4 AQ Plus 100g
B-70A-15-CO-C4AQP-1KG	Bulk Coral 70A 15µm C4 AQ Plus 1kg
B-70A-15-CO-C4AQ-100G	Bulk Coral 70A 15µm C4 AQ 100g
B-70A-15-CO-C4AQ-1KG	Bulk Coral 70A 15µm C4 AQ 1kg
<b>Autres</b>	
B-70A-15-CO-COOH-100G	Bulk Coral 70A 15µm COOH 100g
B-70A-15-CO-COOH-1KG	Bulk Coral 70A 15µm COOH 1kg
B-70A-15-CO-SO3H-100G	Bulk Coral 70A 15µm SO3H 100g
B-70A-15-CO-SO3H-1KG	Bulk Coral 70A 15µm SO3H 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-15-CO-DNH-100G	Bulk Coral 70A 15µm DNH 100g
B-70A-15-CO-DNH-1KG	Bulk Coral 70A 15µm DNH 1kg
B-70A-15-CO-PEI-100G	Bulk Coral 70A 15µm PEI 100g
B-70A-15-CO-PEI-1KG	Bulk Coral 70A 15µm PEI 1kg
B-70A-15-CO-ARG-100G	Bulk Coral 70A 15µm ARG 100g
B-70A-15-CO-ARG-1KG	Bulk Coral 70A 15µm ARG 1kg
B-70A-15-CO-NH-100G	Bulk Coral 70A 15µm NH 100g
B-70A-15-CO-NH-1KG	Bulk Coral 70A 15µm NH 1kg
B-70A-15-CO-NH2-100G	Bulk Coral 70A 15µm NH2 100g
B-70A-15-CO-NH2-1KG	Bulk Coral 70A 15µm NH2 1kg
B-70A-15-CO-CN-100G	Bulk Coral 70A 15µm CN 100g
B-70A-15-CO-CN-1KG	Bulk Coral 70A 15µm CN 1kg
B-70A-15-CO-Phenyl-100G	Bulk Coral 70A 15µm Phenyl 100g
B-70A-15-CO-Phenyl-1KG	Bulk Coral 70A 15µm Phenyl 1kg
B-70A-15-CO-Phenyl HE-100G	Bulk Coral 70A 15µm Phenyl HE 100g
B-70A-15-CO-Phenyl HE-1KG	Bulk Coral 70A 15µm Phenyl HE 1kg
B-70A-15-CO-Diol-100G	Bulk Coral 70A 15µm Diol 100g
B-70A-15-CO-Diol-1KG	Bulk Coral 70A 15µm Diol 1kg
<b>Particule Coral 15, Pore 100A</b>	
B-100A-15-CO-Si-100G	Bulk Coral 100A 15µm 100g
B-100A-15-CO-Si-1KG	Bulk Coral 100A 15µm 1kg
<b>C8</b>	
B-100A-15-CO-C8-100G	Bulk Coral 100A 15µm C8 100g
B-100A-15-CO-C8-1KG	Bulk Coral 100A 15µm C8 1kg
B-100A-15-CO-C8W-100G	Bulk Coral 100A 15µm C8 Wax 100g
B-100A-15-CO-C8W-1KG	Bulk Coral 100A 15µm C8 Wax 1kg
B-100A-15-CO-C8P-100G	Bulk Coral 100A 15µm C8 Polar 100g
B-100A-15-CO-C8P-1KG	Bulk Coral 100A 15µm C8 Polar 1kg
B-100A-15-CO-C8T-100G	Bulk Coral 100A 15µm C8 TriF 100g
B-100A-15-CO-C8T-1KG	Bulk Coral 100A 15µm C8 TriF 1kg
B-100A-15-CO-C8EX-100G	Bulk Coral 100A 15µm C8 Extreme 100g
B-100A-15-CO-C8EX-1KG	Bulk Coral 100A 15µm C8 Extreme 1kg
B-100A-15-CO-C8PT-100G	Bulk Coral 100A 15µm C8 Polar TriF 100g
B-100A-15-CO-C8PT-1KG	Bulk Coral 100A 15µm C8 Polar TriF 1kg
B-100A-15-CO-C8AQP-100G	Bulk Coral 100A 15µm C8 AQ Plus 100g
B-100A-15-CO-C8AQP-1KG	Bulk Coral 100A 15µm C8 AQ Plus 1kg
B-100A-15-CO-C8AQ-100G	Bulk Coral 100A 15µm C8 AQ 100g
B-100A-15-CO-C8AQ-1KG	Bulk Coral 100A 15µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-100A-15-CO-C18-100G	Bulk Coral 100A 15µm C18 100g
B-100A-15-CO-C18-1KG	Bulk Coral 100A 15µm C18 1kg
B-100A-15-CO-C18W-100G	Bulk Coral 100A 15µm C18 Wax 100g
B-100A-15-CO-C18W-1KG	Bulk Coral 100A 15µm C18 Wax 1kg
B-100A-15-CO-C18P-100G	Bulk Coral 100A 15µm C18 Polar 100g
B-100A-15-CO-C18P-1KG	Bulk Coral 100A 15µm C18 Polar 1kg
B-100A-15-CO-C18T-100G	Bulk Coral 100A 15µm C18 TriF 100g
B-100A-15-CO-C18T-1KG	Bulk Coral 100A 15µm C18 TriF 1kg
B-100A-15-CO-C18EX-100G	Bulk Coral 100A 15µm C18 Extreme 100g
B-100A-15-CO-C18EX-1KG	Bulk Coral 100A 15µm C18 Extreme 1kg
B-100A-15-CO-C18PT-100G	Bulk Coral 100A 15µm C18 Polar TriF 100g
B-100A-15-CO-C18PT-1KG	Bulk Coral 100A 15µm C18 Polar TriF 1kg
B-100A-15-CO-C18AQP-100G	Bulk Coral 100A 15µm C18 AQ Plus 100g
B-100A-15-CO-C18AQP-1KG	Bulk Coral 100A 15µm C18 AQ Plus 1kg
B-100A-15-CO-C18AQ-100G	Bulk Coral 100A 15µm C18 AQ 100g
B-100A-15-CO-C18AQ-1KG	Bulk Coral 100A 15µm C18 AQ 1kg
<b>C4</b>	
B-100A-15-CO-C4-100G	Bulk Coral 100A 15µm C4 100g
B-100A-15-CO-C4-1KG	Bulk Coral 100A 15µm C4 1kg
B-100A-15-CO-C4W-100G	Bulk Coral 100A 15µm C4 Wax 100g
B-100A-15-CO-C4W-1KG	Bulk Coral 100A 15µm C4 Wax 1kg
B-100A-15-CO-C4P-100G	Bulk Coral 100A 15µm C4 Polar 100g
B-100A-15-CO-C4P-1KG	Bulk Coral 100A 15µm C4 Polar 1kg
B-100A-15-CO-C4T-100G	Bulk Coral 100A 15µm C4 TriF 100g
B-100A-15-CO-C4T-1KG	Bulk Coral 100A 15µm C4 TriF 1kg
B-100A-15-CO-C4EX-100G	Bulk Coral 100A 15µm C4 Extreme 100g
B-100A-15-CO-C4EX-1KG	Bulk Coral 100A 15µm C4 Extreme 1kg
B-100A-15-CO-C4PT-100G	Bulk Coral 100A 15µm C4 Polar TriF 100g
B-100A-15-CO-C4PT-1KG	Bulk Coral 100A 15µm C4 Polar TriF 1kg
B-100A-15-CO-C4AQP-100G	Bulk Coral 100A 15µm C4 AQ Plus 100g
B-100A-15-CO-C4AQP-1KG	Bulk Coral 100A 15µm C4 AQ Plus 1kg
B-100A-15-CO-C4AQ-100G	Bulk Coral 100A 15µm C4 AQ 100g
B-100A-15-CO-C4AQ-1KG	Bulk Coral 100A 15µm C4 AQ 1kg
<b>Autres</b>	
B-100A-15-CO-COOH-100G	Bulk Coral 100A 15µm COOH 100g
B-100A-15-CO-COOH-1KG	Bulk Coral 100A 15µm COOH 1kg
B-100A-15-CO-SO3H-100G	Bulk Coral 100A 15µm SO3H 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-15-CO-SO3H-1KG	Bulk Coral 100A 15µm SO3H 1kg
B-100A-15-CO-DNH-100G	Bulk Coral 100A 15µm DNH 100g
B-100A-15-CO-DNH-1KG	Bulk Coral 100A 15µm DNH 1kg
B-100A-15-CO-PEI-100G	Bulk Coral 100A 15µm PEI 100g
B-100A-15-CO-PEI-1KG	Bulk Coral 100A 15µm PEI 1kg
B-100A-15-CO-ARG-100G	Bulk Coral 100A 15µm ARG 100g
B-100A-15-CO-ARG-1KG	Bulk Coral 100A 15µm ARG 1kg
B-100A-15-CO-NH-100G	Bulk Coral 100A 15µm NH 100g
B-100A-15-CO-NH-1KG	Bulk Coral 100A 15µm NH 1kg
B-100A-15-CO-NH2-100G	Bulk Coral 100A 15µm NH2 100g
B-100A-15-CO-NH2-1KG	Bulk Coral 100A 15µm NH2 1kg
B-100A-15-CO-CN-100G	Bulk Coral 100A 15µm CN 100g
B-100A-15-CO-CN-1KG	Bulk Coral 100A 15µm CN 1kg
B-100A-15-CO-Phenyl-100G	Bulk Coral 100A 15µm Phenyl 100g
B-100A-15-CO-Phenyl-1KG	Bulk Coral 100A 15µm Phenyl 1kg
B-100A-15-CO-Phenyl HE-100G	Bulk Coral 100A 15µm Phenyl HE 100g
B-100A-15-CO-Phenyl HE-1KG	Bulk Coral 100A 15µm Phenyl HE 1kg
B-100A-15-CO-Diol-100G	Bulk Coral 100A 15µm Diol 100g
B-100A-15-CO-Diol-1KG	Bulk Coral 100A 15µm Diol 1kg
<b>Particule Coral 15, Pore 150A</b>	
B-150A-15-CO-Si-100G	Bulk Coral 150A 15µm 100g
B-150A-15-CO-Si-1KG	Bulk Coral 150A 15µm 1kg
<b>C8</b>	
B-150A-15-CO-C8-100G	Bulk Coral 150A 15µm C8 100g
B-150A-15-CO-C8-1KG	Bulk Coral 150A 15µm C8 1kg
B-150A-15-CO-C8W-100G	Bulk Coral 150A 15µm C8 Wax 100g
B-150A-15-CO-C8W-1KG	Bulk Coral 150A 15µm C8 Wax 1kg
B-150A-15-CO-C8P-100G	Bulk Coral 150A 15µm C8 Polar 100g
B-150A-15-CO-C8P-1KG	Bulk Coral 150A 15µm C8 Polar 1kg
B-150A-15-CO-C8T-100G	Bulk Coral 150A 15µm C8 TriF 100g
B-150A-15-CO-C8T-1KG	Bulk Coral 150A 15µm C8 TriF 1kg
B-150A-15-CO-C8EX-100G	Bulk Coral 150A 15µm C8 Extreme 100g
B-150A-15-CO-C8EX-1KG	Bulk Coral 150A 15µm C8 Extreme 1kg
B-150A-15-CO-C8PT-100G	Bulk Coral 150A 15µm C8 Polar TriF 100g
B-150A-15-CO-C8PT-1KG	Bulk Coral 150A 15µm C8 Polar TriF 1kg
B-150A-15-CO-C8AQP-100G	Bulk Coral 150A 15µm C8 AQ Plus 100g
B-150A-15-CO-C8AQP-1KG	Bulk Coral 150A 15µm C8 AQ Plus 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-15-CO-C8AQ-100G	Bulk Coral 150A 15µm C8 AQ 100g
B-150A-15-CO-C8AQ-1KG	Bulk Coral 150A 15µm C8 AQ 1kg
<b>C18</b>	
B-150A-15-CO-C18-100G	Bulk Coral 150A 15µm C18 100g
B-150A-15-CO-C18-1KG	Bulk Coral 150A 15µm C18 1kg
B-150A-15-CO-C18W-100G	Bulk Coral 150A 15µm C18 Wax 100g
B-150A-15-CO-C18W-1KG	Bulk Coral 150A 15µm C18 Wax 1kg
B-150A-15-CO-C18P-100G	Bulk Coral 150A 15µm C18 Polar 100g
B-150A-15-CO-C18P-1KG	Bulk Coral 150A 15µm C18 Polar 1kg
B-150A-15-CO-C18T-100G	Bulk Coral 150A 15µm C18 TriF 100g
B-150A-15-CO-C18T-1KG	Bulk Coral 150A 15µm C18 TriF 1kg
B-150A-15-CO-C18EX-100G	Bulk Coral 150A 15µm C18 Extreme 100g
B-150A-15-CO-C18EX-1KG	Bulk Coral 150A 15µm C18 Extreme 1kg
B-150A-15-CO-C18PT-100G	Bulk Coral 150A 15µm C18 Polar TriF 100g
B-150A-15-CO-C18PT-1KG	Bulk Coral 150A 15µm C18 Polar TriF 1kg
B-150A-15-CO-C18AQP-100G	Bulk Coral 150A 15µm C18 AQ Plus 100g
B-150A-15-CO-C18AQP-1KG	Bulk Coral 150A 15µm C18 AQ Plus 1kg
B-150A-15-CO-C18AQ-100G	Bulk Coral 150A 15µm C18 AQ 100g
B-150A-15-CO-C18AQ-1KG	Bulk Coral 150A 15µm C18 AQ 1kg
<b>C4</b>	
B-150A-15-CO-C4-100G	Bulk Coral 150A 15µm C4 100g
B-150A-15-CO-C4-1KG	Bulk Coral 150A 15µm C4 1kg
B-150A-15-CO-C4W-100G	Bulk Coral 150A 15µm C4 Wax 100g
B-150A-15-CO-C4W-1KG	Bulk Coral 150A 15µm C4 Wax 1kg
B-150A-15-CO-C4P-100G	Bulk Coral 150A 15µm C4 Polar 100g
B-150A-15-CO-C4P-1KG	Bulk Coral 150A 15µm C4 Polar 1kg
B-150A-15-CO-C4T-100G	Bulk Coral 150A 15µm C4 TriF 100g
B-150A-15-CO-C4T-1KG	Bulk Coral 150A 15µm C4 TriF 1kg
B-150A-15-CO-C4EX-100G	Bulk Coral 150A 15µm C4 Extreme 100g
B-150A-15-CO-C4EX-1KG	Bulk Coral 150A 15µm C4 Extreme 1kg
B-150A-15-CO-C4PT-100G	Bulk Coral 150A 15µm C4 Polar TriF 100g
B-150A-15-CO-C4PT-1KG	Bulk Coral 150A 15µm C4 Polar TriF 1kg
B-150A-15-CO-C4AQP-100G	Bulk Coral 150A 15µm C4 AQ Plus 100g
B-150A-15-CO-C4AQP-1KG	Bulk Coral 150A 15µm C4 AQ Plus 1kg
B-150A-15-CO-C4AQ-100G	Bulk Coral 150A 15µm C4 AQ 100g
B-150A-15-CO-C4AQ-1KG	Bulk Coral 150A 15µm C4 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Autres</b>	
B-150A-15-CO-COOH-100G	Bulk Coral 150A 15µm COOH 100g
B-150A-15-CO-COOH-1KG	Bulk Coral 150A 15µm COOH 1kg
B-150A-15-CO-SO3H-100G	Bulk Coral 150A 15µm SO3H 100g
B-150A-15-CO-SO3H-1KG	Bulk Coral 150A 15µm SO3H 1kg
B-150A-15-CO-DNH-100G	Bulk Coral 150A 15µm DNH 100g
B-150A-15-CO-DNH-1KG	Bulk Coral 150A 15µm DNH 1kg
B-150A-15-CO-PEI-100G	Bulk Coral 150A 15µm PEI 100g
B-150A-15-CO-PEI-1KG	Bulk Coral 150A 15µm PEI 1kg
B-150A-15-CO-ARG-100G	Bulk Coral 150A 15µm ARG 100g
B-150A-15-CO-ARG-1KG	Bulk Coral 150A 15µm ARG 1kg
B-150A-15-CO-NH-100G	Bulk Coral 150A 15µm NH 100g
B-150A-15-CO-NH-1KG	Bulk Coral 150A 15µm NH 1kg
B-150A-15-CO-NH2-100G	Bulk Coral 150A 15µm NH2 100g
B-150A-15-CO-NH2-1KG	Bulk Coral 150A 15µm NH2 1kg
B-150A-15-CO-CN-100G	Bulk Coral 150A 15µm CN 100g
B-150A-15-CO-CN-1KG	Bulk Coral 150A 15µm CN 1kg
B-150A-15-CO-Phenyl-100G	Bulk Coral 150A 15µm Phenyl 100g
B-150A-15-CO-Phenyl-1KG	Bulk Coral 150A 15µm Phenyl 1kg
B-150A-15-CO-Phenyl HE-100G	Bulk Coral 150A 15µm Phenyl HE 100g
B-150A-15-CO-Phenyl HE-1KG	Bulk Coral 150A 15µm Phenyl HE 1kg
B-150A-15-CO-Diol-100G	Bulk Coral 150A 15µm Diol 100g
B-150A-15-CO-Diol-1KG	Bulk Coral 150A 15µm Diol 1kg
<b>Particule Coral 15, Pore 200A</b>	
B-200A-15-CO-Si-100G	Bulk Coral 200A 15µm 100g
B-200A-15-CO-Si-1KG	Bulk Coral 200A 15µm 1kg
<b>C8</b>	
B-200A-15-CO-C8-100G	Bulk Coral 200A 15µm C8 100g
B-200A-15-CO-C8-1KG	Bulk Coral 200A 15µm C8 1kg
B-200A-15-CO-C8W-100G	Bulk Coral 200A 15µm C8 Wax 100g
B-200A-15-CO-C8W-1KG	Bulk Coral 200A 15µm C8 Wax 1kg
B-200A-15-CO-C8P-100G	Bulk Coral 200A 15µm C8 Polar 100g
B-200A-15-CO-C8P-1KG	Bulk Coral 200A 15µm C8 Polar 1kg
B-200A-15-CO-C8T-100G	Bulk Coral 200A 15µm C8 TriF 100g
B-200A-15-CO-C8T-1KG	Bulk Coral 200A 15µm C8 TriF 1kg
B-200A-15-CO-C8EX-100G	Bulk Coral 200A 15µm C8 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-15-CO-C8EX-1KG	Bulk Coral 200A 15µm C8 Extreme 1kg
B-200A-15-CO-C8PT-100G	Bulk Coral 200A 15µm C8 Polar TriF 100g
B-200A-15-CO-C8PT-1KG	Bulk Coral 200A 15µm C8 Polar TriF 1kg
B-200A-15-CO-C8AQP-100G	Bulk Coral 200A 15µm C8 AQ Plus 100g
B-200A-15-CO-C8AQP-1KG	Bulk Coral 200A 15µm C8 AQ Plus 1kg
B-200A-15-CO-C8AQ-100G	Bulk Coral 200A 15µm C8 AQ 100g
B-200A-15-CO-C8AQ-1KG	Bulk Coral 200A 15µm C8 AQ 1kg
<b>C18</b>	
B-200A-15-CO-C18-100G	Bulk Coral 200A 15µm C18 100g
B-200A-15-CO-C18-1KG	Bulk Coral 200A 15µm C18 1kg
B-200A-15-CO-C18W-100G	Bulk Coral 200A 15µm C18 Wax 100g
B-200A-15-CO-C18W-1KG	Bulk Coral 200A 15µm C18 Wax 1kg
B-200A-15-CO-C18P-100G	Bulk Coral 200A 15µm C18 Polar 100g
B-200A-15-CO-C18P-1KG	Bulk Coral 200A 15µm C18 Polar 1kg
B-200A-15-CO-C18T-100G	Bulk Coral 200A 15µm C18 TriF 100g
B-200A-15-CO-C18T-1KG	Bulk Coral 200A 15µm C18 TriF 1kg
B-200A-15-CO-C18EX-100G	Bulk Coral 200A 15µm C18 Extreme 100g
B-200A-15-CO-C18EX-1KG	Bulk Coral 200A 15µm C18 Extreme 1kg
B-200A-15-CO-C18PT-100G	Bulk Coral 200A 15µm C18 Polar TriF 100g
B-200A-15-CO-C18PT-1KG	Bulk Coral 200A 15µm C18 Polar TriF 1kg
B-200A-15-CO-C18AQP-100G	Bulk Coral 200A 15µm C18 AQ Plus 100g
B-200A-15-CO-C18AQP-1KG	Bulk Coral 200A 15µm C18 AQ Plus 1kg
B-200A-15-CO-C18AQ-100G	Bulk Coral 200A 15µm C18 AQ 100g
B-200A-15-CO-C18AQ-1KG	Bulk Coral 200A 15µm C18 AQ 1kg
<b>C4</b>	
B-200A-15-CO-C4-100G	Bulk Coral 200A 15µm C4 100g
B-200A-15-CO-C4-1KG	Bulk Coral 200A 15µm C4 1kg
B-200A-15-CO-C4W-100G	Bulk Coral 200A 15µm C4 Wax 100g
B-200A-15-CO-C4W-1KG	Bulk Coral 200A 15µm C4 Wax 1kg
B-200A-15-CO-C4P-100G	Bulk Coral 200A 15µm C4 Polar 100g
B-200A-15-CO-C4P-1KG	Bulk Coral 200A 15µm C4 Polar 1kg
B-200A-15-CO-C4T-100G	Bulk Coral 200A 15µm C4 TriF 100g
B-200A-15-CO-C4T-1KG	Bulk Coral 200A 15µm C4 TriF 1kg
B-200A-15-CO-C4EX-100G	Bulk Coral 200A 15µm C4 Extreme 100g
B-200A-15-CO-C4EX-1KG	Bulk Coral 200A 15µm C4 Extreme 1kg
B-200A-15-CO-C4PT-100G	Bulk Coral 200A 15µm C4 Polar TriF 100g
B-200A-15-CO-C4PT-1KG	Bulk Coral 200A 15µm C4 Polar TriF 1kg
B-200A-15-CO-C4AQP-100G	Bulk Coral 200A 15µm C4 AQ Plus 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-15-CO-C4AQP-1KG	Bulk Coral 200A 15µm C4 AQ Plus 1kg
B-200A-15-CO-C4AQ-100G	Bulk Coral 200A 15µm C4 AQ 100g
B-200A-15-CO-C4AQ-1KG	Bulk Coral 200A 15µm C4 AQ 1kg
<b>Autres</b>	
B-200A-15-CO-COOH-100G	Bulk Coral 200A 15µm COOH 100g
B-200A-15-CO-COOH-1KG	Bulk Coral 200A 15µm COOH 1kg
B-200A-15-CO-SO3H-100G	Bulk Coral 200A 15µm SO3H 100g
B-200A-15-CO-SO3H-1KG	Bulk Coral 200A 15µm SO3H 1kg
B-200A-15-CO-DNH-100G	Bulk Coral 200A 15µm DNH 100g
B-200A-15-CO-DNH-1KG	Bulk Coral 200A 15µm DNH 1kg
B-200A-15-CO-PEI-100G	Bulk Coral 200A 15µm PEI 100g
B-200A-15-CO-PEI-1KG	Bulk Coral 200A 15µm PEI 1kg
B-200A-15-CO-ARG-100G	Bulk Coral 200A 15µm ARG 100g
B-200A-15-CO-ARG-1KG	Bulk Coral 200A 15µm ARG 1kg
B-200A-15-CO-NH-100G	Bulk Coral 200A 15µm NH 100g
B-200A-15-CO-NH-1KG	Bulk Coral 200A 15µm NH 1kg
B-200A-15-CO-NH2-100G	Bulk Coral 200A 15µm NH2 100g
B-200A-15-CO-NH2-1KG	Bulk Coral 200A 15µm NH2 1kg
B-200A-15-CO-CN-100G	Bulk Coral 200A 15µm CN 100g
B-200A-15-CO-CN-1KG	Bulk Coral 200A 15µm CN 1kg
B-200A-15-CO-Phenyl-100G	Bulk Coral 200A 15µm Phenyl 100g
B-200A-15-CO-Phenyl-1KG	Bulk Coral 200A 15µm Phenyl 1kg
B-200A-15-CO-Phenyl HE-100G	Bulk Coral 200A 15µm Phenyl HE 100g
B-200A-15-CO-Phenyl HE-1KG	Bulk Coral 200A 15µm Phenyl HE 1kg
B-200A-15-CO-Diol-100G	Bulk Coral 200A 15µm Diol 100g
B-200A-15-CO-Diol-1KG	Bulk Coral 200A 15µm Diol 1kg
<b>Particule Coral 15, Pore 300A</b>	
B-300A-15-CO-Si-100G	Bulk Coral 300A 15µm 100g
B-300A-15-CO-Si-1KG	Bulk Coral 300A 15µm 1kg
<b>C8</b>	
B-300A-15-CO-C8-100G	Bulk Coral 300A 15µm C8 100g
B-300A-15-CO-C8-1KG	Bulk Coral 300A 15µm C8 1kg
B-300A-15-CO-C8W-100G	Bulk Coral 300A 15µm C8 Wax 100g
B-300A-15-CO-C8W-1KG	Bulk Coral 300A 15µm C8 Wax 1kg
B-300A-15-CO-C8P-100G	Bulk Coral 300A 15µm C8 Polar 100g
B-300A-15-CO-C8P-1KG	Bulk Coral 300A 15µm C8 Polar 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-15-CO-C8T-100G	Bulk Coral 300A 15µm C8 TriF 100g
B-300A-15-CO-C8T-1KG	Bulk Coral 300A 15µm C8 TriF 1kg
B-300A-15-CO-C8EX-100G	Bulk Coral 300A 15µm C8 Extreme 100g
B-300A-15-CO-C8EX-1KG	Bulk Coral 300A 15µm C8 Extreme 1kg
B-300A-15-CO-C8PT-100G	Bulk Coral 300A 15µm C8 Polar TriF 100g
B-300A-15-CO-C8PT-1KG	Bulk Coral 300A 15µm C8 Polar TriF 1kg
B-300A-15-CO-C8AQP-100G	Bulk Coral 300A 15µm C8 AQ Plus 100g
B-300A-15-CO-C8AQP-1KG	Bulk Coral 300A 15µm C8 AQ Plus 1kg
B-300A-15-CO-C8AQ-100G	Bulk Coral 300A 15µm C8 AQ 100g
B-300A-15-CO-C8AQ-1KG	Bulk Coral 300A 15µm C8 AQ 1kg
<b>C18</b>	
B-300A-15-CO-C18-100G	Bulk Coral 300A 15µm C18 100g
B-300A-15-CO-C18-1KG	Bulk Coral 300A 15µm C18 1kg
B-300A-15-CO-C18W-100G	Bulk Coral 300A 15µm C18 Wax 100g
B-300A-15-CO-C18W-1KG	Bulk Coral 300A 15µm C18 Wax 1kg
B-300A-15-CO-C18P-100G	Bulk Coral 300A 15µm C18 Polar 100g
B-300A-15-CO-C18P-1KG	Bulk Coral 300A 15µm C18 Polar 1kg
B-300A-15-CO-C18T-100G	Bulk Coral 300A 15µm C18 TriF 100g
B-300A-15-CO-C18T-1KG	Bulk Coral 300A 15µm C18 TriF 1kg
B-300A-15-CO-C18EX-100G	Bulk Coral 300A 15µm C18 Extreme 100g
B-300A-15-CO-C18EX-1KG	Bulk Coral 300A 15µm C18 Extreme 1kg
B-300A-15-CO-C18PT-100G	Bulk Coral 300A 15µm C18 Polar TriF 100g
B-300A-15-CO-C18PT-1KG	Bulk Coral 300A 15µm C18 Polar TriF 1kg
B-300A-15-CO-C18AQP-100G	Bulk Coral 300A 15µm C18 AQ Plus 100g
B-300A-15-CO-C18AQP-1KG	Bulk Coral 300A 15µm C18 AQ Plus 1kg
B-300A-15-CO-C18AQ-100G	Bulk Coral 300A 15µm C18 AQ 100g
B-300A-15-CO-C18AQ-1KG	Bulk Coral 300A 15µm C18 AQ 1kg
<b>C4</b>	
B-300A-15-CO-C4-100G	Bulk Coral 300A 15µm C4 100g
B-300A-15-CO-C4-1KG	Bulk Coral 300A 15µm C4 1kg
B-300A-15-CO-C4W-100G	Bulk Coral 300A 15µm C4 Wax 100g
B-300A-15-CO-C4W-1KG	Bulk Coral 300A 15µm C4 Wax 1kg
B-300A-15-CO-C4P-100G	Bulk Coral 300A 15µm C4 Polar 100g
B-300A-15-CO-C4P-1KG	Bulk Coral 300A 15µm C4 Polar 1kg
B-300A-15-CO-C4T-100G	Bulk Coral 300A 15µm C4 TriF 100g
B-300A-15-CO-C4T-1KG	Bulk Coral 300A 15µm C4 TriF 1kg
B-300A-15-CO-C4EX-100G	Bulk Coral 300A 15µm C4 Extreme 100g
B-300A-15-CO-C4EX-1KG	Bulk Coral 300A 15µm C4 Extreme 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-15-CO-C4PT-100G	Bulk Coral 300A 15µm C4 Polar TriF 100g
B-300A-15-CO-C4PT-1KG	Bulk Coral 300A 15µm C4 Polar TriF 1kg
B-300A-15-CO-C4AQP-100G	Bulk Coral 300A 15µm C4 AQ Plus 100g
B-300A-15-CO-C4AQP-1KG	Bulk Coral 300A 15µm C4 AQ Plus 1kg
B-300A-15-CO-C4AQ-100G	Bulk Coral 300A 15µm C4 AQ 100g
B-300A-15-CO-C4AQ-1KG	Bulk Coral 300A 15µm C4 AQ 1kg
<b>Autres</b>	
B-300A-15-CO-COOH-100G	Bulk Coral 300A 15µm COOH 100g
B-300A-15-CO-COOH-1KG	Bulk Coral 300A 15µm COOH 1kg
B-300A-15-CO-SO3H-100G	Bulk Coral 300A 15µm SO3H 100g
B-300A-15-CO-SO3H-1KG	Bulk Coral 300A 15µm SO3H 1kg
B-300A-15-CO-DNH-100G	Bulk Coral 300A 15µm DNH 100g
B-300A-15-CO-DNH-1KG	Bulk Coral 300A 15µm DNH 1kg
B-300A-15-CO-PEI-100G	Bulk Coral 300A 15µm PEI 100g
B-300A-15-CO-PEI-1KG	Bulk Coral 300A 15µm PEI 1kg
B-300A-15-CO-ARG-100G	Bulk Coral 300A 15µm ARG 100g
B-300A-15-CO-ARG-1KG	Bulk Coral 300A 15µm ARG 1kg
B-300A-15-CO-NH-100G	Bulk Coral 300A 15µm NH 100g
B-300A-15-CO-NH-1KG	Bulk Coral 300A 15µm NH 1kg
B-300A-15-CO-NH2-100G	Bulk Coral 300A 15µm NH2 100g
B-300A-15-CO-NH2-1KG	Bulk Coral 300A 15µm NH2 1kg
B-300A-15-CO-CN-100G	Bulk Coral 300A 15µm CN 100g
B-300A-15-CO-CN-1KG	Bulk Coral 300A 15µm CN 1kg
B-300A-15-CO-Phenyl-100G	Bulk Coral 300A 15µm Phenyl 100g
B-300A-15-CO-Phenyl-1KG	Bulk Coral 300A 15µm Phenyl 1kg
B-300A-15-CO-Phenyl HE-100G	Bulk Coral 300A 15µm Phenyl HE 100g
B-300A-15-CO-Phenyl HE-1KG	Bulk Coral 300A 15µm Phenyl HE 1kg
B-300A-15-CO-Diol-100G	Bulk Coral 300A 15µm Diol 100g
B-300A-15-CO-Diol-1KG	Bulk Coral 300A 15µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Coral 15, Pore 500A</b>	
B-500A-15-CO-Si-100G	Bulk Coral 500A 15µm 100g
B-500A-15-CO-Si-1KG	Bulk Coral 500A 15µm 1kg
<b>C8</b>	
B-500A-15-CO-C8-100G	Bulk Coral 500A 15µm C8 100g
B-500A-15-CO-C8-1KG	Bulk Coral 500A 15µm C8 1kg
B-500A-15-CO-C8W-100G	Bulk Coral 500A 15µm C8 Wax 100g
B-500A-15-CO-C8W-1KG	Bulk Coral 500A 15µm C8 Wax 1kg
B-500A-15-CO-C8P-100G	Bulk Coral 500A 15µm C8 Polar 100g
B-500A-15-CO-C8P-1KG	Bulk Coral 500A 15µm C8 Polar 1kg
B-500A-15-CO-C8T-100G	Bulk Coral 500A 15µm C8 TriF 100g
B-500A-15-CO-C8T-1KG	Bulk Coral 500A 15µm C8 TriF 1kg
B-500A-15-CO-C8EX-100G	Bulk Coral 500A 15µm C8 Extreme 100g
B-500A-15-CO-C8EX-1KG	Bulk Coral 500A 15µm C8 Extreme 1kg
B-500A-15-CO-C8PT-100G	Bulk Coral 500A 15µm C8 Polar TriF 100g
B-500A-15-CO-C8PT-1KG	Bulk Coral 500A 15µm C8 Polar TriF 1kg
B-500A-15-CO-C8AQP-100G	Bulk Coral 500A 15µm C8 AQ Plus 100g
B-500A-15-CO-C8AQP-1KG	Bulk Coral 500A 15µm C8 AQ Plus 1kg
B-500A-15-CO-C8AQ-100G	Bulk Coral 500A 15µm C8 AQ 100g
B-500A-15-CO-C8AQ-1KG	Bulk Coral 500A 15µm C8 AQ 1kg
<b>C18</b>	
B-500A-15-CO-C18-100G	Bulk Coral 500A 15µm C18 100g
B-500A-15-CO-C18-1KG	Bulk Coral 500A 15µm C18 1kg
B-500A-15-CO-C18W-100G	Bulk Coral 500A 15µm C18 Wax 100g
B-500A-15-CO-C18W-1KG	Bulk Coral 500A 15µm C18 Wax 1kg
B-500A-15-CO-C18P-100G	Bulk Coral 500A 15µm C18 Polar 100g
B-500A-15-CO-C18P-1KG	Bulk Coral 500A 15µm C18 Polar 1kg
B-500A-15-CO-C18T-100G	Bulk Coral 500A 15µm C18 TriF 100g
B-500A-15-CO-C18T-1KG	Bulk Coral 500A 15µm C18 TriF 1kg
B-500A-15-CO-C18EX-100G	Bulk Coral 500A 15µm C18 Extreme 100g
B-500A-15-CO-C18EX-1KG	Bulk Coral 500A 15µm C18 Extreme 1kg
B-500A-15-CO-C18PT-100G	Bulk Coral 500A 15µm C18 Polar TriF 100g
B-500A-15-CO-C18PT-1KG	Bulk Coral 500A 15µm C18 Polar TriF 1kg
B-500A-15-CO-C18AQP-100G	Bulk Coral 500A 15µm C18 AQ Plus 100g
B-500A-15-CO-C18AQP-1KG	Bulk Coral 500A 15µm C18 AQ Plus 1kg
B-500A-15-CO-C18AQ-100G	Bulk Coral 500A 15µm C18 AQ 100g
B-500A-15-CO-C18AQ-1KG	Bulk Coral 500A 15µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-500A-15-CO-C4-100G	Bulk Coral 500A 15µm C4 100g
B-500A-15-CO-C4-1KG	Bulk Coral 500A 15µm C4 1kg
B-500A-15-CO-C4W-100G	Bulk Coral 500A 15µm C4 Wax 100g
B-500A-15-CO-C4W-1KG	Bulk Coral 500A 15µm C4 Wax 1kg
B-500A-15-CO-C4P-100G	Bulk Coral 500A 15µm C4 Polar 100g
B-500A-15-CO-C4P-1KG	Bulk Coral 500A 15µm C4 Polar 1kg
B-500A-15-CO-C4T-100G	Bulk Coral 500A 15µm C4 TriF 100g
B-500A-15-CO-C4T-1KG	Bulk Coral 500A 15µm C4 TriF 1kg
B-500A-15-CO-C4EX-100G	Bulk Coral 500A 15µm C4 Extreme 100g
B-500A-15-CO-C4EX-1KG	Bulk Coral 500A 15µm C4 Extreme 1kg
B-500A-15-CO-C4PT-100G	Bulk Coral 500A 15µm C4 Polar TriF 100g
B-500A-15-CO-C4PT-1KG	Bulk Coral 500A 15µm C4 Polar TriF 1kg
B-500A-15-CO-C4AQP-100G	Bulk Coral 500A 15µm C4 AQ Plus 100g
B-500A-15-CO-C4AQP-1KG	Bulk Coral 500A 15µm C4 AQ Plus 1kg
B-500A-15-CO-C4AQ-100G	Bulk Coral 500A 15µm C4 AQ 100g
B-500A-15-CO-C4AQ-1KG	Bulk Coral 500A 15µm C4 AQ 1kg
<b>Autres</b>	
B-500A-15-CO-COOH-100G	Bulk Coral 500A 15µm COOH 100g
B-500A-15-CO-COOH-1KG	Bulk Coral 500A 15µm COOH 1kg
B-500A-15-CO-SO3H-100G	Bulk Coral 500A 15µm SO3H 100g
B-500A-15-CO-SO3H-1KG	Bulk Coral 500A 15µm SO3H 1kg
B-500A-15-CO-DNH-100G	Bulk Coral 500A 15µm DNH 100g
B-500A-15-CO-DNH-1KG	Bulk Coral 500A 15µm DNH 1kg
B-500A-15-CO-PEI-100G	Bulk Coral 500A 15µm PEI 100g
B-500A-15-CO-PEI-1KG	Bulk Coral 500A 15µm PEI 1kg
B-500A-15-CO-ARG-100G	Bulk Coral 500A 15µm ARG 100g
B-500A-15-CO-ARG-1KG	Bulk Coral 500A 15µm ARG 1kg
B-500A-15-CO-NH-100G	Bulk Coral 500A 15µm NH 100g
B-500A-15-CO-NH-1KG	Bulk Coral 500A 15µm NH 1kg
B-500A-15-CO-NH2-100G	Bulk Coral 500A 15µm NH2 100g
B-500A-15-CO-NH2-1KG	Bulk Coral 500A 15µm NH2 1kg
B-500A-15-CO-CN-100G	Bulk Coral 500A 15µm CN 100g
B-500A-15-CO-CN-1KG	Bulk Coral 500A 15µm CN 1kg
B-500A-15-CO-Phenyl-100G	Bulk Coral 500A 15µm Phenyl 100g
B-500A-15-CO-Phenyl-1KG	Bulk Coral 500A 15µm Phenyl 1kg
B-500A-15-CO-Phenyl HE-100G	Bulk Coral 500A 15µm Phenyl HE 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-500A-15-CO-Phenyl HE-1KG	Bulk Coral 500A 15µm Phenyl HE 1kg
B-500A-15-CO-Diol-100G	Bulk Coral 500A 15µm Diol 100g
B-500A-15-CO-Diol-1KG	Bulk Coral 500A 15µm Diol 1kg
<b>Particule Coral 15, Pore 800A</b>	
B-800A-15-CO-Si-100G	Bulk Coral 800A 15µm 100g
B-800A-15-CO-Si-1KG	Bulk Coral 800A 15µm 1kg
<b>C8</b>	
B-800A-15-CO-C8-100G	Bulk Coral 800A 15µm C8 100g
B-800A-15-CO-C8-1KG	Bulk Coral 800A 15µm C8 1kg
B-800A-15-CO-C8W-100G	Bulk Coral 800A 15µm C8 Wax 100g
B-800A-15-CO-C8W-1KG	Bulk Coral 800A 15µm C8 Wax 1kg
B-800A-15-CO-C8P-100G	Bulk Coral 800A 15µm C8 Polar 100g
B-800A-15-CO-C8P-1KG	Bulk Coral 800A 15µm C8 Polar 1kg
B-800A-15-CO-C8T-100G	Bulk Coral 800A 15µm C8 TriF 100g
B-800A-15-CO-C8T-1KG	Bulk Coral 800A 15µm C8 TriF 1kg
B-800A-15-CO-C8EX-100G	Bulk Coral 800A 15µm C8 Extreme 100g
B-800A-15-CO-C8EX-1KG	Bulk Coral 800A 15µm C8 Extreme 1kg
B-800A-15-CO-C8PT-100G	Bulk Coral 800A 15µm C8 Polar TriF 100g
B-800A-15-CO-C8PT-1KG	Bulk Coral 800A 15µm C8 Polar TriF 1kg
B-800A-15-CO-C8AQP-100G	Bulk Coral 800A 15µm C8 AQ Plus 100g
B-800A-15-CO-C8AQP-1KG	Bulk Coral 800A 15µm C8 AQ Plus 1kg
B-800A-15-CO-C8AQ-100G	Bulk Coral 800A 15µm C8 AQ 100g
B-800A-15-CO-C8AQ-1KG	Bulk Coral 800A 15µm C8 AQ 1kg
<b>C18</b>	
B-800A-15-CO-C18-100G	Bulk Coral 800A 15µm C18 100g
B-800A-15-CO-C18-1KG	Bulk Coral 800A 15µm C18 1kg
B-800A-15-CO-C18W-100G	Bulk Coral 800A 15µm C18 Wax 100g
B-800A-15-CO-C18W-1KG	Bulk Coral 800A 15µm C18 Wax 1kg
B-800A-15-CO-C18P-100G	Bulk Coral 800A 15µm C18 Polar 100g
B-800A-15-CO-C18P-1KG	Bulk Coral 800A 15µm C18 Polar 1kg
B-800A-15-CO-C18T-100G	Bulk Coral 800A 15µm C18 TriF 100g
B-800A-15-CO-C18T-1KG	Bulk Coral 800A 15µm C18 TriF 1kg
B-800A-15-CO-C18EX-100G	Bulk Coral 800A 15µm C18 Extreme 100g
B-800A-15-CO-C18EX-1KG	Bulk Coral 800A 15µm C18 Extreme 1kg
B-800A-15-CO-C18PT-100G	Bulk Coral 800A 15µm C18 Polar TriF 100g
B-800A-15-CO-C18PT-1KG	Bulk Coral 800A 15µm C18 Polar TriF 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-15-CO-C18AQP-100G	Bulk Coral 800A 15µm C18 AQ Plus 100g
B-800A-15-CO-C18AQP-1KG	Bulk Coral 800A 15µm C18 AQ Plus 1kg
B-800A-15-CO-C18AQ-100G	Bulk Coral 800A 15µm C18 AQ 100g
B-800A-15-CO-C18AQ-1KG	Bulk Coral 800A 15µm C18 AQ 1kg
<b>C4</b>	
B-800A-15-CO-C4-100G	Bulk Coral 800A 15µm C4 100g
B-800A-15-CO-C4-1KG	Bulk Coral 800A 15µm C4 1kg
B-800A-15-CO-C4W-100G	Bulk Coral 800A 15µm C4 Wax 100g
B-800A-15-CO-C4W-1KG	Bulk Coral 800A 15µm C4 Wax 1kg
B-800A-15-CO-C4P-100G	Bulk Coral 800A 15µm C4 Polar 100g
B-800A-15-CO-C4P-1KG	Bulk Coral 800A 15µm C4 Polar 1kg
B-800A-15-CO-C4T-100G	Bulk Coral 800A 15µm C4 TriF 100g
B-800A-15-CO-C4T-1KG	Bulk Coral 800A 15µm C4 TriF 1kg
B-800A-15-CO-C4EX-100G	Bulk Coral 800A 15µm C4 Extreme 100g
B-800A-15-CO-C4EX-1KG	Bulk Coral 800A 15µm C4 Extreme 1kg
B-800A-15-CO-C4PT-100G	Bulk Coral 800A 15µm C4 Polar TriF 100g
B-800A-15-CO-C4PT-1KG	Bulk Coral 800A 15µm C4 Polar TriF 1kg
B-800A-15-CO-C4AQP-100G	Bulk Coral 800A 15µm C4 AQ Plus 100g
B-800A-15-CO-C4AQP-1KG	Bulk Coral 800A 15µm C4 AQ Plus 1kg
B-800A-15-CO-C4AQ-100G	Bulk Coral 800A 15µm C4 AQ 100g
B-800A-15-CO-C4AQ-1KG	Bulk Coral 800A 15µm C4 AQ 1kg
<b>Autres</b>	
B-800A-15-CO-COOH-100G	Bulk Coral 800A 15µm COOH 100g
B-800A-15-CO-COOH-1KG	Bulk Coral 800A 15µm COOH 1kg
B-800A-15-CO-SO3H-100G	Bulk Coral 800A 15µm SO3H 100g
B-800A-15-CO-SO3H-1KG	Bulk Coral 800A 15µm SO3H 1kg
B-800A-15-CO-DNH-100G	Bulk Coral 800A 15µm DNH 100g
B-800A-15-CO-DNH-1KG	Bulk Coral 800A 15µm DNH 1kg
B-800A-15-CO-PEI-100G	Bulk Coral 800A 15µm PEI 100g
B-800A-15-CO-PEI-1KG	Bulk Coral 800A 15µm PEI 1kg
B-800A-15-CO-ARG-100G	Bulk Coral 800A 15µm ARG 100g
B-800A-15-CO-ARG-1KG	Bulk Coral 800A 15µm ARG 1kg
B-800A-15-CO-NH-100G	Bulk Coral 800A 15µm NH 100g
B-800A-15-CO-NH-1KG	Bulk Coral 800A 15µm NH 1kg
B-800A-15-CO-NH2-100G	Bulk Coral 800A 15µm NH2 100g
B-800A-15-CO-NH2-1KG	Bulk Coral 800A 15µm NH2 1kg
B-800A-15-CO-CN-100G	Bulk Coral 800A 15µm CN 100g
B-800A-15-CO-CN-1KG	Bulk Coral 800A 15µm CN 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-15-CO-Phenyl-100G	Bulk Coral 800A 15µm Phenyl 100g
B-800A-15-CO-Phenyl-1KG	Bulk Coral 800A 15µm Phenyl 1kg
B-800A-15-CO-Phenyl HE-100G	Bulk Coral 800A 15µm Phenyl HE 100g
B-800A-15-CO-Phenyl HE-1KG	Bulk Coral 800A 15µm Phenyl HE 1kg
B-800A-15-CO-Diol-100G	Bulk Coral 800A 15µm Diol 100g
B-800A-15-CO-Diol-1KG	Bulk Coral 800A 15µm Diol 1kg
<b>Particule Coral 15, Pore 1000A</b>	
B-1000A-15-CO-Si-100G	Bulk Coral 1000A 15µm 100g
B-1000A-15-CO-Si-1KG	Bulk Coral 1000A 15µm 1kg
<b>C8</b>	
B-1000A-15-CO-C8-100G	Bulk Coral 1000A 15µm C8 100g
B-1000A-15-CO-C8-1KG	Bulk Coral 1000A 15µm C8 1kg
B-1000A-15-CO-C8W-100G	Bulk Coral 1000A 15µm C8 Wax 100g
B-1000A-15-CO-C8W-1KG	Bulk Coral 1000A 15µm C8 Wax 1kg
B-1000A-15-CO-C8P-100G	Bulk Coral 1000A 15µm C8 Polar 100g
B-1000A-15-CO-C8P-1KG	Bulk Coral 1000A 15µm C8 Polar 1kg
B-1000A-15-CO-C8T-100G	Bulk Coral 1000A 15µm C8 TriF 100g
B-1000A-15-CO-C8T-1KG	Bulk Coral 1000A 15µm C8 TriF 1kg
B-1000A-15-CO-C8EX-100G	Bulk Coral 1000A 15µm C8 Extreme 100g
B-1000A-15-CO-C8EX-1KG	Bulk Coral 1000A 15µm C8 Extreme 1kg
B-1000A-15-CO-C8PT-100G	Bulk Coral 1000A 15µm C8 Polar TriF 100g
B-1000A-15-CO-C8PT-1KG	Bulk Coral 1000A 15µm C8 Polar TriF 1kg
B-1000A-15-CO-C8AQP-100G	Bulk Coral 1000A 15µm C8 AQ Plus 100g
B-1000A-15-CO-C8AQP-1KG	Bulk Coral 1000A 15µm C8 AQ Plus 1kg
B-1000A-15-CO-C8AQ-100G	Bulk Coral 1000A 15µm C8 AQ 100g
B-1000A-15-CO-C8AQ-1KG	Bulk Coral 1000A 15µm C8 AQ 1kg
<b>C18</b>	
B-1000A-15-CO-C18-100G	Bulk Coral 1000A 15µm C18 100g
B-1000A-15-CO-C18-1KG	Bulk Coral 1000A 15µm C18 1kg
B-1000A-15-CO-C18W-100G	Bulk Coral 1000A 15µm C18 Wax 100g
B-1000A-15-CO-C18W-1KG	Bulk Coral 1000A 15µm C18 Wax 1kg
B-1000A-15-CO-C18P-100G	Bulk Coral 1000A 15µm C18 Polar 100g
B-1000A-15-CO-C18P-1KG	Bulk Coral 1000A 15µm C18 Polar 1kg
B-1000A-15-CO-C18T-100G	Bulk Coral 1000A 15µm C18 TriF 100g
B-1000A-15-CO-C18T-1KG	Bulk Coral 1000A 15µm C18 TriF 1kg
B-1000A-15-CO-C18EX-100G	Bulk Coral 1000A 15µm C18 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-1000A-15-CO-C18EX-1KG	Bulk Coral 1000A 15µm C18 Extreme 1kg
B-1000A-15-CO-C18PT-100G	Bulk Coral 1000A 15µm C18 Polar TriF 100g
B-1000A-15-CO-C18PT-1KG	Bulk Coral 1000A 15µm C18 Polar TriF 1kg
B-1000A-15-CO-C18AQP-100G	Bulk Coral 1000A 15µm C18 AQ Plus 100g
B-1000A-15-CO-C18AQP-1KG	Bulk Coral 1000A 15µm C18 AQ Plus 1kg
B-1000A-15-CO-C18AQ-100G	Bulk Coral 1000A 15µm C18 AQ 100g
B-1000A-15-CO-C18AQ-1KG	Bulk Coral 1000A 15µm C18 AQ 1kg
<b>C4</b>	
B-1000A-15-CO-C4-100G	Bulk Coral 1000A 15µm C4 100g
B-1000A-15-CO-C4-1KG	Bulk Coral 1000A 15µm C4 1kg
B-1000A-15-CO-C4W-100G	Bulk Coral 1000A 15µm C4 Wax 100g
B-1000A-15-CO-C4W-1KG	Bulk Coral 1000A 15µm C4 Wax 1kg
B-1000A-15-CO-C4P-100G	Bulk Coral 1000A 15µm C4 Polar 100g
B-1000A-15-CO-C4P-1KG	Bulk Coral 1000A 15µm C4 Polar 1kg
B-1000A-15-CO-C4T-100G	Bulk Coral 1000A 15µm C4 TriF 100g
B-1000A-15-CO-C4T-1KG	Bulk Coral 1000A 15µm C4 TriF 1kg
B-1000A-15-CO-C4EX-100G	Bulk Coral 1000A 15µm C4 Extreme 100g
B-1000A-15-CO-C4EX-1KG	Bulk Coral 1000A 15µm C4 Extreme 1kg
B-1000A-15-CO-C4PT-100G	Bulk Coral 1000A 15µm C4 Polar TriF 100g
B-1000A-15-CO-C4PT-1KG	Bulk Coral 1000A 15µm C4 Polar TriF 1kg
B-1000A-15-CO-C4AQP-100G	Bulk Coral 1000A 15µm C4 AQ Plus 100g
B-1000A-15-CO-C4AQP-1KG	Bulk Coral 1000A 15µm C4 AQ Plus 1kg
B-1000A-15-CO-C4AQ-100G	Bulk Coral 1000A 15µm C4 AQ 100g
B-1000A-15-CO-C4AQ-1KG	Bulk Coral 1000A 15µm C4 AQ 1kg
<b>Autres</b>	
B-1000A-15-CO-COOH-100G	Bulk Coral 1000A 15µm COOH 100g
B-1000A-15-CO-COOH-1KG	Bulk Coral 1000A 15µm COOH 1kg
B-1000A-15-CO-SO3H-100G	Bulk Coral 1000A 15µm SO3H 100g
B-1000A-15-CO-SO3H-1KG	Bulk Coral 1000A 15µm SO3H 1kg
B-1000A-15-CO-DNH-100G	Bulk Coral 1000A 15µm DNH 100g
B-1000A-15-CO-DNH-1KG	Bulk Coral 1000A 15µm DNH 1kg
B-1000A-15-CO-PEI-100G	Bulk Coral 1000A 15µm PEI 100g
B-1000A-15-CO-PEI-1KG	Bulk Coral 1000A 15µm PEI 1kg
B-1000A-15-CO-ARG-100G	Bulk Coral 1000A 15µm ARG 100g
B-1000A-15-CO-ARG-1KG	Bulk Coral 1000A 15µm ARG 1kg
B-1000A-15-CO-NH-100G	Bulk Coral 1000A 15µm NH 100g
B-1000A-15-CO-NH-1KG	Bulk Coral 1000A 15µm NH 1kg
B-1000A-15-CO-NH2-100G	Bulk Coral 1000A 15µm NH2 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-1000A-15-CO-NH2-1KG	Bulk Coral 1000A 15µm NH2 1kg
B-1000A-15-CO-CN-100G	Bulk Coral 1000A 15µm CN 100g
B-1000A-15-CO-CN-1KG	Bulk Coral 1000A 15µm CN 1kg
B-1000A-15-CO-Phenyl-100G	Bulk Coral 1000A 15µm Phenyl 100g
B-1000A-15-CO-Phenyl-1KG	Bulk Coral 1000A 15µm Phenyl 1kg
B-1000A-15-CO-Phenyl HE-100G	Bulk Coral 1000A 15µm Phenyl HE 100g
B-1000A-15-CO-Phenyl HE-1KG	Bulk Coral 1000A 15µm Phenyl HE 1kg
B-1000A-15-CO-Diol-100G	Bulk Coral 1000A 15µm Diol 100g
B-1000A-15-CO-Diol-1KG	Bulk Coral 1000A 15µm Diol 1kg
<b>Coral 30µm</b>	
<b>Particule Coral 30, Pore 70A</b>	
B-70A-30-CO-Si-100G	Bulk Coral 70A 30µm 100g
B-70A-30-CO-Si-1KG	Bulk Coral 70A 30µm 1kg
<b>C8</b>	
B-70A-30-CO-C8-100G	Bulk Coral 70A 30µm C8 100g
B-70A-30-CO-C8-1KG	Bulk Coral 70A 30µm C8 1kg
B-70A-30-CO-C8W-100G	Bulk Coral 70A 30µm C8 Wax 100g
B-70A-30-CO-C8W-1KG	Bulk Coral 70A 30µm C8 Wax 1kg
B-70A-30-CO-C8P-100G	Bulk Coral 70A 30µm C8 Polar 100g
B-70A-30-CO-C8P-1KG	Bulk Coral 70A 30µm C8 Polar 1kg
B-70A-30-CO-C8T-100G	Bulk Coral 70A 30µm C8 TriF 100g
B-70A-30-CO-C8T-1KG	Bulk Coral 70A 30µm C8 TriF 1kg
B-70A-30-CO-C8EX-100G	Bulk Coral 70A 30µm C8 Extreme 100g
B-70A-30-CO-C8EX-1KG	Bulk Coral 70A 30µm C8 Extreme 1kg
B-70A-30-CO-C8PT-100G	Bulk Coral 70A 30µm C8 Polar TriF 100g
B-70A-30-CO-C8PT-1KG	Bulk Coral 70A 30µm C8 Polar TriF 1kg
B-70A-30-CO-C8AQP-100G	Bulk Coral 70A 30µm C8 AQ Plus 100g
B-70A-30-CO-C8AQP-1KG	Bulk Coral 70A 30µm C8 AQ Plus 1kg
B-70A-30-CO-C8AQ-100G	Bulk Coral 70A 30µm C8 AQ 100g
B-70A-30-CO-C8AQ-1KG	Bulk Coral 70A 30µm C8 AQ 1kg
<b>C18</b>	
B-70A-30-CO-C18-100G	Bulk Coral 70A 30µm C18 100g
B-70A-30-CO-C18-1KG	Bulk Coral 70A 30µm C18 1kg
B-70A-30-CO-C18W-100G	Bulk Coral 70A 30µm C18 Wax 100g
B-70A-30-CO-C18W-1KG	Bulk Coral 70A 30µm C18 Wax 1kg
B-70A-30-CO-C18P-100G	Bulk Coral 70A 30µm C18 Polar 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-30-CO-C18P-1KG	Bulk Coral 70A 30µm C18 Polar 1kg
B-70A-30-CO-C18T-100G	Bulk Coral 70A 30µm C18 TriF 100g
B-70A-30-CO-C18T-1KG	Bulk Coral 70A 30µm C18 TriF 1kg
B-70A-30-CO-C18EX-100G	Bulk Coral 70A 30µm C18 Extreme 100g
B-70A-30-CO-C18EX-1KG	Bulk Coral 70A 30µm C18 Extreme 1kg
B-70A-30-CO-C18PT-100G	Bulk Coral 70A 30µm C18 Polar TriF 100g
B-70A-30-CO-C18PT-1KG	Bulk Coral 70A 30µm C18 Polar TriF 1kg
B-70A-30-CO-C18AQP-100G	Bulk Coral 70A 30µm C18 AQ Plus 100g
B-70A-30-CO-C18AQP-1KG	Bulk Coral 70A 30µm C18 AQ Plus 1kg
B-70A-30-CO-C18AQ-100G	Bulk Coral 70A 30µm C18 AQ 100g
B-70A-30-CO-C18AQ-1KG	Bulk Coral 70A 30µm C18 AQ 1kg
<b>C4</b>	
B-70A-30-CO-C4-100G	Bulk Coral 70A 30µm C4 100g
B-70A-30-CO-C4-1KG	Bulk Coral 70A 30µm C4 1kg
B-70A-30-CO-C4W-100G	Bulk Coral 70A 30µm C4 Wax 100g
B-70A-30-CO-C4W-1KG	Bulk Coral 70A 30µm C4 Wax 1kg
B-70A-30-CO-C4P-100G	Bulk Coral 70A 30µm C4 Polar 100g
B-70A-30-CO-C4P-1KG	Bulk Coral 70A 30µm C4 Polar 1kg
B-70A-30-CO-C4T-100G	Bulk Coral 70A 30µm C4 TriF 100g
B-70A-30-CO-C4T-1KG	Bulk Coral 70A 30µm C4 TriF 1kg
B-70A-30-CO-C4EX-100G	Bulk Coral 70A 30µm C4 Extreme 100g
B-70A-30-CO-C4EX-1KG	Bulk Coral 70A 30µm C4 Extreme 1kg
B-70A-30-CO-C4PT-100G	Bulk Coral 70A 30µm C4 Polar TriF 100g
B-70A-30-CO-C4PT-1KG	Bulk Coral 70A 30µm C4 Polar TriF 1kg
B-70A-30-CO-C4AQP-100G	Bulk Coral 70A 30µm C4 AQ Plus 100g
B-70A-30-CO-C4AQP-1KG	Bulk Coral 70A 30µm C4 AQ Plus 1kg
B-70A-30-CO-C4AQ-100G	Bulk Coral 70A 30µm C4 AQ 100g
B-70A-30-CO-C4AQ-1KG	Bulk Coral 70A 30µm C4 AQ 1kg
<b>Autres</b>	
B-70A-30-CO-COOH-100G	Bulk Coral 70A 30µm COOH 100g
B-70A-30-CO-COOH-1KG	Bulk Coral 70A 30µm COOH 1kg
B-70A-30-CO-SO3H-100G	Bulk Coral 70A 30µm SO3H 100g
B-70A-30-CO-SO3H-1KG	Bulk Coral 70A 30µm SO3H 1kg
B-70A-30-CO-DNH-100G	Bulk Coral 70A 30µm DNH 100g
B-70A-30-CO-DNH-1KG	Bulk Coral 70A 30µm DNH 1kg
B-70A-30-CO-PEI-100G	Bulk Coral 70A 30µm PEI 100g
B-70A-30-CO-PEI-1KG	Bulk Coral 70A 30µm PEI 1kg
B-70A-30-CO-ARG-100G	Bulk Coral 70A 30µm ARG 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70A-30-CO-ARG-1KG	Bulk Coral 70A 30µm ARG 1kg
B-70A-30-CO-NH-100G	Bulk Coral 70A 30µm NH 100g
B-70A-30-CO-NH-1KG	Bulk Coral 70A 30µm NH 1kg
B-70A-30-CO-CO-NH2-100G	Bulk Coral 70A 30µm NH2 100g
B-70A-30-CO-CO-NH2-1KG	Bulk Coral 70A 30µm NH2 1kg
B-70A-30-CO-CN-100G	Bulk Coral 70A 30µm CN 100g
B-70A-30-CO-CN-1KG	Bulk Coral 70A 30µm CN 1kg
B-70A-30-CO-Phenyl-100G	Bulk Coral 70A 30µm Phenyl 100g
B-70A-30-CO-Phenyl-1KG	Bulk Coral 70A 30µm Phenyl 1kg
B-70A-30-CO-Phenyl HE-100G	Bulk Coral 70A 30µm Phenyl HE 100g
B-70A-30-CO-Phenyl HE-1KG	Bulk Coral 70A 30µm Phenyl HE 1kg
B-70A-30-CO-Diol-100G	Bulk Coral 70A 30µm Diol 100g
B-70A-30-CO-Diol-1KG	Bulk Coral 70A 30µm Diol 1kg
<b>Particule Coral 30, Pore 100A</b>	
B-100A-30-CO-Si-100G	Bulk Coral 100A 30µm 100g
B-100A-30-CO-Si-1KG	Bulk Coral 100A 30µm 1kg
<b>C8</b>	
B-100A-30-CO-C8-100G	Bulk Coral 100A 30µm C8 100g
B-100A-30-CO-C8-1KG	Bulk Coral 100A 30µm C8 1kg
B-100A-30-CO-C8W-100G	Bulk Coral 100A 30µm C8 Wax 100g
B-100A-30-CO-C8W-1KG	Bulk Coral 100A 30µm C8 Wax 1kg
B-100A-30-CO-C8P-100G	Bulk Coral 100A 30µm C8 Polar 100g
B-100A-30-CO-C8P-1KG	Bulk Coral 100A 30µm C8 Polar 1kg
B-100A-30-CO-C8T-100G	Bulk Coral 100A 30µm C8 TriF 100g
B-100A-30-CO-C8T-1KG	Bulk Coral 100A 30µm C8 TriF 1kg
B-100A-30-CO-C8EX-100G	Bulk Coral 100A 30µm C8 Extreme 100g
B-100A-30-CO-C8EX-1KG	Bulk Coral 100A 30µm C8 Extreme 1kg
B-100A-30-CO-C8PT-100G	Bulk Coral 100A 30µm C8 Polar TriF 100g
B-100A-30-CO-C8PT-1KG	Bulk Coral 100A 30µm C8 Polar TriF 1kg
B-100A-30-CO-C8AQP-100G	Bulk Coral 100A 30µm C8 AQ Plus 100g
B-100A-30-CO-C8AQP-1KG	Bulk Coral 100A 30µm C8 AQ Plus 1kg
B-100A-30-CO-C8AQ-100G	Bulk Coral 100A 30µm C8 AQ 100g
B-100A-30-CO-C8AQ-1KG	Bulk Coral 100A 30µm C8 AQ 1kg
<b>C18</b>	
B-100A-30-CO-C18-100G	Bulk Coral 100A 30µm C18 100g
B-100A-30-CO-C18-1KG	Bulk Coral 100A 30µm C18 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-30-CO-C18W-100G	Bulk Coral 100A 30µm C18 Wax 100g
B-100A-30-CO-C18W-1KG	Bulk Coral 100A 30µm C18 Wax 1kg
B-100A-30-CO-C18P-100G	Bulk Coral 100A 30µm C18 Polar 100g
B-100A-30-CO-C18P-1KG	Bulk Coral 100A 30µm C18 Polar 1kg
B-100A-30-CO-C18T-100G	Bulk Coral 100A 30µm C18 TriF 100g
B-100A-30-CO-C18T-1KG	Bulk Coral 100A 30µm C18 TriF 1kg
B-100A-30-CO-C18EX-100G	Bulk Coral 100A 30µm C18 Extreme 100g
B-100A-30-CO-C18EX-1KG	Bulk Coral 100A 30µm C18 Extreme 1kg
B-100A-30-CO-C18PT-100G	Bulk Coral 100A 30µm C18 Polar TriF 100g
B-100A-30-CO-C18PT-1KG	Bulk Coral 100A 30µm C18 Polar TriF 1kg
B-100A-30-CO-C18AQP-100G	Bulk Coral 100A 30µm C18 AQ Plus 100g
B-100A-30-CO-C18AQP-1KG	Bulk Coral 100A 30µm C18 AQ Plus 1kg
B-100A-30-CO-C18AQ-100G	Bulk Coral 100A 30µm C18 AQ 100g
B-100A-30-CO-C18AQ-1KG	Bulk Coral 100A 30µm C18 AQ 1kg
<b>C4</b>	
B-100A-30-CO-C4-100G	Bulk Coral 100A 30µm C4 100g
B-100A-30-CO-C4-1KG	Bulk Coral 100A 30µm C4 1kg
B-100A-30-CO-C4W-100G	Bulk Coral 100A 30µm C4 Wax 100g
B-100A-30-CO-C4W-1KG	Bulk Coral 100A 30µm C4 Wax 1kg
B-100A-30-CO-C4P-100G	Bulk Coral 100A 30µm C4 Polar 100g
B-100A-30-CO-C4P-1KG	Bulk Coral 100A 30µm C4 Polar 1kg
B-100A-30-CO-C4T-100G	Bulk Coral 100A 30µm C4 TriF 100g
B-100A-30-CO-C4T-1KG	Bulk Coral 100A 30µm C4 TriF 1kg
B-100A-30-CO-C4EX-100G	Bulk Coral 100A 30µm C4 Extreme 100g
B-100A-30-CO-C4EX-1KG	Bulk Coral 100A 30µm C4 Extreme 1kg
B-100A-30-CO-C4PT-100G	Bulk Coral 100A 30µm C4 Polar TriF 100g
B-100A-30-CO-C4PT-1KG	Bulk Coral 100A 30µm C4 Polar TriF 1kg
B-100A-30-CO-C4AQP-100G	Bulk Coral 100A 30µm C4 AQ Plus 100g
B-100A-30-CO-C4AQP-1KG	Bulk Coral 100A 30µm C4 AQ Plus 1kg
B-100A-30-CO-C4AQ-100G	Bulk Coral 100A 30µm C4 AQ 100g
B-100A-30-CO-C4AQ-1KG	Bulk Coral 100A 30µm C4 AQ 1kg
<b>Autres</b>	
B-100A-30-CO-COOH-100G	Bulk Coral 100A 30µm COOH 100g
B-100A-30-CO-COOH-1KG	Bulk Coral 100A 30µm COOH 1kg
B-100A-30-CO-SO3H-100G	Bulk Coral 100A 30µm SO3H 100g
B-100A-30-CO-SO3H-1KG	Bulk Coral 100A 30µm SO3H 1kg
B-100A-30-CO-DNH-100G	Bulk Coral 100A 30µm DNH 100g
B-100A-30-CO-DNH-1KG	Bulk Coral 100A 30µm DNH 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100A-30-CO-PEI-100G	Bulk Coral 100A 30µm PEI 100g
B-100A-30-CO-PEI-1KG	Bulk Coral 100A 30µm PEI 1kg
B-100A-30-CO-ARG-100G	Bulk Coral 100A 30µm ARG 100g
B-100A-30-CO-ARG-1KG	Bulk Coral 100A 30µm ARG 1kg
B-100A-30-CO-NH-100G	Bulk Coral 100A 30µm NH 100g
B-100A-30-CO-NH-1KG	Bulk Coral 100A 30µm NH 1kg
B-100A-30-CO-CO-NH2-100G	Bulk Coral 100A 30µm NH2 100g
B-100A-30-CO-CO-NH2-1KG	Bulk Coral 100A 30µm NH2 1kg
B-100A-30-CO-CN-100G	Bulk Coral 100A 30µm CN 100g
B-100A-30-CO-CN-1KG	Bulk Coral 100A 30µm CN 1kg
B-100A-30-CO-Phenyl-100G	Bulk Coral 100A 30µm Phenyl 100g
B-100A-30-CO-Phenyl-1KG	Bulk Coral 100A 30µm Phenyl 1kg
B-100A-30-CO-Phenyl HE-100G	Bulk Coral 100A 30µm Phenyl HE 100g
B-100A-30-CO-Phenyl HE-1KG	Bulk Coral 100A 30µm Phenyl HE 1kg
B-100A-30-CO-Diol-100G	Bulk Coral 100A 30µm Diol 100g
B-100A-30-CO-Diol-1KG	Bulk Coral 100A 30µm Diol 1kg
<b>Particule Coral 30, Pore 150A</b>	
B-150A-30-CO-Si-100G	Bulk Coral 150A 30µm 100g
B-150A-30-CO-Si-1KG	Bulk Coral 150A 30µm 1kg
<b>C8</b>	
B-150A-30-CO-C8-100G	Bulk Coral 150A 30µm C8 100g
B-150A-30-CO-C8-1KG	Bulk Coral 150A 30µm C8 1kg
B-150A-30-CO-C8W-100G	Bulk Coral 150A 30µm C8 Wax 100g
B-150A-30-CO-C8W-1KG	Bulk Coral 150A 30µm C8 Wax 1kg
B-150A-30-CO-C8P-100G	Bulk Coral 150A 30µm C8 Polar 100g
B-150A-30-CO-C8P-1KG	Bulk Coral 150A 30µm C8 Polar 1kg
B-150A-30-CO-C8T-100G	Bulk Coral 150A 30µm C8 TriF 100g
B-150A-30-CO-C8T-1KG	Bulk Coral 150A 30µm C8 TriF 1kg
B-150A-30-CO-C8EX-100G	Bulk Coral 150A 30µm C8 Extreme 100g
B-150A-30-CO-C8EX-1KG	Bulk Coral 150A 30µm C8 Extreme 1kg
B-150A-30-CO-C8PT-100G	Bulk Coral 150A 30µm C8 Polar TriF 100g
B-150A-30-CO-C8PT-1KG	Bulk Coral 150A 30µm C8 Polar TriF 1kg
B-150A-30-CO-C8AQP-100G	Bulk Coral 150A 30µm C8 AQ Plus 100g
B-150A-30-CO-C8AQP-1KG	Bulk Coral 150A 30µm C8 AQ Plus 1kg
B-150A-30-CO-C8AQ-100G	Bulk Coral 150A 30µm C8 AQ 100g
B-150A-30-CO-C8AQ-1KG	Bulk Coral 150A 30µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-150A-30-CO-C18-100G	Bulk Coral 150A 30µm C18 100g
B-150A-30-CO-C18-1KG	Bulk Coral 150A 30µm C18 1kg
B-150A-30-CO-C18W-100G	Bulk Coral 150A 30µm C18 Wax 100g
B-150A-30-CO-C18W-1KG	Bulk Coral 150A 30µm C18 Wax 1kg
B-150A-30-CO-C18P-100G	Bulk Coral 150A 30µm C18 Polar 100g
B-150A-30-CO-C18P-1KG	Bulk Coral 150A 30µm C18 Polar 1kg
B-150A-30-CO-C18T-100G	Bulk Coral 150A 30µm C18 TriF 100g
B-150A-30-CO-C18T-1KG	Bulk Coral 150A 30µm C18 TriF 1kg
B-150A-30-CO-C18EX-100G	Bulk Coral 150A 30µm C18 Extreme 100g
B-150A-30-CO-C18EX-1KG	Bulk Coral 150A 30µm C18 Extreme 1kg
B-150A-30-CO-C18PT-100G	Bulk Coral 150A 30µm C18 Polar TriF 100g
B-150A-30-CO-C18PT-1KG	Bulk Coral 150A 30µm C18 Polar TriF 1kg
B-150A-30-CO-C18AQP-100G	Bulk Coral 150A 30µm C18 AQ Plus 100g
B-150A-30-CO-C18AQP-1KG	Bulk Coral 150A 30µm C18 AQ Plus 1kg
B-150A-30-CO-C18AQ-100G	Bulk Coral 150A 30µm C18 AQ 100g
B-150A-30-CO-C18AQ-1KG	Bulk Coral 150A 30µm C18 AQ 1kg
<b>C4</b>	
B-150A-30-CO-C4-100G	Bulk Coral 150A 30µm C4 100g
B-150A-30-CO-C4-1KG	Bulk Coral 150A 30µm C4 1kg
B-150A-30-CO-C4W-100G	Bulk Coral 150A 30µm C4 Wax 100g
B-150A-30-CO-C4W-1KG	Bulk Coral 150A 30µm C4 Wax 1kg
B-150A-30-CO-C4P-100G	Bulk Coral 150A 30µm C4 Polar 100g
B-150A-30-CO-C4P-1KG	Bulk Coral 150A 30µm C4 Polar 1kg
B-150A-30-CO-C4T-100G	Bulk Coral 150A 30µm C4 TriF 100g
B-150A-30-CO-C4T-1KG	Bulk Coral 150A 30µm C4 TriF 1kg
B-150A-30-CO-C4EX-100G	Bulk Coral 150A 30µm C4 Extreme 100g
B-150A-30-CO-C4EX-1KG	Bulk Coral 150A 30µm C4 Extreme 1kg
B-150A-30-CO-C4PT-100G	Bulk Coral 150A 30µm C4 Polar TriF 100g
B-150A-30-CO-C4PT-1KG	Bulk Coral 150A 30µm C4 Polar TriF 1kg
B-150A-30-CO-C4AQP-100G	Bulk Coral 150A 30µm C4 AQ Plus 100g
B-150A-30-CO-C4AQP-1KG	Bulk Coral 150A 30µm C4 AQ Plus 1kg
B-150A-30-CO-C4AQ-100G	Bulk Coral 150A 30µm C4 AQ 100g
B-150A-30-CO-C4AQ-1KG	Bulk Coral 150A 30µm C4 AQ 1kg
<b>Autres</b>	
B-150A-30-CO-COOH-100G	Bulk Coral 150A 30µm COOH 100g
B-150A-30-CO-COOH-1KG	Bulk Coral 150A 30µm COOH 1kg
B-150A-30-CO-SO3H-100G	Bulk Coral 150A 30µm SO3H 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-150A-30-CO-S03H-1KG	Bulk Coral 150A 30µm S03H 1kg
B-150A-30-CO-DNH-100G	Bulk Coral 150A 30µm DNH 100g
B-150A-30-CO-DNH-1KG	Bulk Coral 150A 30µm DNH 1kg
B-150A-30-CO-PEI-100G	Bulk Coral 150A 30µm PEI 100g
B-150A-30-CO-PEI-1KG	Bulk Coral 150A 30µm PEI 1kg
B-150A-30-CO-ARG-100G	Bulk Coral 150A 30µm ARG 100g
B-150A-30-CO-ARG-1KG	Bulk Coral 150A 30µm ARG 1kg
B-150A-30-CO-NH-100G	Bulk Coral 150A 30µm NH 100g
B-150A-30-CO-NH-1KG	Bulk Coral 150A 30µm NH 1kg
B-150A-30-CO-CO-NH2-100G	Bulk Coral 150A 30µm NH2 100g
B-150A-30-CO-CO-NH2-1KG	Bulk Coral 150A 30µm NH2 1kg
B-150A-30-CO-CN-100G	Bulk Coral 150A 30µm CN 100g
B-150A-30-CO-CN-1KG	Bulk Coral 150A 30µm CN 1kg
B-150A-30-CO-Phenyl-100G	Bulk Coral 150A 30µm Phenyl 100g
B-150A-30-CO-Phenyl-1KG	Bulk Coral 150A 30µm Phenyl 1kg
B-150A-30-CO-Phenyl HE-100G	Bulk Coral 150A 30µm Phenyl HE 100g
B-150A-30-CO-Phenyl HE-1KG	Bulk Coral 150A 30µm Phenyl HE 1kg
B-150A-30-CO-Diol-100G	Bulk Coral 150A 30µm Diol 100g
B-150A-30-CO-Diol-1KG	Bulk Coral 150A 30µm Diol 1kg
<b>Particule Coral 30, Pore 200A</b>	
B-200A-30-CO-Si-100G	Bulk Coral 200A 30µm 100g
B-200A-30-CO-Si-1KG	Bulk Coral 200A 30µm 1kg
<b>C8</b>	
B-200A-30-CO-C8-100G	Bulk Coral 200A 30µm C8 100g
B-200A-30-CO-C8-1KG	Bulk Coral 200A 30µm C8 1kg
B-200A-30-CO-C8W-100G	Bulk Coral 200A 30µm C8 Wax 100g
B-200A-30-CO-C8W-1KG	Bulk Coral 200A 30µm C8 Wax 1kg
B-200A-30-CO-C8P-100G	Bulk Coral 200A 30µm C8 Polar 100g
B-200A-30-CO-C8P-1KG	Bulk Coral 200A 30µm C8 Polar 1kg
B-200A-30-CO-C8T-100G	Bulk Coral 200A 30µm C8 TriF 100g
B-200A-30-CO-C8T-1KG	Bulk Coral 200A 30µm C8 TriF 1kg
B-200A-30-CO-C8EX-100G	Bulk Coral 200A 30µm C8 Extreme 100g
B-200A-30-CO-C8EX-1KG	Bulk Coral 200A 30µm C8 Extreme 1kg
B-200A-30-CO-C8PT-100G	Bulk Coral 200A 30µm C8 Polar TriF 100g
B-200A-30-CO-C8PT-1KG	Bulk Coral 200A 30µm C8 Polar TriF 1kg
B-200A-30-CO-C8AQP-100G	Bulk Coral 200A 30µm C8 AQ Plus 100g
B-200A-30-CO-C8AQP-1KG	Bulk Coral 200A 30µm C8 AQ Plus 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-200A-30-CO-C8AQ-100G	Bulk Coral 200A 30µm C8 AQ 100g
B-200A-30-CO-C8AQ-1KG	Bulk Coral 200A 30µm C8 AQ 1kg
<b>C18</b>	
B-200A-30-CO-C18-100G	Bulk Coral 200A 30µm C18 100g
B-200A-30-CO-C18-1KG	Bulk Coral 200A 30µm C18 1kg
B-200A-30-CO-C18W-100G	Bulk Coral 200A 30µm C18 Wax 100g
B-200A-30-CO-C18W-1KG	Bulk Coral 200A 30µm C18 Wax 1kg
B-200A-30-CO-C18P-100G	Bulk Coral 200A 30µm C18 Polar 100g
B-200A-30-CO-C18P-1KG	Bulk Coral 200A 30µm C18 Polar 1kg
B-200A-30-CO-C18T-100G	Bulk Coral 200A 30µm C18 TriF 100g
B-200A-30-CO-C18T-1KG	Bulk Coral 200A 30µm C18 TriF 1kg
B-200A-30-CO-C18EX-100G	Bulk Coral 200A 30µm C18 Extreme 100g
B-200A-30-CO-C18EX-1KG	Bulk Coral 200A 30µm C18 Extreme 1kg
B-200A-30-CO-C18PT-100G	Bulk Coral 200A 30µm C18 Polar TriF 100g
B-200A-30-CO-C18PT-1KG	Bulk Coral 200A 30µm C18 Polar TriF 1kg
B-200A-30-CO-C18AQP-100G	Bulk Coral 200A 30µm C18 AQ Plus 100g
B-200A-30-CO-C18AQP-1KG	Bulk Coral 200A 30µm C18 AQ Plus 1kg
B-200A-30-CO-C18AQ-100G	Bulk Coral 200A 30µm C18 AQ 100g
B-200A-30-CO-C18AQ-1KG	Bulk Coral 200A 30µm C18 AQ 1kg
<b>C4</b>	
B-200A-30-CO-C4-100G	Bulk Coral 200A 30µm C4 100g
B-200A-30-CO-C4-1KG	Bulk Coral 200A 30µm C4 1kg
B-200A-30-CO-C4W-100G	Bulk Coral 200A 30µm C4 Wax 100g
B-200A-30-CO-C4W-1KG	Bulk Coral 200A 30µm C4 Wax 1kg
B-200A-30-CO-C4P-100G	Bulk Coral 200A 30µm C4 Polar 100g
B-200A-30-CO-C4P-1KG	Bulk Coral 200A 30µm C4 Polar 1kg
B-200A-30-CO-C4T-100G	Bulk Coral 200A 30µm C4 TriF 100g
B-200A-30-CO-C4T-1KG	Bulk Coral 200A 30µm C4 TriF 1kg
B-200A-30-CO-C4EX-100G	Bulk Coral 200A 30µm C4 Extreme 100g
B-200A-30-CO-C4EX-1KG	Bulk Coral 200A 30µm C4 Extreme 1kg
B-200A-30-CO-C4PT-100G	Bulk Coral 200A 30µm C4 Polar TriF 100g
B-200A-30-CO-C4PT-1KG	Bulk Coral 200A 30µm C4 Polar TriF 1kg
B-200A-30-CO-C4AQP-100G	Bulk Coral 200A 30µm C4 AQ Plus 100g
B-200A-30-CO-C4AQP-1KG	Bulk Coral 200A 30µm C4 AQ Plus 1kg
B-200A-30-CO-C4AQ-100G	Bulk Coral 200A 30µm C4 AQ 100g
B-200A-30-CO-C4AQ-1KG	Bulk Coral 200A 30µm C4 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Autres</b>	
B-200A-30-CO-COOH-100G	Bulk Coral 200A 30µm COOH 100g
B-200A-30-CO-COOH-1KG	Bulk Coral 200A 30µm COOH 1kg
B-200A-30-CO-SO3H-100G	Bulk Coral 200A 30µm SO3H 100g
B-200A-30-CO-SO3H-1KG	Bulk Coral 200A 30µm SO3H 1kg
B-200A-30-CO-DNH-100G	Bulk Coral 200A 30µm DNH 100g
B-200A-30-CO-DNH-1KG	Bulk Coral 200A 30µm DNH 1kg
B-200A-30-CO-PEI-100G	Bulk Coral 200A 30µm PEI 100g
B-200A-30-CO-PEI-1KG	Bulk Coral 200A 30µm PEI 1kg
B-200A-30-CO-ARG-100G	Bulk Coral 200A 30µm ARG 100g
B-200A-30-CO-ARG-1KG	Bulk Coral 200A 30µm ARG 1kg
B-200A-30-CO-NH-100G	Bulk Coral 200A 30µm NH 100g
B-200A-30-CO-NH-1KG	Bulk Coral 200A 30µm NH 1kg
B-200A-30-CO-CO-NH2-100G	Bulk Coral 200A 30µm NH2 100g
B-200A-30-CO-CO-NH2-1KG	Bulk Coral 200A 30µm NH2 1kg
B-200A-30-CO-CN-100G	Bulk Coral 200A 30µm CN 100g
B-200A-30-CO-CN-1KG	Bulk Coral 200A 30µm CN 1kg
B-200A-30-CO-Phenyl-100G	Bulk Coral 200A 30µm Phenyl 100g
B-200A-30-CO-Phenyl-1KG	Bulk Coral 200A 30µm Phenyl 1kg
B-200A-30-CO-Phenyl HE-100G	Bulk Coral 200A 30µm Phenyl HE 100g
B-200A-30-CO-Phenyl HE-1KG	Bulk Coral 200A 30µm Phenyl HE 1kg
B-200A-30-CO-Diol-100G	Bulk Coral 200A 30µm Diol 100g
B-200A-30-CO-Diol-1KG	Bulk Coral 200A 30µm Diol 1kg
<b>Particule Coral 30, Pore 300A</b>	
B-300A-30-CO-Si-100G	Bulk Coral 300A 30µm 100g
B-300A-30-CO-Si-1KG	Bulk Coral 300A 30µm 1kg
<b>C8</b>	
B-300A-30-CO-C8-100G	Bulk Coral 300A 30µm C8 100g
B-300A-30-CO-C8-1KG	Bulk Coral 300A 30µm C8 1kg
B-300A-30-CO-C8W-100G	Bulk Coral 300A 30µm C8 Wax 100g
B-300A-30-CO-C8W-1KG	Bulk Coral 300A 30µm C8 Wax 1kg
B-300A-30-CO-C8P-100G	Bulk Coral 300A 30µm C8 Polar 100g
B-300A-30-CO-C8P-1KG	Bulk Coral 300A 30µm C8 Polar 1kg
B-300A-30-CO-C8T-100G	Bulk Coral 300A 30µm C8 TriF 100g
B-300A-30-CO-C8T-1KG	Bulk Coral 300A 30µm C8 TriF 1kg
B-300A-30-CO-C8EX-100G	Bulk Coral 300A 30µm C8 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-30-CO-C8EX-1KG	Bulk Coral 300A 30µm C8 Extreme 1kg
B-300A-30-CO-C8PT-100G	Bulk Coral 300A 30µm C8 Polar TriF 100g
B-300A-30-CO-C8PT-1KG	Bulk Coral 300A 30µm C8 Polar TriF 1kg
B-300A-30-CO-C8AQP-100G	Bulk Coral 300A 30µm C8 AQ Plus 100g
B-300A-30-CO-C8AQP-1KG	Bulk Coral 300A 30µm C8 AQ Plus 1kg
B-300A-30-CO-C8AQ-100G	Bulk Coral 300A 30µm C8 AQ 100g
B-300A-30-CO-C8AQ-1KG	Bulk Coral 300A 30µm C8 AQ 1kg
<b>C18</b>	
B-300A-30-CO-C18-100G	Bulk Coral 300A 30µm C18 100g
B-300A-30-CO-C18-1KG	Bulk Coral 300A 30µm C18 1kg
B-300A-30-CO-C18W-100G	Bulk Coral 300A 30µm C18 Wax 100g
B-300A-30-CO-C18W-1KG	Bulk Coral 300A 30µm C18 Wax 1kg
B-300A-30-CO-C18P-100G	Bulk Coral 300A 30µm C18 Polar 100g
B-300A-30-CO-C18P-1KG	Bulk Coral 300A 30µm C18 Polar 1kg
B-300A-30-CO-C18T-100G	Bulk Coral 300A 30µm C18 TriF 100g
B-300A-30-CO-C18T-1KG	Bulk Coral 300A 30µm C18 TriF 1kg
B-300A-30-CO-C18EX-100G	Bulk Coral 300A 30µm C18 Extreme 100g
B-300A-30-CO-C18EX-1KG	Bulk Coral 300A 30µm C18 Extreme 1kg
B-300A-30-CO-C18PT-100G	Bulk Coral 300A 30µm C18 Polar TriF 100g
B-300A-30-CO-C18PT-1KG	Bulk Coral 300A 30µm C18 Polar TriF 1kg
B-300A-30-CO-C18AQP-100G	Bulk Coral 300A 30µm C18 AQ Plus 100g
B-300A-30-CO-C18AQP-1KG	Bulk Coral 300A 30µm C18 AQ Plus 1kg
B-300A-30-CO-C18AQ-100G	Bulk Coral 300A 30µm C18 AQ 100g
B-300A-30-CO-C18AQ-1KG	Bulk Coral 300A 30µm C18 AQ 1kg
<b>C4</b>	
B-300A-30-CO-C4-100G	Bulk Coral 300A 30µm C4 100g
B-300A-30-CO-C4-1KG	Bulk Coral 300A 30µm C4 1kg
B-300A-30-CO-C4W-100G	Bulk Coral 300A 30µm C4 Wax 100g
B-300A-30-CO-C4W-1KG	Bulk Coral 300A 30µm C4 Wax 1kg
B-300A-30-CO-C4P-100G	Bulk Coral 300A 30µm C4 Polar 100g
B-300A-30-CO-C4P-1KG	Bulk Coral 300A 30µm C4 Polar 1kg
B-300A-30-CO-C4T-100G	Bulk Coral 300A 30µm C4 TriF 100g
B-300A-30-CO-C4T-1KG	Bulk Coral 300A 30µm C4 TriF 1kg
B-300A-30-CO-C4EX-100G	Bulk Coral 300A 30µm C4 Extreme 100g
B-300A-30-CO-C4EX-1KG	Bulk Coral 300A 30µm C4 Extreme 1kg
B-300A-30-CO-C4PT-100G	Bulk Coral 300A 30µm C4 Polar TriF 100g
B-300A-30-CO-C4PT-1KG	Bulk Coral 300A 30µm C4 Polar TriF 1kg
B-300A-30-CO-C4AQP-100G	Bulk Coral 300A 30µm C4 AQ Plus 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300A-30-CO-C4AQP-1KG	Bulk Coral 300A 30µm C4 AQ Plus 1kg
B-300A-30-CO-C4AQ-100G	Bulk Coral 300A 30µm C4 AQ 100g
B-300A-30-CO-C4AQ-1KG	Bulk Coral 300A 30µm C4 AQ 1kg
<b>Autres</b>	
B-300A-30-CO-COOH-100G	Bulk Coral 300A 30µm COOH 100g
B-300A-30-CO-COOH-1KG	Bulk Coral 300A 30µm COOH 1kg
B-300A-30-CO-SO3H-100G	Bulk Coral 300A 30µm SO3H 100g
B-300A-30-CO-SO3H-1KG	Bulk Coral 300A 30µm SO3H 1kg
B-300A-30-CO-DNH-100G	Bulk Coral 300A 30µm DNH 100g
B-300A-30-CO-DNH-1KG	Bulk Coral 300A 30µm DNH 1kg
B-300A-30-CO-PEI-100G	Bulk Coral 300A 30µm PEI 100g
B-300A-30-CO-PEI-1KG	Bulk Coral 300A 30µm PEI 1kg
B-300A-30-CO-ARG-100G	Bulk Coral 300A 30µm ARG 100g
B-300A-30-CO-ARG-1KG	Bulk Coral 300A 30µm ARG 1kg
B-300A-30-CO-NH-100G	Bulk Coral 300A 30µm NH 100g
B-300A-30-CO-NH-1KG	Bulk Coral 300A 30µm NH 1kg
B-300A-30-CO-CO-NH2-100G	Bulk Coral 300A 30µm NH2 100g
B-300A-30-CO-CO-NH2-1KG	Bulk Coral 300A 30µm NH2 1kg
B-300A-30-CO-CN-100G	Bulk Coral 300A 30µm CN 100g
B-300A-30-CO-CN-1KG	Bulk Coral 300A 30µm CN 1kg
B-300A-30-CO-Phenyl-100G	Bulk Coral 300A 30µm Phenyl 100g
B-300A-30-CO-Phenyl-1KG	Bulk Coral 300A 30µm Phenyl 1kg
B-300A-30-CO-Phenyl HE-100G	Bulk Coral 300A 30µm Phenyl HE 100g
B-300A-30-CO-Phenyl HE-1KG	Bulk Coral 300A 30µm Phenyl HE 1kg
B-300A-30-CO-Diol-100G	Bulk Coral 300A 30µm Diol 100g
B-300A-30-CO-Diol-1KG	Bulk Coral 300A 30µm Diol 1kg
<b>Particule Coral 30, Pore 500A</b>	
B-500A-30-CO-Si-100G	Bulk Coral 500A 30µm 100g
B-500A-30-CO-Si-1KG	Bulk Coral 500A 30µm 1kg
<b>C8</b>	
B-500A-30-CO-C8-100G	Bulk Coral 500A 30µm C8 100g
B-500A-30-CO-C8-1KG	Bulk Coral 500A 30µm C8 1kg
B-500A-30-CO-C8W-100G	Bulk Coral 500A 30µm C8 Wax 100g
B-500A-30-CO-C8W-1KG	Bulk Coral 500A 30µm C8 Wax 1kg
B-500A-30-CO-C8P-100G	Bulk Coral 500A 30µm C8 Polar 100g
B-500A-30-CO-C8P-1KG	Bulk Coral 500A 30µm C8 Polar 1kg
B-500A-30-CO-C8T-100G	Bulk Coral 500A 30µm C8 TriF 100g
B-500A-30-CO-C8T-1KG	Bulk Coral 500A 30µm C8 TriF 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-500A-30-CO-C8EX-100G	Bulk Coral 500A 30µm C8 Extreme 100g
B-500A-30-CO-C8EX-1KG	Bulk Coral 500A 30µm C8 Extreme 1kg
B-500A-30-CO-C8PT-100G	Bulk Coral 500A 30µm C8 Polar TriF 100g
B-500A-30-CO-C8PT-1KG	Bulk Coral 500A 30µm C8 Polar TriF 1kg
B-500A-30-CO-C8AQP-100G	Bulk Coral 500A 30µm C8 AQ Plus 100g
B-500A-30-CO-C8AQP-1KG	Bulk Coral 500A 30µm C8 AQ Plus 1kg
B-500A-30-CO-C8AQ-100G	Bulk Coral 500A 30µm C8 AQ 100g
B-500A-30-CO-C8AQ-1KG	Bulk Coral 500A 30µm C8 AQ 1kg
<b>C18</b>	
B-500A-30-CO-C18-100G	Bulk Coral 500A 30µm C18 100g
B-500A-30-CO-C18-1KG	Bulk Coral 500A 30µm C18 1kg
B-500A-30-CO-C18W-100G	Bulk Coral 500A 30µm C18 Wax 100g
B-500A-30-CO-C18W-1KG	Bulk Coral 500A 30µm C18 Wax 1kg
B-500A-30-CO-C18P-100G	Bulk Coral 500A 30µm C18 Polar 100g
B-500A-30-CO-C18P-1KG	Bulk Coral 500A 30µm C18 Polar 1kg
B-500A-30-CO-C18T-100G	Bulk Coral 500A 30µm C18 TriF 100g
B-500A-30-CO-C18T-1KG	Bulk Coral 500A 30µm C18 TriF 1kg
B-500A-30-CO-C18EX-100G	Bulk Coral 500A 30µm C18 Extreme 100g
B-500A-30-CO-C18EX-1KG	Bulk Coral 500A 30µm C18 Extreme 1kg
B-500A-30-CO-C18PT-100G	Bulk Coral 500A 30µm C18 Polar TriF 100g
B-500A-30-CO-C18PT-1KG	Bulk Coral 500A 30µm C18 Polar TriF 1kg
B-500A-30-CO-C18AQP-100G	Bulk Coral 500A 30µm C18 AQ Plus 100g
B-500A-30-CO-C18AQP-1KG	Bulk Coral 500A 30µm C18 AQ Plus 1kg
B-500A-30-CO-C18AQ-100G	Bulk Coral 500A 30µm C18 AQ 100g
B-500A-30-CO-C18AQ-1KG	Bulk Coral 500A 30µm C18 AQ 1kg
<b>C4</b>	
B-500A-30-CO-C4-100G	Bulk Coral 500A 30µm C4 100g
B-500A-30-CO-C4-1KG	Bulk Coral 500A 30µm C4 1kg
B-500A-30-CO-C4W-100G	Bulk Coral 500A 30µm C4 Wax 100g
B-500A-30-CO-C4W-1KG	Bulk Coral 500A 30µm C4 Wax 1kg
B-500A-30-CO-C4P-100G	Bulk Coral 500A 30µm C4 Polar 100g
B-500A-30-CO-C4P-1KG	Bulk Coral 500A 30µm C4 Polar 1kg
B-500A-30-CO-C4T-100G	Bulk Coral 500A 30µm C4 TriF 100g
B-500A-30-CO-C4T-1KG	Bulk Coral 500A 30µm C4 TriF 1kg
B-500A-30-CO-C4EX-100G	Bulk Coral 500A 30µm C4 Extreme 100g
B-500A-30-CO-C4EX-1KG	Bulk Coral 500A 30µm C4 Extreme 1kg
B-500A-30-CO-C4PT-100G	Bulk Coral 500A 30µm C4 Polar TriF 100g
B-500A-30-CO-C4PT-1KG	Bulk Coral 500A 30µm C4 Polar TriF 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-500A-30-CO-C4AQP-100G	Bulk Coral 500A 30µm C4 AQ Plus 100g
B-500A-30-CO-C4AQP-1KG	Bulk Coral 500A 30µm C4 AQ Plus 1kg
B-500A-30-CO-C4AQ-100G	Bulk Coral 500A 30µm C4 AQ 100g
B-500A-30-CO-C4AQ-1KG	Bulk Coral 500A 30µm C4 AQ 1kg
<b>Autres</b>	
B-500A-30-CO-COOH-100G	Bulk Coral 500A 30µm COOH 100g
B-500A-30-CO-COOH-1KG	Bulk Coral 500A 30µm COOH 1kg
B-500A-30-CO-SO3H-100G	Bulk Coral 500A 30µm SO3H 100g
B-500A-30-CO-SO3H-1KG	Bulk Coral 500A 30µm SO3H 1kg
B-500A-30-CO-DNH-100G	Bulk Coral 500A 30µm DNH 100g
B-500A-30-CO-DNH-1KG	Bulk Coral 500A 30µm DNH 1kg
B-500A-30-CO-PEI-100G	Bulk Coral 500A 30µm PEI 100g
B-500A-30-CO-PEI-1KG	Bulk Coral 500A 30µm PEI 1kg
B-500A-30-CO-ARG-100G	Bulk Coral 500A 30µm ARG 100g
B-500A-30-CO-ARG-1KG	Bulk Coral 500A 30µm ARG 1kg
B-500A-30-CO-NH-100G	Bulk Coral 500A 30µm NH 100g
B-500A-30-CO-NH-1KG	Bulk Coral 500A 30µm NH 1kg
B-500A-30-CO-CO-NH2-100G	Bulk Coral 500A 30µm NH2 100g
B-500A-30-CO-CO-NH2-1KG	Bulk Coral 500A 30µm NH2 1kg
B-500A-30-CO-CN-100G	Bulk Coral 500A 30µm CN 100g
B-500A-30-CO-CN-1KG	Bulk Coral 500A 30µm CN 1kg
B-500A-30-CO-Phenyl-100G	Bulk Coral 500A 30µm Phenyl 100g
B-500A-30-CO-Phenyl-1KG	Bulk Coral 500A 30µm Phenyl 1kg
B-500A-30-CO-Phenyl HE-100G	Bulk Coral 500A 30µm Phenyl HE 100g
B-500A-30-CO-Phenyl HE-1KG	Bulk Coral 500A 30µm Phenyl HE 1kg
B-500A-30-CO-Diol-100G	Bulk Coral 500A 30µm Diol 100g
B-500A-30-CO-Diol-1KG	Bulk Coral 500A 30µm Diol 1kg
<b>Particule Coral 30, Pore 800A</b>	
B-800A-30-CO-Si-100G	Bulk Coral 800A 30µm 100g
B-800A-30-CO-Si-1KG	Bulk Coral 800A 30µm 1kg
<b>C8</b>	
B-800A-30-CO-C8-100G	Bulk Coral 800A 30µm C8 100g
B-800A-30-CO-C8-1KG	Bulk Coral 800A 30µm C8 1kg
B-800A-30-CO-C8W-100G	Bulk Coral 800A 30µm C8 Wax 100g
B-800A-30-CO-C8W-1KG	Bulk Coral 800A 30µm C8 Wax 1kg
B-800A-30-CO-C8P-100G	Bulk Coral 800A 30µm C8 Polar 100g
B-800A-30-CO-C8P-1KG	Bulk Coral 800A 30µm C8 Polar 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-30-CO-C8P-1KG	Bulk Coral 800A 30µm C8 Polar 1kg
B-800A-30-CO-C8T-100G	Bulk Coral 800A 30µm C8 TriF 100g
B-800A-30-CO-C8T-1KG	Bulk Coral 800A 30µm C8 TriF 1kg
B-800A-30-CO-C8EX-100G	Bulk Coral 800A 30µm C8 Extreme 100g
B-800A-30-CO-C8EX-1KG	Bulk Coral 800A 30µm C8 Extreme 1kg
B-800A-30-CO-C8PT-100G	Bulk Coral 800A 30µm C8 Polar TriF 100g
B-800A-30-CO-C8PT-1KG	Bulk Coral 800A 30µm C8 Polar TriF 1kg
B-800A-30-CO-C8AQP-100G	Bulk Coral 800A 30µm C8 AQ Plus 100g
B-800A-30-CO-C8AQP-1KG	Bulk Coral 800A 30µm C8 AQ Plus 1kg
B-800A-30-CO-C8AQ-100G	Bulk Coral 800A 30µm C8 AQ 100g
B-800A-30-CO-C8AQ-1KG	Bulk Coral 800A 30µm C8 AQ 1kg
<b>C18</b>	
B-800A-30-CO-C18-100G	Bulk Coral 800A 30µm C18 100g
B-800A-30-CO-C18-1KG	Bulk Coral 800A 30µm C18 1kg
B-800A-30-CO-C18W-100G	Bulk Coral 800A 30µm C18 Wax 100g
B-800A-30-CO-C18W-1KG	Bulk Coral 800A 30µm C18 Wax 1kg
B-800A-30-CO-C18P-100G	Bulk Coral 800A 30µm C18 Polar 100g
B-800A-30-CO-C18P-1KG	Bulk Coral 800A 30µm C18 Polar 1kg
B-800A-30-CO-C18T-100G	Bulk Coral 800A 30µm C18 TriF 100g
B-800A-30-CO-C18T-1KG	Bulk Coral 800A 30µm C18 TriF 1kg
B-800A-30-CO-C18EX-100G	Bulk Coral 800A 30µm C18 Extreme 100g
B-800A-30-CO-C18EX-1KG	Bulk Coral 800A 30µm C18 Extreme 1kg
B-800A-30-CO-C18PT-100G	Bulk Coral 800A 30µm C18 Polar TriF 100g
B-800A-30-CO-C18PT-1KG	Bulk Coral 800A 30µm C18 Polar TriF 1kg
B-800A-30-CO-C18AQP-100G	Bulk Coral 800A 30µm C18 AQ Plus 100g
B-800A-30-CO-C18AQP-1KG	Bulk Coral 800A 30µm C18 AQ Plus 1kg
B-800A-30-CO-C18AQ-100G	Bulk Coral 800A 30µm C18 AQ 100g
B-800A-30-CO-C18AQ-1KG	Bulk Coral 800A 30µm C18 AQ 1kg
<b>C4</b>	
B-800A-30-CO-C4-100G	Bulk Coral 800A 30µm C4 100g
B-800A-30-CO-C4-1KG	Bulk Coral 800A 30µm C4 1kg
B-800A-30-CO-C4W-100G	Bulk Coral 800A 30µm C4 Wax 100g
B-800A-30-CO-C4W-1KG	Bulk Coral 800A 30µm C4 Wax 1kg
B-800A-30-CO-C4P-100G	Bulk Coral 800A 30µm C4 Polar 100g
B-800A-30-CO-C4P-1KG	Bulk Coral 800A 30µm C4 Polar 1kg
B-800A-30-CO-C4T-100G	Bulk Coral 800A 30µm C4 TriF 100g
B-800A-30-CO-C4T-1KG	Bulk Coral 800A 30µm C4 TriF 1kg
B-800A-30-CO-C4EX-100G	Bulk Coral 800A 30µm C4 Extreme 100g



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800A-30-CO-C4EX-1KG	Bulk Coral 800A 30µm C4 Extreme 1kg
B-800A-30-CO-C4PT-100G	Bulk Coral 800A 30µm C4 Polar TriF 100g
B-800A-30-CO-C4PT-1KG	Bulk Coral 800A 30µm C4 Polar TriF 1kg
B-800A-30-CO-C4AQP-100G	Bulk Coral 800A 30µm C4 AQ Plus 100g
B-800A-30-CO-C4AQP-1KG	Bulk Coral 800A 30µm C4 AQ Plus 1kg
B-800A-30-CO-C4AQ-100G	Bulk Coral 800A 30µm C4 AQ 100g
B-800A-30-CO-C4AQ-1KG	Bulk Coral 800A 30µm C4 AQ 1kg
<b>Autres</b>	
B-800A-30-CO-COOH-100G	Bulk Coral 800A 30µm COOH 100g
B-800A-30-CO-COOH-1KG	Bulk Coral 800A 30µm COOH 1kg
B-800A-30-CO-S03H-100G	Bulk Coral 800A 30µm S03H 100g
B-800A-30-CO-S03H-1KG	Bulk Coral 800A 30µm S03H 1kg
B-800A-30-CO-DNH-100G	Bulk Coral 800A 30µm DNH 100g
B-800A-30-CO-DNH-1KG	Bulk Coral 800A 30µm DNH 1kg
B-800A-30-CO-PEI-100G	Bulk Coral 800A 30µm PEI 100g
B-800A-30-CO-PEI-1KG	Bulk Coral 800A 30µm PEI 1kg
B-800A-30-CO-ARG-100G	Bulk Coral 800A 30µm ARG 100g
B-800A-30-CO-ARG-1KG	Bulk Coral 800A 30µm ARG 1kg
B-800A-30-CO-NH-100G	Bulk Coral 800A 30µm NH 100g
B-800A-30-CO-NH-1KG	Bulk Coral 800A 30µm NH 1kg
B-800A-30-CO-CO-NH2-100G	Bulk Coral 800A 30µm NH2 100g
B-800A-30-CO-CO-NH2-1KG	Bulk Coral 800A 30µm NH2 1kg
B-800A-30-CO-CN-100G	Bulk Coral 800A 30µm CN 100g
B-800A-30-CO-CN-1KG	Bulk Coral 800A 30µm CN 1kg
B-800A-30-CO-Phenyl-100G	Bulk Coral 800A 30µm Phenyl 100g
B-800A-30-CO-Phenyl-1KG	Bulk Coral 800A 30µm Phenyl 1kg
B-800A-30-CO-Phenyl HE-100G	Bulk Coral 800A 30µm Phenyl HE 100g
B-800A-30-CO-Phenyl HE-1KG	Bulk Coral 800A 30µm Phenyl HE 1kg
B-800A-30-CO-Diol-100G	Bulk Coral 800A 30µm Diol 100g
B-800A-30-CO-Diol-1KG	Bulk Coral 800A 30µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Coral 30, Pore 1000A</b>	
B-1000A-30-CO-Si-100G	Bulk Coral 1000A 30µm 100g
B-1000A-30-CO-Si-1KG	Bulk Coral 1000A 30µm 1kg
<b>C8</b>	
B-1000A-30-CO-C8-100G	Bulk Coral 1000A 30µm C8 100g
B-1000A-30-CO-C8-1KG	Bulk Coral 1000A 30µm C8 1kg
B-1000A-30-CO-C8W-100G	Bulk Coral 1000A 30µm C8 Wax 100g
B-1000A-30-CO-C8W-1KG	Bulk Coral 1000A 30µm C8 Wax 1kg
B-1000A-30-CO-C8P-100G	Bulk Coral 1000A 30µm C8 Polar 100g
B-1000A-30-CO-C8P-1KG	Bulk Coral 1000A 30µm C8 Polar 1kg
B-1000A-30-CO-C8T-100G	Bulk Coral 1000A 30µm C8 TriF 100g
B-1000A-30-CO-C8T-1KG	Bulk Coral 1000A 30µm C8 TriF 1kg
B-1000A-30-CO-C8EX-100G	Bulk Coral 1000A 30µm C8 Extreme 100g
B-1000A-30-CO-C8EX-1KG	Bulk Coral 1000A 30µm C8 Extreme 1kg
B-1000A-30-CO-C8PT-100G	Bulk Coral 1000A 30µm C8 Polar TriF 100g
B-1000A-30-CO-C8PT-1KG	Bulk Coral 1000A 30µm C8 Polar TriF 1kg
B-1000A-30-CO-C8AQP-100G	Bulk Coral 1000A 30µm C8 AQ Plus 100g
B-1000A-30-CO-C8AQP-1KG	Bulk Coral 1000A 30µm C8 AQ Plus 1kg
B-1000A-30-CO-C8AQ-100G	Bulk Coral 1000A 30µm C8 AQ 100g
B-1000A-30-CO-C8AQ-1KG	Bulk Coral 1000A 30µm C8 AQ 1kg
<b>C18</b>	
B-1000A-30-CO-C18-100G	Bulk Coral 1000A 30µm C18 100g
B-1000A-30-CO-C18-1KG	Bulk Coral 1000A 30µm C18 1kg
B-1000A-30-CO-C18W-100G	Bulk Coral 1000A 30µm C18 Wax 100g
B-1000A-30-CO-C18W-1KG	Bulk Coral 1000A 30µm C18 Wax 1kg
B-1000A-30-CO-C18P-100G	Bulk Coral 1000A 30µm C18 Polar 100g
B-1000A-30-CO-C18P-1KG	Bulk Coral 1000A 30µm C18 Polar 1kg
B-1000A-30-CO-C18T-100G	Bulk Coral 1000A 30µm C18 TriF 100g
B-1000A-30-CO-C18T-1KG	Bulk Coral 1000A 30µm C18 TriF 1kg
B-1000A-30-CO-C18EX-100G	Bulk Coral 1000A 30µm C18 Extreme 100g
B-1000A-30-CO-C18EX-1KG	Bulk Coral 1000A 30µm C18 Extreme 1kg
B-1000A-30-CO-C18PT-100G	Bulk Coral 1000A 30µm C18 Polar TriF 100g
B-1000A-30-CO-C18PT-1KG	Bulk Coral 1000A 30µm C18 Polar TriF 1kg
B-1000A-30-CO-C18AQP-100G	Bulk Coral 1000A 30µm C18 AQ Plus 100g
B-1000A-30-CO-C18AQP-1KG	Bulk Coral 1000A 30µm C18 AQ Plus 1kg
B-1000A-30-CO-C18AQ-100G	Bulk Coral 1000A 30µm C18 AQ 100g
B-1000A-30-CO-C18AQ-1KG	Bulk Coral 1000A 30µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-1000A-30-CO-C4-100G	Bulk Coral 1000A 30µm C4 100g
B-1000A-30-CO-C4-1KG	Bulk Coral 1000A 30µm C4 1kg
B-1000A-30-CO-C4W-100G	Bulk Coral 1000A 30µm C4 Wax 100g
B-1000A-30-CO-C4W-1KG	Bulk Coral 1000A 30µm C4 Wax 1kg
B-1000A-30-CO-C4P-100G	Bulk Coral 1000A 30µm C4 Polar 100g
B-1000A-30-CO-C4P-1KG	Bulk Coral 1000A 30µm C4 Polar 1kg
B-1000A-30-CO-C4T-100G	Bulk Coral 1000A 30µm C4 TriF 100g
B-1000A-30-CO-C4T-1KG	Bulk Coral 1000A 30µm C4 TriF 1kg
B-1000A-30-CO-C4EX-100G	Bulk Coral 1000A 30µm C4 Extreme 100g
B-1000A-30-CO-C4EX-1KG	Bulk Coral 1000A 30µm C4 Extreme 1kg
B-1000A-30-CO-C4PT-100G	Bulk Coral 1000A 30µm C4 Polar TriF 100g
B-1000A-30-CO-C4PT-1KG	Bulk Coral 1000A 30µm C4 Polar TriF 1kg
B-1000A-30-CO-C4AQP-100G	Bulk Coral 1000A 30µm C4 AQ Plus 100g
B-1000A-30-CO-C4AQP-1KG	Bulk Coral 1000A 30µm C4 AQ Plus 1kg
B-1000A-30-CO-C4AQ-100G	Bulk Coral 1000A 30µm C4 AQ 100g
B-1000A-30-CO-C4AQ-1KG	Bulk Coral 1000A 30µm C4 AQ 1kg
<b>Autres</b>	
B-1000A-30-CO-COOH-100G	Bulk Coral 1000A 30µm COOH 100g
B-1000A-30-CO-COOH-1KG	Bulk Coral 1000A 30µm COOH 1kg
B-1000A-30-CO-SO3H-100G	Bulk Coral 1000A 30µm SO3H 100g
B-1000A-30-CO-SO3H-1KG	Bulk Coral 1000A 30µm SO3H 1kg
B-1000A-30-CO-DNH-100G	Bulk Coral 1000A 30µm DNH 100g
B-1000A-30-CO-DNH-1KG	Bulk Coral 1000A 30µm DNH 1kg
B-1000A-30-CO-PEI-100G	Bulk Coral 1000A 30µm PEI 100g
B-1000A-30-CO-PEI-1KG	Bulk Coral 1000A 30µm PEI 1kg
B-1000A-30-CO-ARG-100G	Bulk Coral 1000A 30µm ARG 100g
B-1000A-30-CO-ARG-1KG	Bulk Coral 1000A 30µm ARG 1kg
B-1000A-30-CO-NH-100G	Bulk Coral 1000A 30µm NH 100g
B-1000A-30-CO-NH-1KG	Bulk Coral 1000A 30µm NH 1kg
B-1000A-30-CO-CO-NH2-100G	Bulk Coral 1000A 30µm NH2 100g
B-1000A-30-CO-CO-NH2-1KG	Bulk Coral 1000A 30µm NH2 1kg
B-1000A-30-CO-CN-100G	Bulk Coral 1000A 30µm CN 100g
B-1000A-30-CO-CN-1KG	Bulk Coral 1000A 30µm CN 1kg
B-1000A-30-CO-Phenyl-100G	Bulk Coral 1000A 30µm Phenyl 100g
B-1000A-30-CO-Phenyl-1KG	Bulk Coral 1000A 30µm Phenyl 1kg
B-1000A-30-CO-Phenyl HE-100G	Bulk Coral 1000A 30µm Phenyl HE 100g
B-1000A-30-CO-Phenyl HE-1KG	Bulk Coral 1000A 30µm Phenyl HE 1kg
B-1000A-30-CO-Diol-100G	Bulk Coral 1000A 30µm Diol 100g
B-1000A-30-CO-Diol-1KG	Bulk Coral 1000A 30µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Coral 60µm</b>	
<b>Particule Coral 60, Pore 55A</b>	
B-55-60-CO-Si-100G	Bulk Coral 55A 60µm 100g
B-55-60-CO-Si-1KG	Bulk Coral 55A 60µm 1kg
<b>C8</b>	
B-55-60-CO-C8T-100G	Bulk Coral 55A 60µm C8 TriF 100g
B-55-60-CO-C8T-1KG	Bulk Coral 55A 60µm C8 TriF 1kg
B-55-60-CO-C8EX-100G	Bulk Coral 55A 60µm C8 Extreme 100g
B-55-60-CO-C8EX-1KG	Bulk Coral 55A 60µm C8 Extreme 1kg
B-55-60-CO-C8PT-100G	Bulk Coral 55A 60µm C8 Polar TriF 100g
B-55-60-CO-C8PT-1KG	Bulk Coral 55A 60µm C8 Polar TriF 1kg
B-55-60-CO-C8AQP-100G	Bulk Coral 55A 60µm C8 AQ Plus 100g
B-55-60-CO-C8AQP-1KG	Bulk Coral 55A 60µm C8 AQ Plus 1kg
B-55-60-CO-C8AQ-100G	Bulk Coral 55A 60µm C8 AQ 100g
B-55-60-CO-C8AQ-1KG	Bulk Coral 55A 60µm C8 AQ 1kg
<b>C18</b>	
B-55-60-CO-C18T-100G	Bulk Coral 55A 60µm C18 TriF 100g
B-55-60-CO-C18T-1KG	Bulk Coral 55A 60µm C18 TriF 1kg
B-55-60-CO-C18EX-100G	Bulk Coral 55A 60µm C18 Extreme 100g
B-55-60-CO-C18EX-1KG	Bulk Coral 55A 60µm C18 Extreme 1kg
B-55-60-CO-C18PT-100G	Bulk Coral 55A 60µm C18 Polar TriF 100g
B-55-60-CO-C18PT-1KG	Bulk Coral 55A 60µm C18 Polar TriF 1kg
B-55-60-CO-C18AQP-100G	Bulk Coral 55A 60µm C18 AQ Plus 100g
B-55-60-CO-C18AQP-1KG	Bulk Coral 55A 60µm C18 AQ Plus 1kg
B-55-60-CO-C18AQ-100G	Bulk Coral 55A 60µm C18 AQ 100g
B-55-60-CO-C18AQ-1KG	Bulk Coral 55A 60µm C18 AQ 1kg
<b>C4</b>	
B-55-60-CO-C4T-100G	Bulk Coral 55A 60µm C4 TriF 100g
B-55-60-CO-C4T-1KG	Bulk Coral 55A 60µm C4 TriF 1kg
B-55-60-CO-C4EX-100G	Bulk Coral 55A 60µm C4 Extreme 100g
B-55-60-CO-C4EX-1KG	Bulk Coral 55A 60µm C4 Extreme 1kg
B-55-60-CO-C4PT-100G	Bulk Coral 55A 60µm C4 Polar TriF 100g
B-55-60-CO-C4PT-1KG	Bulk Coral 55A 60µm C4 Polar TriF 1kg
B-55-60-CO-C4AQP-100G	Bulk Coral 55A 60µm C4 AQ Plus 100g
B-55-60-CO-C4AQP-1KG	Bulk Coral 55A 60µm C4 AQ Plus 1kg
B-55-60-CO-C4AQ-100G	Bulk Coral 55A 60µm C4 AQ 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-55-60-CO-C4AQ-1KG	Bulk Coral 55A 60µm C4 AQ 1kg
<b>Autres</b>	
B-55-60-CO-COOH-100G	Bulk Coral 55A 60µm COOH 100g
B-55-60-CO-COOH-1KG	Bulk Coral 55A 60µm COOH 1kg
B-55-60-CO-SO3H-100G	Bulk Coral 55A 60µm SO3H 100g
B-55-60-CO-SO3H-1KG	Bulk Coral 55A 60µm SO3H 1kg
B-55-60-CO-DNH-100G	Bulk Coral 55A 60µm DNH 100g
B-55-60-CO-DNH-1KG	Bulk Coral 55A 60µm DNH 1kg
B-55-60-CO-PEI-100G	Bulk Coral 55A 60µm PEI 100g
B-55-60-CO-PEI-1KG	Bulk Coral 55A 60µm PEI 1kg
B-55-60-CO-ARG-100G	Bulk Coral 55A 60µm ARG 100g
B-55-60-CO-ARG-1KG	Bulk Coral 55A 60µm ARG 1kg
B-55-60-CO-NH-100G	Bulk Coral 55A 60µm NH 100g
B-55-60-CO-NH-1KG	Bulk Coral 55A 60µm NH 1kg
B-55-60-CO-NH2-100G	Bulk Coral 55A 60µm NH2 100g
B-55-60-CO-NH2-1KG	Bulk Coral 55A 60µm NH2 1kg
B-55-60-CO-CO-CN-100G	Bulk Coral 55A 60µm CN 100g
B-55-60-CO-CO-CN-1KG	Bulk Coral 55A 60µm CN 1kg
B-55-60-CO-Phenyl-100G	Bulk Coral 55A 60µm Phenyl 100g
B-55-60-CO-Phenyl-1KG	Bulk Coral 55A 60µm Phenyl 1kg
B-55-60-CO-Diol-100G	Bulk Coral 55A 60µm Diol 100g
B-55-60-CO-Diol-1KG	Bulk Coral 55A 60µm Diol 1kg
<b>Particule Coral 60, Pore 70A</b>	
B-70-60-100G	Bulk Coral 70A 60µm 100g
B-70-60-1KG	Bulk Coral 70A 60µm 1kg
<b>C8</b>	
B-70-60-CO-C8T-100G	Bulk Coral 70A 60µm C8 TriF 100g
B-70-60-CO-C8T-1KG	Bulk Coral 70A 60µm C8 TriF 1kg
B-70-60-CO-C8EX-100G	Bulk Coral 70A 60µm C8 Extreme 100g
B-70-60-CO-C8EX-1KG	Bulk Coral 70A 60µm C8 Extreme 1kg
B-70-60-CO-C8PT-100G	Bulk Coral 70A 60µm C8 Polar TriF 100g
B-70-60-CO-C8PT-1KG	Bulk Coral 70A 60µm C8 Polar TriF 1kg
B-70-60-CO-C8AQP-100G	Bulk Coral 70A 60µm C8 AQ Plus 100g
B-70-60-CO-C8AQP-1KG	Bulk Coral 70A 60µm C8 AQ Plus 1kg
B-70-60-CO-C8AQ-100G	Bulk Coral 70A 60µm C8 AQ 100g
B-70-60-CO-C8AQ-1KG	Bulk Coral 70A 60µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-70-60-CO-C18T-100G	Bulk Coral 70A 60µm C18 TriF 100g
B-70-60-CO-C18T-1KG	Bulk Coral 70A 60µm C18 TriF 1kg
B-70-60-CO-C18EX-100G	Bulk Coral 70A 60µm C18 Extreme 100g
B-70-60-CO-C18EX-1KG	Bulk Coral 70A 60µm C18 Extreme 1kg
B-70-60-CO-C18PT-100G	Bulk Coral 70A 60µm C18 Polar TriF 100g
B-70-60-CO-C18PT-1KG	Bulk Coral 70A 60µm C18 Polar TriF 1kg
B-70-60-CO-C18AQP-100G	Bulk Coral 70A 60µm C18 AQ Plus 100g
B-70-60-CO-C18AQP-1KG	Bulk Coral 70A 60µm C18 AQ Plus 1kg
B-70-60-CO-C18AQ-100G	Bulk Coral 70A 60µm C18 AQ 100g
B-70-60-CO-C18AQ-1KG	Bulk Coral 70A 60µm C18 AQ 1kg
<b>C4</b>	
B-70-60-CO-C4T-100G	Bulk Coral 70A 60µm C4 TriF 100g
B-70-60-CO-C4T-1KG	Bulk Coral 70A 60µm C4 TriF 1kg
B-70-60-CO-C4EX-100G	Bulk Coral 70A 60µm C4 Extreme 100g
B-70-60-CO-C4EX-1KG	Bulk Coral 70A 60µm C4 Extreme 1kg
B-70-60-CO-C4PT-100G	Bulk Coral 70A 60µm C4 Polar TriF 100g
B-70-60-CO-C4PT-1KG	Bulk Coral 70A 60µm C4 Polar TriF 1kg
B-70-60-CO-C4AQP-100G	Bulk Coral 70A 60µm C4 AQ Plus 100g
B-70-60-CO-C4AQP-1KG	Bulk Coral 70A 60µm C4 AQ Plus 1kg
B-70-60-CO-C4AQ-100G	Bulk Coral 70A 60µm C4 AQ 100g
B-70-60-CO-C4AQ-1KG	Bulk Coral 70A 60µm C4 AQ 1kg
<b>Autres</b>	
B-70-60-CO-COOH-100G	Bulk Coral 70A 60µm COOH 100g
B-70-60-CO-COOH-1KG	Bulk Coral 70A 60µm COOH 1kg
B-70-60-CO-SO3H-100G	Bulk Coral 70A 60µm SO3H 100g
B-70-60-CO-SO3H-1KG	Bulk Coral 70A 60µm SO3H 1kg
B-70-60-CO-DNH-100G	Bulk Coral 70A 60µm DNH 100g
B-70-60-CO-DNH-1KG	Bulk Coral 70A 60µm DNH 1kg
B-70-60-CO-PEI-100G	Bulk Coral 70A 60µm PEI 100g
B-70-60-CO-PEI-1KG	Bulk Coral 70A 60µm PEI 1kg
B-70-60-CO-ARG-100G	Bulk Coral 70A 60µm ARG 100g
B-70-60-CO-ARG-1KG	Bulk Coral 70A 60µm ARG 1kg
B-70-60-CO-NH-100G	Bulk Coral 70A 60µm NH 100g
B-70-60-CO-NH-1KG	Bulk Coral 70A 60µm NH 1kg
B-70-60-CO-NH2-100G	Bulk Coral 70A 60µm NH2 100g
B-70-60-CO-NH2-1KG	Bulk Coral 70A 60µm NH2 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70-60-CO-CO-CN-100G	Bulk Coral 70A 60µm CN 100g
B-70-60-CO-CO-CN-1KG	Bulk Coral 70A 60µm CN 1kg
B-70-60-CO-Phenyl-100G	Bulk Coral 70A 60µm Phenyl 100g
B-70-60-CO-Phenyl-1KG	Bulk Coral 70A 60µm Phenyl 1kg
B-70-60-CO-Diol-100G	Bulk Coral 70A 60µm Diol 100g
B-70-60-CO-Diol-1KG	Bulk Coral 70A 60µm Diol 1kg
<b>Particule Coral 60, Pore 100A</b>	
B-100-60-CO-Si-100G	Bulk Coral 100A 60µm 100g
B-100-60-CO-Si-1KG	Bulk Coral 100A 60µm 1kg
<b>C8</b>	
B-100-60-CO-C8T-100G	Bulk Coral 100A 60µm C8 TriF 100g
B-100-60-CO-C8T-1KG	Bulk Coral 100A 60µm C8 TriF 1kg
B-100-60-CO-C8EX-100G	Bulk Coral 100A 60µm C8 Extreme 100g
B-100-60-CO-C8EX-1KG	Bulk Coral 100A 60µm C8 Extreme 1kg
B-100-60-CO-C8PT-100G	Bulk Coral 100A 60µm C8 Polar TriF 100g
B-100-60-CO-C8PT-1KG	Bulk Coral 100A 60µm C8 Polar TriF 1kg
B-100-60-CO-C8AQP-100G	Bulk Coral 100A 60µm C8 AQ Plus 100g
B-100-60-CO-C8AQP-1KG	Bulk Coral 100A 60µm C8 AQ Plus 1kg
B-100-60-CO-C8AQ-100G	Bulk Coral 100A 60µm C8 AQ 100g
B-100-60-CO-C8AQ-1KG	Bulk Coral 100A 60µm C8 AQ 1kg
<b>C18</b>	
B-100-60-CO-C18T-100G	Bulk Coral 100A 60µm C18 TriF 100g
B-100-60-CO-C18T-1KG	Bulk Coral 100A 60µm C18 TriF 1kg
B-100-60-CO-C18EX-100G	Bulk Coral 100A 60µm C18 Extreme 100g
B-100-60-CO-C18EX-1KG	Bulk Coral 100A 60µm C18 Extreme 1kg
B-100-60-CO-C18PT-100G	Bulk Coral 100A 60µm C18 Polar TriF 100g
B-100-60-CO-C18PT-1KG	Bulk Coral 100A 60µm C18 Polar TriF 1kg
B-100-60-CO-C18AQP-100G	Bulk Coral 100A 60µm C18 AQ Plus 100g
B-100-60-CO-C18AQP-1KG	Bulk Coral 100A 60µm C18 AQ Plus 1kg
B-100-60-CO-C18AQ-100G	Bulk Coral 100A 60µm C18 AQ 100g
B-100-60-CO-C18AQ-1KG	Bulk Coral 100A 60µm C18 AQ 1kg
<b>C4</b>	
B-100-60-CO-C4T-100G	Bulk Coral 100A 60µm C4 TriF 100g
B-100-60-CO-C4T-1KG	Bulk Coral 100A 60µm C4 TriF 1kg
B-100-60-CO-C4EX-100G	Bulk Coral 100A 60µm C4 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-100-60-CO-C4EX-1KG	Bulk Coral 100A 60µm C4 Extreme 1kg
B-100-60-CO-C4PT-100G	Bulk Coral 100A 60µm C4 Polar TriF 100g
B-100-60-CO-C4PT-1KG	Bulk Coral 100A 60µm C4 Polar TriF 1kg
B-100-60-CO-C4AQP-100G	Bulk Coral 100A 60µm C4 AQ Plus 100g
B-100-60-CO-C4AQP-1KG	Bulk Coral 100A 60µm C4 AQ Plus 1kg
B-100-60-CO-C4AQ-100G	Bulk Coral 100A 60µm C4 AQ 100g
B-100-60-CO-C4AQ-1KG	Bulk Coral 100A 60µm C4 AQ 1kg
<b>Autres</b>	
B-100-60-CO-COOH-100G	Bulk Coral 100A 60µm COOH 100g
B-100-60-CO-COOH-1KG	Bulk Coral 100A 60µm COOH 1kg
B-100-60-CO-SO3H-100G	Bulk Coral 100A 60µm SO3H 100g
B-100-60-CO-SO3H-1KG	Bulk Coral 100A 60µm SO3H 1kg
B-100-60-CO-DNH-100G	Bulk Coral 100A 60µm DNH 100g
B-100-60-CO-DNH-1KG	Bulk Coral 100A 60µm DNH 1kg
B-100-60-CO-PEI-100G	Bulk Coral 100A 60µm PEI 100g
B-100-60-CO-PEI-1KG	Bulk Coral 100A 60µm PEI 1kg
B-100-60-CO-ARG-100G	Bulk Coral 100A 60µm ARG 100g
B-100-60-CO-ARG-1KG	Bulk Coral 100A 60µm ARG 1kg
B-100-60-CO-NH-100G	Bulk Coral 100A 60µm NH 100g
B-100-60-CO-NH-1KG	Bulk Coral 100A 60µm NH 1kg
B-100-60-CO-NH2-100G	Bulk Coral 100A 60µm NH2 100g
B-100-60-CO-NH2-1KG	Bulk Coral 100A 60µm NH2 1kg
B-100-60-CO-CO-CN-100G	Bulk Coral 100A 60µm CN 100g
B-100-60-CO-CO-CN-1KG	Bulk Coral 100A 60µm CN 1kg
B-100-60-CO-Phenyl-100G	Bulk Coral 100A 60µm Phenyl 100g
B-100-60-CO-Phenyl-1KG	Bulk Coral 100A 60µm Phenyl 1kg
B-100-60-CO-Diol-100G	Bulk Coral 100A 60µm Diol 100g
B-100-60-CO-Diol-1KG	Bulk Coral 100A 60µm Diol 1kg
<b>Particule Coral 60, Pore 300A</b>	
B-300-60-CO-Si-100G	Bulk Coral 300A 60µm 100g
B-300-60-CO-Si-1KG	Bulk Coral 300A 60µm 1kg
<b>C8</b>	
B-300-60-CO-C8T-100G	Bulk Coral 300A 60µm C8 TriF 100g
B-300-60-CO-C8T-1KG	Bulk Coral 300A 60µm C8 TriF 1kg
B-300-60-CO-C8EX-100G	Bulk Coral 300A 60µm C8 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300-60-CO-C8EX-1KG	Bulk Coral 300A 60µm C8 Extreme 1kg
B-300-60-CO-C8PT-100G	Bulk Coral 300A 60µm C8 Polar TriF 100g
B-300-60-CO-C8PT-1KG	Bulk Coral 300A 60µm C8 Polar TriF 1kg
B-300-60-CO-C8AQP-100G	Bulk Coral 300A 60µm C8 AQ Plus 100g
B-300-60-CO-C8AQP-1KG	Bulk Coral 300A 60µm C8 AQ Plus 1kg
B-300-60-CO-C8AQ-100G	Bulk Coral 300A 60µm C8 AQ 100g
B-300-60-CO-C8AQ-1KG	Bulk Coral 300A 60µm C8 AQ 1kg
<b>C18</b>	
B-300-60-CO-C18T-100G	Bulk Coral 300A 60µm C18 TriF 100g
B-300-60-CO-C18T-1KG	Bulk Coral 300A 60µm C18 TriF 1kg
B-300-60-CO-C18EX-100G	Bulk Coral 300A 60µm C18 Extreme 100g
B-300-60-CO-C18EX-1KG	Bulk Coral 300A 60µm C18 Extreme 1kg
B-300-60-CO-C18PT-100G	Bulk Coral 300A 60µm C18 Polar TriF 100g
B-300-60-CO-C18PT-1KG	Bulk Coral 300A 60µm C18 Polar TriF 1kg
B-300-60-CO-C18AQP-100G	Bulk Coral 300A 60µm C18 AQ Plus 100g
B-300-60-CO-C18AQP-1KG	Bulk Coral 300A 60µm C18 AQ Plus 1kg
B-300-60-CO-C18AQ-100G	Bulk Coral 300A 60µm C18 AQ 100g
B-300-60-CO-C18AQ-1KG	Bulk Coral 300A 60µm C18 AQ 1kg
<b>C4</b>	
B-300-60-CO-C4T-100G	Bulk Coral 300A 60µm C4 TriF 100g
B-300-60-CO-C4T-1KG	Bulk Coral 300A 60µm C4 TriF 1kg
B-300-60-CO-C4EX-100G	Bulk Coral 300A 60µm C4 Extreme 100g
B-300-60-CO-C4EX-1KG	Bulk Coral 300A 60µm C4 Extreme 1kg
B-300-60-CO-C4PT-100G	Bulk Coral 300A 60µm C4 Polar TriF 100g
B-300-60-CO-C4PT-1KG	Bulk Coral 300A 60µm C4 Polar TriF 1kg
B-300-60-CO-C4AQP-100G	Bulk Coral 300A 60µm C4 AQ Plus 100g
B-300-60-CO-C4AQP-1KG	Bulk Coral 300A 60µm C4 AQ Plus 1kg
B-300-60-CO-C4AQ-100G	Bulk Coral 300A 60µm C4 AQ 100g
B-300-60-CO-C4AQ-1KG	Bulk Coral 300A 60µm C4 AQ 1kg
<b>Autres</b>	
B-300-60-CO-COOH-100G	Bulk Coral 300A 60µm COOH 100g
B-300-60-CO-COOH-1KG	Bulk Coral 300A 60µm COOH 1kg
B-300-60-CO-SO3H-100G	Bulk Coral 300A 60µm SO3H 100g
B-300-60-CO-SO3H-1KG	Bulk Coral 300A 60µm SO3H 1kg
B-300-60-CO-DNH-100G	Bulk Coral 300A 60µm DNH 100g
B-300-60-CO-DNH-1KG	Bulk Coral 300A 60µm DNH 1kg
B-300-60-CO-PEI-100G	Bulk Coral 300A 60µm PEI 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300-60-CO-PEI-1KG	Bulk Coral 300A 60µm PEI 1kg
B-300-60-CO-ARG-100G	Bulk Coral 300A 60µm ARG 100g
B-300-60-CO-ARG-1KG	Bulk Coral 300A 60µm ARG 1kg
B-300-60-CO-NH-100G	Bulk Coral 300A 60µm NH 100g
B-300-60-CO-NH-1KG	Bulk Coral 300A 60µm NH 1kg
B-300-60-CO-NH2-100G	Bulk Coral 300A 60µm NH2 100g
B-300-60-CO-NH2-1KG	Bulk Coral 300A 60µm NH2 1kg
B-300-60-CO-CO-CN-100G	Bulk Coral 300A 60µm CN 100g
B-300-60-CO-CO-CN-1KG	Bulk Coral 300A 60µm CN 1kg
B-300-60-CO-Phenyl-100G	Bulk Coral 300A 60µm Phenyl 100g
B-300-60-CO-Phenyl-1KG	Bulk Coral 300A 60µm Phenyl 1kg
B-300-60-CO-Diol-100G	Bulk Coral 300A 60µm Diol 100g
B-300-60-CO-Diol-1KG	Bulk Coral 300A 60µm Diol 1kg
<b>Particule Coral 60, Pore 500A</b>	
B-500-60-CO-Si-100G	Bulk Coral 500A 60µm 100g
B-500-60-CO-Si-1KG	Bulk Coral 500A 60µm 1kg
<b>C8</b>	
B-500-60-CO-C8T-100G	Bulk Coral 500A 60µm C8 TriF 100g
B-500-60-CO-C8T-1KG	Bulk Coral 500A 60µm C8 TriF 1kg
B-500-60-CO-C8EX-100G	Bulk Coral 500A 60µm C8 Extreme 100g
B-500-60-CO-C8EX-1KG	Bulk Coral 500A 60µm C8 Extreme 1kg
B-500-60-CO-C8PT-100G	Bulk Coral 500A 60µm C8 Polar TriF 100g
B-500-60-CO-C8PT-1KG	Bulk Coral 500A 60µm C8 Polar TriF 1kg
B-500-60-CO-C8AQP-100G	Bulk Coral 500A 60µm C8 AQ Plus 100g
B-500-60-CO-C8AQP-1KG	Bulk Coral 500A 60µm C8 AQ Plus 1kg
B-500-60-CO-C8AQ-100G	Bulk Coral 500A 60µm C8 AQ 100g
B-500-60-CO-C8AQ-1KG	Bulk Coral 500A 60µm C8 AQ 1kg
<b>C18</b>	
B-500-60-CO-C18T-100G	Bulk Coral 500A 60µm C18 TriF 100g
B-500-60-CO-C18T-1KG	Bulk Coral 500A 60µm C18 TriF 1kg
B-500-60-CO-C18EX-100G	Bulk Coral 500A 60µm C18 Extreme 100g
B-500-60-CO-C18EX-1KG	Bulk Coral 500A 60µm C18 Extreme 1kg
B-500-60-CO-C18PT-100G	Bulk Coral 500A 60µm C18 Polar TriF 100g
B-500-60-CO-C18PT-1KG	Bulk Coral 500A 60µm C18 Polar TriF 1kg
B-500-60-CO-C18AQP-100G	Bulk Coral 500A 60µm C18 AQ Plus 100g
B-500-60-CO-C18AQP-1KG	Bulk Coral 500A 60µm C18 AQ Plus 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-500-60-CO-C18AQ-100G	Bulk Coral 500A 60µm C18 AQ 100g
B-500-60-CO-C18AQ-1KG	Bulk Coral 500A 60µm C18 AQ 1kg
<b>C4</b>	
B-500-60-CO-C4T-100G	Bulk Coral 500A 60µm C4 TriF 100g
B-500-60-CO-C4T-1KG	Bulk Coral 500A 60µm C4 TriF 1kg
B-500-60-CO-C4EX-100G	Bulk Coral 500A 60µm C4 Extreme 100g
B-500-60-CO-C4EX-1KG	Bulk Coral 500A 60µm C4 Extreme 1kg
B-500-60-CO-C4PT-100G	Bulk Coral 500A 60µm C4 Polar TriF 100g
B-500-60-CO-C4PT-1KG	Bulk Coral 500A 60µm C4 Polar TriF 1kg
B-500-60-CO-C4AQP-100G	Bulk Coral 500A 60µm C4 AQ Plus 100g
B-500-60-CO-C4AQP-1KG	Bulk Coral 500A 60µm C4 AQ Plus 1kg
B-500-60-CO-C4AQ-100G	Bulk Coral 500A 60µm C4 AQ 100g
B-500-60-CO-C4AQ-1KG	Bulk Coral 500A 60µm C4 AQ 1kg
<b>Autres</b>	
B-500-60-CO-COOH-100G	Bulk Coral 500A 60µm COOH 100g
B-500-60-CO-COOH-1KG	Bulk Coral 500A 60µm COOH 1kg
B-500-60-CO-SO3H-100G	Bulk Coral 500A 60µm SO3H 100g
B-500-60-CO-SO3H-1KG	Bulk Coral 500A 60µm SO3H 1kg
B-500-60-CO-DNH-100G	Bulk Coral 500A 60µm DNH 100g
B-500-60-CO-DNH-1KG	Bulk Coral 500A 60µm DNH 1kg
B-500-60-CO-PEI-100G	Bulk Coral 500A 60µm PEI 100g
B-500-60-CO-PEI-1KG	Bulk Coral 500A 60µm PEI 1kg
B-500-60-CO-ARG-100G	Bulk Coral 500A 60µm ARG 100g
B-500-60-CO-ARG-1KG	Bulk Coral 500A 60µm ARG 1kg
B-500-60-CO-NH-100G	Bulk Coral 500A 60µm NH 100g
B-500-60-CO-NH-1KG	Bulk Coral 500A 60µm NH 1kg
B-500-60-CO-NH2-100G	Bulk Coral 500A 60µm NH2 100g
B-500-60-CO-NH2-1KG	Bulk Coral 500A 60µm NH2 1kg
B-500-60-CO-CO-CN-100G	Bulk Coral 500A 60µm CN 100g
B-500-60-CO-CO-CN-1KG	Bulk Coral 500A 60µm CN 1kg
B-500-60-CO-Phenyl-100G	Bulk Coral 500A 60µm Phenyl 100g
B-500-60-CO-Phenyl-1KG	Bulk Coral 500A 60µm Phenyl 1kg
B-500-60-CO-Diol-100G	Bulk Coral 500A 60µm Diol 100g
B-500-60-CO-Diol-1KG	Bulk Coral 500A 60µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Coral 60, Pore 800A</b>	
B-800-60-CO-Si-100G	Bulk Coral 800A 60µm 100g
B-800-60-CO-Si-1KG	Bulk Coral 800A 60µm 1kg
<b>C8</b>	
B-800-60-CO-C8T-100G	Bulk Coral 800A 60µm C8 TriF 100g
B-800-60-CO-C8T-1KG	Bulk Coral 800A 60µm C8 TriF 1kg
B-800-60-CO-C8EX-100G	Bulk Coral 800A 60µm C8 Extreme 100g
B-800-60-CO-C8EX-1KG	Bulk Coral 800A 60µm C8 Extreme 1kg
B-800-60-CO-C8PT-100G	Bulk Coral 800A 60µm C8 Polar TriF 100g
B-800-60-CO-C8PT-1KG	Bulk Coral 800A 60µm C8 Polar TriF 1kg
B-800-60-CO-C8AQP-100G	Bulk Coral 800A 60µm C8 AQ Plus 100g
B-800-60-CO-C8AQP-1KG	Bulk Coral 800A 60µm C8 AQ Plus 1kg
B-800-60-CO-C8AQ-100G	Bulk Coral 800A 60µm C8 AQ 100g
B-800-60-CO-C8AQ-1KG	Bulk Coral 800A 60µm C8 AQ 1kg
<b>C18</b>	
B-800-60-CO-C18T-100G	Bulk Coral 800A 60µm C18 TriF 100g
B-800-60-CO-C18T-1KG	Bulk Coral 800A 60µm C18 TriF 1kg
B-800-60-CO-C18EX-100G	Bulk Coral 800A 60µm C18 Extreme 100g
B-800-60-CO-C18EX-1KG	Bulk Coral 800A 60µm C18 Extreme 1kg
B-800-60-CO-C18PT-100G	Bulk Coral 800A 60µm C18 Polar TriF 100g
B-800-60-CO-C18PT-1KG	Bulk Coral 800A 60µm C18 Polar TriF 1kg
B-800-60-CO-C18AQP-100G	Bulk Coral 800A 60µm C18 AQ Plus 100g
B-800-60-CO-C18AQP-1KG	Bulk Coral 800A 60µm C18 AQ Plus 1kg
B-800-60-CO-C18AQ-100G	Bulk Coral 800A 60µm C18 AQ 100g
B-800-60-CO-C18AQ-1KG	Bulk Coral 800A 60µm C18 AQ 1kg
<b>C4</b>	
B-800-60-CO-C4T-100G	Bulk Coral 800A 60µm C4 TriF 100g
B-800-60-CO-C4T-1KG	Bulk Coral 800A 60µm C4 TriF 1kg
B-800-60-CO-C4EX-100G	Bulk Coral 800A 60µm C4 Extreme 100g
B-800-60-CO-C4EX-1KG	Bulk Coral 800A 60µm C4 Extreme 1kg
B-800-60-CO-C4PT-100G	Bulk Coral 800A 60µm C4 Polar TriF 100g
B-800-60-CO-C4PT-1KG	Bulk Coral 800A 60µm C4 Polar TriF 1kg
B-800-60-CO-C4AQP-100G	Bulk Coral 800A 60µm C4 AQ Plus 100g
B-800-60-CO-C4AQP-1KG	Bulk Coral 800A 60µm C4 AQ Plus 1kg
B-800-60-CO-C4AQ-100G	Bulk Coral 800A 60µm C4 AQ 100g
B-800-60-CO-C4AQ-1KG	Bulk Coral 800A 60µm C4 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Autres</b>	
B-800-60-CO-COOH-100G	Bulk Coral 800A 60µm COOH 100g
B-800-60-CO-COOH-1KG	Bulk Coral 800A 60µm COOH 1kg
B-800-60-CO-SO3H-100G	Bulk Coral 800A 60µm SO3H 100g
B-800-60-CO-SO3H-1KG	Bulk Coral 800A 60µm SO3H 1kg
B-800-60-CO-DNH-100G	Bulk Coral 800A 60µm DNH 100g
B-800-60-CO-DNH-1KG	Bulk Coral 800A 60µm DNH 1kg
B-800-60-CO-PEI-100G	Bulk Coral 800A 60µm PEI 100g
B-800-60-CO-PEI-1KG	Bulk Coral 800A 60µm PEI 1kg
B-800-60-CO-ARG-100G	Bulk Coral 800A 60µm ARG 100g
B-800-60-CO-ARG-1KG	Bulk Coral 800A 60µm ARG 1kg
B-800-60-CO-NH-100G	Bulk Coral 800A 60µm NH 100g
B-800-60-CO-NH-1KG	Bulk Coral 800A 60µm NH 1kg
B-800-60-CO-NH2-100G	Bulk Coral 800A 60µm NH2 100g
B-800-60-CO-NH2-1KG	Bulk Coral 800A 60µm NH2 1kg
B-800-60-CO-CO-CN-100G	Bulk Coral 800A 60µm CN 100g
B-800-60-CO-CO-CN-1KG	Bulk Coral 800A 60µm CN 1kg
B-800-60-CO-Phenyl-100G	Bulk Coral 800A 60µm Phenyl 100g
B-800-60-CO-Phenyl-1KG	Bulk Coral 800A 60µm Phenyl 1kg
B-800-60-CO-Diol-100G	Bulk Coral 800A 60µm Diol 100g
B-800-60-CO-Diol-1KG	Bulk Coral 800A 60µm Diol 1kg
<b>Particule Coral 60, Pore 1000A</b>	
B-1000-60-CO-Si-100G	Bulk Coral 1000A 60µm 100g
B-1000-60-CO-Si-1KG	Bulk Coral 1000A 60µm 1kg
<b>C8</b>	
B-1000-60-CO-C8T-100G	Bulk Coral 1000A 60µm C8 TriF 100g
B-1000-60-CO-C8T-1KG	Bulk Coral 1000A 60µm C8 TriF 1kg
B-1000-60-CO-C8EX-100G	Bulk Coral 1000A 60µm C8 Extreme 100g
B-1000-60-CO-C8EX-1KG	Bulk Coral 1000A 60µm C8 Extreme 1kg
B-1000-60-CO-C8PT-100G	Bulk Coral 1000A 60µm C8 Polar TriF 100g
B-1000-60-CO-C8PT-1KG	Bulk Coral 1000A 60µm C8 Polar TriF 1kg
B-1000-60-CO-C8AQP-100G	Bulk Coral 1000A 60µm C8 AQ Plus 100g
B-1000-60-CO-C8AQP-1KG	Bulk Coral 1000A 60µm C8 AQ Plus 1kg
B-1000-60-CO-C8AQ-100G	Bulk Coral 1000A 60µm C8 AQ 100g
B-1000-60-CO-C8AQ-1KG	Bulk Coral 1000A 60µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-1000-60-CO-C18T-100G	Bulk Coral 1000A 60µm C18 TriF 100g
B-1000-60-CO-C18T-1KG	Bulk Coral 1000A 60µm C18 TriF 1kg
B-1000-60-CO-C18EX-100G	Bulk Coral 1000A 60µm C18 Extreme 100g
B-1000-60-CO-C18EX-1KG	Bulk Coral 1000A 60µm C18 Extreme 1kg
B-1000-60-CO-C18PT-100G	Bulk Coral 1000A 60µm C18 Polar TriF 100g
B-1000-60-CO-C18PT-1KG	Bulk Coral 1000A 60µm C18 Polar TriF 1kg
B-1000-60-CO-C18AQP-100G	Bulk Coral 1000A 60µm C18 AQ Plus 100g
B-1000-60-CO-C18AQP-1KG	Bulk Coral 1000A 60µm C18 AQ Plus 1kg
B-1000-60-CO-C18AQ-100G	Bulk Coral 1000A 60µm C18 AQ 100g
B-1000-60-CO-C18AQ-1KG	Bulk Coral 1000A 60µm C18 AQ 1kg
<b>C4</b>	
B-1000-60-CO-C4T-100G	Bulk Coral 1000A 60µm C4 TriF 100g
B-1000-60-CO-C4T-1KG	Bulk Coral 1000A 60µm C4 TriF 1kg
B-1000-60-CO-C4EX-100G	Bulk Coral 1000A 60µm C4 Extreme 100g
B-1000-60-CO-C4EX-1KG	Bulk Coral 1000A 60µm C4 Extreme 1kg
B-1000-60-CO-C4PT-100G	Bulk Coral 1000A 60µm C4 Polar TriF 100g
B-1000-60-CO-C4PT-1KG	Bulk Coral 1000A 60µm C4 Polar TriF 1kg
B-1000-60-CO-C4AQP-100G	Bulk Coral 1000A 60µm C4 AQ Plus 100g
B-1000-60-CO-C4AQP-1KG	Bulk Coral 1000A 60µm C4 AQ Plus 1kg
B-1000-60-CO-C4AQ-100G	Bulk Coral 1000A 60µm C4 AQ 100g
B-1000-60-CO-C4AQ-1KG	Bulk Coral 1000A 60µm C4 AQ 1kg
<b>Autres</b>	
B-1000-60-CO-COOH-100G	Bulk Coral 1000A 60µm COOH 100g
B-1000-60-CO-COOH-1KG	Bulk Coral 1000A 60µm COOH 1kg
B-1000-60-CO-SO3H-100G	Bulk Coral 1000A 60µm SO3H 100g
B-1000-60-CO-SO3H-1KG	Bulk Coral 1000A 60µm SO3H 1kg
B-1000-60-CO-DNH-100G	Bulk Coral 1000A 60µm DNH 100g
B-1000-60-CO-DNH-1KG	Bulk Coral 1000A 60µm DNH 1kg
B-1000-60-CO-PEI-100G	Bulk Coral 1000A 60µm PEI 100g
B-1000-60-CO-PEI-1KG	Bulk Coral 1000A 60µm PEI 1kg
B-1000-60-CO-ARG-100G	Bulk Coral 1000A 60µm ARG 100g
B-1000-60-CO-ARG-1KG	Bulk Coral 1000A 60µm ARG 1kg
B-1000-60-CO-NH-100G	Bulk Coral 1000A 60µm NH 100g
B-1000-60-CO-NH-1KG	Bulk Coral 1000A 60µm NH 1kg
B-1000-60-CO-NH2-100G	Bulk Coral 1000A 60µm NH2 100g
B-1000-60-CO-NH2-1KG	Bulk Coral 1000A 60µm NH2 1kg
B-1000-60-CO-CO-CN-100G	Bulk Coral 1000A 60µm CN 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-1000-60-CO-CO-CN-1KG	Bulk Coral 1000A 60µm CN 1kg
B-1000-60-CO-Phenyl-100G	Bulk Coral 1000A 60µm Phenyl 100g
B-1000-60-CO-Phenyl-1KG	Bulk Coral 1000A 60µm Phenyl 1kg
B-1000-60-CO-Diol-100G	Bulk Coral 1000A 60µm Diol 100g
B-1000-60-CO-Diol-1KG	Bulk Coral 1000A 60µm Diol 1kg
<b>Coral 140µm</b>	
<b>Particule Coral 140, Pore 55A</b>	
B-55-140-CO-Si-100G	Bulk Coral 55A 140µm 100g
B-55-140-CO-Si-1KG	Bulk Coral 55A 140µm 1kg
<b>C8</b>	
B-55-140-CO-C8T-100G	Bulk Coral 55A 140µm C8 TriF 100g
B-55-140-CO-C8T-1KG	Bulk Coral 55A 140µm C8 TriF 1kg
B-55-140-CO-C8EX-100G	Bulk Coral 55A 140µm C8 Extreme 100g
B-55-140-CO-C8EX-1KG	Bulk Coral 55A 140µm C8 Extreme 1kg
B-55-140-CO-C8PT-100G	Bulk Coral 55A 140µm C8 Polar TriF 100g
B-55-140-CO-C8PT-1KG	Bulk Coral 55A 140µm C8 Polar TriF 1kg
B-55-140-CO-C8AQP-100G	Bulk Coral 55A 140µm C8 AQ Plus 100g
B-55-140-CO-C8AQP-1KG	Bulk Coral 55A 140µm C8 AQ Plus 1kg
B-55-140-CO-C8AQ-100G	Bulk Coral 55A 140µm C8 AQ 100g
B-55-140-CO-C8AQ-1KG	Bulk Coral 55A 140µm C8 AQ 1kg
<b>C18</b>	
B-55-140-CO-C18T-100G	Bulk Coral 55A 140µm C18 TriF 100g
B-55-140-CO-C18T-1KG	Bulk Coral 55A 140µm C18 TriF 1kg
B-55-140-CO-C18EX-100G	Bulk Coral 55A 140µm C18 Extreme 100g
B-55-140-CO-C18EX-1KG	Bulk Coral 55A 140µm C18 Extreme 1kg
B-55-140-CO-C18PT-100G	Bulk Coral 55A 140µm C18 Polar TriF 100g
B-55-140-CO-C18PT-1KG	Bulk Coral 55A 140µm C18 Polar TriF 1kg
B-55-140-CO-C18AQP-100G	Bulk Coral 55A 140µm C18 AQ Plus 100g
B-55-140-CO-C18AQP-1KG	Bulk Coral 55A 140µm C18 AQ Plus 1kg
B-55-140-CO-C18AQ-100G	Bulk Coral 55A 140µm C18 AQ 100g
B-55-140-CO-C18AQ-1KG	Bulk Coral 55A 140µm C18 AQ 1kg
<b>C4</b>	
B-55-140-CO-C4T-100G	Bulk Coral 55A 140µm C4 TriF 100g
B-55-140-CO-C4T-1KG	Bulk Coral 55A 140µm C4 TriF 1kg
B-55-140-CO-C4EX-100G	Bulk Coral 55A 140µm C4 Extreme 100g
B-55-140-CO-C4EX-1KG	Bulk Coral 55A 140µm C4 Extreme 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-55-140-CO-C4PT-100G	Bulk Coral 55A 140µm C4 Polar TriF 100g
B-55-140-CO-C4PT-1KG	Bulk Coral 55A 140µm C4 Polar TriF 1kg
B-55-140-CO-C4AQP-100G	Bulk Coral 55A 140µm C4 AQ Plus 100g
B-55-140-CO-C4AQP-1KG	Bulk Coral 55A 140µm C4 AQ Plus 1kg
B-55-140-CO-C4AQ-100G	Bulk Coral 55A 140µm C4 AQ 100g
B-55-140-CO-C4AQ-1KG	Bulk Coral 55A 140µm C4 AQ 1kg
<b>Autres</b>	
B-55-140-CO-COOH-100G	Bulk Coral 55A 140µm COOH 100g
B-55-140-CO-COOH-1KG	Bulk Coral 55A 140µm COOH 1kg
B-55-140-CO-SO3H-100G	Bulk Coral 55A 140µm SO3H 100g
B-55-140-CO-SO3H-1KG	Bulk Coral 55A 140µm SO3H 1kg
B-55-140-CO-DNH-100G	Bulk Coral 55A 140µm DNH 100g
B-55-140-CO-DNH-1KG	Bulk Coral 55A 140µm DNH 1kg
B-55-140-CO-PEI-100G	Bulk Coral 55A 140µm PEI 100g
B-55-140-CO-PEI-1KG	Bulk Coral 55A 140µm PEI 1kg
B-55-140-CO-ARG-100G	Bulk Coral 55A 140µm ARG 100g
B-55-140-CO-ARG-1KG	Bulk Coral 55A 140µm ARG 1kg
B-55-140-CO-NH-100G	Bulk Coral 55A 140µm NH 100g
B-55-140-CO-NH-1KG	Bulk Coral 55A 140µm NH 1kg
B-55-140-CO-NH2-100G	Bulk Coral 55A 140µm NH2 100g
B-55-140-CO-NH2-1KG	Bulk Coral 55A 140µm NH2 1kg
B-55-140-CO-CN-100G	Bulk Coral 55A 140µm CN 100g
B-55-140-CO-CN-1KG	Bulk Coral 55A 140µm CN 1kg
B-55-140-CO-Phenyl-100G	Bulk Coral 55A 140µm Phenyl 100g
B-55-140-CO-Phenyl-1KG	Bulk Coral 55A 140µm Phenyl 1kg
B-55-140-CO-Diol-100G	Bulk Coral 55A 140µm Diol 100g
B-55-140-CO-Diol-1KG	Bulk Coral 55A 140µm Diol 1kg
<b>Particule Coral 140, Pore 55A</b>	
B-70-140-CO-Si-100G	Bulk Coral 70A 140µm 100g
B-70-140-CO-Si-1KG	Bulk Coral 70A 140µm 1kg
<b>C8</b>	
B-70-140-CO-C8T-100G	Bulk Coral 70A 140µm C8 TriF 100g
B-70-140-CO-C8T-1KG	Bulk Coral 70A 140µm C8 TriF 1kg
B-70-140-CO-C8EX-100G	Bulk Coral 70A 140µm C8 Extreme 100g
B-70-140-CO-C8EX-1KG	Bulk Coral 70A 140µm C8 Extreme 1kg



# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70-140-CO-C8PT-100G	Bulk Coral 70A 140µm C8 Polar TriF 100g
B-70-140-CO-C8PT-1KG	Bulk Coral 70A 140µm C8 Polar TriF 1kg
B-70-140-CO-C8AQP-100G	Bulk Coral 70A 140µm C8 AQ Plus 100g
B-70-140-CO-C8AQP-1KG	Bulk Coral 70A 140µm C8 AQ Plus 1kg
B-70-140-CO-C8AQ-100G	Bulk Coral 70A 140µm C8 AQ 100g
B-70-140-CO-C8AQ-1KG	Bulk Coral 70A 140µm C8 AQ 1kg
<b>C18</b>	
B-70-140-CO-C18T-100G	Bulk Coral 70A 140µm C18 TriF 100g
B-70-140-CO-C18T-1KG	Bulk Coral 70A 140µm C18 TriF 1kg
B-70-140-CO-C18EX-100G	Bulk Coral 70A 140µm C18 Extreme 100g
B-70-140-CO-C18EX-1KG	Bulk Coral 70A 140µm C18 Extreme 1kg
B-70-140-CO-C18PT-100G	Bulk Coral 70A 140µm C18 Polar TriF 100g
B-70-140-CO-C18PT-1KG	Bulk Coral 70A 140µm C18 Polar TriF 1kg
B-70-140-CO-C18AQP-100G	Bulk Coral 70A 140µm C18 AQ Plus 100g
B-70-140-CO-C18AQP-1KG	Bulk Coral 70A 140µm C18 AQ Plus 1kg
B-70-140-CO-C18AQ-100G	Bulk Coral 70A 140µm C18 AQ 100g
B-70-140-CO-C18AQ-1KG	Bulk Coral 70A 140µm C18 AQ 1kg
<b>C4</b>	
B-70-140-CO-C4T-100G	Bulk Coral 70A 140µm C4 TriF 100g
B-70-140-CO-C4T-1KG	Bulk Coral 70A 140µm C4 TriF 1kg
B-70-140-CO-C4EX-100G	Bulk Coral 70A 140µm C4 Extreme 100g
B-70-140-CO-C4EX-1KG	Bulk Coral 70A 140µm C4 Extreme 1kg
B-70-140-CO-C4PT-100G	Bulk Coral 70A 140µm C4 Polar TriF 100g
B-70-140-CO-C4PT-1KG	Bulk Coral 70A 140µm C4 Polar TriF 1kg
B-70-140-CO-C4AQP-100G	Bulk Coral 70A 140µm C4 AQ Plus 100g
B-70-140-CO-C4AQP-1KG	Bulk Coral 70A 140µm C4 AQ Plus 1kg
B-70-140-CO-C4AQ-100G	Bulk Coral 70A 140µm C4 AQ 100g
B-70-140-CO-C4AQ-1KG	Bulk Coral 70A 140µm C4 AQ 1kg
<b>Autres</b>	
B-70-140-CO-COOH-100G	Bulk Coral 70A 140µm COOH 100g
B-70-140-CO-COOH-1KG	Bulk Coral 70A 140µm COOH 1kg
B-70-140-CO-SO3H-100G	Bulk Coral 70A 140µm SO3H 100g
B-70-140-CO-SO3H-1KG	Bulk Coral 70A 140µm SO3H 1kg
B-70-140-CO-DNH-100G	Bulk Coral 70A 140µm DNH 100g
B-70-140-CO-DNH-1KG	Bulk Coral 70A 140µm DNH 1kg
B-70-140-CO-PEI-100G	Bulk Coral 70A 140µm PEI 100g
B-70-140-CO-PEI-1KG	Bulk Coral 70A 140µm PEI 1kg
B-70-140-CO-ARG-100G	Bulk Coral 70A 140µm ARG 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-70-140-CO-ARG-1KG	Bulk Coral 70A 140µm ARG 1kg
B-70-140-CO-NH-100G	Bulk Coral 70A 140µm NH 100g
B-70-140-CO-NH-1KG	Bulk Coral 70A 140µm NH 1kg
B-70-140-CO-NH2-100G	Bulk Coral 70A 140µm NH2 100g
B-70-140-CO-NH2-1KG	Bulk Coral 70A 140µm NH2 1kg
B-70-140-CO-CN-100G	Bulk Coral 70A 140µm CN 100g
B-70-140-CO-CN-1KG	Bulk Coral 70A 140µm CN 1kg
B-70-140-CO-Phenyl-100G	Bulk Coral 70A 140µm Phenyl 100g
B-70-140-CO-Phenyl-1KG	Bulk Coral 70A 140µm Phenyl 1kg
B-70-140-CO-Diol-100G	Bulk Coral 70A 140µm Diol 100g
B-70-140-CO-Diol-1KG	Bulk Coral 70A 140µm Diol 1kg
<b>Particule Coral 140, Pore 100A</b>	
B-100-140-CO-Si-100G	Bulk Coral 100A 140µm 100g
B-100-140-CO-Si-1KG	Bulk Coral 100A 140µm 1kg
<b>C8</b>	
B-100-140-CO-C8T-100G	Bulk Coral 100A 140µm C8 TriF 100g
B-100-140-CO-C8T-1KG	Bulk Coral 100A 140µm C8 TriF 1kg
B-100-140-CO-C8EX-100G	Bulk Coral 100A 140µm C8 Extreme 100g
B-100-140-CO-C8EX-1KG	Bulk Coral 100A 140µm C8 Extreme 1kg
B-100-140-CO-C8PT-100G	Bulk Coral 100A 140µm C8 Polar TriF 100g
B-100-140-CO-C8PT-1KG	Bulk Coral 100A 140µm C8 Polar TriF 1kg
B-100-140-CO-C8AQP-100G	Bulk Coral 100A 140µm C8 AQ Plus 100g
B-100-140-CO-C8AQP-1KG	Bulk Coral 100A 140µm C8 AQ Plus 1kg
B-100-140-CO-C8AQ-100G	Bulk Coral 100A 140µm C8 AQ 100g
B-100-140-CO-C8AQ-1KG	Bulk Coral 100A 140µm C8 AQ 1kg
<b>C18</b>	
B-100-140-CO-C18T-100G	Bulk Coral 100A 140µm C18 TriF 100g
B-100-140-CO-C18T-1KG	Bulk Coral 100A 140µm C18 TriF 1kg
B-100-140-CO-C18EX-100G	Bulk Coral 100A 140µm C18 Extreme 100g
B-100-140-CO-C18EX-1KG	Bulk Coral 100A 140µm C18 Extreme 1kg
B-100-140-CO-C18PT-100G	Bulk Coral 100A 140µm C18 Polar TriF 100g
B-100-140-CO-C18PT-1KG	Bulk Coral 100A 140µm C18 Polar TriF 1kg
B-100-140-CO-C18AQP-100G	Bulk Coral 100A 140µm C18 AQ Plus 100g
B-100-140-CO-C18AQP-1KG	Bulk Coral 100A 140µm C18 AQ Plus 1kg
B-100-140-CO-C18AQ-100G	Bulk Coral 100A 140µm C18 AQ 100g
B-100-140-CO-C18AQ-1KG	Bulk Coral 100A 140µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-100-140-CO-C4T-100G	Bulk Coral 100A 140µm C4 TriF 100g
B-100-140-CO-C4T-1KG	Bulk Coral 100A 140µm C4 TriF 1kg
B-100-140-CO-C4EX-100G	Bulk Coral 100A 140µm C4 Extreme 100g
B-100-140-CO-C4EX-1KG	Bulk Coral 100A 140µm C4 Extreme 1kg
B-100-140-CO-C4PT-100G	Bulk Coral 100A 140µm C4 Polar TriF 100g
B-100-140-CO-C4PT-1KG	Bulk Coral 100A 140µm C4 Polar TriF 1kg
B-100-140-CO-C4AQP-100G	Bulk Coral 100A 140µm C4 AQ Plus 100g
B-100-140-CO-C4AQP-1KG	Bulk Coral 100A 140µm C4 AQ Plus 1kg
B-100-140-CO-C4AQ-100G	Bulk Coral 100A 140µm C4 AQ 100g
B-100-140-CO-C4AQ-1KG	Bulk Coral 100A 140µm C4 AQ 1kg
<b>Autres</b>	
B-100-140-CO-COOH-100G	Bulk Coral 100A 140µm COOH 100g
B-100-140-CO-COOH-1KG	Bulk Coral 100A 140µm COOH 1kg
B-100-140-CO-SO3H-100G	Bulk Coral 100A 140µm SO3H 100g
B-100-140-CO-SO3H-1KG	Bulk Coral 100A 140µm SO3H 1kg
B-100-140-CO-DNH-100G	Bulk Coral 100A 140µm DNH 100g
B-100-140-CO-DNH-1KG	Bulk Coral 100A 140µm DNH 1kg
B-100-140-CO-PEI-100G	Bulk Coral 100A 140µm PEI 100g
B-100-140-CO-PEI-1KG	Bulk Coral 100A 140µm PEI 1kg
B-100-140-CO-ARG-100G	Bulk Coral 100A 140µm ARG 100g
B-100-140-CO-ARG-1KG	Bulk Coral 100A 140µm ARG 1kg
B-100-140-CO-NH-100G	Bulk Coral 100A 140µm NH 100g
B-100-140-CO-NH-1KG	Bulk Coral 100A 140µm NH 1kg
B-100-140-CO-NH2-100G	Bulk Coral 100A 140µm NH2 100g
B-100-140-CO-NH2-1KG	Bulk Coral 100A 140µm NH2 1kg
B-100-140-CO-CN-100G	Bulk Coral 100A 140µm CN 100g
B-100-140-CO-CN-1KG	Bulk Coral 100A 140µm CN 1kg
B-100-140-CO-Phenyl-100G	Bulk Coral 100A 140µm Phenyl 100g
B-100-140-CO-Phenyl-1KG	Bulk Coral 100A 140µm Phenyl 1kg
B-100-140-CO-Diol-100G	Bulk Coral 100A 140µm Diol 100g
B-100-140-CO-Diol-1KG	Bulk Coral 100A 140µm Diol 1kg
<b>Particule Coral 140, Pore 300A</b>	
B-300-140-CO-Si-100G	Bulk Coral 300A 140µm 100g
B-300-140-CO-Si-1KG	Bulk Coral 300A 140µm 1kg
<b>C8</b>	
B-300-140-CO-C8T-100G	Bulk Coral 300A 140µm C8 TriF 100g
B-300-140-CO-C8T-1KG	Bulk Coral 300A 140µm C8 TriF 1kg
B-300-140-CO-C8EX-100G	Bulk Coral 300A 140µm C8 Extreme 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300-140-CO-C8EX-1KG	Bulk Coral 300A 140µm C8 Extreme 1kg
B-300-140-CO-C8PT-100G	Bulk Coral 300A 140µm C8 Polar TriF 100g
B-300-140-CO-C8PT-1KG	Bulk Coral 300A 140µm C8 Polar TriF 1kg
B-300-140-CO-C8AQP-100G	Bulk Coral 300A 140µm C8 AQ Plus 100g
B-300-140-CO-C8AQP-1KG	Bulk Coral 300A 140µm C8 AQ Plus 1kg
B-300-140-CO-C8AQ-100G	Bulk Coral 300A 140µm C8 AQ 100g
B-300-140-CO-C8AQ-1KG	Bulk Coral 300A 140µm C8 AQ 1kg
<b>C18</b>	
B-300-140-CO-C18T-100G	Bulk Coral 300A 140µm C18 TriF 100g
B-300-140-CO-C18T-1KG	Bulk Coral 300A 140µm C18 TriF 1kg
B-300-140-CO-C18EX-100G	Bulk Coral 300A 140µm C18 Extreme 100g
B-300-140-CO-C18EX-1KG	Bulk Coral 300A 140µm C18 Extreme 1kg
B-300-140-CO-C18PT-100G	Bulk Coral 300A 140µm C18 Polar TriF 100g
B-300-140-CO-C18PT-1KG	Bulk Coral 300A 140µm C18 Polar TriF 1kg
B-300-140-CO-C18AQP-100G	Bulk Coral 300A 140µm C18 AQ Plus 100g
B-300-140-CO-C18AQP-1KG	Bulk Coral 300A 140µm C18 AQ Plus 1kg
B-300-140-CO-C18AQ-100G	Bulk Coral 300A 140µm C18 AQ 100g
B-300-140-CO-C18AQ-1KG	Bulk Coral 300A 140µm C18 AQ 1kg
<b>C4</b>	
B-300-140-CO-C4T-100G	Bulk Coral 300A 140µm C4 TriF 100g
B-300-140-CO-C4T-1KG	Bulk Coral 300A 140µm C4 TriF 1kg
B-300-140-CO-C4EX-100G	Bulk Coral 300A 140µm C4 Extreme 100g
B-300-140-CO-C4EX-1KG	Bulk Coral 300A 140µm C4 Extreme 1kg
B-300-140-CO-C4PT-100G	Bulk Coral 300A 140µm C4 Polar TriF 100g
B-300-140-CO-C4PT-1KG	Bulk Coral 300A 140µm C4 Polar TriF 1kg
B-300-140-CO-C4AQP-100G	Bulk Coral 300A 140µm C4 AQ Plus 100g
B-300-140-CO-C4AQP-1KG	Bulk Coral 300A 140µm C4 AQ Plus 1kg
B-300-140-CO-C4AQ-100G	Bulk Coral 300A 140µm C4 AQ 100g
B-300-140-CO-C4AQ-1KG	Bulk Coral 300A 140µm C4 AQ 1kg
<b>Autres</b>	
B-300-140-CO-COOH-100G	Bulk Coral 300A 140µm COOH 100g
B-300-140-CO-COOH-1KG	Bulk Coral 300A 140µm COOH 1kg
B-300-140-CO-SO3H-100G	Bulk Coral 300A 140µm SO3H 100g
B-300-140-CO-SO3H-1KG	Bulk Coral 300A 140µm SO3H 1kg
B-300-140-CO-DNH-100G	Bulk Coral 300A 140µm DNH 100g
B-300-140-CO-DNH-1KG	Bulk Coral 300A 140µm DNH 1kg
B-300-140-CO-PEI-100G	Bulk Coral 300A 140µm PEI 100g
B-300-140-CO-PEI-1KG	Bulk Coral 300A 140µm PEI 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-300-140-CO-ARG-100G	Bulk Coral 300A 140µm ARG 100g
B-300-140-CO-ARG-1KG	Bulk Coral 300A 140µm ARG 1kg
B-300-140-CO-NH-100G	Bulk Coral 300A 140µm NH 100g
B-300-140-CO-NH-1KG	Bulk Coral 300A 140µm NH 1kg
B-300-140-CO-NH2-100G	Bulk Coral 300A 140µm NH2 100g
B-300-140-CO-NH2-1KG	Bulk Coral 300A 140µm NH2 1kg
B-300-140-CO-CN-100G	Bulk Coral 300A 140µm CN 100g
B-300-140-CO-CN-1KG	Bulk Coral 300A 140µm CN 1kg
B-300-140-CO-Phenyl-100G	Bulk Coral 300A 140µm Phenyl 100g
B-300-140-CO-Phenyl-1KG	Bulk Coral 300A 140µm Phenyl 1kg
B-300-140-CO-Diol-100G	Bulk Coral 300A 140µm Diol 100g
B-300-140-CO-Diol-1KG	Bulk Coral 300A 140µm Diol 1kg
<b>Particule Coral 140, Pore 500A</b>	
B-500-140-CO-Si-100G	Bulk Coral 500A 140µm 100g
B-500-140-CO-Si-1KG	Bulk Coral 500A 140µm 1kg
<b>C8</b>	
B-500-140-CO-C8T-100G	Bulk Coral 500A 140µm C8 TriF 100g
B-500-140-CO-C8T-1KG	Bulk Coral 500A 140µm C8 TriF 1kg
B-500-140-CO-C8EX-100G	Bulk Coral 500A 140µm C8 Extreme 100g
B-500-140-CO-C8EX-1KG	Bulk Coral 500A 140µm C8 Extreme 1kg
B-500-140-CO-C8PT-100G	Bulk Coral 500A 140µm C8 Polar TriF 100g
B-500-140-CO-C8PT-1KG	Bulk Coral 500A 140µm C8 Polar TriF 1kg
B-500-140-CO-C8AQP-100G	Bulk Coral 500A 140µm C8 AQ Plus 100g
B-500-140-CO-C8AQP-1KG	Bulk Coral 500A 140µm C8 AQ Plus 1kg
B-500-140-CO-C8AQ-100G	Bulk Coral 500A 140µm C8 AQ 100g
B-500-140-CO-C8AQ-1KG	Bulk Coral 500A 140µm C8 AQ 1kg
<b>C18</b>	
B-500-140-CO-C18T-100G	Bulk Coral 500A 140µm C18 TriF 100g
B-500-140-CO-C18T-1KG	Bulk Coral 500A 140µm C18 TriF 1kg
B-500-140-CO-C18EX-100G	Bulk Coral 500A 140µm C18 Extreme 100g
B-500-140-CO-C18EX-1KG	Bulk Coral 500A 140µm C18 Extreme 1kg
B-500-140-CO-C18PT-100G	Bulk Coral 500A 140µm C18 Polar TriF 100g
B-500-140-CO-C18PT-1KG	Bulk Coral 500A 140µm C18 Polar TriF 1kg
B-500-140-CO-C18AQP-100G	Bulk Coral 500A 140µm C18 AQ Plus 100g
B-500-140-CO-C18AQP-1KG	Bulk Coral 500A 140µm C18 AQ Plus 1kg
B-500-140-CO-C18AQ-100G	Bulk Coral 500A 140µm C18 AQ 100g
B-500-140-CO-C18AQ-1KG	Bulk Coral 500A 140µm C18 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C4</b>	
B-500-140-CO-C4T-100G	Bulk Coral 500A 140µm C4 TriF 100g
B-500-140-CO-C4T-1KG	Bulk Coral 500A 140µm C4 TriF 1kg
B-500-140-CO-C4EX-100G	Bulk Coral 500A 140µm C4 Extreme 100g
B-500-140-CO-C4EX-1KG	Bulk Coral 500A 140µm C4 Extreme 1kg
B-500-140-CO-C4PT-100G	Bulk Coral 500A 140µm C4 Polar TriF 100g
B-500-140-CO-C4PT-1KG	Bulk Coral 500A 140µm C4 Polar TriF 1kg
B-500-140-CO-C4AQP-100G	Bulk Coral 500A 140µm C4 AQ Plus 100g
B-500-140-CO-C4AQP-1KG	Bulk Coral 500A 140µm C4 AQ Plus 1kg
B-500-140-CO-C4AQ-100G	Bulk Coral 500A 140µm C4 AQ 100g
B-500-140-CO-C4AQ-1KG	Bulk Coral 500A 140µm C4 AQ 1kg
<b>Autres</b>	
B-500-140-CO-COOH-100G	Bulk Coral 500A 140µm COOH 100g
B-500-140-CO-COOH-1KG	Bulk Coral 500A 140µm COOH 1kg
B-500-140-CO-SO3H-100G	Bulk Coral 500A 140µm SO3H 100g
B-500-140-CO-SO3H-1KG	Bulk Coral 500A 140µm SO3H 1kg
B-500-140-CO-DNH-100G	Bulk Coral 500A 140µm DNH 100g
B-500-140-CO-DNH-1KG	Bulk Coral 500A 140µm DNH 1kg
B-500-140-CO-PEI-100G	Bulk Coral 500A 140µm PEI 100g
B-500-140-CO-PEI-1KG	Bulk Coral 500A 140µm PEI 1kg
B-500-140-CO-ARG-100G	Bulk Coral 500A 140µm ARG 100g
B-500-140-CO-ARG-1KG	Bulk Coral 500A 140µm ARG 1kg
B-500-140-CO-NH-100G	Bulk Coral 500A 140µm NH 100g
B-500-140-CO-NH-1KG	Bulk Coral 500A 140µm NH 1kg
B-500-140-CO-NH2-100G	Bulk Coral 500A 140µm NH2 100g
B-500-140-CO-NH2-1KG	Bulk Coral 500A 140µm NH2 1kg
B-500-140-CO-CN-100G	Bulk Coral 500A 140µm CN 100g
B-500-140-CO-CN-1KG	Bulk Coral 500A 140µm CN 1kg
B-500-140-CO-Phenyl-100G	Bulk Coral 500A 140µm Phenyl 100g
B-500-140-CO-Phenyl-1KG	Bulk Coral 500A 140µm Phenyl 1kg
B-500-140-CO-Diol-100G	Bulk Coral 500A 140µm Diol 100g
B-500-140-CO-Diol-1KG	Bulk Coral 500A 140µm Diol 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>Particule Coral 140, Pore 800A</b>	
B-800-140-CO-Si-100G	Bulk Coral 800A 140µm 100g
B-800-140-CO-Si-1KG	Bulk Coral 800A 140µm 1kg
<b>C8</b>	
B-800-140-CO-C8T-100G	Bulk Coral 800A 140µm C8 TriF 100g
B-800-140-CO-C8T-1KG	Bulk Coral 800A 140µm C8 TriF 1kg
B-800-140-CO-C8EX-100G	Bulk Coral 800A 140µm C8 Extreme 100g
B-800-140-CO-C8EX-1KG	Bulk Coral 800A 140µm C8 Extreme 1kg
B-800-140-CO-C8PT-100G	Bulk Coral 800A 140µm C8 Polar TriF 100g
B-800-140-CO-C8PT-1KG	Bulk Coral 800A 140µm C8 Polar TriF 1kg
B-800-140-CO-C8AQP-100G	Bulk Coral 800A 140µm C8 AQ Plus 100g
B-800-140-CO-C8AQP-1KG	Bulk Coral 800A 140µm C8 AQ Plus 1kg
B-800-140-CO-C8AQ-100G	Bulk Coral 800A 140µm C8 AQ 100g
B-800-140-CO-C8AQ-1KG	Bulk Coral 800A 140µm C8 AQ 1kg
<b>C18</b>	
B-800-140-CO-C18T-100G	Bulk Coral 800A 140µm C18 TriF 100g
B-800-140-CO-C18T-1KG	Bulk Coral 800A 140µm C18 TriF 1kg
B-800-140-CO-C18EX-100G	Bulk Coral 800A 140µm C18 Extreme 100g
B-800-140-CO-C18EX-1KG	Bulk Coral 800A 140µm C18 Extreme 1kg
B-800-140-CO-C18PT-100G	Bulk Coral 800A 140µm C18 Polar TriF 100g
B-800-140-CO-C18PT-1KG	Bulk Coral 800A 140µm C18 Polar TriF 1kg
B-800-140-CO-C18AQP-100G	Bulk Coral 800A 140µm C18 AQ Plus 100g
B-800-140-CO-C18AQP-1KG	Bulk Coral 800A 140µm C18 AQ Plus 1kg
B-800-140-CO-C18AQ-100G	Bulk Coral 800A 140µm C18 AQ 100g
B-800-140-CO-C18AQ-1KG	Bulk Coral 800A 140µm C18 AQ 1kg
<b>C4</b>	
B-800-140-CO-C4T-100G	Bulk Coral 800A 140µm C4 TriF 100g
B-800-140-CO-C4T-1KG	Bulk Coral 800A 140µm C4 TriF 1kg
B-800-140-CO-C4EX-100G	Bulk Coral 800A 140µm C4 Extreme 100g
B-800-140-CO-C4EX-1KG	Bulk Coral 800A 140µm C4 Extreme 1kg
B-800-140-CO-C4PT-100G	Bulk Coral 800A 140µm C4 Polar TriF 100g
B-800-140-CO-C4PT-1KG	Bulk Coral 800A 140µm C4 Polar TriF 1kg
B-800-140-CO-C4AQP-100G	Bulk Coral 800A 140µm C4 AQ Plus 100g
B-800-140-CO-C4AQP-1KG	Bulk Coral 800A 140µm C4 AQ Plus 1kg
B-800-140-CO-C4AQ-100G	Bulk Coral 800A 140µm C4 AQ 100g
B-800-140-CO-C4AQ-1KG	Bulk Coral 800A 140µm C4 AQ 1kg
<b>Autres</b>	

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-800-140-CO-COOH-100G	Bulk Coral 800A 140µm COOH 100g
B-800-140-CO-COOH-1KG	Bulk Coral 800A 140µm COOH 1kg
B-800-140-CO-SO3H-100G	Bulk Coral 800A 140µm SO3H 100g
B-800-140-CO-SO3H-1KG	Bulk Coral 800A 140µm SO3H 1kg
B-800-140-CO-DNH-100G	Bulk Coral 800A 140µm DNH 100g
B-800-140-CO-DNH-1KG	Bulk Coral 800A 140µm DNH 1kg
B-800-140-CO-PEI-100G	Bulk Coral 800A 140µm PEI 100g
B-800-140-CO-PEI-1KG	Bulk Coral 800A 140µm PEI 1kg
B-800-140-CO-ARG-100G	Bulk Coral 800A 140µm ARG 100g
B-800-140-CO-ARG-1KG	Bulk Coral 800A 140µm ARG 1kg
B-800-140-CO-NH-100G	Bulk Coral 800A 140µm NH 100g
B-800-140-CO-NH-1KG	Bulk Coral 800A 140µm NH 1kg
B-800-140-CO-NH2-100G	Bulk Coral 800A 140µm NH2 100g
B-800-140-CO-NH2-1KG	Bulk Coral 800A 140µm NH2 1kg
B-800-140-CO-CN-100G	Bulk Coral 800A 140µm CN 100g
B-800-140-CO-CN-1KG	Bulk Coral 800A 140µm CN 1kg
B-800-140-CO-Phenyl-100G	Bulk Coral 800A 140µm Phenyl 100g
B-800-140-CO-Phenyl-1KG	Bulk Coral 800A 140µm Phenyl 1kg
B-800-140-CO-Diol-100G	Bulk Coral 800A 140µm Diol 100g
B-800-140-CO-Diol-1KG	Bulk Coral 800A 140µm Diol 1kg
<b>Particule Coral 140, Pore 1000A</b>	
B-1000-140-CO-Si-100G	Bulk Coral 1000A 140µm 100g
B-1000-140-CO-Si-1KG	Bulk Coral 1000A 140µm 1kg
<b>C8</b>	
B-1000-140-CO-C8T-100G	Bulk Coral 1000A 140µm C8 TriF 100g
B-1000-140-CO-C8T-1KG	Bulk Coral 1000A 140µm C8 TriF 1kg
B-1000-140-CO-C8EX-100G	Bulk Coral 1000A 140µm C8 Extreme 100g
B-1000-140-CO-C8EX-1KG	Bulk Coral 1000A 140µm C8 Extreme 1kg
B-1000-140-CO-C8PT-100G	Bulk Coral 1000A 140µm C8 Polar TriF 100g
B-1000-140-CO-C8PT-1KG	Bulk Coral 1000A 140µm C8 Polar TriF 1kg
B-1000-140-CO-C8AQP-100G	Bulk Coral 1000A 140µm C8 AQ Plus 100g
B-1000-140-CO-C8AQP-1KG	Bulk Coral 1000A 140µm C8 AQ Plus 1kg
B-1000-140-CO-C8AQ-100G	Bulk Coral 1000A 140µm C8 AQ 100g
B-1000-140-CO-C8AQ-1KG	Bulk Coral 1000A 140µm C8 AQ 1kg

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
<b>C18</b>	
B-1000-140-CO-C18T-100G	Bulk Coral 1000A 140µm C18 TriF 100g
B-1000-140-CO-C18T-1KG	Bulk Coral 1000A 140µm C18 TriF 1kg
B-1000-140-CO-C18EX-100G	Bulk Coral 1000A 140µm C18 Extreme 100g
B-1000-140-CO-C18EX-1KG	Bulk Coral 1000A 140µm C18 Extreme 1kg
B-1000-140-CO-C18PT-100G	Bulk Coral 1000A 140µm C18 Polar TriF 100g
B-1000-140-CO-C18PT-1KG	Bulk Coral 1000A 140µm C18 Polar TriF 1kg
B-1000-140-CO-C18AQP-100G	Bulk Coral 1000A 140µm C18 AQ Plus 100g
B-1000-140-CO-C18AQP-1KG	Bulk Coral 1000A 140µm C18 AQ Plus 1kg
B-1000-140-CO-C18AQ-100G	Bulk Coral 1000A 140µm C18 AQ 100g
B-1000-140-CO-C18AQ-1KG	Bulk Coral 1000A 140µm C18 AQ 1kg
<b>C4</b>	
B-1000-140-CO-C4T-100G	Bulk Coral 1000A 140µm C4 TriF 100g
B-1000-140-CO-C4T-1KG	Bulk Coral 1000A 140µm C4 TriF 1kg
B-1000-140-CO-C4EX-100G	Bulk Coral 1000A 140µm C4 Extreme 100g
B-1000-140-CO-C4EX-1KG	Bulk Coral 1000A 140µm C4 Extreme 1kg
B-1000-140-CO-C4PT-100G	Bulk Coral 1000A 140µm C4 Polar TriF 100g
B-1000-140-CO-C4PT-1KG	Bulk Coral 1000A 140µm C4 Polar TriF 1kg
B-1000-140-CO-C4AQP-100G	Bulk Coral 1000A 140µm C4 AQ Plus 100g
B-1000-140-CO-C4AQP-1KG	Bulk Coral 1000A 140µm C4 AQ Plus 1kg
B-1000-140-CO-C4AQ-100G	Bulk Coral 1000A 140µm C4 AQ 100g
B-1000-140-CO-C4AQ-1KG	Bulk Coral 1000A 140µm C4 AQ 1kg
<b>Autres</b>	
B-1000-140-CO-COOH-100G	Bulk Coral 1000A 140µm COOH 100g
B-1000-140-CO-COOH-1KG	Bulk Coral 1000A 140µm COOH 1kg
B-1000-140-CO-SO3H-100G	Bulk Coral 1000A 140µm SO3H 100g
B-1000-140-CO-SO3H-1KG	Bulk Coral 1000A 140µm SO3H 1kg
B-1000-140-CO-DNH-100G	Bulk Coral 1000A 140µm DNH 100g
B-1000-140-CO-DNH-1KG	Bulk Coral 1000A 140µm DNH 1kg
B-1000-140-CO-PEI-100G	Bulk Coral 1000A 140µm PEI 100g
B-1000-140-CO-PEI-1KG	Bulk Coral 1000A 140µm PEI 1kg
B-1000-140-CO-ARG-100G	Bulk Coral 1000A 140µm ARG 100g
B-1000-140-CO-ARG-1KG	Bulk Coral 1000A 140µm ARG 1kg
B-1000-140-CO-NH-100G	Bulk Coral 1000A 140µm NH 100g
B-1000-140-CO-NH-1KG	Bulk Coral 1000A 140µm NH 1kg
B-1000-140-CO-NH2-100G	Bulk Coral 1000A 140µm NH2 100g
B-1000-140-CO-NH2-1KG	Bulk Coral 1000A 140µm NH2 1kg
B-1000-140-CO-CN-100G	Bulk Coral 1000A 140µm CN 100g

# Phases Bulk

## POUR COMMANDER

Référence	Désignation
B-1000-140-CO-CN-1KG	Bulk Coral 1000A 140µm CN 1kg
B-1000-140-CO-Phenyl-100G	Bulk Coral 1000A 140µm Phenyl 100g
B-1000-140-CO-Phenyl-1KG	Bulk Coral 1000A 140µm Phenyl 1kg
B-1000-140-CO-Diol-100G	Bulk Coral 1000A 140µm Diol 100g
B-1000-140-CO-Diol-1KG	Bulk Coral 1000A 140µm Diol 1kg

# Silice en vrac pour HPLC

Notes

# Silice en vrac pour HPLC

Notes





*im*Chem

164 Avenue Joseph Kessel - Bâtiment 7 - 78960 Voisins le Bretonneux - France  
Tel. : +33 9 72 32 10 17 06 - Fax : +33 1 83 64 60 45 - info@imchem.fr - www.imchem.fr