

## Application Data Alphabetical List Page1

Compound	Compound No.	$\alpha$	e.o.	CSP	Mobile Phase
Abscisic acid	H14	1.06		OA-4700	1040
Acebutolol	B36	4.25		CBH	22101
N-Acetylalanine	G22	1.12		OA-3200	2
N-Acetylalanine	G22	1.12		OA-3300	2
N-Acetylalanine	G22	1.17		OA-4700	1041
N-Acetylalanine	G22	1.17		OA-5000	131
N-Acetylasnartic acid	G49	1.11		OA-4700	1041
N-Acetylleucine	G24	1.04		OA-3100	1
N-Acetylleucine	G24	1.10		OA-3200	2
N-Acetylleucine	G24	1.39		OA-4700	1041
N-Acetylleucine	G24	1.39		OA-5000	130
O-Acetylmandelic acid	D12	1.06	S,R	OA-4700	1040
N-Acetylmethionine	G25	1.06		OA-3200	2
N-Acetylmethionirie	G25	1.05		OA-3300	2
N-Acetylmethionine	G25	1.36		OA-4700	1041
N-Acetylmethionine	G25	1.27		OA-5000	130
N-Acetylphenylalanine	A34	1.08		OA-4000	1041
N-Acetylphenylalanine	A34	1.18		OA-4100	1041
N-Acetylphenylalanine	A34	1.08		OA-4500	1041
N-Acetylphenylalanine	A34	1.09	D,L	OA-4600	3
N-Acetylphenylalanine	A34	1.15		OA-4600	1041
N-Acetylphenylalanine	A34	1.24	D,L	OA-4700	3
N-Acetylphenylalanine	A34	1.38	D,L	OA-4700	1041
N-Acetylserine	G50	1.18		OA-4700	1041
N-Acetyltryptophan	A36	1.25	D,L	OA-4700	1041
N-Acetylvaline	G23	1.11		OA-3300	2
N-Acetylvaline	G23	1.11		OA-4000	1041
N-Acetylvaline	G23	1.10		OA-4100	1041
N-Acetylvaline	G23	1.11	D,L	OA-4600	3
N-Acetylvaline	G23	1.18		OA-4600	1041
N-Acetylvaline	G23	1.18	D,L	OA-4700	3
N-Acetylvaline	G23	1.24	D,L	OA-4700	1041
N-Acetylvaline	G23	1.50		OA-5000	130
Alanine	G01	1.66	L,D	OA-5000	21
Alanine	G01	2.28	D,L	OA-6000	21
Alanine	G01	2.23	D,L	OA-6100	21
alanine- $\beta$ -naphthylamide	C71	1.17		OA-8000	60001
Alanylglycine	G36	2.26		OA-5000	21
Alanylglycyl-glycine	G38	1.62		OA-5000	21
Albuterol	B35	1.18	+,-	OA-4900	1054
Allethrin	I03	I03		OA-2000I	56
Allethrolone	I06	1.09	S,R	OA-4100	58
Allethrolone	I06	1.08		OA-4500	58
Allethrolone	I06	1.08		OA-4700	58
Alloisoleucine	G15	1.67	L,D	OA-5000	138
Alloisoleucine	G15	1.20	L,D	OA-5500	21
Alloisoleucine	G15	2.79	D,L	OA-6000	130

Allothreonine	G18	1.19	L,D	OA-5000	11
Allothreonine	G18	1.72	D,L	OA-5500	21
Allothreonine	G18	3.93	D,L	OA-6000	21
Allothreonine	G18	6.34	D,L	OA-6100	21
erythro-2-Amino-1,2-diphenylethanol	B27	1.14	+,-	OA-4500	1047
erythro-2-Amino-1,2-diphenylethanol	B27	1.20	+,-	OA-4900	1047
erythro-2-Amino-1,2-diphenviethanol	B27	1.13		OA-5500	130
erythro-2-Amino-1,2-diphenviethanol	B27	1.10		OA-7000	43060
2-Amino-1-phenylethanol	B15	1.19		OA-5500	21
3-Amino-2-methylpropionic acid	G34	1.08		OA-5500	21
2-Amino-4-phosphonobutylic acid	G45	1.31		OA-5000	131
2-Amino-4-phosphonobutylic acid	G45	1.86		OA-5500	130

PAGE : Ab ~ Am Am ~ Be Be ~ Ci Cl ~ Di Di ~ Fl Fl ~ Ib Is ~ Mo Na ~ Ph Ph ~ Pr Pr ~ Te Te ~ Un Va ~

CLOSE