



Vials & Closures



Syringe Filters




SPE, QuEChERS



GC Columns



HPLC Columns



Flash Chromatography



Standards



Reagents



Spare parts
& Lamps

Centrifuge tubes

Designed for QuEChERS workflow

Features

- **Medical-Grade PP Material** : High-purity polypropylene with excellent chemical resistance.
- **Resistant to e-beam & gamma sterilization**
- **DNase/RNase and pyrogen-free.**
- **Leak-Proof Cap Design**, prevents leakage & sample loss during centrifugation or storage.
- **Withstands temperatures from - 80°C to 120°C.**
- **Centrifugal speed up to 10,000 rpm.**
- **Clear & Accurate** : Easy-to-read graduations & writable frosted area for quick labeling.
- **Cleanroom Manufacturing** : - Produced with strict quality control.

Validated for Common Analytical Tests

- Extractable Test
- Leakage Test
- Pesticide & Veterinary Drug Residue Test

Application

- Ideal for food safety testing, environmental analysis, pharmaceutical R&D, & general laboratory sample preparation.



IC015CT-1



IC015SC



ICBAK050



IC050-1

Ordering information

imChem item code	Designation	Qty
IC015CT-1	Cryopreservation tubes, sterile 2mL	500pcs/bag
IC015SC	Centrifuge Tubes, 15 mL, with separate cap, non-sterile	500 Pcs/Carton
ICBAK050	Centrifuge Tubes, 50 mL, with assembled cap, non-sterile	25 Pcs/Bag - 20 bags/Carton
IC050-1	Centrifuge Tubes, 50 mL, with assembled cap, Self-standing, non-sterile	50 Pcs/Bag - 10 bags/Carton



Your expert source for chromatography consumables

164 Av. Joseph Kessel – 78960 Voisins le Bretonneux – France
Tel +33 9 72 32 10 17 - info@imchem.fr

www.imchem.fr